GE Healthcare

MR Quiet Upgrade

OPERATING DOCUMENTATION

5670077
Revision 3.0

WARNING

STRONG MAGNETIC FIELD

NO PACEMAKERS*
NO NEUROSTIMULATORS*
NO CONDUCTIVE/METALLIC IMPLANTS*

**Persons with pacemakers,
neurostimulators or metallic
implants must not enter this area.**
Serious injury may result.

* In general, patients with conductive (e.g. metallic) implants are contraindicated for MR scans. For patients with implants that are labeled as 'MR Safe' or 'MR Conditional', consult the implant device manufacturer's documentation.

* WARNING: Only use quadrature transmit for 'MR Conditional' devices.

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**WARNING
(EN)**

This service manual is available in English only.

- If a customer's service provider requires a language other than English, it is the customer's responsibility to provide translation services.
- Do not attempt to service the equipment unless this service manual has been consulted and is understood.
- Failure to heed this warning may result in injury to the service provider, operator or patient from electric shock, mechanical or other hazards.

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This page left intentionally blank.

Revision History

Revision	Date	Description
1.0	May 1, 2013	Initial release.
2.0	July 23, 2013	Updated revision to assign part classification – no content change
3.0	July 23, 2013	Updated revision to synchronize revision history with MyWorkshop – no content change

This page left intentionally blank.

Table of Contents

Chapter 1 OVERVIEW	15
1 Getting Started	15
1.1 Upgrade Introduction	15
1.2 Site Ready Check	15
1.3 Pre-Delivery FDO Checks	15
1.4 Upgrade Schedule Considerations	15
1.5 Product Delivery Instructions and Packing	15
1.6 Damage in Transportation	16
1.7 Product Locator	16
1.8 Required Equipment and Tools	17
1.9 Discarded Material Return Policy	17
1.9.1 Material Policy Overview	18
1.9.1.1 Returned Material	18
1.9.1.2 Disposed Material	18
1.9.2 Material to be Returned	18
1.9.3 Return Kit	18
1.9.4 Deinstallation Procedures	18
1.9.5 Region Specific Return Information	18
1.9.5.1 Sites in Americas	18
1.9.5.2 Sites in Asia	19
1.9.5.3 Sites in Europe	19
1.9.6 Product Locator Information	19
1.10 Upgrade Components	19
1.10.1 Components Supplied for MR450w GEM Systems (M7001SB)	20
1.10.2 Components Supplied for MR750w GEM Systems (M7001SD)	21
Chapter 2 UPGRADE PROCEDURE	23
1 Upgrade Procedures	23
1.1 Pre-Mechanical Upgrade Procedures	23
1.2 Mechanical Upgrade Procedures	23

1.2.1 Power Down CAM Chassis.....	23
1.2.2 Board Replacement Identification.....	25
1.2.2.1 MR450w Board Positions.....	25
1.2.2.2 MR750w Board Positions.....	26
1.2.3 Circuit Board Removal.....	26
1.2.4 Circuit Board Replacement.....	27
1.2.5 Return Removed Modules.....	28
1.3 Post-Mechanical Upgrade Procedures.....	28
1.4 Completion.....	30
Chapter 3 APPENDIX.....	31
1 Upgrade Test Record.....	31

Chapter 1 Overview

1 Getting Started

1.1 Upgrade Introduction

This publication documents the mechanical installation information for upgrading a MR450w GEM or MR750w GEM to MR Quiet.

WARNING

FERROUS MATERIAL HAZARD!

DELIVERY EQUIPMENT AND OTHER TOOLS AND PARTS REQUIRED FOR THIS INSTALLATION CONTAIN FERROUS MATERIAL AND WILL BE STRONGLY ATTRACTED TO MAGNET AND MAY BECOME DANGEROUS PROJECTILES.

KEEP ALL FERROUS TOOLS AT LEAST 10 FEET AWAY FROM THE MAGNET.

1.2 Site Ready Check

Pre-installation work must be completed before equipment is delivered to avoid delays and confusion.

1.3 Pre-Delivery FDO Checks

Before the shipment is made, Be sure that FDO review has been completed by consultation with the Upgrade Specialist. Also verify the Surface coil compatibility issues have also been resolved.

1.4 Upgrade Schedule Considerations

FMI, periodic maintenance, installation of prerequisite upgrades and/or options will directly add to the man-hour and system down-time total. Be sure that all planned activity is considered when predicting the return of system to customer. Also verify Application Specialist is scheduled at the end of the installation to instruct users. Additional time must also be scheduled to instruct users following completion of upgrade.

1.5 Product Delivery Instructions and Packing

The "Shipping Document" lists catalog numbers delivered. **Review and confirm that order is delivered complete.** Check impact on installation schedule if Catalogs and/or packing boxes are missing and/or noted as shipped short.

Labels that are attached to the outside of packing boxes summarize box contents. PDI (Product Delivery Instructions) specify box contents, part numbers, and shipping procedure. The PDI is numbered according to the catalog number. Lists of items included with each box are detailed by separate checklists, or a separate sheet that provides a summary of that box contents. Refer to PDI and packing lists for information specific to your shipment.

A set of service and operator manuals is delivered with the Upgrade Kit. Refer to checklists packed with "Technical Publication" boxes for a list of delivered documents.

1.6 Damage in Transportation

All packages should be closely examined at time of delivery. If damage is apparent, have notation **"damage in shipment"** written on **all** copies of the freight or express bill **before** delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier **immediately** upon discovery, or in any event, within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this **14** day period. File a report with:

- Call 1-800-548-3366. Select "Install Support Services for FOA and MIS"
- Contact your local service coordinator for more information on this process.

1.7 Product Locator

The Global Install Base Database (also known as Product Locator System) tracks shipment, trans-shipment, and field location of the serialized models. There are now two methods for submitting Product Locator information.

1. At this time, for **U.S. ONLY**, the preferred method of submitting information is the FE Site Verification Web Site at:

<http://gein2.med.ge.com/gib>

The FE Site Verification consists of three components that are available on the web from the main menu. They are:

- Install/deinstall product locator model and serial numbers
- Add/modify ship to address information
- Update CARES FE data for primary/secondary FEs

To obtain a copy of the FE training tutorial for using this website, a downloadable copy is available at the following: Product Locator Support Central Page - <http://supportcentral.ge.com/15563>

2. One "Shipping Card" is filled out and submitted when shipped (extra cards are supplied for trans-shipments between storage and distribution points), and the "Installation Card" and extra shipment cards are attached.

Verify that serial and model number on each rating plate matches installation card numbers before removing installation card. Note that there may be one or more shipment cards and bar code labels with the installation card. These shipment cards are used to trace the transfer of serialized units between various inventory storage and distribution points until the product reaches its final installation destination. Process just the installation card and discard any extra shipment cards and labels.

1.8 Required Equipment and Tools

NOTICE

Follow all required safety and PPE procedures customary for your organization when working on this product.

Table 1-1: General Safety Requirements

Item	Description
1	Wear composite safety shoes for foot protection.
2	Wear gloves for hand protection.
3	Wear safety glasses for eye protection.

DANGER

FERROUS MATERIAL HAZARD!

IF MAGNET IS ENERGIZED, TOOLS, COMPONENTS, AND PARTS REQUIRED FOR THIS INSTALLATION THAT CONTAIN FERROUS MATERIAL WILL BE STRONGLY ATTRACTED TO MAGNET AND MAY BECOME DANGEROUS PROJECTILES.

KEEP ALL FERROUS TOOLS AS FAR FROM THE MAGNET AS POSSIBLE.

Tools to upgrade the system are listed in following table:

Table 1-2: Installation Equipment

Item	GE Part Number	Description
1		Screwdriver, slotted (long & narrow)
2	46-271138G1	Restricted Access Control Kit. Contains two plastic warning signs for posting at site during installation and service activity.
3		Service Methods DVD
4	46-198094P1	Wrist grounding strap
5	---	Plastic or aluminum flashlight
6	5307511-3 or (5434817-20 & 5434817 set)	RF Power Measurement Kit
7	46-194427P284	Digital Multi-meter (DMM)

1.9 Discarded Material Return Policy

The equipment or material being removed by the upgrade have a potential value for service of the installed base. The disposition of such must be controlled by GE. Disposition of the removed material is per GE Medical Systems Policy and Procedures.

1.9.1 Material Policy Overview

1.9.1.1 Returned Material

Per GE Healthcare Policy 3.9: all requested material is to be returned, as directed, intact and promptly to Renewable Resources. Do not remove or swap any parts from the material to be returned.

1.9.1.2 Disposed Material

Regulations governing the disposal of this equipment or material vary from location to location and from country to country. In order to insure compliance with all local, state, or national environment, health, and safety regulations, dispose of all items per local GEMS policy and procedures and other applicable published guidance. Contact GEMS Recycling Center for assistance.

1.9.2 Material to be Returned

All deinstalled UPM processor boards must be returned to GE Renewable Resources for recycling. If a 3.0T system, also return the GEM HNU coil.

- Mark each box: MR QUIET UPGRADE
- Send the boxes to the following address:

GE Renewable Resources
Attn: Asset Recovery Leader, Harvest
7624 South 10th Street
Oak Creek, WI 53154

1.9.3 Return Kit

If needed, packaging and/or transportation details can be obtained through the local Installation Specialist.

Reuse appropriate containers that new hardware and modules were shipped to site for return of above material.

1.9.4 Deinstallation Procedures

Refer to instructions in this manual for deinstallation procedures of site.

1.9.5 Region Specific Return Information

1.9.5.1 Sites in Americas

All deinstalled UPM processor boards must be returned to GE Renewable Resources for recycling. If a 3.0T system, also return the GEM HNU coil.

- Mark each box: MR QUIET UPGRADE
- Send the boxes to the following address:

GE Renewable Resources
Attn: Asset Recovery Leader, Harvest
7624 South 10th Street
Oak Creek, WI 53154

1.9.5.2 Sites in Asia

Return procedure for GEMS–Asia countries, contact:

GE HEALTHCARE JAPAN CORPORATION
ATTN: Yoshiaki Tsujii (SWaP/HARVEST)
67-4, Takakura-machi, Hachioji-shi-Tokyo
192-0033 Japan
Contact at Destination: 81-426-48-2595
FAX: 81-426-48-2902

1.9.5.3 Sites in Europe

For non-Gold Seal returns in Europe:



European Renewable Resources Center
GEODIS LOGISTICS - CD 26 - Pièce de la Remise - Bâtiment EVL2
Quai 45 à 49 - 91090 LISSES - FRANCE
TEL = (+33)-1-64-97-55-41 FAX = (+33)-1-60-86-12-31

1.9.6 Product Locator Information

Process product locator information for all serialized components removed from the system according to policy. Refer to [Section 1.7](#).

1.10 Upgrade Components

Dependent on system being upgraded, the following components are provided for upgrading a MR450w GEM or MR750w GEM sites.

Illustration 1-1: Location of Replacement Boards in PGR Cabinet**Illustration 1-2: 3.0T Quiet GEM HNU Coil**

1.10.1 Components Supplied for MR450w GEM Systems (M7001SB)

- Quantity 2 of RF Detector boards
- Quantity 2 of UPM Processor boards
- DV24 software collector
- Silent MR option key
- MR Quiet Upgrade manual

1.10.2 Components Supplied for MR750w GEMSystems (M7001SD)

- Quantity 4 of RF Detector boards
- Quantity 2 of UPM Processor boards
- DV24 software collector
- Silent MR option key
- 3.0T Quiet GEM HNU coil
- MR Quiet Upgrade manual

This page left intentionally blank.

Chapter 2 Upgrade Procedure

1 Upgrade Procedures

The general sequence of this upgrade procedure is:

[Section 1.1](#) Pre-Mechanical Upgrade Procedures

[Section 1.2](#) Mechanical Upgrade Procedures

[Section 1.3](#) Post-Mechanical Upgrade Procedures

[Section 1.4](#) Completion

1.1 Pre-Mechanical Upgrade Procedures

Refer to the applicable MR Service Methods DVD to complete the following procedures.

1. Verify that all required upgrades have arrived complete and site is ready to shut down.
2. Do check scan to verify proper system operation – see service methods Check Scan (Good-Bye)
3. To confirm the UPM Processor boards are ready to harvest, run the following two UPM diagnostics. If the diagnostic completes, the UPM Processor boards are good to return. Fill in the attached Upgrade Test Procedure page located in the **Chapter 3, Appendix** of this manual. Make sure the test procedure is dated, signed, and all the * entries are completed.
 - Analog Data
 - I/O Data

NOTE: The diagnostic location is:

- Diagnostics >> System Function >> RF >> UPM Diagnostics

OR

- Diagnostics >> Hardware Location >> PGR Cabinet >> CAM >> UPM Diagnostics

4. Perform Save Info process – see service methods Save and Restore Tool


1.2 Mechanical Upgrade Procedures

1.2.1 Power Down CAM Chassis

Prior to powering down the entire CAM Chassis, power down the ICNs from the VRE power control. The VRE System, which resides in the PGR Cabinet, can become corrupted if it loses power when the applications software is active.


1. Use DMM to check unused outlet on terminal strip for presence of 115 VAC.

Illustration 2-1: Service Power Strip

2. At the host computer, open the service browser and select **Utilities > Tool Box > VRE > Power Control**
3. Next, select [Power Off VRE]
4. At the front panel of the PDU, switch the RECON circuit breaker to the OFF (down) position.

Illustration 2-2: RECON Circuit Breaker on PDU

5. At the front of the PDU, switch the PEN Cabinet or PDM-RF breaker to the OFF (down) position.
6. At the front panel of the PDU covering the circuit breakers, apply a personal LOTO device for each FE working on the equipment.

Illustration 2-3: PDU LOTO Device

7. Confirm that the power is OFF. Check the service power strip for 115 VAC power with the DMM. If the power is still ON, confirm that the RECON CB on the PDU is in the OFF position.


1.2.2 Board Replacement Identification

Refer to one of the following two sections, dependent on system being upgraded (MR450w or MR750w), for the applicable board configuration,

1.2.2.1 MR450w Board Positions

For slots 20U and 20L, two new RF Detector boards (6250264-2) are to be installed.

For slots 21U and 21L, two new UPM Processor boards (6250266) are to be installed.


Illustration 2-4: MR450w GEM Boards

1.2.2.2 MR750w Board Positions

For slots 19U, 19L, 20U and 20L, four new RF Detector boards (6250264) are to be installed.

For slots 21U and 21L, two new UPM Processor boards (6250266) are to be installed.


Illustration 2-5: MR750w GEM Boards

1.2.3 Circuit Board Removal

1. Attach an ESD grounding wrist strap.
2. Where applicable, remove the cables (Sub D, SMB, etc.) from the face of the board.
NOTE: When removing the cables, verify that the cable labels are accurate.
3. Loosen the faceplate holding screw. (These screws are captive to the faceplate, and do not need to be completely removed.)
4. Using your thumb, push in the gray release tab and gently pull the extractor tab located on the bottom of the circuit board. The lever action of the ejector disengages the circuit board connections from the chassis.

Illustration 2-6: Board Extractor Tabs

5. Being careful not to bend or bind the card in any way, slide the circuit board all the way out of the card cage.

NOTE: Handle the circuit boards using only the faceplate, ejectors, and the edges of the board.

6. Place the UPM Processor boards, with the completed test procedure, in an appropriate shipping package to return for recycling. The Detector boards should be disposed of per local regulations. Refer to the *Disposed Material* section in [Chapter 1, Overview](#).

1.2.4 Circuit Board Replacement

1. Remove the packing material from the applicable replacement board. Position board in assigned slot and gently slide the circuit board along the slot guides into the chassis.
2. Gently engage the circuit board connectors with the backplane.
 - a. If the circuit board is not properly engaging the backplane, remove the board and check the connectors for bent or broken pins.
 - b. Once the board ejectors contact the CAM chassis faceplate, gently push the ejector into the center of the board. This causes the ejector to insert the board into the chassis. (Do **not** attempt to use anything other than the ejector to insert or remove the board from the chassis.)

3. Tighten, but do not overtighten, the faceplate holding screws to secure the circuit board to the card rack.
4. Reattach the cables to the face of the board.

1.2.5 Return Removed Modules

All discarded material (i.e. removed/discarded hardware, packing material) must be returned in accordance with GE recycling policy. Refer to [Chapter 1, Discarded Material Return Policy](#)

1.3 Post-Mechanical Upgrade Procedures

1. Remove LOTO. Restore power.
2. Load DV24 software — see service methods Loading Host System Software.
 - When in Guided Install the second time (when restoring all SaveInfo files), go to the Hardware Config menu and select “QUPM 1 Mhz Sampling Rate” for UPM type


Illustration 2-7: Guided Install Hardware Configuration Menu

- When asked to load all applicable options, load the Silent MR option key. (This is not the Silent MRA option key)

- If upgrading from a pre-DV24 software version, select the “Install Scantool” button from the Option Key page of the “Guided Install” to get the full list of scan tools loaded onto the system.

Illustration 2-8: Guided Install Option Key page

3. Perform UPM Calibrations and UPM Functional Checks. Refer to Service Methods for procedure details. (Do not use the service key)
4. If upgrading from a pre-DV24 software version, insert service key and perform the HOEC Calibration. Refer to Service Methods for procedure details, otherwise skip this step.
5. If MR750w 3T GEM system, replace old GEM HNU coil with new 3.0T Quiet GEM HNU coil. Harvest the old HNU coil by using the packaging that the new Quiet GEM HNU coil came in for return shipping.
6. Run MCQA for any new coils.
7. For 1.5T/3T: Run a Quiet scan.
 - Select GE protocols tab
 - Filter list by: Quiet
 - Select adult_head_GEM Brain — Silent MR protocol
8. Perform SaveInfo — see service methods Save and Restore Tool

1.4 Completion

The following steps need to be completed by the GE Field Engineer.

1. Record and enter applicable data into applicable site configuration files and records.
2. Complete Product Locator information for all installed serialized components, new or updated, via either:

- (U.S. Only) FE Site Verification Web Site

or

- Process and return product locator installation cards for all serialized components to:

Product Locator File, P.O. Box 414, W-523, Milwaukee, WI 53201-0414

Refer to [Chapter 1, Product Locator](#) for details on submitting Product Locator information.

NOTE: Failure to fill out and return Product Locator Cards may result in failure of your site to receive future FMI's.

3. Make sure all replaced components are packaged and returned for recycling. Include the Test Procedure results.
4. Make sure all Service and Operator documentation is installed on the host computer.

Chapter 3 Appendix

1 Upgrade Test Record

Table 3-1: MR Quiet Upgrade Master Harvest Test Procedure

MR Quiet Upgrade					
System ID *		Dispatch Number *			
Cabinet Serial Number *		FE Name *			
		SSO *			
		Signature *			
		Test Date *			
Test Site Name *		Country *			
After test completion, send this document with harvested parts for review of test results. Completed test record should be signed and dated. Ensure all cells with * requiring information are completed or marked with "N/A" if not applicable.					
Manual Documentation Number		Revision Number	FE Confirm Rev.	If revision number is different, what rev. was used?	
5670014: Discovery MR450, MR750, Optima MR450w BASE & GEM, MR750w GEM Service Methods		2	*		
Test Number	Test & Procedure Documentation	Expected Results	Measured Results	Test Result Pass / Fail	Specific Failures
Test					
UPM Diagnostic Procedure Discovery MR450 & MR750, Optima MR450w BASE & GEM, MR750w GEM - <i>Troubleshooting->Diagnostics->System Function->RF->UPM Diagnostics</i>					
1	UPM Diagnostics I/O Data	Test completes	Test completes	*Circle One Pass / Fail	
2	UPM Diagnostics Analog Data	Test completes	Test completes	*Circle One Pass / Fail	
Patient Data Removal Verification					
3	Confirm Patient Data has been removed from the hard disc drives	N/A - System does not contain components which hold patient data		N/A	

This page left intentionally blank.

© 2013 General Electric Company.
General Electric Company, doing business as GE Healthcare.

www.gehealthcare.com