

Mobile unit # _____

MRI Mobile Unit Checklist

Name: _____

FE Name: _____

Date: _____

Checklist

EXTERIOR CONDITION AND STORAGE AREA

	OK	Cleaned	Adjusted	Needs Correction
Front				
Paint and body clean, no dents and in good condition.....	☐	☐	☐	☐
Welding and riveting on frame of exterior all in place no holes, all covers in place for AC and generator.....	☐	☐	☐	☐
Air break and electrical connectors in good condition, all running lights and covers on.....	☐	☐	☐	☐

Paint and body clean, no dents and in good condition.....	☐	☐	☐	☐
Welding and riveting on frame of exterior all in place no holes, all covers in place for AC and generator.....	☐	☐	☐	☐
Rear door operational and clean with all hardware and lock.....	☐	☐	☐	☐
Corner and gutter extrusions in place.....	☐	☐	☐	☐
Lower bay storage - check for corrosion, leakage and door locks working, all doors clean and look good.....	☐	☐	☐	☐
Lower bay storage - everything cleaned out, only necessary items and a few pieces of wood for leveling.....	☐	☐	☐	☐
Awning over lift gate operational and clean, awning handle and emergency wrench in storage place.....	☐	☐	☐	☐
Patient lift and roll door working properly with both remote and panel controls.....	☐	☐	☐	☐
Lift gate railings secure with all hardware and working properly, everything clean with caution markings.....	☐	☐	☐	☐
Aluminum stair assembly stable with rails and all bolts or hardware.....	☐	☐	☐	☐
Main entry door and Emergency exit door hardware working (exit signs and path marked).....	☐	☐	☐	☐
Exterior light and rain deflector drip rail at entry door in place and working.....	☐	☐	☐	☐

INTERIOR CONDITION

Floor in tech area clean and in good condition, computers or cabinets mounted with nothing else on floors..	☐	☐	☐	☐
Counter tops clean and in good condition, everything on them mounted or secure for moving.....	☐	☐	☐	☐
Condition of storage cabinets, everything cleaned with only necessary items.....	☐	☐	☐	☐
Current save Info with MOD options and manuals, folder with all data sheets and log book in cabinet.....	☐	☐	☐	☐
Floor in magnet room clean and in good condition, nothing on floors.....	☐	☐	☐	☐
Magnet covers, table, coils, pads and phantoms all cleaned and stored in a secure location.....	☐	☐	☐	☐
Floor in equipment room clean and in good condition, nothing on floors.....	☐	☐	☐	☐
All equipment has screws and secure, everything clean with all cabinet cover on and doors locked.....	☐	☐	☐	☐
Ceiling material is acoustic tile, carpet, Other_____. clean and in good condition.....	☐	☐	☐	☐
Conditions of wall coverings, holes patched and walls clean in good condition.....	☐	☐	☐	☐
All electrical outlets, light switches and e-stop cover plates in place and clean.....	☐	☐	☐	☐
Fire detection system, current non ferrous fire extinguisher installed in at least 2 locations.....	☐	☐	☐	☐
Optional equipment mounted and secure, Injector system, stereo system, camera, ect.....	☐	☐	☐	☐

Checklist

MOBILE SYSTEM CHECK

	OK	Cleaned	Adjusted	Needs Correction
Condition of RussellStoll power connector, cable secure and lugs tight.....	☐	☐	☐	☐
Check all power and ground connections, all covers and breaker panels labeled and secure.....	☐	☐	☐	☐
Check diesel fuel level and generator operation, verify condition.....	☐	☐	☐	☐
Check shore/generator manual transfer switch.....	☐	☐	☐	☐
Check 12V DC battery and charger/inverter with 12V circuit breaker panel and master switch.....	☐	☐	☐	☐
Check all lighting inside, including all high intensity LED lights, make sure all lights and switches work.....	☐	☐	☐	☐
Check working condition of AM/FM stereo CD player in Control Room(if applicable).....	☐	☐	☐	☐
Working condition of HVAC units, coils cleaned and system PM complete everything in spec.....	☐	☐	☐	☐
Check working condition of A/C and heater in equipment, operator and magnet rooms.....	☐	☐	☐	☐
Twin I-beam, landing legs, and support stands on chassis (check for cracks or damage).....	☐	☐	☐	☐
Verify tire tread, wear, and pressure check that lug nuts are all there.....	☐	☐	☐	☐
Anti-lock brakes(if available).....	☐	☐	☐	☐
Condition of roof, roof hatch and magnet exit, all seals and RF mesh in place.....	☐	☐	☐	☐

FINAL BEFORE MOVING

Do PM on Unit and Chillers (if it has not been done recently).....	☐	☐	☐	☐
Secure Keyboard, monitor and mouse.....	☐	☐	☐	☐
Bungee magnet covers.....	☐	☐	☐	☐
Lower table, lock wheels and make sure it is secure.....	☐	☐	☐	☐
Check license plate, tags, registration and HUD are on vehicle.....	☐	☐	☐	☐
Check to make sure shim lead is up.....	☐	☐	☐	☐
Check Helium Level and Pressure.....	☐	☐	☐	☐

Compressor hrs. _____

Helium _____

Pres. _____

Notes