Technical Publication

Direction 2296444-100

Revision 03

GE Medical Systems LightSpeed 3.X Functional Interconnects

The information in this service manual applies to the following
LightSpeed 3.X CT systems:

- LightSpeed Ultra (8-slice, MDAS)
- LightSpeed Plus (4-slice, MDAS)
- LightSpeed QX/i (4-slice, MDAS)

ADVANCED SERVICE DOCUMENTATION
PROPERTY OF GE
FOR GE SERVICE PERSONNEL ONLY
NO RIGHTS LICENSED
DO NOT USE OR COPY
DISCLOSURE TO THIRD PARTIES PROHIBITED

Copyright © 2001-2005 by General Electric Company.
All rights reserved.

PROPRIETARY TO GE MEDICAL SYSTEMS

Permission to use this Advanced Service Software and related documentation (herein called the material) by persons other than GE Medical Systems employees is provided only under an Advanced Service Package License relating specifically to this Proprietary Material. This is a different agreement from the one under which operating and basic service software is lessened. A license to use operating or basic service software does not extend to or cover this software or related documentation.

If you are a GE Medical Systems employee or a customer who has entered into such a license agreement with GE Medical Systems to use this proprietary software, you are authorized to use this Material according to the conditions stated in your license agreement. However, you do not have the permission of GE Medical Systems to alter, decompose or reverse-assemble the software, and unless you are a GE employee, you may not copy the Material. The Material is protected by Copyright and Trade Secret laws; the violation of which can result in civil damages and criminal prosecution.

If you are not party to such a license agreement or a GE Medical Systems Employee, you must exit this Material now.

IMPORTANT PRECAUTIONS

LANGUAGE

WARNING

- THIS SERVICE MANUAL IS AVAILABLE IN ENGLISH ONLY.
- IF A CUSTOMER’S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER’S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR OR PATIENT FROM ELECTRIC SHOCK, MECHANICAL OR OTHER HAZARDS.

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DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation “damage in shipment” written on all copies of the freight or express bill before delivery is accepted or “signed for” by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report:

- Call 1-800-548-3366 and use option 8.
- Fill out a report on <http://us44hdd21/sctq/InstallFulfill/InstalFulfillment.htm>
- Contact your local service coordinator for more information on this process.

Rev. Jan. 5, 2005

CERTIFIED ELECTRICAL CONTRACTOR STATEMENT

All electrical Installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE Medical personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE’s electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE’s field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

IMPORTANT...X-RAY PROTECTION

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The General Electric Company, Medical Systems Group, will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the General Electric Company, Medical Systems Group, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective materials and devices are available. It is urged that such materials or devices be used.

LITHIUM BATTERY CAUTIONARY STATEMENTS

CAUTION
Risk of
Explosion

Danger of explosion if battery is incorrectly replaced. Replace only with the same or equivalent type recommended by the manufacturer. Discard used batteries according to the manufacturer’s instructions.

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OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the GEMS PQR Process to report all omissions, errors, and defects in this publication.

Revision History

Rev	Date	Reason for change
0	11/13/03	Removed information from the LightSpeed 3.X System Service Manual - Advanced, and created a stand-alone manual.
1	12/18/03	CQA 1033884: Added Section 1.0 — Patient Alignment Lights - 2126249, to Chapter 4
2	08/17/04	PQR 13012916: Added Section 5.0 — Table Cradle Drive (Vertical & Longitudinal Drives) - 2123540 to Chapter 4
3	05/10/05	PQR 13024644: Updated drawing 226770SCH, Added drawing 2267713SCH

Table of Contents

Chapter 1		Chapter 5	
System Communications.....	7	Axial Control	39
Section 1.0 — Gantry Slip Ring Communications - 2267732	8	Section 1.0 — Axial Servo Control Loop - 2267757.....	40
Section 2.0 — Rotating Controller Interface Bus (RCIB) - 2305060.....	9	Section 2.0 — Axial Drive Power Control - 2267784.....	41
Section 3.0 — Final Scan Control / Intercom / AutoVoice - 2151953.....	10		
Chapter 2		Chapter 6	
X-Ray Generation	13	Operator I/O.....	43
Section 1.0 — Rotor Control - 2267317	14	Section 1.0 — Gantry Operator Interface - 2267760	44
Section 2.0 — Exposure Enable - 2267735	15	Section 2.0 — Site X-Ray ON Light - 2267761	45
Section 3.0 — Rotating HVON Sense - 2267736	16	Section 3.0 — Site System ON Light	47
Section 4.0 — Filament Select and mA Loop - 2267737	17	Section 4.0 — RX Tilt & Gantry Console Reset - 2267765	48
Section 5.0 — Tube ID, Cooling and Pressure Sense - 2267738.....	18		
Section 6.0 — Back-Up Timer - 2267739	19	Chapter 7	
Section 7.0 — kV Loop - 2267740	20	System Monitoring	49
Section 8.0 — BU Contactor / Exposure Interlocks - 2267741	22	Section 1.0 — Mains Undervoltage - 2267731 (Not applicable to NGPDU)	50
Section 9.0 — X-Ray Collimation and Filtration - 2267744	25	Section 2.0 — Gantry Temp. Sensors & Power Supply Monitors - 2267766....	51
Section 10.0 — Final Exposure Command - 2267746	26		
Chapter 3		Chapter 8	
Data Acquisition	27	System Power Control	53
Section 1.0 — DAS Power Supplies - 2302185	28	Section 1.0 — HVDC Bus Control - 2267769	54
Section 2.0 — Detector Heater Control Board - 2302184.....	29	Section 2.0 — Gantry 120VAC Control - 2267770.....	57
		Section 3.0 — E-Stop Drives Reset - 2267773	59
Chapter 4		Section 4.0 — Gantry Stationary Power and Ground - 2267713	61
Patient Positioning	31		
Section 1.0 — Patient Alignment Lights - 2126249	32	Chapter 9	
Section 2.0 — Table Gantry Sync - 2267745	33	Image Generation (Global Console - Linux).....	65
Section 3.0 — Table Elevation Drive - 2267754	34	Section 1.0 — Scan Data Save	66
Section 4.0 — Gantry Tilt Function - 2267756	35	Section 2.0 — Data Restore.....	67
Section 5.0 — Table Cradle Drive (Vertical & Longitudinal Drives) - 2123540 .	36	Section 3.0 — Pre-Processing, Filtered Back Projection & Post Processing..	68
		Section 4.0 — Image Transfer	69
		Chapter 10	
		Image Generation (Global Console - Octane2).....	71
		Section 1.0 — Scan Data Save	72

Section 2.0 — Data Restore	73
Section 3.0 — Pre-Processing, Filtered Back Projection & Post Processing..	74
Section 4.0 — Image Transfer	75
Chapter 11	
Image Generation (H3 Console).....	77
Section 1.0 — Scan Data Save	78
Section 2.0 — Data Restore	79
Section 3.0 — Pre-Processing, Filtered Back Projection & Post Processing..	80
Section 4.0 — Image Transfer	81

Chapter 1

System Communications

Section 1.0 — Gantry Slip Ring Communications - 2267732


Figure 1-1 2267732 Sheet 1, Rev. 0

Section 3.0 — Final Scan Control / Intercom / AutoVoice - 2151953


Figure 1-3 2151953 Sheet 1, Rev 1


Figure 1-4 2151953 Sheet 2, Rev 1

Section 1.0 — Rotor Control - 2267317




Figure 2-1 2267317 Sheet 1, Rev. 0

Section 2.0 — Exposure Enable - 2267735

Section 3.0 — Rotating HVON Sense - 2267736

2267736SCH

BLOCK PATHNAME \$GEMS_CT_FUNCT_INTCNT/XRGEN/HVON_SENSE/H2_2267736SCH/H2_HVON_SENSE					BLOCK SHT 1 OF 1				
REV	PCN	DESCRIPTION	APPR BY	DATE	MADE	RYAN TOONEN	DATE	15/MAR/2000	DESIGN TITLE ROTATING HVON SENSE
MADE BY		DATE			CHECK	R. VANDEN BOOGARD	DATE	10/AUG/2000	FIRST MADE FOR H2 SYS FUNC
REV	PCN	DESCRIPTION	APPR BY	DATE	APPR	J. R. SCHMIDT	DATE	10/AUG/2000	DRAWING NO. 2267736SCH
MADE BY		DATE			OTHER		DATE		REV 0
					GE MEDICAL SYSTEMS MILWAUKEE, WI USA				
					SIZE B				
					SHT/CONT. ON 1 / -				

Figure 2-3 2267736 Sheet 1, Rev. 0

Section 4.0 — Filament Select and mA Loop - 2267737


Figure 2-4 2267737 Sheet 1, Rev. 0

Section 5.0 — Tube ID, Cooling and Pressure Sense - 2267738




Figure 2-5 2267738 Sheet 1, Rev. 0

Section 6.0 — Back-Up Timer - 2267739

Section 7.0 — kV Loop - 2267740

NOTE
FIBER OPTIC CABLES SHOWN
ON SHEET 2 ARE IN
THE SAME ORDER AS THEY
APPEAR ON THIS PAGE.

REFER TO
FINAL EXP COMMAND
AND EXP ENABLE
FUNCTIONAL
DRAWINGS
FOR DETAILS

Figure 2-7 2267740 Sheet 1, Rev. 0




Figure 2-8 2267740 Sheet 2, Rev. 0

Section 8.0 —BU Contactor / Exposure Interlocks - 2267741



Figure 2-11 BU Contactor / Exposure Interlocks (NGPDU shown; see Figure 2-10 for CPDU)

Section 10.0 — Final Exposure Command - 2267746

2267746SCH

BLOCK PATHNAME \$GEMS_CT_FUNCT_INTCONT/XRGEN/EXP_CMD/H2_2267746SCH/H2_EXP_CMD										BLOCK SHT 1 OF 1			
REV	PCN	DESCRIPTION			MADE	R. VANDEN BOOGARD	DATE	28/JUL/2000	GE MEDICAL SYSTEMS MILWAUKEE, WI USA	DESIGN TITLE FINAL EXPOSURE COMMAND			
MADE BY		DATE			CHECK		DATE			FIRST MADE FOR H2 SYS FUNC			
REV	PCN	DESCRIPTION			APPR	J. R. SCHMIDT	DATE	10/AUG/2000	DRAWING NO. 2267746SCH	REV 0	SIZE B	SHT/CONT. ON 1 / -	
MADE BY		DATE			OTHER		DATE						

Figure 2-13 2267746 Sheet 1, Rev. 0

Chapter 3

Data Acquisition

Section 1.0 — DAS Power Supplies - 2302185


Figure 3-1 2302185 Sheet 1, Rev. 0

Section 2.0 — Detector Heater Control Board - 2302184



Figure 3-2 2302184 Sheet 1, Rev. 0

Chapter 4

Patient Positioning

Figure 4-1 2126249 Sheet 1, Rev 0

Section 2.0 — Table Gantry Sync - 2267745


Figure 4-2 2267745 Sheet 1, Rev. 0

Section 3.0 — Table Elevation Drive - 2267754



Figure 4-3 2267754 Sheet 1, Rev. 0


Section 4.0 — Gantry Tilt Function - 2267756

BLOCK PATHNAME \$GEMS_CT_FUNCNT/PATIENT_POS/TILT_DRIVE/H2_2267756SCH/H2_TILT_DRIVE				BLOCK SHT 1 OF 2			
REV	PCN	DESCRIPTION	DATE	MADE	MICHAEL PEETE	DATE	12/SEP/2000
MADE BY				CHECK	MICHAEL PEETE	DATE	10/OCT/2000
REV	PCN	DESCRIPTION	DATE	APPR	PHIL PEARSON	DATE	10/OCT/2000
MADE BY				OTHER		DATE	
				GE MEDICAL SYSTEMS MILWAUKEE, WI USA		DESIGN TITLE H-2 GANTRY TILT FUNCTION FIRST MADE FOR H2 SYS FUNC	
				DRAWING NO. 2267756SCH		REV 0	SIZE B
						SHT/CONT. ON 1 / 2	

Figure 4-4 2267756 Sheet 1, Rev. 0

Section 5.0—Table Cradle Drive (Vertical & Longitudinal Drives) - 2123540


Figure 4-6 2123540 Sheet 2, Rev 2

Chapter 5

Axial Control

Section 1.0 — Axial Servo Control Loop - 2267757


Figure 5-1 2267757 Sheet 1 Rev. 0

Chapter 6

Operator I/O

COMPONENTS:

- DISPLAY (FRONT)
- PB SWITCH PANEL (LEFT)
- PB SWITCH PANEL (RIGHT)
- FOOT SWITCH ASSEMBLY (LEFT)
- FOOT SWITCH ASSEMBLY (RIGHT)
- GANTRY TOUCH SENSORS
- CHASSIS GND STUD
- POWER REGULATION CRT
- GANTRY OPERATOR INTERFACE

TABLE 1: COMPONENTS

REV	PCN	DESCRIPTION	DATE
1	1	1	1

TABLE 2: REVISIONS

REV	PCN	DESCRIPTION	DATE
1	1	1	1

Section 1.0 —Gantry Operator Interface - 2267760

Section 2.0 — Site X-Ray ON Light - 2267761



Figure 6-2 2267761 Sheet 1, Rev. 0 (CPDU shown; see Figure 6-3 for NGPDU)

NOTE: * CONNECTIONS TO TS6 ARE
SUPPLIED BY THE CUSTOMER

Figure 6-3 Site X-Ray ON Light (NGPDU shown; see Figure 6-2 for CPDU)

Section 3.0 — Site System ON Light


Figure 6-4 Site System ON Light (NGPDU shown)

Section 4.0 — RX Tilt & Gantry Console Reset - 2267765


Figure 6-5 2267765 Sheet 1, Rev. 0

Section 1.0 — Mains Undervoltage - 2267731 (Not applicable to NGPDU)


Figure 7-1 2267731 Sheet 1, Rev. 0

Section 2.0 — Gantry Temp. Sensors & Power Supply Monitors - 2267766



Figure 7-2 2267766 Sheet 1, Rev. 0

Chapter 8

System Power Control

Section 1.0 — HVDC Bus Control - 2267769

2267769SCH

TO SHEET 2

THE REVISION NUMBER OF THIS SHEET IS THE CONTROL NUMBER OF THIS ENTIRE DOCUMENT. THIS REVISION NUMBER MUST INCREASE WHENEVER ANY SHEET IS CHANGED.					
PCN #	DATE	18APR00	11AUG00		
SHT	2	A	0		

BLOCK PATHNAME \$GEMS_CT_FUNCT_INTCNT/SYS_PWR_CTRL/HVDC_BUS/H2_2267769SCH/H2_HVDC_BUS										BLOCK SHT 1 OF 2							
REV	PCN	DESCRIPTION			MADE	RYAN TOONEN	DATE	17/APR/2000	GE MEDICAL SYSTEMS MILWAUKEE, WI USA	DESIGN TITLE HVDC BUS CONTROL							
MADE BY		DATE			APPR BY		DATE	9/AUG/2000		FIRST MADE FOR H2 SYS FUNC							
REV	PCN	DESCRIPTION			APPR	J. R. SCHMIDT	DATE	10/AUG/2000		DRAWING NO.		2267769SCH	REV	0	SIZE	B	SHT/CONT. ON
MADE BY		DATE			APPR BY		DATE										

Figure 8-1 2267769 Sheet 1, Rev. 0 (CPDU shown; see Figure 8-2 for NGPDU)



Figure 8-2 HVDC Bus Control (NGPDU shown; see Figure 8-1 for CPDU)

2267769SCH

SYSTEMS WITH PERFORMIX TUBE

BLOCK PATHNAME \$GEMS_CT_FUNCT_INTCNT/SYS_PWR_CTRL/HVDC_BUS/H2_2267769SCH/H2_HVDC_BUS						BLOCK SHT 2 OF 2								
REV	PCN	DESCRIPTION		MADE	RYAN TOONEN	DATE	17/APR/2000	GE MEDICAL SYSTEMS		DESIGN TITLE HVDC BUS CONTROL				
MADE BY		DATE	APPR BY	DATE	CHECK	R. VANDEN BOOGARD	DATE			10/AUG/2000	FIRST MADE FOR H2 SYS FUNC			
REV	PCN	DESCRIPTION		APPR	J. R. SCHMIDT	DATE	10/AUG/2000	MILWAUKEE, WI USA		DRAWING NO.		REV	SIZE	SHT/CONT. ON
MADE BY		DATE	APPR BY	DATE	OTHER	DATE					2267769SCH	0	B	2 / -

Figure 8-3 2267769 Sheet 2, Rev. 0


Figure 8-10 2267713 Sheet 3, Rev. 1


Section 4.0 —Gantry Stationary Power and Ground - 2267713

Chapter 9


Image Generation (Global Console - Linux)

Section 1.0 — Scan Data Save

rev. 03/11/02
Figure 9-1 Scan Data Save


Section 3.0 — Pre-Processing, Filtered Back Projection & Post Processing

rev. 03/11/02

Figure 9-3 Pre-Processing, Filtered Back Projection & Post Processing

Section 4.0 — Image Transfer


rev. 03/11/02

Figure 9-4 Image Transfer

Chapter 10

Image Generation (Global Console - Octane2)

Section 1.0 — Scan Data Save

Section 2.0 — Data Restore


Figure 10-2 Data Restore

Section 3.0 — Pre-Processing, Filtered Back Projection & Post Processing


Figure 10-3 Pre-Processing, Filtered Back Projection & Post Processing

Section 4.0 — Image Transfer


Figure 10-4 Image Transfer

Chapter 11

Image Generation (H3 Console)

Section 1.0 — Scan Data Save


Figure 11-1 Scan Data Save

Section 2.0 — Data Restore


Figure 11-2 Data Restore

Section 3.0 — Pre-Processing, Filtered Back Projection & Post Processing


Figure 11-3 Pre-Processing, Filtered Back Projection & Post Processing

Section 4.0 — Image Transfer


Figure 11-4 Image Transfer

GE MEDICAL SYSTEMS

GE MEDICAL SYSTEMS-AMERICAS: FAX 262.312.7434
3000 N. GRANDVIEW BLVD., WAUKESHA, WI 53188 U.S.A.
GE MEDICAL SYSTEMS-EUROPE: FAX 33.1.40.93.33.33
PARIS, FRANCE
GE MEDICAL SYSTEMS-ASIA: FAX 65.291.7006
SINGAPORE