

**ADMINISTRATIVE INFORMATION****FMI No: 25308****FMI Title: Grease Inspection and Protocol Check and Change**

FMI TYPE	CHARGE CODES	IMPLEMENTATION DATES
☐ Mandatory Safety	GEMS - AM:	Issue Date: 04/10/02
☒ Mandatory	GEMS - E:	
☐ Mandatory on Request	GEMS - A:	Close Date: 08/10/02

☒ Staged FMI	Customer List provided by FMI Administration. If you identify additional units within effectivity range, but they are not shown on the customer list, contact FMI Administration in your Pole.
☐	Parts and any special tools are furnished automatically. To address shortages and ordering exceptions, contact FMI Administration in your Pole.
☐	Parts Ordering: GEMS - Americas, Contact your Region logistics coordinator GEMS - Europe, Order FMI kits from GPO/EPL Buc GEMS - Asia, Contact your region logistics coordinator
☒	No Parts are required to complete this FMI.
☐ Batch with	_____
	To minimize costs, combine implementation of this FMI with the FMI's listed for all customers where the effectivity is the same (Please see your customer lists).
☐ Non-Staged	No Customer List is provided. Each service location must identify affected units, and place part orders as indicated in the Special Instructions.

Technical Support

GEMS - Americas	GEMS - Europe	GEMS - Asia
GE Medical Systems Insite Center (414) 524-5300 800-321-7973 Milwaukee, WI 53201, USA	GEMS - E / ESC BP34 F 78533 BUC CEDEX, FRANCE Tel: (33 1) 39 20 00 07 or ESC Fax: (33 1) 30 70 99 70	Asia Support Center Phone: 81-426-56-0033 Fax: 81-426-48-2905

Special Instructions:**Unit Tracked: Model # 2180482 and 2260619****FMI Instructions Part No.: 2338752-100****Do NOT leave unsecured on site**

FMI Completion Certification Report

Mail To

GEMS - Americas

G.E. Medical Systems
Attn: FMI Admin. W-523
P.O. Box 414
Milwaukee, WI 53201

GEMS - Europe

GEMS - Europe
FMI Admin. 322
BP34
F78533 Buc Cedex, France

Send to your
Country FMI or Service
Administrator

FMI No.: 25308

FMI Title: Grease Inspection and Protocol Check and Change

Customer Full Name:

Number and Street Address:

City and Country:

Country Code

Site Number

System I.D.

Unit Tracked:

Model Number	Serial Number	Compl. Date DD/MM/YY	FMI Comp. Code	Service Hours	
				Labor	Travel

FMI Completion Code

1. Modification Complete or previously modified.
2. FMI not required per FMI instruction.
3. Not modified. Unit in storage. FMI kit returned to stock.
4. Modification refused or equipment altered. Customer signature and title required.
5. Unit location unknown or unit not in area. Product Locator notified.

Customer Refusal or Equipment Altered - Modification Pending.

Installer
Comments

Effectivity

The following CT Equipment is affected: All LightSpeed QXi and LightSpeed Plus Systems. Model number 2180482 for QXi and Model Number 2260619 for the Plus systems.

Purpose

The purpose of this FMI is to close two FPRs. The first FPR is FPR1016957-grease deposits from a collimator, which can smudge and cause an artifact. Image Artifacts have been generated and reported on some LightSpeed QXI and LightSpeed Plus systems due to the contamination of the bowtie and the primary copper filter. This contamination is from the lubricating grease used on the filter positioning drive screw assembly. A manufacturing change was made in April 2001, to reduce the amount of grease used.

The second FPR is FPR 10112115-Sinus supine scan protocol contains prone patient orientation. A discrepancy in the protocol name and the patient position indicated exists in the first protocol in the Adult Orbit area on some LightSpeed CT systems. The issue relates to an incorrect patient position setting. This issue may or may not be present on your system currently, depending if this standard default protocol has already been user modified on your system.

Furnished Materials

FMI Kit Part Number 2339937 includes the items listed below.

Qty	Description	Part #
1	FMI Document	2338752-100
1	Directions to Check & Change Protocol	2339144-100
1	Artifacts Caused by Collimator Grease	2339305-100

Tools Required

FOR GREASE INSPECTION:

Bright Flashlight

FOR PROTOCOL CHECK AND CHANGE:

None Needed

People Power

Labor Estimates: 1 hour

Travel Estimates: 1 hour

Procedure 1

Grease Inspection - FPR 1016957

The inspection takes less than 5 minutes. It consists of simply examining the Copper primary filter.

Section 1 – Required Tools

- Bright Flashlight

Section 2 – General Inspection Procedure

- 1) Remove Mylar Scan Window
- 2) Launch Diagnostic Data Collection, DDC, from the Service Desktop.
 - a) Select DIAGNOSTICS.
 - b) Select DIAGNOSTIC DATA COLLECTION.
 - c) Select POSITION TUBE.
 - d) Enter 180 and execute.
 - e) Select STATIC X-RAY OFF.
 - 1) Filter **AIR**.
 - 2) Slice Collimation **Largest Aperture**, 4 x 5.00 for example.
 - f) Select **ACCEPT RX**.
- 3) Using your flashlight, inspect the primary copper filter by looking down into the collimator output port.

Comment: The copper filter should be clean, dent and scratch free. Discoloration is acceptable.

If you see contamination and need to proceed with cleaning, please refer to Service Note T-1449 (P/N 2339305-100) for directions (the Service Note is included in this FMI kit).

Procedure 2

Protocol Check and Change Procedure - FPR 10112115

To determine if your system currently has this issue please check the first protocol in the Adult Orbit area to insure that the name of the protocol and the patient position indicated in this protocol are consistent:

1. Click on 'Protocol Management' at the bottom of the left screen.
2. Click on 'Protocol Management' on the pop up box that comes up.
3. On the list of protocol groups that comes up, click on the second group (2.x)
4. Check the first protocol name that comes up in the resulting list. If the name of the first protocol in this 2nd group list contains the word 'Supine' then continue on with the next step. If the first protocol name here does not contain the word 'Supine', then you can stop, as your system does not have this issue.
5. If the first protocol name in the 2nd protocol group list contains the word 'Supine' then click on this protocol to highlight it and then click the 'Edit' button at the bottom of the list.
6. On the first scan series screen that comes up, check the Patient Position setting. If it is set to Supine then continue to the next step. If it is not then click above the patient icon at the top left until the Patient Position changes to Supine.
7. Click the 'Next Series' button near the bottom of the screen.
8. On the second scan series screen that comes up, check the Patient Position setting. If it is set to Supine then continue to the next step. If it is not then click above the patient icon at the top left until the Patient Position changes to Supine.
9. Click the 'Accept' button at the bottom right of the screen to save the protocol settings.
10. Exit out of the Protocol Management function and locate your system save state MOD or other MOD that is used to backup your scan protocols and insert it into the MOD drive.
11. Select the 'Service' desktop icon on the right screen.
12. Select 'PM' on the service desktop.
13. Select 'System State'.
14. Click on 'Protocols', then click on 'Save'.
15. After the save function completes, click on the 'Dismiss' button.
16. You can then eject the system state MOD and switch out of the service desktop.

FMI Completion

Fill out and return the FMI Completion Certification Report paperwork (or electronically close your FMI dispatch) in accordance with your local procedure.

This Page Intentionally Left Blank.