



| FMI TYPE                                      | CHARGE CODES | IMPLEMENTATION DATES |
|-----------------------------------------------|--------------|----------------------|
| <input type="checkbox"/> Mandatory Safety     | GEMS - AM:   | Issue Date:          |
| <input checked="" type="checkbox"/> Mandatory | GEMS - E:    |                      |
| <input type="checkbox"/> Mandatory on Request | GEMS - A:    | Close Date:          |

☒ **Staged FMI** **Customer List** provided by FMI Administration. If you identify additional units within effectivity range, but they are not shown on the customer list, contact FMI Administration in your Pole.

☒ **Parts** and any special tools are furnished automatically. To address shortages and ordering exceptions, contact FMI Administration in your Pole.

☐ **Parts Ordering:** GEMS - Americas, Contact your Region logistics coordinator  
GEMS - Europe, Order FMI kits from GPO/EPL Buc  
GEMS - Asia, Contact your region logistics coordinator

☐ **No Parts** are required to complete this FMI.

☐ **Batch with**

---

To minimize costs, combine implementation of this FMI with the FMI's listed for all customers where the effectivity is the same (Please see your customer lists).

☐ **Non-Staged** **No Customer List is provided.** Each service location must identify affected units, and place part orders as indicated in the Special Instructions.

| <b>GEMS - Americas</b>                                                                                                               | <b>GEMS - Europe</b>                                                                                                                                | <b>GEMS - Asia</b>                                                                       |
|--------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------|
| <b>GE Medical Systems</b><br><b>Insite Center</b><br><b>(414) 524-5300</b><br><b>800-321-7973</b><br><b>Milwaukee, WI 53201, USA</b> | <b>GEMS - E / ESC</b><br><b>BP34</b><br><b>F 78533 BUC CEDEX, FRANCE</b><br><b>Tel: (33 1) 39 20 00 07 or ESC</b><br><b>Fax: (33 1) 30 70 99 70</b> | <b>Asia Support Center</b><br><b>Phone: 81-426-56-0033</b><br><b>Fax: 81-426-48-2905</b> |

**Do NOT leave unsecured on site**

# **FMI** Completion Certification Report

## Mail To

**GEMS - Americas**

**GEMS - Europe**

G.E. Medical Systems  
Attn: FMI Admin. W-523  
P.O. Box 414  
Milwaukee, WI 53201

GEMS - Europe  
FMI Admin. 322  
BP34  
F78533 Buc Cedex, France

Send to your  
Country FMI or Service  
Administrator

FMI No.: 25350

FMI Title: LightSpeed 5.x w/ Performix Pro (MCS\_7079) Tube -  
Pressure Switch & Pump Replacement

Customer Full Name:

Number and Street Address:

City and Country:

Country Code

Site Number

System I.D.

### Unit Tracked:

| Model Number | Serial Number | Compl. Date<br>DD/MM/YY | FMI Comp.<br>Code | Service Hours |        |
|--------------|---------------|-------------------------|-------------------|---------------|--------|
|              |               |                         |                   | Labor         | Travel |
|              |               |                         |                   |               |        |
|              |               |                         |                   |               |        |
|              |               |                         |                   |               |        |
|              |               |                         |                   |               |        |

### FMI Completion Code

1. Modification Complete or previously modified.
2. FMI not required per FMI instruction.
3. Not modified. Unit in storage. FMI kit returned to stock.
4. Modification refused or equipment altered. Customer signature and title required.
5. Unit location unknown or unit not in area. Product Locator notified.

### Customer Refusal or Equipment Altered - Modification Pending.

Installer

Comments

## Effectivity

| Costumer Name                                               | Scanner Type         | Sys ID       |
|-------------------------------------------------------------|----------------------|--------------|
| Molinette Italy                                             | HP80                 | A5177358     |
| HIA Val De Grace                                            | HP80                 | M4434860     |
| CED                                                         | HP80                 | A5269709     |
| Centro Cardiologico Monzino (Milano)                        | <a href="#">HP80</a> | A5624105     |
| Pontificia Universidad Catolic (Santiago, Chile)            | HP80                 | 26329CT4     |
| Ctr de Lutte Contre Cancer Rene Huguenin (St-Cloud, France) | HP80                 | M4143946     |
| MIA. MATER HOSPITAL                                         | RT                   | 0910 21 2110 |
| AUSTIN & REPAT. HOSP.                                       | RT                   | 0910 21 3079 |
| MIA. JOHN FAUKNER HOSP.                                     | RT                   | 0910 21 3072 |

## Purpose

The purpose of this FMI is to upgrade the Heat Exchanger and Coolant Pump. After these changes, an additional kit must be used to purge all air that may be present in the Heat Exchanger Coolant system.

## Prerequisite FMIs

none

## Related FMIs

none

## Furnished Materials

| Qty | Description                 | Part #      |
|-----|-----------------------------|-------------|
| 1   | 10 psi switch               | 5122148     |
| 1   | Thumbscrew Clamp            | 5122103     |
| 1   | Heat Exchanger Coolant Pump | 5112969     |
| 1   | Teflon Tape Roll            | 46-208982P1 |

## ***Tool Required***

| Qty | Description                           | SW Version | Part #      |
|-----|---------------------------------------|------------|-------------|
| 1   | FMI Document                          | N/a        | 5126334-100 |
| 1   | Performix Pro (MCS_7079)<br>Purge Kit | N/a        | 5122149     |
| 1   | Permanent Marker - Black              | N/a        | N/a         |

## ***People Power***

"4" hours labor and "1.5" hours of travel

## ***FMI Process Flow***




Figure 1: FMI Process Flow

## Procedure

### Preparation or Before you Begin

1. Read this entire document before you begin this procedure
2. Ensure the contents of the package are complete
3. Ensure that your system has a Performix Pro (MCS\_7079) X-Ray tube. There are several ways to do this:
  - a. Go to Reconfig and click on Hardware to look at tube type. The table below lists the possibilities

| Software Version | LightSpeed 5x HP80 Tube Type | LightSpeed 5x RT Tube Type |
|------------------|------------------------------|----------------------------|
| M3               | MX_240MCT                    | MCS_7079                   |
| 04MV13.7         | OEM                          | OEM                        |

**Table 1: Tube Configuration per system type as shown in Reconfig**


Only the OEM and MCS\_7079 selections correspond to the MCS\_7079 tube

- b. Open the right side gantry cover and verify that the hoses are blue and are labeled "Water-Based Coolant"

### Heat Exchanger Pressure Switch Replacement


1. Remove the Side, Top, and Front gantry covers.
2. Disable the Axial Drive power by turning off the AXIAL ENABLE switch on the Service Switch Board.

3. Position the heat exchanger so that the pressure switch is at the top of the heat exchanger. The surface that the pressure switch is threaded into should be as close to horizontal as possible.



**Figure 2: Oil Pressure Switch**

4. Turn off all three service switches on the Service Switch Board (120VAC, Axial Drive Enable, HVDC Enable).
5. Engage gantry rotational lock. Verify that the gantry is locked by attempting to rotate the gantry by hand.
6. With the screw clamp provided, clamp the black hose so that no fluid can pass (see figure).



**Figure 3: Black hose clamp in use**

7. Prepare the new pressure switch with Teflon tape as shown the figure below.

Note: Do not cover the opening of the switch or the system will not work.



**Figure 4: How to prepare the new oil pressure switch**

8. Disconnect the wires connecting the heat exchanger to the pressure switch.
9. Wrap absorbent pad around the threads of the pressure switch on the heat exchanger.
10. Unscrew the pressure switch from the heat exchanger. Just before the switch is free go slowly to absorb any fluid.

Note: There should be no release of fluid from if the clamp is set properly on the black hose.

11. Install the new pressure switch finger tight and then turn an additional 1.5 turns with a wrench.
12. Connect the wires to the pressure switch at the quick disconnect.
13. Remove the clamp from the black hose.
14. Put any tie-wraps back in place to hold wires to the black hose. See figure 3.

**NOTE: If you do not remove the clamp – the system will report overpressure errors via the JEDI generator.**

15. Disengage gantry rotational lock.
16. With a marker, add a “-2” to the part number listed on the heat exchanger, such that the heat exchanger now indicates that it is part number 2383952-2.
17. Next you will need to replace the Pump, and then purge any air from the system.
18. Replace the Heat Exchanger Pump.
  - a. See Appendix A or the System Service Manual for instructions if needed.
19. Purge Air from the Heat Exchanger Subsystem.
  - a. See Appendix B for instructions if needed.
20. Run Tube Warmup

**25350 LightSpeed 5.x w/ Performix Pro (MCS\_7079) Tube - Pressure Switch & Pump Replacement Pg 6/20**

---

21. Perform scout scan
22. If no issues arise, safely restore the system.



# Appendix A

## Heat Exchanger Pump Replacement (Pro16)

### 1 Personnel Requirements

| Required Persons | Preliminary Reqs      | Procedure | Finalization          |
|------------------|-----------------------|-----------|-----------------------|
| 1                | Timing Not Applicable | 30 mins   | Timing Not Applicable |

### 2 Overview

To replace the X-Ray heat exchanger pump:

- Remove old pump
- Mount new pump

|               |            |
|---------------|------------|
| Last Revised: | March 2004 |
|---------------|------------|

### 3 Preliminary Requirements

#### 3.1 Tools and Test Equipment

No tools or test equipment required.

#### 3.2 Consumables

No consumables required.

#### 3.3 Replacement Parts

No replacement parts required.

#### 3.4 Safety



POTENTIAL FOR SHOCK.

VOLTAGE MAY BE PRESENT.

USE PROPER LOCKOUT/TAGOUT PROCEDURES BEFORE REPLACING COMPONENTS.

#### 3.5 Required Conditions

No required conditions.

## 4 Procedure

### 4.1 HX Pump Removal

#### **WARNING**




POTENTIAL FOR SHOCK.

VOLTAGE MAY BE PRESENT.

USE PROPER LOCKOUT/TAGOUT PROCEDURES BEFORE REPLACING COMPONENTS.

1. Remove gantry power. Be certain to follow appropriate Lockout/Tagout procedures.
2. Remove side and front gantry covers.
3. Turn OFF all three switches (**AXIAL DRIVE ENABLE, HVDC ENABLE, 120 VAC**) on the Service Switch Panel.
4. Rotate the gantry so the tube is at the 2 o'clock position.
5. Remove the 4 M6 screws ([Illustration 1](#)) that mount the pump to the 107 weight stack bracket.

**Illustration 1:** HX pump mounting M6 screws




6. Disconnect HX Pump 120 VAC at quick disconnect (see [Illustration 2](#)).
7. Remove two M6 nuts (shown in [Illustration 2](#) behind 120 VAC quick disconnect) for the pump's front mounting plate.

**Illustration 2:** 120VAC connector and front M6 nuts



8. Disconnect pump from heat exchanger at quick disconnect. Line the slot up with the pin and push. The pump should easily disengage (see [Illustration 3](#)).

**Illustration 3:** Align Pin and Slot




9. Disconnect oil hose from X-Ray tube.
10. Remove power interface board cover.
11. Remove pump and hose from gantry.

## 4.2 HX Pump Installation

1. Rotate gantry so that X-Ray tube is at the 2 o'clock position.
2. Connect pump to Heat Exchanger at quick disconnect. Make sure you align the pin and slot, then connect. See [Illustration 3](#) in [HX Pump Removal](#).
3. Route the hose to the tube. It routes underneath the Heat Exchanger shroud, but between the shroud and the clear plastic power interface board cover.
4. Attach quick disconnect and safety mechanism at tube.
5. Reattach the Power Interface board cover.
6. Connect pump 120 VAC connector.
7. Unmount the ORP assembly by removing the two M12 bolts.
8. Install the 4 M6 HX pump mounting screws. If you just replaced the heat exchanger, it may require some force to realign the pump. Using a lever helps but be sure not to damage the hour meter. See [Illustration 4](#).

**Illustration 4:** Using lever to push pump around to align holes



9. Attach the ORP assembly.
  - a. Apply "pre-load" torque to the two M12 bolts:

|                |            |                |               |
|----------------|------------|----------------|---------------|
| 25 ft.-<br>lbs | 33 N-<br>m | 295 in-<br>lbs | 340 kg-<br>cm |
|----------------|------------|----------------|---------------|

- b. Apply final torque to the two M12 bolts:

|                |            |                |               |
|----------------|------------|----------------|---------------|
| 49 ft.-<br>lbs | 66 N-<br>m | 590 in-<br>lbs | 680 kg-<br>cm |
|----------------|------------|----------------|---------------|

10. Attach the two M6 nuts for the front pump.

## 5 Finalization

1. Ensure that proper torque specifications (see [Torque Wrench Information](#)) are followed for all fasteners.

## ***Appendix B – Heat Exchanger Air Purging***

### **Objective of procedure**

The objective of this procedure is to purge air from the liquid tube cooling system. High amounts of air present in the system can lead to a decrease in flow which can then result in potential tube overheating / over pressure problems and / or image quality problems.

### **Estimated Time to Complete**

45 minutes

### **Procedure summary**

Power down gantry

Remove air from purge kit hoses

Attach purge kit between heat exchanger and x-ray tube

Turn on gantry 120VAC to enable heat exchanger pump

Keep power on until there is no milky colored liquid in purge kit

Remove purge kit

Safely restore system

### **Special Tools**

Performix Pro Purge Kit (Part number 5122149)

### **Procedure Details**


#### **GANTRY PREP:**

1. Remove side, and top gantry covers.
2. Turn off all 3 service switched on the Service Switch Board (120 VAC, Axial Drive Enable, HVDC Enable)
3. Remove front Gantry Cover. Cable disconnection is optional based on amount of space desired for procedure.



**25350 LightSpeed 5.x w/ Performix Pro (MCS\_7079) Tube - Pressure Switch & Pump Replacement Pg 13/20**

4. Position Gantry such that quick disconnect on blue hose (between tube and Heat Exchanger) is easily accessible.
5. Loosen Jam Nuts at bottom of quick disconnect bolts and remove the bolts. Do NOT disconnect the hose at this point.



6.




7. Rotate the gantry such that the tube is in between the 5 and 6 o'clock positions and place absorbent pads in base of cover under the hose (shown in figure).

**CAUTION: Pinch Point**

The following work to be performed could create inadvertent gantry rotation causing injury

Engage gantry rotational lock to prevent gantry rotation

8. Engage gantry rotational lock. Ensure that the gantry is locked by attempting to rotate the gantry by hand.



### PURGE BOTTLE PREP:

1. Unpack the purge kit






**25350 LightSpeed 5.x w/ Performix Pro (MCS\_7079) Tube - Pressure Switch & Pump Replacement Pg 15/20**

2. Connect the quick disconnects on the purge kit together.
3. Make sure the cap is tight.
4. Turn the bottle upside down until all air in the hoses is removed. Verify that all air is removed from hoses by noting lack of bubbles flowing into the bottle.
5. With the bottle still inverted, disconnect the hoses from each other.
6. Return the bottle to the upright position.

**ATTACH PURGE KIT TO GANTRY:**

1. Disconnect the hose between the heat exchanger and the tube.
2. Connect the purge kit between the tube and the heat exchanger and set it on absorbent pads on the floor.



3. Loosen the cap on the bottle 1 complete turn to allow for venting during the purge process.

**PURGE PROCEDURE:**

1. Turn on the Heat Exchanger pump by turning on the Gantry 120 VAC switch on the service switch board.
2. Watch for the following:
  - a. LEAKS – if any leaks are detected, turn off the 120VAC switch and repair the leak before continuing
  - b. FLUID COLOR CHANGE – If there is any air in the system, the fluid in the purge kit will become a milky color. If there is no air in the system, the fluid will remain a clear lime green color

---

**25350 LightSpeed 5.x w/ Performix Pro (MCS\_7079) Tube - Pressure Switch & Pump Replacement Pg 16/20**

---

- c. FLUID LEVEL – During purging the fluid level in the bottle may rise and fall approximately 5cm. If the fluid level rises to the top of the bottle, remove enough fluid from the bottle to ensure that the bottle does not overflow.

3. Once the fluid is a clear color (not milky colored) AND the fluid level in the bottle has stabilized turn off the gantry 120VAC to remove power from the heat exchanger pump.

Note: This will usually happen either immediately (indicating no air was present in system before purging), or after some period of time (fluid turns milky color then becomes clear again after approximately 10-15 minutes)

4. Tighten the cap on the purge bottle
5. Elevate the purge bottle so that it is higher than the hose connections.
6. At this height above the hose connections, hold the purge bottle upside-down so that the fluid can drain into the heat exchanger's accumulator.
7. Hold purge bottle in this inverted position for approximately 15-45 seconds.
8. Return purge bottle to upright position on the floor.
9. Disconnect the purge kit from the heat exchanger and tube.
10. Reconnect the tube to the heat exchanger
11. Disengage gantry rotational lock.
12. Rotate the gantry such that the quick disconnects for the hose are easily accessible.
13. Tighten bolts on quick disconnect.
14. Snug the jam nuts on the quick disconnect.
15. Safely restore system for customer use.

## Finalization

Run Tube warmup and verify no errors.