

**ADMINISTRATIVE INFORMATION****FMI No: 25388****FMI Title: GraphicRx File Transfer**

FMI TYPE	CHARGE CODES	IMPLEMENTATION DATES
☒ Mandatory Safety	GEMS - AM:	Issue Date:
☐ Mandatory	GEMS - E:	
☐ Mandatory on Request	GEMS - A:	Close Date:

☒ Staged FMI	Customer List provided by FMI Administration. If you identify additional units within effectivity range, but they are not shown on the customer list, contact FMI Administration in your Pole.
☒	Parts and any special tools are furnished automatically. To address shortages and ordering exceptions, contact FMI Administration in your Pole.
☐	Parts Ordering: GEMS - Americas, Contact your Region logistics coordinator GEMS - Europe, Order FMI kits from GPO/EPL Buc GEMS - Asia, Contact your region logistics coordinator
☐	No Parts are required to complete this FMI.
☐ Batch with	_____
	To minimize costs, combine implementation of this FMI with the FMIs listed for all customers where the effectivity is the same (Please see your customer lists).
☐ Non-Staged	No Customer List is provided. Each service location must identify affected units, and place part orders as indicated in the Special Instructions.

Technical Support

GEMS - Americas	GEMS - Europe	GEMS - Asia
GE Medical Systems Insite Center 800-321-7937 Milwaukee, WI 53201, USA	GEMS - E / ESC BP34 F 78533 BUC CEDEX, FRANCE Tel: (33 1) 39 20 00 07 or ESC Fax: (33 1) 30 70 99 70	Asia Support Center Phone: 81-426-56-0033 Fax: 81-426-48-2905

Special Instructions:**Note: This FMI is not applicable to CT/PET Systems.****This FMI can be installed remotely. For sites with connectivity, no site visit is necessary!****Unit Tracked: Refer to Effectivity list on page 2.****FMI Instructions Part No.: 5191712-100****Do NOT leave unsecured on site**

FMI Completion Certification Report

Mail To

GEMS - Americas

GEMS - Europe

G.E. Medical Systems
Attn: FMI Admin. W-523
P.O. Box 414
Milwaukee, WI 53201

GEMS - Europe
FMI Admin. 322
BP34
F78533 Buc Cedex, France

Send to your
Country FMI or Service
Administrator

FMI No.: 25388

FMI Title: GraphicRx File Transfer

Customer Full Name:

Number and Street Address:

City and Country:

Country Code

Site Number

System I.D.

Unit Tracked:

Model Number	Serial Number	Compl. Date DD/MM/YY	FMI Comp. Code	Service Hours	
				Labor	Travel

FMI Completion Code

1. Modification Complete or previously modified.
2. FMI not required per FMI instruction.
3. Not modified. Unit in storage. FMI kit returned to stock.
4. Modification refused or equipment altered. Customer signature and title required.
5. Unit location unknown or unit not in area. Product Locator notified.

Customer Refusal or Equipment Altered - Modification Pending.

Installer
Comments

Table of Contents

Table of Contents 1

Revision History2

Effectivity2

Purpose2

 SAFETY (PSR/FPR) FIXES INCLUDE:.....2

 CUSTOMER IMPACT2

Prerequisite FMIs3

Related FMIs3

Furnished Materials.....3

Tool Requirements3

Personnel Requirements.....3

Before You Begin4

 VERIFY PREREQUISITE FMI INSTALLED4

 VERIFY FMI 25388 IS NOT ALREADY COMPLETED.4

Procedure5

Revision History

Revision	Date	Reason for Change
1	10/11/06	Initial release
2	10/26/06	Updated Verify FMI 25388 is not Already Complete Section.

Effectivity

The effectivity for this FMI applies to all LightSpeed 1X – 4X and HiSpeed QXI systems with IRIX based GOC1 consoles.

CT systems using Irix/Pegasus. Console numbers to be tracked:

- 2266832, 2266832-2, 2266832-3, 2266832-4
- 2304732, 2304732-2, 2304732-3, 2304732-4
- 2341104, 2341104-2, 2341104-3, 2341104-4
- 2341923, 2341923-2

Note: This FMI is not applicable to CT/PET Systems.

Purpose

The purpose of this FMI is to upgrade IRIX-based CT scanners with the latest safety fixes.

Safety (PSR/FPR) fixes include:

1.1 H3.1M5: Scan location issue on GraphicRx , additional Slices (PSR 13067856)

2.1 H3.1M5: Additional slices added (PSR 13070654)

Customer Impact

This issue involved a few sites that had experienced the addition of one or two extra slices during tilted- multi group head studies. This issue had been addressed by patch F, but the program file misplacement in patch G has prevented the fix from being active. One program file, graphicRx was installed into the wrong directory; therefore, leaving the previous version (patch C1) still on the system. This FMI provides instructions to correct this fault.

Prerequisite FMIs

Kit P/N	FMI Number	Title
5136941	25363	SmartTube Phase 2 SW Patch for IRIX/Pegasus (Patch G.6)

Related FMIs

None.

Furnished Materials

FMI Kit 5191712 consists of the following materials:

Item	Part No	Qty	Description
1	5191712-100	1	FMI 25388 Document
2	2348269-100	1	LS2002 Software Installation Procedure

Tool Requirements

No special tools required for this FMI.

Personnel Requirements

One (1) Service Engineer for .5 hours, plus 1 hour travel time.

NOTE: This patch may be done remotely.

Before You Begin

NOTE: This FMI does not require the field engineer to be on site. This FMI can be done remotely. If these steps are performed remotely, verify the scanner is not in use by the customer. Applications may be up while making the changes. The fix will be implemented after their next system reboot.

Verify Prerequisite FMI Installed

1. Check the CT Software Version

1. As ctuser:

- a. Open a Unix Shell from the Tool Chest menu on the desktop.
- b. Type: **more /usr/g/patch/inst_patches <enter>**

Within the output you should see a line similar to the following:

```
Is2002_M5_patch_G_6
```

NOTE: If Patch G.6 is not installed, STOP! Patch G.6 must be present for this fix to be installed. Patch G.6 can be installed via FMI 25363, SW P/N 5142334. Or, Patch G.6 can be ordered via Coakley Tech.

<http://coakleytech.med.ge.com/cgi-win/default.aspx>

Verify FMI 25388 is not Already Completed.

1. As ctuser, check the file in this directory:

- a. Type: **cd /usr/g/ctuser/display/arch/irix65 <enter>**
- b. Type: **ls -al bin/graphicRx <enter>**

2. If the FMI has already been done, the shell should display the following (Ignore the dates):

```
-r-xr--xr--x  1 ctuser  ctuser  1315368 Apr  8 2006 bin/graphicRx*
```

If the size of the bin/graphicRx file is 1315368, this FMI has already been completed on this site and there is no further action required. Close the dispatch as Code 1 as FMI done. If not proceed with the FMI.

Procedure

1. Change directory to location of fix:
 - a. As ctuser, type: **cd /usr/g/ctuser/display/arch/irix65 <enter>**
2. Save the existing graphicRx file that is in the bin directory:
 - a. As ctuser, type: **mv bin/graphicRx bin/graphicRx.OLD <enter>**
3. Copy the correct graphicRx file to the bin directory:
 - a. As ctuser, type: **cp graphicRx bin/graphicRx <enter>**
4. Set the file permissions of the graphicRx file that is now in the bin directory
 - a. As ctuser, type: **chmod 555 bin/graphicRx <enter>**
5. Verify that the file changes have been correctly made
 - a. As ctuser, type: **ls -al bin/graphicRx* <enter>**

The shell should display the following (ignoring the date and order):

```
-r-xr-xr-x 1 ctuser ctuser 1315232 Apr 8 2006 bin/graphicRx.OLD*  
-r-xr-xr-x 1 ctuser ctuser 1315368 Sep 18 2006 bin/graphicRx*
```

The size of the bin/graphicRx file should be 1315368 .

This verifies that the copy worked correctly.
6. This fix will become active the next time the system restarts.
 - If done onsite, or the system is not in use, reboot at this time
 - or
 - wait until the next reboot by the customer.
7. Close the dispatch for this FMI to indicate completion. Do not leave this document on site.