GE Healthcare

Prospective Gating LightSpeed Plus

OPERATING DOCUMENTATION

2288566-100

Rev 4

© 2000, General Electric Company,
All Rights Reserved

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WARNING

(EN)

- IF A CUSTOMER'S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER'S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR, OR PATIENT FROM ELECTRIC SHOCK, OR FROM MECHANICAL OR OTHER HAZARDS.

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DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent write "Damage In Shipment" on ALL copies of the freight or express bill BEFORE delivery is accepted or "signed for" by a General Electric representative or hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report:

- Call 1-800-548-3366 and use option 8.
- Fill out a report on <http://egems.med.ge.com/edq/home.jsp>
- Contact your local service coordinator for more information on this process.

CERTIFIED ELECTRICAL CONTRACTOR STATEMENT

All electrical Installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE Medical personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

IMPORTANT...X-RAY PROTECTION

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The General Electric Company, Healthcare Group, will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP

Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the General Electric Company, Healthcare Group, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective materials and devices are available. It is urged that such materials or devices be used.

IMPORTANT...RADIOACTIVE MATERIAL HANDLING

Only employees formally trained in radioactive materials handling and this equipment are authorized by the GE Healthcare Radiation Safety Officer to use radioactive materials to service this equipment.

GE Healthcare Services is required to notify the applicable U.S. state agency PRIOR to any source service event involving pin source handling. See NUC/PET Radioactive material guides for specific instruction or contact your EHS Specialist.

A radiation survey must be performed when a pin source has been removed and replaced. See Radiation Survey Form Instructions or contact your EHS Specialist.

Rev 2 (July 21, 2005)

OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the GE Healthcare iTrak Process to report all omissions, errors, and defects in this publication.

REVISION HISTORY

Date	Rev	Change Description
12/22/00	0	Initial Release of Document
04/23/01	1	CQA 1013197: Corrected STC backplane connector
07/24/01	2	Removed references to cables not included with Cardiac option
02/21/03	3	Removed information on unsupported features (RPM device)
4/24/2008	4	Clarified the installation of the RJ45 connector. Updated the Language front matter and use of current branding template. Removed reference information not related to the installation of this kit.

TABLE OF CONTENTS

IMPORTANT PRECAUTIONS	1
REVISION HISTORY	11
TABLE OF CONTENTS	12
1 - LIGHTSPEED PROSPECTIVE GATING OPTION INSTALLATION	13
1.0 Overview	13
1.1 Installation Procedure	13
1-2-1 Preparation - Cover Removal	13
1-2-2 Installation	15
1-2-3 Operators Console RJ 45 Connection	19
1-2-4 Finishing Up - Testing & Cover Re-Installation.....	21
2 - CARDIAC HARDWARE FUNCTION TEST	21
2-1 - Setup Cardiac Monitor -Test EKG/ECG Waveform.....	21
2-1-1 Systems with IVY monitor- EKG/ECG Viewer	21

1 - LIGHTSPEED PROSPECTIVE GATING OPTION INSTALLATION

1.0 Overview

The CT Cardiac (Prospective) Gating Option hardware is shipped ready for immediate installation. The right side is where the Stationary Controller (STC) subsystem is located. However, the Gantry Options Board (GO Board) can be installed on the left side of the Gantry, if desired. There is sufficient cable length in the connection between the STC and the GO board to allow the board to be mounted on either side.

Note

Cardiac can only be operated from one side of the gantry-the side where the Gantry Options board is located. If the customer wishes to operate Cardiac from the left side, then the GO board must be mounted in the left base cover (as opposed to the right, as is described in this manual).

1.1 Installation Procedure

Table 1-1 Part List

NUMBER	QTY	DESCRIPTION/NAME
2228295	1.0	RJ45 FEEDTHRU CONNECTOR
2290289	2.0	PLASTIC RECEPTACLE CAP FOR S & B SERIES LEMO CONNECTORS SIZE 1B
5198566	1.0	cable-gantry - exg pro
5199332	1.0	cable-cardiac - pulmonary
2280889	1.0	GANTRY OPTION BOARD
2286932	1.0	CARDIAC/SV
2288882-3	1.0	11 HOLE BULKHEAD FOR SMARTVIEW AND CARDIAC OR PULMONARY
2288882-4	1.0	8 HOLE BULKHEAD FOR CARDIAC AND/OR PULMINARY
46-328417P74	1.0	SCREW HEXAGON SOCKET 3 MM 10 MM
46-328425P1	4.0	NUT HEXAGON 4 MM
2288566-100	1.0	PERSPECTIVE GATING OPTION - LIGHTSP
2288879	1.0	BRACKET, GANTRY OPTIONS BD MTG
46-328417P4	6.0	SCREW HEXAGON SOCKET 4 MM 12 MM ZINC PLATED, STEEL
46-328432P2	10.0	WASHER LOCK - SPRING 4.1 MM 7.6 MM ZINC PLATED, SPRING STEEL
46-328430P2	6.0	WASHER PLAIN - NORMAL 4.3 MM 9 MM ZINC PLATED, STEEL

The following procedure assumes a right side setup. That is, it assumes Customer access to the cardiac feature via the right side of the gantry when viewing gantry from the rear.

1-2-1 Preparation - Cover Removal

1. Lower table to home (lowest) position.
2. Use proper Lockout/Tagout procedures as required.
3. Use an 8mm Hex wrench to unlatch the side cover from the front cover. See Figure 1-1.


FIGURE 1-1 SIDE COVER LATCHES

4. Remove the right side cover by lifting it upward to release the two (2) latches, located on the top edge of the cover. Once removed, the STC backplane and switches should be exposed.

FIGURE 1-2 SIDES AND TOP COVER CLASP

5. Turn OFF all 3 switches (Axial Drive, HVDC, 120 VAC) on the STC Backplane.
6. If installed, remove both lower left/right/front rear base covers (2, 3 and 4) respectively, to gain access to the Option mount Bracket 1, as shown in Figure 1-3).


FIGURE 1-3 GANTRY BASE COVERS, EXPLODED VIEW

7. Determine which side to mount the gantry option on, and remove the existing blank bulkhead cover from the mounting bracket.

Note:

If there is no Customer preference install this option on the left side of the gantry when viewing the gantry from the back. (Opposite side of Gantry Drive System).

1-2-2 Installation

There are two bulkhead covers included with the kit with different graphics/icons on them. Their physical installation to bracket 1 is the same. Bracket 1 is not included with this kit. It is located on gantry base behind the base covers. Determine which option the customer has and install the most appropriate cover. (In most cases you will be installing cover 2288882-3) which is used in this procedure as the example for bulkhead cover installation.

- 11 Hole Bulkhead Cover- (2288882-3)
 - o Cardiac (ECG/EKG)
 - o Respiratory
 - o SmartView
 - o Ethernet/Remote Cardiac Monitor
- 8 Hole Bulkhead Cover (2288882-4)
 - o Cardiac (ECG/EKG)
 - o Respiratory

1. Attach the appropriate bulkhead cover to the base cover-mounting bracket (Bracket 1), using two (2) nuts (Figure 1-4).


FIGURE 1-4 MOUNTING BRACKET 1 AND COVER PLATES

2. Attach the “Y” cable (PN# 5199332) connectors to the bulkhead cover. (See Figure 1-5).

Note:

The red mark on both cable ends should be either both up or both down.

- a. Remove the outer conical nut from the Lemo cable ends.
 - b. Insert the Lemo connector marked “Cardiac” through the appropriate hole (Heart Icon) from the back of the bulkhead cover. If necessary, adjust the inner nut to accommodate for the thickness of the plate.
 - c. Tighten the outside conical nut first, so that the connector and nut are flush. If necessary, gently adjust the inside nut to ensure that the connector and external conical nut are flush with each other.
 - d. Gently tighten the inner nut.
3. Repeat step 2 for the end of the cable marked “Pulmonary”.
 4. Secure the Cable- Cardiac/Pulmonary Bulkhead cable (PN# 5199332) to J3 of the Gantry Options Board.


FIGURE 1-5 RIGHT SIDE MOUNTING BRACKET, WITH GO BOARD, BULKHEAD, AND CABLES

5. Attach the Gantry Options Board mounting bracket (Bracket 2, shown in Figure 1-6) to the base cover mounting bracket 1. (See Figure 1-4).


FIGURE 1-6 GANTRY OPTIONS BOARD MOUNTING BRACKET (BRACKET 2)

6. Attach the Gantry Options (GO) Board to the mounting bracket (Bracket 2). Make sure that J2 is at the top of the board. (See Figure 1-4 for proper orientation.)


FIGURE 1-7 FUNCTIONAL INTERCONNECT DIAGRAM CARDIAC CABLE LAYOUT

7. On the STC Backplane, connect cable (PN# 2286932) to the STC connector J16.


FIGURE 1-8 LOCATING J16 ON STC BACKPLANE

8. Secure the other end of the cable (PN# 2286932) to J2 of the Gantry Options Board. Route the cable with any cable bundles across the base of the gantry and use cable ties to secure the cable to the existing cable bundles.
9. Attach the two plastic receptacle caps for the Lemo connectors. (PN# 2290289) using the hardware provided.

FIGURE 1-9 ATTACH RECEPTICAL CAPS

1-2-3 Operators Console RJ 45 Connection

1. At the back of the Operators Console locate the RJ45 Cardiac Monitor cable pre-wired at system Installation (PN# 2373436-4 or equivalent. This cable does not come with this kit). Plug this cable into the EKG/ECG output. Refer to figure 1-10.


FIGURE 1-10 CONSOLE CONNECTORS USED FOR EKG/ECG VIEWER

2. Make sure the Cardiac Monitor is setup to output EKG/ECG Data to the Ethernet port. (For IVY monitors this switch is located on the back of the monitor).
3. Locate the RJ45 Feed through (PN# 2228295) located in the kit. There are two types of RJ45 feed through connectors. A metal cased type and a plastic type. To proceed:

FIGURE 1-11: TWO TYPES OF RJ45 FEED THROUGH CONNECTORS

- a. If you have the metal cased version, insert it into the square hole from the front side of the Gating Options panel and snap it in place with the flexible metal snap at the **BOTTOM** of the hole.
- b. If you have the plastic type, refer to Figure 1-12. Push the RJ 45 connector into the square hole from the **FRONT** of the Option Panel by inserting the top catch behind the panel and then swinging the connector down forcing the Flexible Snap to engage into the bottom of the square hole, locking it into place. The Block must be installed so the longest part of the block is facing **OUT**.

WRONG OREINTATION

CORRECT OREINTATION

FIGURE 1-12 CORRECT OREINTATION AND INSTALLATION OF THE RJ45 FEEDTHROUGH

IMPORTANT NOTE!

This plastic type of RJ45 Feed through connector (PN# 2228295) has a hard stop that should be at the top of the connector and a flexible snap that should be at the BOTTOM of the connector when properly installed.

4. At the back of the Gantry Options panel at the Gantry, attach the other end of the RJ45 Cardiac Monitor cable pre-wired at system Installation (PN# 2373436-4 or equivalent) to the RJ45 Feed through block just installed. (Behind the Gating Options panel).

1-2-4 Finishing Up - Testing & Cover Re-Installation

1. Restore Power to the system. The Power LED on the Gantry Options Board should be lit.
2. Install all covers.

2 - CARDIAC HARDWARE FUNCTION TEST

2-1 - Setup Cardiac Monitor -Test EKG/ECG Waveform

Because the Cardiac Monitor types may vary over time, always refer to the operator's manual that came with the Cardiac Monitor for the specific procedures on how to correctly attach the version of monitor present, and prep it for a test waveform. The procedure in section 2-1-1 describes a typical setup.

2-1-1 Systems with IVY monitor- EKG/ECG Viewer

Systems using the IVY Cardiac Monitor or equivalent cardiac monitor design, may route the cardiac information back to the console. Systems with older versions of cardiac monitors may not have this feature. If the system has an IVY cardiac monitor or equivalent:

1. At the back of the IVY Monitor there is a slide switch marked "Disk Storage/ Ethernet. Make sure this switch is in the Ethernet position to enable the data to be viewed at the Operators Console.
2. Connect the data cable (RJ 45) that came with the cardiac monitor from the back of the IVY cardiac monitor to the Cardiac Option panel RJ45 connector. (This will route the cardiac information to the console).
3. Connect the EKG/ECG Leads to the test output of the Cardiac Monitor.
4. Open the Cardiac Viewer at the Host display, to view the cardiac information.

Note:

Any remaining cable length (if extra) should be wrapped in a figure 8 to avoid creation of any loops.


FIGURE 2-1: SAMPLE SETUP

1. Refer to the scanner Operators documentation/Learning & Reference Guide to prepare the scanner for a prospective or cardiac scan.
2. Ensure waveform appears on the Cardiac Viewer on the Host Display. If not:
 - a. Check to make sure the PWR LED is lit on the Gating Options Board.
 - b. Check the cable connection between the console and the Cardiac Option panel (RJ 45 connector) at the gantry.
 - c. For IVY Cardiac monitors with an Ethernet output, check to make sure the switch at the rear of the IVY cardiac monitor is set in the "Ethernet" position.
 - d. Check the ECG/EKG Leads are set up to the test output on the cardiac monitor and there is a test waveform on the cardiac monitor display.
3. Check to ensure the GO Board is receiving the signal by selecting a cardiac-gated protocol and ensuring the gating button on the console reads the (BPM) Beats Per Minute.
4. If HR (Cardiac Trace) does not appear, signal may not be reaching the GOB. Check the cabling. Continue troubleshooting the configuration until the gating option is functioning, as it should.

Note:

To test further functionality, requires the aid of the CT Scanner Operator/ Technologist or CT Applications Specialist depending on the options purchased and version of scanner this option has been installed on. Cardiac Protocols may vary. Refer to the scanner Operators documentation/Learning & Reference Guide for details on how to setup and run tests for different applications functions of the Cardiac feature.

5. This completes the service functional test of the hardware installed.

This Page intentionally blank

GE Medical Systems, a General Electric Company, going to market as GE Healthcare.
3000 N. Grandview Boulevard
Waukesha, Wisconsin 53188
USA

www.gehealthcare.com

Imagination at work