

FCT SERVICE NOTE

SN06030809A

Class A Service information enclosed.

Date: 8-Mar-06

Subject: **Memory Leakage Due to Detailed or Fast Calibrations May Cause Host Console Slowdown or Halt**

Effectivity: **VCT M4 SW RIs 06MW03.4 & GRE non-VCT SW RIs 06MW03.5**

Details: When running "Detailed Calibrations" under "Scanner Utilities" or "Fast Calibrations" under "Daily Prep" on:

- VCT systems running software release 06MW03.4, or
- GRE-based, non-VCT systems running software release 06MW03.5,

these software tools allocate memory, but does not return it for use by other software routines. The host console becomes memory-starved and operation slows, or may even halt/lockup.

Resolution: Reboot the host console before and after running "Detailed Calibrations" under "Scanner Utilities."

Rebooting the host console before running "Detailed Calibrations" will make the optimal amount of memory available for this tool to complete and optimize system-speed.

Rebooting the host console after running "Detailed Calibrations," and before returning the system for customer use, will make the optimal amount of memory available, optimize host console speed, and make the host console most stable for patient scanning.

Recommended daily reboot of the host console will address the memory leakage issue associated with the "Fast Calibrations" tool found under "Daily Prep."

NOTE: 1½ hrs is the typical runtime for a full set of non-VCT detailed calibrations to complete. 2 hrs is the typical runtime for a full set of VCT detail calibrations to complete. Given these runtimes, it is important these tools complete without system halt.

Mike Dejewski
LSI, Molecular Imaging & Computed Tomography