

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.
Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

MI&CT SERVICE NOTE

SN07041216a

Class C information is Property of GE.

This short-term informational note is invalid 6 months after the issue date.

Date: April 12th, 2007

Subject: Pre-mature Fuse holder failure is occurring on some MDAS fuse box assemblies due to overheating of the fuse holder.

Effectivity: Lightspeed CT Systems with MDAS.

Details:

The current design of the Fuse holder (P/N46-170042P26 -20amp) allows the force due to CT gantry rotation to push the fuse against the spring within the fuse holder, possibly allowing separation or a reduced electrical contact to occur. The resultant effect seen by the Fuse holder is an increase in the contact resistance of the Fuse holder, which can create a localized overheating to occur at the Fuse holder and fuse interface. Some of the field sites have reported either a melted solder at the Fuse holder/cable interface or have had early life failures on the Fuse Holder.

A new Fuse holder fix on failure kit, GE part number 5191271, has been tested and verified to show a decrease in this overall temperature rise in the Fuse holder. If the current MDAS Fuse holder shows any problems due to temperature related effects, the fix on failure kit containing the new Fuse Holder should be ordered and the replacement performed.

Resolution:

The recommended corrective action is as follows, upon failure:

1. Order a fix on failure Fuse holder kit, GE P/N 5191271.
2. Perform this fix on failure upgrade per the service instructions internal to this fix on failure kit, GE part number 5191271-100.

Tim Rose,
MI&CT Development Engineer

MI&CT Engineering Manager Signature _____ Mark Nyman _____