

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.  
Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

**MI&CT ESERVICE NOTE  
SN07062211a**

**Class C information is Property of GE.**

This short-term informational note is invalid 6 months after the issue date.

**Date:** June 22, 2007

**Subject:** **load\_darc\_subsystem -restart** script required for all DARC Node software reloads.

**Effectivity:** GOC3 or GOC4 Operator Consoles with **07MW11.x** software.

**Details:** The 07MW11.x LFC procedure and DARC Replacement procedures do Not include the correct replacement script. The normal **darc\_load\_subsystem** script works properly during a LFC or immediately (1<sup>st</sup> time) after the Host Computer has successfully loaded but during a DARC replacement the **darc\_load\_subsystem -restart** switch (flag) command is needed.

When the **darc\_load\_subsystem** script is invoked without any switch it will attempt to load both the OS & Apps in the DARC. If OS has already been loaded once successfully (during the LFC) then the **darc\_load\_subsystem** this script will only load the DARC Application software.

**Failure Example:** The DARC Node OS has loaded successfully, however, the DARC Apps will not load. Next the DARC Node is replaced. If the user invokes the **darc\_load\_subsystem** script after the DARC Node is replaced, the **darc\_load\_subsystem** script sees that the OS was previously loaded and will not reload the OS again but will load the DARC Applications resident on the OS/System Disk. The Operator Console will halt because no OS is present on the DARC Node.

A switch '-restart' was introduced with 07MW11.x software for GOC3 & GOC4 consoles. In the case mentioned above, the script **darc\_load\_subsystem - restart** is required to reload both the OS and Application software resident on the Host during the Replacement DARC Node process.

Failure Scenarios:

- 1.DARC Node Reload
- 2.DARC Node Replacement

**Resolution:** When loading DARC Node software on an Operator Console that already has been loaded with 07MW11.x DARC Node software, the following script will need to be performed: See the latest (attached) DARC LFC procedure details listed below.

[root@hostname scripts] **./darc\_load\_subsystem -restart**

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.  
Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

#### New 07MW11.x LFC Subsection:

### 8.4 OS / Applications Installation for DARC from Host Computer

**No DVD is required during a LFC.** OS will be loaded onto the DARC Node from the Host Computer as long as the DARC BIOS Order is configured properly for your DARC Node hardware type. The Apps will be loaded at this time also. This process takes approximately 35 minutes.

#### **DVD is required for Reloading / Replacement of a DARC Node.**

The user must perform the `[root@hostname scripts] ./load_darc_subsystem --restart` for all DARC Node replacements or any DARC Node reloads.

**NOTE:** The DARC OS was created during the initial Host Applications Software Installation during the LFC. The DARC Apps software was loaded during the initial Host Applications Software Installation. Basically, the DARC OS and Application software is stored on the Host computer.

1. Insert the OS DVD into the Host Computer for a DARC Node reload or replacement.
2. Open a Unix Shell and type the following:
  - a. `{ctuser@hostname} su -`
  - b. Password: `<password>`

**NOTE:** If the Password has not been modified by the site and is not known, the correct course of action is to contact Local GE Service.

- c. `[root@hostname] cd /usr/g/scripts`

**NOTE:** The user must perform the `[root@hostname scripts] ./load_darc_subsystem --restart` for all DARC Node replacements or any DARC Node reloads. If the DARC Application software is corrupt and requires a reload then a LFC must be performed.

- d. Perform one of the following commands as required:
  - **For all full LFC's:**

`[root@hostname scripts] ./load_darc_subsystem`

- **For DARC Node Reloads or Replacements:**

`[root@hostname scripts] ./load_darc_subsystem -restart`

An Attention window appears.

Follow the instructions in the Attention window.

Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

3. Select **[OK]**.
4. The Unix Shell will display output during the DARC load. Example output is as follows:

```
[root@hostnamescripts]# ./load_darc_subsystem -restart
```

```
[load_darc_subsystem] [Start time] Loading the darc subsystem(OS + Apps) takes 30 minutes
```

```
[load_darc_subsystem] [Start time] Loading the darc OS
```

```
[load_darc_subsystem] [Time] Loading the darc OS is complete
```

```
[load_darc_subsystem] [Time] Loading the darc apps
```

```
[load_darc_subsystem] [End time] Loading the darc is complete
```

```
[load_darc_subsystem] [End time] Load started: (Date) (Start time)
```

```
[load_darc_subsystem] [End time] Load ended: (Date) (End time)
```

```
[load_darc_subsystem] [End time] Finished loading darc subsystem.
```

```
[root@hostname]#
```

5. Additionally, an *sh* window displays the OS and Apps load messages.

The DARC Load will stop prematurely ONLY if an error is encountered. An Error Dialog Box will appear with a brief description of the error. Select **[OK]** after performing the instructions displayed on the Monitor.

**NOTE:** If issues occur during the DARC OS load, a pop-up may appear stating that the DARC BIOS must be checked. Jarrell DARC's will perform the Restore Settings - DARC BIOS Boot Order process in [Section 8.6](#). For Jarrell there is only one mode of BIOS operation/settings. Westville DARC's perform the BIOS setup for load and then restore (for normal operation).

**NOTE:** A %done status is displayed during the DARC OS load. This is a timed countdown sequence and does not provide status of the actual load process.

7. When the DARC load has completed, the *sh* window will close. The Unix Shell provides the status of the load process.

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.

Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

**Affected Service Publication:**

5193575-800 - LS 7.X Limited Access Service Info (Adv, Rev 2

5178786-800 - Discovery VCT PET-LS7.X Service Methods (Adv), Rev 2

Associated Class A Docs

**PQR/SPR #:** FCTge 29294

Resa Lambert

CT Console - Service Methods

MI&CT Engineering Manager Reviewer: \_\_\_\_\_