

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare. Disclosure to persons not employed by GE Healthcare or not in possession of an Advanced Service Agreement with GE Healthcare is prohibited. Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

**MI&CT eSERVICE NOTE
SN08091610a**

Class C information. Property of GE Healthcare

This short-term informational note is invalid 6 months after the issue date.

Date: September 16, 2008

Subject: VCT Detector Module FRU kit contents

Effectivity: VCT (64 slice only) detector module FRU kits for HALO 5196401 and 5196402

Details:

Due to an error at the vendor, there were 25 VCT HALO detector module FRU kits that were shipped with 2 hex driver bits of the same size. Two 3mm driver bits are shipped instead of a 2.5mm and 3mm driver bit. The bits are delivered in part numbered bags but there were 3mm bits incorrectly placed in some of the 2.5mm bit bags. The 2.5mm driver bit is used to torque the digital cable screws.

This has been isolated to only 25 total kits but those kits have already been distributed globally and have started showing up on site.

Resolution:

The driver bits are also available in the DAS/Detector service kit (5133871) shipped with every VCT system to be kept on site. If you receive a kit with the wrong bit size, use the bits in the DAS/Detector service kit delivered with the system. This issue has been resolved at the vendor. All tools delivered with the replacement detector modules must be returned per normal parts return process, the bits from these 25 kits will be corrected upon return before reuse of the tool set.

Note: any 2.5mm hex driver bit can be used.

Affected Service Publication:

None

Tracking Number (PQR/SPR/ECR):

PQR 13206947

Steve Merritt
MICT Service Engineering