

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare. Disclosure to persons not employed by GE Healthcare or not in possession of an Advanced Service Agreement with GE Healthcare is prohibited. Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

MI&CT eSERVICE NOTE
SN08102716a

Class C information. Property of GE Healthcare

Date: October 27, 2008

Subject: Flashing LS Gantry Pushbutton and Gantry Display no longer required

Effectivity: CT Lightspeed 2x, 3x and 4x Gantrys

Details: DIR 2385487-100 and the associated software gcflash_package.exe, and cable (PN#2381519) are no longer required. The procedure is obsolete and should not be used.

Resolution:

Any parts remaining in stock for the Gantry Display part numbers 2256730, 2256730-2 and Gantry Control Panels part # 2256728, 2256729 have been re-worked flashing these units is no longer required. If one of these parts fail the correct course of action is to use up remaining stock and follow the FE Assist instructions to acquire the correct replacement parts, which do not require any flashing in the Field.

FE Assist notes have been added to these parts with replacement part instructions.

Need to replace Gantry Display	2256730	Order	2256730-3
Need to replace Gantry Display	2256730-2	Order	2256730-3
Need to replace Gantry Control panel	2256728	Order	5139591
Need to replace Gantry Control panel	2256729	Order	5139590

Service publication DIR 2385487-100 has been obsoleted and has been removed from affected Service publications.

Affected Service Publication:

DIR 2385487-100

DIR 2360027-300

Tracking Number (Complaint # or SPR #): 13201982

Douglas R. Hofstetter
MICT Service Engineering