



GE Medical Systems

Technical Publication

Direction 2237373-100

Revision 0

**Service Note T-1421, Field Engineer's
using Service CDRom**

RESTRICTED TO GENERAL ELECTRIC MEDICAL SYSTEMS (GEMS)

**THIS MATERIAL IS THE PROPERTY OF GEMS AND CONTAINS
INFORMATION WHICH IS RESTRICTED TO USE BY GEMS
MANUFACTURING, ENGINEERING AND SERVICE PERSONNEL.
ANY REPRODUCTION IN WHOLE OR IN PART IS STRICTLY
PROHIBITED. NEITHER THIS MATERIAL NOR ANY INFORMATION
CONTAINED HEREIN SHALL BE DISCLOSED TO ANYONE
OUTSIDE OF GEMS MANUFACTURING, ENGINEERING AND
SERVICE PERSONNEL.**

Copyrighted © by GE Medical Systems, 1999



Page Left Intentionally Blank

LEGAL NOTES

PROPRIETARY TO GE MEDICAL SYSTEMS

Permission to use this Advanced Service Software and related documentation (herein called the material) by persons other than GE Medical Systems employees is provided only under an Advanced Service Package License relating specifically to this Proprietary Material. This is a different agreement from the one under which operating and basic service software is lessened. A license to use operating or basic service software does not extend to or cover this software or related documentation.

If you are a GE Medical Systems employee or a customer who has entered into such a license agreement with GE Medical Systems to use this proprietary software, you are authorized to use this Material according to the conditions stated in your license agreement. However, you do not have the permission of GE Medical Systems to alter, decompose or reverse-assemble the software, and unless you are a GE employee, you may not copy the Material. The Material is protected by Copyright and Trade Secret laws; the violation of which can result in civil damages and criminal prosecution.

If you are not party to such a license agreement or a GE Medical Systems Employee, you must exit this Material now.

GE EXCLUSIVE USE ONLY

Restricted To General Electric Medical Systems (GEMS)

This Material is the property of GEMS and contains information which is restricted to use by GEMS Manufacturing, Engineering and Service personnel. Any reproduction in whole or in part is strictly prohibited. Neither this material nor any information contained herein shall be disclosed to anyone outside of GWMS Manufacturing, Engineering and Service personnel.

If you are not a GE Medical Systems Employee, you must exit this material now.

TRADEMARKS

Adobe, the Adobe logo, Acrobat, the Acrobat logo, Exchange, and PostScript are trademarks of Adobe Systems Incorporated or its subsidiaries and may be registered in certain jurisdictions.

Microsoft is a registered trademark and Windows is a trademark of Microsoft Cooperation.

All other products and their name brands are trademarks of their respective holders.

COPYRIGHTS

All Material Copyrighted(c) 1999 by the GE Company, All rights reserved.

OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the GEMS CQA Process to report all omissions, errors, and defects in this publication.

REVISION HISTORY

Revision	Date	Reason for change
0	02/09/99	Initial Release

TABLE OF CONTENTS

PROPRIETARY TO GE MEDICAL SYSTEMS	3
GE EXCLUSIVE USE ONLY	3
TRADEMARKS	3
COPYRIGHTS	3
OMISSIONS & ERRORS.....	3
PURPOSE	5
PROBLEM STATEMENT	5
TOOL.SOLUTION STATEMENT	5
TOOLS REQUIRED	5
PROCEDURE:	5

SETUP AND USING CD-ROM 2152913

PURPOSE

Procedure for loading Adobe Acrobat Search 3.01. Also included, is the procedure for setting up a shortcut icon to the Table of Contents (TOC) for the CDROM. Procedure for setting up indexes for advanced searches is also included.

PROBLEM STATEMENT

Field Engineer's have not been trained properly to utilize their laptop's as an effective service

TOOL.SOLUTION STATEMENT

Provide a procedure to load Adobe Acrobat, create a short-cut icon to the Table of Contents, and select individual indexes for advanced search capabilities.

TOOLS REQUIRED

Laptop Computer and 2152913 Rev. x CT Service Information CDROM, distributed by Service Engineering (X being the most recent rev).

PROCEDURE:

- 1.) Boot your laptop to the desktop.
- 2.) Insert Vol. 2 of the CT Service Information CDROM, 2152913 Rev. X (*X being the most recent rev*).
- 3.) Double click and open the **MY COMPUTER** icon.
- 4.) From **MY COMPUTER** double click and open the **CDROM** icon (*normally designated as drive (D)*).
- 5.) Once you have opened the CDROM directory, double click on the **READERS** folder. In readers, double click on the **WINDOWS** folder. In **WINDOWS**, double click on the **ENGLISH** folder. In the english folder, double click on **WIN95**. In the Win95 folder, double click on the **SETUP.EXE** icon. For those Field Engineer's who are comfortable using the laptop, simply use the following path from the **RUN** command:
D:\readers\windows\english\Win95\setup.exe
- 6.) Click on **YES** in the InstallShield Self-extracting EXE window
- 7.) After the initial setup is complete, you will have a Adobe Acrobat 3.01 Setup window. If you wish to install Adobe, select the **NEXT** button.
- 8.) The next window will be the Software License Agreement window, select **YES**.
- 9.) At the Choose Destination Location window, select **NEXT**. Adobe is being loaded.

- 10.) At the Setup Complete window, select **FINISH**. When the load is finished you will get the read me notes for Adobe. You can close the read me notes by single clicking on the **X** at the top right corner of your screen. Then click on the **OK** on the setup is complete window.
- 11.) Adobe Acrobat search 3.01 is now loaded, please close all open programs. We will now create a short-cut icon to link you directly to your CDROM Table of Contents (TOC).
- 12.) Replace Vol. 2 of Service CDROM with Vol. 1.
- 13.) Right click once on an open area of your desktop. This will produce a box of selections.
- 14.) Drag your pointer to the new selection. This will activate a new box of selections. Highlight **shortcut** and single click.
- 15.) Choose the **BROWSE** button on the Create Shortcut window.
- 16.) In the browse window, select **DRIVE D** (CDROM) from the **Look in** box. If **drive D** is not selected, simply place your pointer on the down arrow in the **LOOK IN** box and select the **D DRIVE**.
- 17.) Go down to the **Files of Type** box and select all files. You will now see a **cd_toc.pdf** folder.
- 18.) Double click on the **cd_toc.pdf** folder, and you will now see it in the command line of the Create Shortcut window. Select the **NEXT** button and you will be at the Select a Title for the Program window. If you would like to rename your icon, you may do it here, or simply press the **FINISH** button.
- 19.) You now have an icon that will link you directly to Adobe Acrobat, and take you directly to the Table of Contents(TOC). Double click on the Table of Contents(TOC) icon, this will take you directly to the CDROM TOC. Use your page down key and you will find at the bottom of page two, the hot buttons for the different directories (i.e., AW, HSA-Cti, Cti, HSA-Z, HSA-RP, PM, HLA, CT9800, XTRAS, YMS, SRVNOTES, and FMI's. The hot buttons are an excellent way to start an initial search. You will see which documents are included in each directory.
- 20.) We will now setup your indexes for more advanced Adobe searches. With Adobe running, single click on the button with the sheet behind the binoculars.
- 21.) You should now have the Adobe Acrobat Search window on your desktop. Single click on the **INDEXES** button. You now have the index selection window.
- 22.) Single click on the **ADD** button. The add index window is now opened. Be sure you are on the (D) directory in the "**Look in**" box.
- 23.) In the D directory you will see the different choices you have for search indexes. (i.e, **HLA,AW, ...**).
- 24.) Double click the index you want, it will then appear in the "**Look in**" box of the **Add Index** window. You will now also notice a choice called **index.pdx** in the selection window.

25.) Double click on `index.pdx` and you will see your selection in the Index Selection window. Now single click on the button.

26.) You can now setup your search indexes any way you like, by simply repeating steps 20 through step 25.

Having the Adobe Acrobat search software loaded, along with Table of Contents shortcut established, will improve the ability of a Field Engineer to interact smoothly with the OLC and improve Customer Satisfaction and Uptime.

Dale Borowiak CT Online Engineer



GE Medical Systems

CT SERVICE ENGINEERING, NB-918

P.O. BOX 414

MILWAUKEE WISCONSIN 53201-0414 U.S.A.

Fax: (414) 785-5874