

Page Description – Network for BB

View Data

- [View Alarms](#)
- [View Event Log](#)
- [View Minute Data](#)
- [View Insite2 Log](#)
- [View Insite2 Status](#)
- [View Analog Inputs](#)
- [View Digital I/O](#)

Configuration

- [Set Time/Date](#)
- [Magnet Type](#)
- [Controls](#)
- [System ID](#)
- [Network](#)
- [Insite2](#)

Help Windows

- [Error Codes](#)
- [Minute Data Help](#)

GE Healthcare

Configure Network IP

IP Configuration

☐ No network connection when not in Web Mode.

☐ IP addresses assigned by DHCP server.

☒ Static IP Address

IP Address: 3.231.48.61

Subnet Mask: 255.255.252.0

Gateway: 3.231.48.254

Note: Changes made on this page will not take effect until after you exit Web Mode.

DNS Configuration

☐ DNS addresses assigned by DHCP.

☐ Static DNS Addresses

1st DNS: []

2nd DNS: [] (optional)

3rd DNS: [] (optional)

[Help](#) [Submit Changes](#)

Only for direct internet connection. Click in **DHCP** if site go to the internet using DHCP server , consult IT site admin

Direct Internet no DHCP & Insite1: please fill Static IP address of the MM3 along with Subnet mask.

Gateway: Direct connection: Ask Hospital IT manager if needed ; Insite1 BB types: Check e-checkout web page

Under Insite1 type of connection, don't select any box bellow

Only for direct internet: Select only If site has DNS Server , consult IT site admin

Only used for direct internet, do not select DNS configuration box, if you are connecting via Insite1 pipeline. In case you select DNS by mistake, uncheck the same and reboot MM3. Select only If site has no DNS server, consult IT site admin to get DNS IP address (universal DNS IP could be used 194.2.0.50 ; 195.68.0.1).

Note 4: User must log out of Web mode for the device to start the communication process

Page Description – Insite 2 for BB

[Leave Message](#)

View Data
[View Alarms](#)
[View Event Log](#)
[View Minute Data](#)
[View Insite2 Log](#)
[View Insite2 Status](#)
[View Analog Inputs](#)
[View Digital I/O](#)

Configuration
[Set Time/Date](#)
[Magnet Type](#)
[Controls](#)
[System ID](#)
[Network](#)
[Insite2](#)

Help Windows

GE Healthcare

Insite2 Configuration

[Help](#)[Submit Changes](#)

System ID and Magnet Type

WARNING: Do not click the submit button on this page until after the correct System ID and Magnet Type have been selected in previous configuration pages. On Monitor has registered with the back office, it is difficult to change that System ID and Magnet Type combination in the back office.

System ID	EDUCTRMR9C	Change
Magnet Type	HFO	Change

InSite Connection Type

☒ Broadband	Broadband Polling Time 5 minutes
☐ Dialup Modem	
☐ Mobile Satellite Modem	
☐ No InSite Connection	

If on a ISDN line this number will need to be changed to reduce telephone line cost. This field exists only in 1.5.2 soft (will do in the checkout phase)

Broadband selected

Note: User must select the correct magnet type prior to connecting to InSite2. If not another site visit will be required

Note 2: User must log out of Web mode for the device to start the communication process

Page Description – Insite 2 for BB

The screenshot shows the 'Remote Service Configuration' and 'HTTP Proxy Server Configuration' sections of the Insite 2 for BB web interface. A sidebar on the left contains navigation links for 'View Data', 'Configuration', and 'Help Windows'. The 'Insite2' link under 'Configuration' is highlighted with a red arrow. Callouts provide additional context for the configuration fields.

View Data

- [View Alarms](#)
- [View Event Log](#)
- [View Minute Data](#)
- [View Insite2 Log](#)
- [View Insite2 Status](#)
- [View Analog Inputs](#)
- [View Digital I/O](#)

Configuration

- [Set Time/Date](#)
- [Magnet Type](#)
- [Controls](#)
- [System ID](#)
- [Network](#)
- [Insite2](#) (highlighted with a red arrow)

Help Windows

- [Error Codes](#)
- [Minute Data Help](#)
- [Patch Software](#)

Remote Service Configuration

Enterprise URL:

Protocol:

HTTP Proxy Server Configuration

☒ Proxy Name/IP Port

☐ Account Username Password

Modem Parameters

Phone Number:

User Name:

Password:

Country/Region:

AT%T19,0,34

ATS0=0

AT&F0 X4 &D3 E0TM

Insite2 Retry Number:

Callouts:

- Enterprise url from installation procedure: <https://us1-ws.service.gehealthcare.com:443>
- Direct Connection: Proxy and port information, only if customer uses proxy server, can be obtained from Hospital IT Manager. Insite1 type of connection always use this proxies, according to the connection type:
 - VPN/ADSL: proxy setting to 150.2.1.251, port 8002
 - ISDN: 150.2.1.251, port 8002
 - NHS: 150.2.1.251, port 8002
 - NHS N2: 62.130.238.140, port 8002
 - NHS N3: 10.190.64.89, port 8002
 - SJU: 82.136.152.78, port 8002
- Direct Connection: Check with Hospital IT manager if authentication is needed
- Modem information is grayed out once broadband is selected

Note: User enters data obtain from the InSite2 Assessment form or troubleshooting flowchart to set up correct communication configuration

Note 2: User must log out of Web mode for the device to start the communication process

Page Description – Insite 2 for BB

Note: ISDN will work only in the following countries:

AUSTRIA
BELGIUM
CZECK REPUBLIC
DENMARK
FINLAND
FRANCE
GERMANY
GREECE
GUADELOUPE
HUNGARY
ISRAEL
ITALY

LUXEMBOURG
MARTINIQUE
NETHERLANDS
NORWAY
POLAND
PORTUGAL
SOUTH AFRICA
SPAIN
SWEDEN
SWITZERLAND
TURKEY
UK