SERVICE MANUAL

GE Healthcare Coils USA Instruments

Document 2381683-9

Revision 3

GE Signa[®] Excite[™] 3.0T 8-CHANNEL TORSO ARRAY

GE Catalog Part Number: M1385AY

2381683-9
Rev 3
© 2010, General Electric Company,
All Rights Reserved

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation “**damage in shipment**” written on **all** copies of the freight or express bill **before** delivery is accepted or “signed for” by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier **immediately** upon discovery, or in any event, within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this **14**-day period.

File a report with

- Call 1–800–548–3366 and use option 8.
- Fill out a report on <http://egems.med.ge.com/edq/home.jsp>
- Contact your local service coordinator for more information on this process.

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TABLE OF CONTENTS

SECTION 1 – INTRODUCTION	8
1-1 Product Identification and Shipping List	8
1-2 Compatibility	9
1-3 Related Documentation	9
1-4 Environmental Requirements	9
1-5 Theory of Operation	9
section 2 – setup and calibration	11
section 2 – setup and calibration	12
2-1 Coil Installation	12
2-1-1 Special Install Notes	12
2-1-2 Configuration	12
2-2 Installation Functional Checks	12
2-3 Periodic Quality Assurance Check	12
Section 3 – Functional Checks	13
3-1 Scanner Verification	13
3-2 Coil Imaging Performance Verification	13
3-2-1 Tools Required	13
3-2-2 Data Collection	13
3-2-3 Test Explanation and Results Interpretation	25
3-3 External Cable Check	25
3-4 PIN Diodes Check	25
3-5 Mechanical Hardware Check	25
3-6 Troubleshooting Tips	26
Section 4 – Maintenance	28
4-1 Coil Care	28
4-2 Special Care Requirements	28
Section 5 – Replacement	29
SECTION 6 – Renewal Parts	30
6-1 Field Replaceable Units	30
6-2 Other Replaceable Accessories	30
SECTION 7 – Appendix	31
7-1 Schematics	31
7-2 Coil Configuration	32
REVISION HISTORY	33

SECTION 1 – INTRODUCTION

1-1 Product Identification and Shipping List

To identify the 8-Channel Torso Array, refer to the coil labels (see *Figure 1*).

Figure 1: Coil label

Figure 2 shows a picture of the 8-Channel Torso Array.

Figure 2: The 8-Channel Torso Array.

SHIPPING LIST – TABLE 1-1

Box #	Part Name	GE Part #	USAI Part #	Qty
1	Coil	2381683-2	100401	1
1	Operator Manual Packet	2422146-500	2422146-500	1
1	Service Manual	2381683-9	780034	1
2	Patient Liner	2381683-5	502548	1
2	Table Pad	2381683-6	150485	1
2	Spacer Pad	2381683-7	150486	1

1-2 Compatibility

This coil is compatible with the Signa[®] Excite[™] 3.0T MR system.

1-3 Related Documentation

Operator's Guide, GE Part Number 2422146.

Signa[®] Excite[™] 3.0T Service Methods CD, 2389160-1.

1-4 Environmental Requirements

Storage Requirements

The coil should be stored in the scanner room.

Dimensions

43.3" W x 15.0" D x 13.0" H (110.0 cm x 38.1 cm x 101.9 cm)

Weight

14.0 lbs (6.4 kg)

1-5 Theory of Operation

The block diagram of the 8-Channel Torso Array is shown in *Figure 3*. *Figure 4* shows the block diagram of the circuit layout. Also refer to the schematic in the Appendix. The receive-only coil is designed to give optimum signal to noise ratio and uniform coverage of the abdominal and pelvic regions, including the liver, spleen and renal arteries. The coil package consists of anterior and posterior sections, which are positioned above and below the region of interest.

The coil is actively decoupled from the RF transmit coil during transmit by means of multiple RF choking circuits or RF chokes. Each coil element has three active decoupling circuit elements and two passive decoupling circuit elements, which can be switched on actively or passively. The passive RF chokes are switched on using small-signal diodes. The small-signal diodes are turned on by the induced RF voltage in the chokes, coupled from the transmit field (passive coupling). Each small-signal diode is turned on when the induced RF voltage reaches about 0.5 volt. The active RF chokes are switched on using PIN diodes. The PIN diodes are turned on by a forward biasing DC current source supplied by the MRI system (active decoupling). When either the PIN diodes or the small-signal diodes are turned on, the RF choking circuits become very high impedance blocks (typically above 3 k-ohm), compared to the other circuit elements (less than 70 ohm). These high impedance elements segregate the coil circuitry into several isolated electrical segments, preventing any current flow in the coil circuit that might be induced by the body coil transmission. In the event that one of the diodes fails, the signal diode will form a short circuit and permanently de-resonate the coil.


Figure 3: Block diagram of the coil.


Figure 4: Block diagram of the circuit layout.

SECTION 2 – SETUP AND CALIBRATION

2-1 Coil Installation

2-1-1 Special Install Notes

None.

2-1-2 Configuration

The name for this coil is 8US TORSOPA. Add the coil using the Configuration File Manager. Refer to: Service Methods CD; Software; System Level Software Utilities; Configuration File Manager. For reference, the configuration information is included in the Appendix; however, that configuration is only current as of the printing of this manual.

2-2 Installation Functional Checks

1. Perform system level Signal to Noise Check. Refer to Service Methods CD; Functional Checks; System Level Functional Checks; Signal to Noise Check.
2. Perform Section 3-2 Coil Imaging Performance Verification.

2-3 Periodic Quality Assurance Check

On a periodic basis, such as during planned maintenance, perform the quality assurance checks as outlined below to ensure the coil is operating properly.

1. Check external cable for cracks or cuts.
2. Perform Section 3-2 Coil Imaging Performance Verification.

SECTION 3 – FUNCTIONAL CHECKS

3-1 Scanner Verification

Perform system level Signal to Noise Check. Refer to Service Methods CD; Functional Checks; RF Path; Signal to Noise Check.

3-2 Coil Imaging Performance Verification

NOTE: Legacy coils were shipped with a coil specified phantom and positioner. New coils do not ship with a phantom; phantoms come in a unified phantom set with the MR system. New coils do ship with a phantom positioner.

3-2-1 Tools Required

TOOLS REQUIRED – TABLE 3-1

Description	Qty	Effectivity	GE Part #
Phantom	1	Legacy	2381683-3
Phantom Positioner	1	Legacy	2381683-4
TL Unified Phantom, SiOil	1	Unified	5343347-2

3-2-2 Data Collection

3-2-2-1 Legacy Phantom Procedure

Step 1 – Place the 8-Channel Torso Array on the patient cradle (see *Figure 5*). Note the serial number of the coil.

Step 2 – Position the phantom positioner on the coil as shown in *Figure 6*.

Figure 5 8 Channel Torso Coil Setup

Figure 6 Phantom Positioner Setup

Step 3 – Position the phantom on the phantom positioner as shown in *Figure 7*.

Step 4 – Place the anterior section on the phantom and fasten the Velcro straps on both sides (see *Figure 8*).


Figure 7 Phantom Setup

Figure Anterior Coil Setup

Step 5 – Establish an axial landmark on the line on top of the coil. Refer to Section 3-2-2-3 to run the MCQA test with the Legacy Phantom.

3-2-2-2 Unified Phantom Procedure

Figure 3-2-2-2-1

Place the posterior section of the coil on the table. (Figure 3-2-2-2-1)


Figure 3-2-2-2-2

Place the positioner of Torso into gaps present on the posterior section of the coil as shown in Figure 3-2-2-2-2

Figure 3-2-2-2-3

Place the TL Unified phantom as shown in figure 3-2-2-2-3


Figure 3-2-2-2-4

Place the TL unified phantom such that the phantom is well centered on both S/I direction as shown in Figure 3-2-2-2-4 and R/L direction as shown in arrows Figure 3-2-2-2-5.


Figure 3-2-2-2-5

Figure 3-2-2-2-6


Figure 3-2-2-2-7

Figure 3-2-2-2-8

Place the anterior section (figure 3-2-2-2-6) in such a way that the anterior of the coil is well set in to the positioner as shown in figure 2.3-7 and fasten the Velcro straps so that the coil firmly snugs the phantom (Figure 3-2-2-2-8).

Figure 3-2-2-2-9

Connect the coil to respective port of the system LPCA. Landmark the coil on the cross mark as shown in Figure 3-2-2-2-9 and press advance to scan.

Refer to Section 3-2-2-3 to run the MCQA test with the Unified Phantom Set.

3-2-2-3 Steps to run MCQA

SETUP AND STARTING Multi-Coil Quality Assurance (MCQA) tool USING LEGACY PHANTOM SET

WARNING: All RF coil related tests must be run on a system that is well-calibrated and passing all system tests, especially white pixel, correlated noise, and MCR (Multi-Coil-Receive) tool.

1. If possible, run the Multi-Coil Chain Receive (MCRv) tool to validate receive chain function before running the MCQA tool. If the MCRv tool has missing or failed results, the MCQA tool will detect it and a pop-up will appear, asking if you want to continue.
2. From the Common Service Desktop (CSD) select [Image Quality], *Multi-Coil QA Tool* and Click here to start this tool. The Coil Log Viewer will open as shown in figure 4.1.3-1.


Figure 4.1.3-1

MCQA Starting window

3. The current coil field will be automatically filled in, based on the Coil ID of the coil connected to the LPCA. The names listed will be the physical coil names for each coil.

NOTE: Coil names are specific to the coil connector. Legacy coils should not have reference to the connector, but any new coil with the Bendix connector will have an HDx/HDe identifier in the coil name. During download, the tool should detect if the wrong coil was selected, and give an error. The user will then have the ability to pick again.

4. Enter the *Coil Serial #*: into the provided field.
5. From the MCQA Tool GUI, press the [Start] button to begin the test.
6. A warning dialog (figure 4.1.3-2) will be shown, asking the user to verify the coil/phantom set-up and position. Ensure the phantom and coil setup is proper and select [Yes] to continue.
 - a. The status window of the MCQA Tool GUI will continuously update to give information on what the tool is doing at any point in time.
 - b. A time bar will appear, showing approximate total test time, elapsed time and percent complete, as shown in figure 4.1.3-3.


Figure 4.1.3-2

Phantom/Coil, Landmark Warning Popup


Figure 4.1.3-3

Status and Time Bar

7. Once the test has completed, the name of the Results file will be listed, along with the Pass/Fail status of the test.
8. Record the Measured ISNR and the LSL for each channel displayed in the Status Window on a separate sheet. (See Table 4.1.3-1)

Table 4.1.3. 1: Data Collection Table Example

Si No	Test No.	Element No.	LSL	ISNR

9. Select the [QUIT] button to exit out of the MCQA Tool.
10. Remove the coil and phantom from the system.

Setup and Starting Multi-coil QA (MCQA) Tool USING unified phantom KIT

All RF coil related tests must be run on a system that is well-calibrated and passing all system tests. The system should have passed "Install In Spec" (IIS), especially white pixel, correlated noise, and MCR (Multi-Coil-Receive) tool.

1. From Common Service Desktop (CSD) go to Service Browser and select [Image Quality], Multi-Coil QA Tool and then Click here to start this tool as shown in Figure 4.2.3-1


Figure 4.2.3-1

2. Upon start-up, the coil that is connected to the system gets selected in the *Current Coil* field as shown in Figure 4.2.3-2


Figure 4.2.3-2

3. A pop-up displays requesting you to, **Please select correct Phantom Set** as shown in Figure 4.2.3-3 below. Click on [OK] to continue.

Figure 4.2.3-3

NOTE: For Optima MR450w, the Unified Phantom Kit is the **ONLY** available phantom set and will be automatically selected in the [Phantom Set] field of the tool.

4. Select the Unified Phantom Set from the drop-down menu as shown in figure 4.2.3-4


Figure 4.2.3-4

NOTE: For the Legacy Phantom Kit, please select LEGACY (old phantom).

5. The current coil field will be automatically filled in, based on the CoilID of the coil connected to the LPCA. Enter the serial number of the coil being tested in the *Coil Serial Number* field.
6. Click on [Start] to begin the automated test as shown in figure 4.2.3-5. Depending on the number of test locations (complexity of the coil) the test may take from 3 to 20 minutes.


Figure 4.2.3-5

7. A confirmation pop-up displays, asking to continue. Click on [Yes]. See figure 4.2.3-6.

Figure 4.2.3-6

8. Upon start-up, a message stating, **Phantom placement and coil landmarking are critical for repeatable results** will appear. If the landmark has been set correctly and there are no air bubbles in the phantom, click [Yes] to continue. See figure 4.2.3-7.

Figure 4.2.3-7

NOTE: The Status window of the MCQA Tool GUI will continuously update to give information on what the tool is doing at any point in time. A time bar (figure 4.2.3-8) will appear, showing approximate total test time, elapsed time and percent complete. Once the test has completed, the name of the Results file will be listed, along with the Pass/Fail status and the results of the test in the MCQA Tool GUI.


Figure 4.2.3-8

9. When the test is complete, test results display on the screen (figure 4.2.3-9) The PASS/FAIL status shows PASS if all coil elements are functioning properly. If any coil elements display FAIL, call the GE Service Representative for coil repair.


Figure 4.2.3-9

10. The MCQA Tool GUI displays **Fail** for reasons including, but not limited to:
 - a. Bad coil element
 - b. Incorrect phantom used for the test
 - c. Incorrect positioning/placement of the phantom

11. Click on [Quit] button to exit MCQA Tool

Viewing Report Results:

After running the Multi-Coil QA test, select [View Report Details] to review the results. To review past results, select *File, Open Results File* and select the desired coil results file. The Results Viewer will open as shown in *figure 4.3-1*. The results file name and Pass/Fail Results shown on the tool GUI will also be listed across the top of the viewer.


Figure 4.3-1

3-2-3 Test Explanation and Results Interpretation

The 8-Channel Torso Array is a moderately complex phased array. The Quality Assurance tool acquires data looking at the composite image (like the customer sees). Individual element performance is evaluated using the composite image.

Signal images are acquired using a FSE-XL sequence with CV saveinter set to 1 to save intermediate images into the image database. The scan is taken twice for a total of 18 images. Next, a no-RF noise scan is taken, also with saveinter on for a total of 9 noise images. Individual element signal and noise images no longer yield results that are displayed to the user.

The test employs two types of image analyses. Separate signal and noise images are evaluated for individual element signal, noise and SNR. In addition, a filtered version of the NEMA SNR measurement method is also used. The filter takes out adverse effects of subtle ghosting on noise measured in difference.

3-3 External Cable Check

If the signal on any channel is low, refer to *Figure 15* to ensure that none of the B, C or D row pins are bent.

3-4 PIN Diodes Check

No diode checks can be made from the connector. All PIN diode checks must be made individually with the coil open.


Figure 15: 48-pin Hypertronics connector.

3-5 Mechanical Hardware Check

None.

3-6 Troubleshooting Tips

Symptom #1: The system reports a coil fault during prescan or does not recognize the coil connection to the system when selected in the software.

Probable Cause	Suggested Actions	Resolution
The coil connector has become disconnected from the system interface.	Check to make sure the coil connector is fully engaged.	Engage connector and try the scan again. Make sure the indicator light on the patient table turns green with the proper coil selected on the console.
One or more PIN diodes have failed.	From inside the coil, check the PIN diodes using a DVM. Use diode mode for forward-bias measurements and resistance mode for reverse-bias measurements. Note that there should be no reverse-bias leakage. Bad PIN diodes are usually shorted or exhibit reverse-bias leakage current.	Replace the coil.
There is a DC short or open somewhere within the DC path inside the coil.	Check the coil and associated DC paths with a DVM, in both diode check and resistance modes.	Replace the coil.
There is a problem with the DC bias from the system interface.	Check to ensure the correct DC bias is being provided to the coil during transmit and receive per GEMS System Service Manual. Try a similar 2-channel coil to see if the same problem occurs	Correct MRI system problem.
Coil ID issues.	From CSD, select Diagnostics, Manual Tests, SRI Functional Checks, and run the SRI Coil ID test. Compare the coil ID value found to that listed in Section 7-2 of this manual.	If the value does not match, or no coil ID is found, replace the coil.

Symptom #2: The coil does not pass SNR tests or exhibits poor image quality on patient scans.

Probable Cause	Suggested Actions	Resolution
One or more of the coil channels has a high noise level.	Check QA results from Data Sheet	Replace the coil.
One or more of the coil channels has low signal.	Look at the uncombined images to determine which channel has the problem. Compare signal mean measurements between channels and system performance logs.	Replace the coil.
Excessive ghosting is causing the noise std. dev. measurements to be artificially high.	Window and level the images down to look at the background for signs of ghosting. Try padding the phantom in the holder. Try running the scan without the phantom holder to see if the ghosting diminishes. Try turning off the cold heads to see if the ghosting is vibration-related.	If padding or changing the phantom position can minimize ghosting, then the problem is caused by excessive vibration elsewhere in the system. If the ghosting is not symptomatic of phantom positioning or padding, then replace the coil.
There is a problem with the DC bias from the system interface.	Check to ensure the correct DC bias is being provided to the coils during transmit and receive per GEMS System Service Manual. Check pins for preamp bias.	Correct MRI system DC bias problem – otherwise replace the coil.

SECTION 4 – MAINTENANCE

4-1 Coil Care

WARNING!

Detach the coil connector from the scanner before attempting to clean the coil. Do not touch the connectors with bare fingers. Never press a sharp object against the surface of the connector. Do not reattach the connector after cleaning the coil until the coil has dried completely. Electric shock may result if the coil is attached to the system during cleaning or when it is wet.

CAUTION

Do not spray or pour cleaning solution directly on the coil. Do not submerge the coil in the solution. The coil contains sensitive electronic components that could be damaged by the solution. The coil cannot be sterilized and should be cleaned only according to the procedure outlined in this section.

The following solutions are recommended for the coil and pad surfaces: a ten percent bleach solution (some discoloration may occur) or warm water. Apply cleaning solution to a soft cotton cloth and proceed to clean. To prevent soiling of the coil, the user should place a cotton sheet over the coil before positioning the patient. If the coil is soiled, clean the coil as described above.

4-2 Special Care Requirements

Prior to returning a coil for service, use a ten percent bleach solution (as described above) to eliminate risk of exposure to potentially infectious materials.

SECTION 5 – REPLACEMENT

Simple removals that are clearly obvious are not described here.

Unless otherwise noted, the steps for re-assembly are simply the reverse order of the steps described for disassembly.

SECTION 6 – RENEWAL PARTS

6-1 Field Replaceable Units

FIELD REPLACEABLE UNITS – TABLE 6-1

Part Name	GE Part #	USAI Part #
Coil	2381683-2	165149
Phantom	2381683-3	155190
Phantom Positioner	2381683-4	155191

6-2 Other Replaceable Accessories

OTHER REPLACEABLE ACCESSORIES LIST – TABLE 6-2

Part Name	GE Part #	USAI Part #
Patient Liner	2381683-5	155192
Table Pad	2381683-6	155193
Spacer Pad	2381683-7	155194
Lower Pad	-	2406572




Figure 16: Parts of the coil.

SECTION 7 – APPENDIX

7-1 Schematics

Schematic of Channel 1 and Channel 2

Schematic of Channel 3 and Channel 4

7-2 Coil Configuration

Parameter	Mode
coilName	8US TORSOPA
coilCode	mb156q7appvweb2is03talpsaw63fsn5
physicalCoilDescription	GE by USAI 8 ch torso array
coilType	3
multiCoil	yes
xmitAtten	0
fastTGstartRG	12
fastTGstartTA	90
cableLoss	1.163
coilLoss	0.89
reconScale	0
autoshimRcvr	-1
headDefaultFreqDir	0
mCoilSeperationDir	7
maxFOV	480
assetThresh	10000
nuclide	1H
scic_a	1.5 0.2 0.2 2.0 128 0.4 1.25
scic_s	1.5 0.2 0.2 2.0 128 0.4 1.25
scic_c	1.5 0.2 0.2 2.0 128 0.4 1.25
mcSwitchSel	255
mcBiasEnable	255
mcVoltConfig0	255
mcVoltConfig 1	255
mcOpenCktConfig0	255
mcOpenCktConfig1	0
mcOpenCktConfig2	255
mcOpenCktConfig3	0
mcShortCktConfig0	255
mcShortCktConfig1	0
mcShortCktConfig2	255
mcShortCktConfig3	0
mcChannelSelect	255
startRec	0
endRec	7
numRec	8
userInterfaceGroup	4
prescanWeight0	1.0
prescanWeight1	1.0
prescanWeight2	1.0
prescanWeight3	1.0
prescanWeight4	1.0
prescanWeight5	1.0
prescanWeight6	1.0
prescanWeight7	1.0
aps1_mod	0
aps1_mod_PStloc	0
aps1_mod_PSrloc	0
aps1_mod_fov	200
aps1_plane	1
opthickPS_mod	5
coilPercentExposed	100
optxpos	3
extremity	no
mcReconEnable	0
mcPortEnable	6
mcErrorEnable	6
startFastRec	4
numFastRec	0
endFastRec	4
linearQuad	1

REVISION HISTORY

Rev	Date	Author	Primary Reason for Change
1	10/03	Teresa DeMarco	First Issue
2	07/04	C. Cohara	Sec. 7, Added FRU Lower Pad
3	11Jan10	D. Jackson	Updated damage in transportation and language policy. Updated shipping contents. Updated title page and all headers and footers. Updated operator manual packet number in shipping list and operator manual number in section 1-3. Added new MCQA procedure for Unified phantom set per the Phantom Elimination Program.