SERVICE MANUAL

Invivo Corporation
3545 SW 47th Avenue
Gainesville, FL 32608

Document TR0273S

Revision 3

GE 1.5T HD TR Knee Array

GE Part Number: M3087JF

MRIDC Part Number: 800273

© Invivo Corporation, 2005

Damage in Transportation

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation “**damage in shipment**” written on **all** copies of the freight or express bill **before** delivery is accepted or “signed for” by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier **immediately** upon discovery, or in any event, within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this **14** day period.

File a report with

- Call 1–800–548–3366 and use option 8.
- Fill out a report on <http://egems.med.ge.com/edq/home.jsp>
- Contact your local service coordinator for more information on this process.

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TABLE OF CONTENTS

Section 1 – Introduction	9
1-1 Product Identification and Shipping List	9
1-2 Compatibility	9
1-3 Related Documentation	9
1-4 Environmental Requirements	10
1-5 Theory of Operation	10
Section 2 – Setup and Calibration	13
2-1 Coil Installation	13
2-1-1 Reverse Polarity Configuration	13
2-1-2 Installing the Coil	14
2-2 Quality Assurance Check	15
Section 3 – Functional Checks	16
3-1 Scanner Verification	16
3-2 Coil Imaging Performance Verification	16
3-2-1 Tools Required	16
3-2-2 Data Collection Setup with Legacy Phantom	16
3-2-3 Data Collection Setup with Unified Phantom	19
3-2-4 MCQA / SNR Test	21
3-2-5 Data Analysis 11.x and 12.x only	22
3-3 External Cable Check	26
3-3-1 RF Signal Channels Check	27
3-3-2 PIN Diode Control Channels Check	27
3-4 PIN Diodes Check	28
3-5 Mechanical Hardware Check	28
3-6 Troubleshooting Tips	28
Section 4 – Maintenance	30
4-1 Coil Care	30
4-2 Special Care Requirements	30
Section 5 – Replacement	31
5-1 Disassembly of Coil	31
5-2 External Cable Replacement	31
5-2-1 Tools Required	31
5-2-2 Replacement Time	31
5-2-3 Replacement Procedure	31
5-3 Mechanical Hardware Replacement	32
5-3-1 Baseplate Replacement	32
5-3-2 Latch Replacement	32
Section 6 – Renewal Parts	34
6-1 Field Replaceable Units	34
6-2 Other Replaceable Accessories	34
Section 7 – Appendix	35
7-1 Schematic	35
7-2 Coil Configuration	35

SECTION 1 – INTRODUCTION

1-1 Product Identification and Shipping List

This is a service manual for the GE 1.5T HD TR Knee Array.

Note: Coils do not ship with phantoms. Phantoms come in a unified phantom set with the MR System.

SHIPPING LIST – TABLE 1-1

Description	GE Part #	Qty
1.5T HD TR Knee Array	5114356-2	1
Operator Manual	5114356-5	1
Service Manual	5114356-6	1
Thin Knee Pad	5114356-7	1
Thick Knee Pad	5114356-8	1
Ramp Pad	5114356-9	1
Ankle Pad	5114356-10	1
Phantom Positioner	5114356-12	1

1-2 Compatibility

This coil is compatible with the following 1.5T hardware configurations:

Signa Excite HD, HDx, HDe, Signa HDxT, Discovery MR450, and Optima MR450w systems

1-3 Related Documentation

- 5133315 - Signa EXCITE HD 1.5T 11.x Service Methods
- 5155500 – Signa HDx 14.x Service Methods
- 5323683-2EN – Excite HDe Service Methods
- 5440000 – Signa 15.x HDxT Service Methods

- 5500110 - Discovery 20.x MR450 1.5T Service Methods
- 5670000 – Optima MR450w 1.5T Service Methods
- 1.5T HD TR Knee Array Operator's Manual, 5114356-5

1-4 Environmental Requirements

Storage

Store the coil in an air-conditioned scan room or equipment room. To store the coil, a storage space of 45.7 cm/18.0 in (width) x 72.4 cm/28.5 in (depth) x 42.7 cm/16.8 in (height) is required.

Dimensions

Coil:	45.7 cm x 72.4 cm x 42.7 cm	(18.0 in x 28.5 in x 16.8 in)
Phantom Positioner:	17.8 cm x 17.8 cm x 19.9 cm	(7.0 in x 7.0 in x 7.83 in)


Weight


Coil:	8.56 Kg	(18.9 lb)
Phantom Positioner:	1.5 Kg	(3.3 lb)

1-5 Theory of Operation

Refer to Figure 1 below and the Block Diagram on the following page.

The 1.5T HD TR Knee Array consists of 8 receive elements, and a quad birdcage transmit coil. Each receive element is comprised of a loop arranged azimuthally around the knee. Each loop is approximately 45° in azimuthal extent. The center of Loop #1 is located 45° rotated in the clockwise direction (as viewed from the head) with respect to the center of the patella. Each of the next seven loops is rotated clockwise by 45° successively. During transmit, each receive loop is decoupled with both passive and active decoupling, with an additional active preamp protection network. During receive, the transmit bird-cage is decoupled and signals from elements 1 – 8 are directed to separate ultra-low-input impedance preamp modules placed directly at the output of the loop. These preamp modules provide additional element decoupling. In addition, the 1.5T HD TR Knee Array allows receiving with the bird-cage coil to obtain a uniform signal from the imaging volume.

SECTION 2 – SETUP AND CALIBRATION

2-1 Coil Installation

2-1-1 Reverse Polarity Configuration

DURING REVERSE POLARITY CONFIGURATION, THE 1.5T HD TR KNEE ARRAY COMBINER BOARD IS EXPOSED. **USE ESD SAFETY PRECAUTIONS**

2-1-1-1 Tools Required

1ea Philips Screwdriver

2-1-1-2 Configuration Time


15 to 20 minutes

2-1-1-3 Configuration Procedure

If required, follow these steps to configure the 1.5T HD TR Knee Array for reverse polarity.


1. Remove the 1.5T HD TR Knee Array coil cable as detailed in section **5-2 External Cable Replacement**.
2. Separate the anterior and posterior coil sections.
3. Using the phillips screw driver, remove 4 screws from the anterior coil section as illustrated in Figure 2-1-1-3-1.
4. Remove the anterior coil section cover as illustrated in Figure 2-1-1-3-2.

Figure 2-1-1-3-1

Remove 4 screws


Figure 2-1-1-3-2

Remove Cover


5. Locate the two *Trap/Combiner Coaxial Cables*; they are marked 1 and 2, and are connected to the *Combiner Printed Circuit Board* at locations *Port 1* and *Port 2* respectively. Reference Figure 2-1-1-3-3.
6. Disconnect both *Trap/Combiner Coaxial Cables*. Connect cable 1 to *Port 2* and connect cable 2 to *Port 1*. Reference Figure 2-1-1-3-4 and Table 2-1-1-3.

Figure 2-1-1-3-3

Trap/Combiner Cable Connectors

Figure 2-1-1-3-4

Exchange Connectors

POLARITY CABLE CONFIGURATION - TABLE 2-1-1-3

Trap/Combiner Coaxial Cable	For Normal Polarity	For Reverse Polarity
Number 1	Connect to Port 1	Connect to Port 2
Number 2	Connect to Port 2	Connect to Port 1

7. Reinstall the anterior coil section cover using the 4 associated phillips head screws.
8. Reinstall the 1.5T HD TR Knee Array coil cable as detailed in section **5-2 External Cable Replacement**.
9. Perform the **Coil Imaging Performance Verification, Section 3-2**.

2-1-2 Installing the Coil

There are no installation steps other than plugging in the coil and performing functional checks. Coil configurations will already be on the system software. For Discovery MR450 and Optima MR450w the coil configuration is HD TRKnee PA. There are 3 configurations available for HDx, HDxt and HDe:

HD TRKnee PA
HDTRKneeService
HDTRK1chService

HD TRKnee PA is the only configuration visible in clinical mode.

HDTRKneeService and **HDTRK1chService** are visible in service mode (geservice).
HDTRKneeService and **HDTRK1chService** must be installed for the SNR tool to run.

2-2 Quality Assurance Check

QA checks should be performed during installation and during planned maintenance. Perform the quality assurance checks as outlined below to ensure the coil is operating properly.

1. Check external cable for cracks or cuts.
2. Perform Section 3 - Coil Imaging Performance Verification.

SECTION 3 – FUNCTIONAL CHECKS

3-1 Scanner Verification

Perform system level Signal to Noise Check. Refer to Service Methods CD; System Level Procedures; Functional Checks; Signal to Noise Check.

3-2 Coil Imaging Performance Verification

3-2-1 Tools Required

TOOLS REQUIRED USING LEGACY PHANTOM – TABLE 3-2-1

Description	GE Part #	MRIDC #	Qty
Phantom	5114356-11	870005	1
Phantom Positioner	5114356-12	870006	1


TOOLS REQUIRED USING UNIFIED PHANTOM – TABLE 3-2-2

Description	GE Part #	MRIDC #	Qty
Large cylindrical unified phantom	5342679	N/a	1
Phantom Positioner	5114356-12	870006	1

3-2-2 Data Collection Setup with Legacy Phantom


1. Place the 1.5T HD TR Knee Array, the phantom positioner, and the phantom on the patient cradle as shown in Figure 3-2-2-1. Note the serial number of the coil.
2. Ensure that the phantom positioner is aligned with the posterior coil section as indicated by the marks on the positioner. (Figure 3-2-2-2)

Figure 3-2-2-1

Coil and Phantom Setup


Figure 3-2-2-2

Phantom Positioner Alignment Marks


3. Being careful to align the connector pins, place the anterior coil section onto the posterior coil section and latch into place. (Figure 3-2-2-3) The scanner will not operate if the coil sections are not latched correctly.
4. Ensure that the coil is aligned and centered with respect to the baseplate as indicated by the arrows on the baseplate offset scale. (Figure 3-2-2-4)

Figure 3-2-2-3

Anterior Coil Section in Place


Figure 3-2-2-4

Align and Center the Coil

5. Establish an axial landmark using the reference mark on top of the coil. (Figure 3-2-2-5)
Advance to scan.

Figure 3-2-2-5

Establish Landmark

3-2-3 Data Collection Setup with Unified Phantom

Visual	Description
 Figure 3.2.3.1	1. Place the lower section of the coil on the patient table as shown in Figure 3.2.3.1.
 Figure 3.2.3.2	2. Align the alignment mark on the coil to Zero position on the scale as shown in Figure 3.2.3.2
 Figure 3.2.3.3	3. Place the large cylindrical unified phantom into the Positioner as shown in Figure 3.2.3.3.

Figure 3.2.3.4


Figure 3.2.3.5

4. Place this Phantom-positioner setup on the lower part of the coil as shown in Figure 3.2.3.4.
5. Ensure that the phantom positioner is aligned with the posterior coil section as indicated by the marks on the positioner and the coil (Figure 3.2.3.5)

Figure 3.2.3.6

6. Carefully align the connector pins on the lower part of the coil with the upper part of the coil and latch both the coil sections into place as in Figure 3.2.3.6. (The scanner will not operate if the coil sections are not latched correctly).


Figure 3.2.3.7

7. Connect the coil to Port A of the system LPCA and Landmark the coil at the cross mark available near the handle of the coil (Figure 3.2.3.7) and advance to scan.

3-2-4 MCQA / SNR Test

1. Move back to the operator workspace and end the previous exam, if required.
2. For 11x/12x systems, go to Common Service Desktop, select Troubleshooting, Image Quality Troubleshooting Tools, and then coil Surface Coils SNR Test. Continue with Step 3 of this section.

For 14x or later systems, go to Common Service Desktop, select Image Quality, and then select MCQA. Steps 3-6 do not apply.

If more information is needed for MCQA, it can be found on 14x/15x service CD or on web at: Troubleshooting/Image Quality/multicoil Quality Assurance. For MR450/450w the path is: Troubleshooting/ System / Multi-Coil Quality Assurance Tool

3. When the Multicoil SNR Tool Control window opens, select the **8 Channel TRKnee** Coil from the Coil Name selection menu. (Figure 3-2-4-1)

Figure 3-2-4-1

The screenshot shows a software window titled "Multicoil SNR Tool Control". At the top is a "File" menu. Below it is a "Scan Panel" with three input fields: "Coil Name" containing "8 Channel TRKnee", "Coil serial #" containing "e1277", and "Passes" containing "1". There are two buttons: "Start Scan" and "Abort Scan". To the right of the "Scan Panel" is a "View Archive Data" button. Below the "Scan Panel" is a "Results" section, which is currently empty and displays the word "EXAMPLE" in large, blue, outlined letters.

Multicoil SNR Tool Control Window

4. Enter the Coil serial number.
5. Enter the number of repetitions in the Passes field. Typically 1 pass is sufficient.
6. Click on the Start Scan button. The system will acquire the images and automatically start data analysis.


3-2-5 Data Analysis 11.x and 12.x only

For information on MCQA for systems later than 12.x, it can be found on 14x/15x service CD or on web at: Troubleshooting/Image Quality/multicoil Quality Assurance. For MR450/450w the path is: Troubleshooting/ System / Multi-Coil Quality Assurance Tool

The data analysis is automated.

1. When the scan completes, the data will be displayed in the Multicoil SNR Tool Control window. (Figure 3-2-5-1) Each measurement is assigned a pass/fail value. The coil passes the SNR test if all measurements pass.


Figure 3-2-5-1

Results in the Multicoil SNR
Tool Control Window

2. A plot window with SNR, indicating the measured value and specification limit, will open. (Figure 3-2-5-2)

Figure 3-2-5-2


Results in the SNR Plot Window

To view archived data files:

- a) Click on the View Data button in the Multicoil SNR Tool Control window.
(Figure 3-2-2-6)
- b) A new window, Multicoil SNR Display Control, will open. (Figure 3-2-5-3)


Figure 3-2-5-3

Multicoil SNR Display Control Window

- c) Select the File Tab and turn Expert Mode on.
- d) Select the **8 Channel TRKnee** coil in the coil selection menu.
- e) Click the Plot Data button.
- f) The Data Plot Control window opens. (Figure 3-2-5-4)


Figure 3-2-5-4

Data Plot Control Window

- g) Select the desired file and click Plot Data.
- h) A plot window with SNR, indicating the measured value and specification limit, will open. (Figure 3-2-5-5)


Figure 3-2-5-5

Results in the SNR Plot Window

Use Figure 3-3 for pin identification while performing the procedures in sections 3-3 and 3-4. Referring to Section 5-2, External Cable Replacement, disconnect the coil cable from the coil.

System Connector

3-3-1 RF Signal Channels Check

NOTE: The center pins of the System Connector *receive coaxial connectors* are extremely fragile [Figure 3-3, C1-C4 and D1-D4]. **Do not attempt to make any measurements at these center pins.**

With the coil cable disconnected, use a multimeter to assure an open circuit condition exists between the designated pin and GND, for each signal channel in Table 3-3-1.

EXTERNAL CABLE EXPECTED READINGS – TABLE 3-3-1

Signal Channel	Coil Cable Connector	
MC-1	A1	GND
MC-2	A15	GND
MC-3	A13	GND
MC-4	A11	GND
MC-5	A9	GND
MC-6	A7	GND
MC-7	A5	GND
MC-8	A3	GND
Transmit Signal	A20	GND
Reflected Power	A17	GND
GND	NOTE: The following points are common to GND A4, A6, A8, A10, A12, A16, A18, A22, E1 Shield, E6 Shield	

If one or more of the readings indicate a short, replace the cable. If all of the above readings indicate an open circuit condition, proceed to the **PIN Diode Control Channels Check**.

3-3-2 PIN Diode Control Channels Check

With the coil cable disconnected, use a multimeter to determine continuity between the designated pins for each control channel in Table 3-3-2. Also verify that the control pins are not shorted to GND.

EXTERNAL CABLE EXPECTED READINGS – TABLE 3-3-2

Control Channel	Coil Cable Connector	System Connector
MC bias 1	A24	F2
MC bias 2	B24	F1
MC bias 3	B23	F3
MC bias 4	A23	F4
MC bias 5	B6	F5
MC bias 6	B8	F6
MC bias 7	B10	F7
MC bias 8	B4	F8
GND	A4, A6, A8, A10, A12, A16, A18, A22	G1-G8, I6, MC Shield

If one or more of the readings are > 3 ohms, replace the coil cable. If all readings are < 3 ohms, proceed to the **PIN Diodes Check**.

3-4 PIN Diodes Check

With the coil cable connected to the coil, use a multimeter to measure the diode drop across the coil connector pins as detailed in Table 3-4. Use Figure 3-3 for pin identification while performing this procedure.

PIN DIODE EXPECTED READINGS – TABLE 3-4

System Connector		Reading (volts)
Positive Lead	Negative Lead	
F1	GND [G1]	2.2 $\pm$ 0.3
GND [G1]	F1	Open
F2	GND [G2]	2.2 $\pm$ 0.3
GND [G2]	F2	Open
F3	GND [G3]	1.1 $\pm$ 0.3
GND [G3]	F3	Open
F4	GND [G4]	1.1 $\pm$ 0.3
GND [G4]	F4	Open
F5	GND [G5]	1.1 $\pm$ 0.3
GND [G5]	F5	Open
F6	GND [G6]	2.2 $\pm$ 0.3
GND [G6]	F6	Open
F7	GND [G7]	2.2 $\pm$ 0.3
GND [G7]	F7	Open
F8	GND [G8]	1.1 $\pm$ 0.3
GND [G8]	F8	Open
E1	E1 Shield	1.5 $\pm$ 0.3
E1 Shield	E1	1.7 $\pm$ 0.3
E6	E6 Shield	Open
E6 Shield	E6	Open

If one or more of the readings indicate a short, replace the coil. If all of the above readings are correct, redo 3-2 *Coil Imaging Performance Verification*. If the problem still exists, replace the 1.5T HD TR Knee Array.

3-5 Mechanical Hardware Check

Not Applicable

3-6 Troubleshooting Tips

If the following error is displayed, there may be a problem with the Coil ID coil or system hardware: *The prescribed RF coil and detected RF coil do not match.*

Troubleshooting instructions are available through the service browser or service methods CDs:

11.x/12.x: Instructions not available in service methods. Access tool via common service browser path: **Diagnostics -> Peripherals -> SRI Functional Tests -> Coil ID**

Start the diagnostic by selecting **SRI Coil ID** and **Run**. The diagnostic will display the contents of the coil ID chip whenever the connector is inserted. Follow the trouble-shooting instructions if the correct ID is not displayed. When the trouble-shooting is complete, and if the correct ID is still not displayed, replace the cable.

HDx(14.x)/HDxT(15.x) Service Methods: **Troubleshooting -> System -> Common Service Desktop Help Screens -> SRI Coil ID – 16-channel switch**

450/450w Service Methods: **Troubleshooting -> Diagnostics -> Hardware Location -> Magnet Room -> Coil ID Functional Diagnostic**

SECTION 4 – MAINTENANCE

4-1 Coil Care

WARNING!

Detach coil connector from scanner before attempting to clean. Do not reattach after cleaning until coil has dried completely. Having the coil attached to the system during cleaning or when it is wet may result in electrical shock.

CAUTION

Do not spray or pour cleaning solution directly on coil. Do not submerge coil in solution. The coil contains sensitive electronics components that could be damaged by the solution.

The 1.5T HD TR Knee Array, baseplate and pads must be cleaned and stored using the following procedure:

- Wipe with a cloth that has been dampened in a solution of 10% bleach and 90% tap water, or 30% isopropyl alcohol and 70% tap water.
- **Do not pour any cleaning solution directly on the coil!**
- Let the coil housing and pads dry before use.
- **Under no circumstances should the coil be placed into any type of sterilizer.**

4-2 Special Care Requirements

None

SECTION 5 – REPLACEMENT

Simple removals that are clearly obvious are not described here.

Unless otherwise noted, the steps for re-assembly are simply the reverse order of the steps described for disassembly.

5-1 Disassembly of Coil

Not Applicable

5-2 External Cable Replacement

DURING EXTERNAL CABLE REPLACEMENT, THE 1.5T HD TR KNEE ARRAY COMBINER BOARD IS EXPOSED. **USE ESD SAFETY PRECAUTIONS**

5-2-1 Tools Required

- 1ea Philips Screwdriver

5-2-2 Replacement Time


- 10 to 15 minutes

5-2-3 Replacement Procedure

The 2 coil connector mounting screws are located at the inferior end of the coil.


1. Remove the 2 coil connector mounting screws as detailed in Figure 5-2-3-1.
2. Remove the coil connector by grasping the coil connector housing and pulling away from the coil in the inferior direction. [Figure 5-2-3-2]
3. Using the FRU cable, reverse the above process for installation.

Figure 5-2-3-1

Remove 2 Screws

Figure 5-2-3-2

Disconnect Cable

5-3 Mechanical Hardware Replacement

5-3-1 Baseplate Replacement

5-3-1-1 Tools Required

- 1ea Philips Screwdriver


5-3-1-2 Replacement Time

- 20 to 30 minutes

5-3-1-3 Replacement Procedure

1. Open the 1.5T HD TR Knee Array Baseplate FRU kit and find the FRU Replacement Instructions for GE 1.5T HD TR Knee Array Baseplate Kit, as illustrated in Figure 5-3-1-3.
2. Follow the FRU Replacement Instructions while performing the baseplate replacement.

Figure 5-3-1-3

FRU Replacement Instructions for
GE 1.5T HD TR Knee Array Baseplate Kit

5-3-2 Latch Replacement

5-3-2-1 Tools Required

- Locking Bail Insert Removal Tool [included in kit]


5-3-2-2 Replacement Time

- 20 to 30 minutes

5-3-2-3 Replacement Procedure

3. Open the 1.5T HD TR Knee Array Latch FRU kit and find the FRU Replacement Instructions for GE 1.5T HD TR Knee Array Latch Kit, as illustrated in Figure 5-3-2-3.
4. Follow the FRU Replacement Instructions while performing the latch replacement.

Figure 5-3-2-3

**FRU Replacement Instructions for
GE 1.5T HD TR Knee Array Latch Kit**

SECTION 6 – RENEWAL PARTS

6-1 Field Replaceable Units

FIELD REPLACEABLE UNITS LIST – TABLE 6-1

Description	GE Part #	MRIDC #
1.5T HD TR Knee Array	5114356-2	870001
1.5T HD TR Knee Array Cable	5114356-3	870002
1.5T HD TR Knee Array Baseplate	5114356-4	870003
Latch Kit	5114356-13	870004
Legacy Phantom	5114356-11	870005
Phantom Positioner	5114356-12	870006
Large cylindrical unified phantom	5342679	N/A

6-2 Other Replaceable Accessories

The following table lists replaceable accessories that may be purchased from MRI Devices Corporation.

OTHER REPLACEABLE ACCESSORIES LIST – TABLE 6-2

Description	MRIDC #
Thin Knee Pad	870007
Thick Knee Pad	870008
Ramp Pad	870009
Ankle Pad	870010

SECTION 7 – APPENDIX

7-1 Schematic

Refer to Block Diagram under Section 1-5 Theory of Operation.

7-2 Coil Configuration

For Discovery MR450 and Optima MR450w the coil configuration is HDTRKneePA

There are 3 configurations available for HDx, HDxt, and HDe:

- HD TRKnee PA
- HDTRKneeService
- HDTRK1chService

HD TRKnee PA is the only configuration visible in clinical mode.

HDTRKneeService and **HDTRK1chService** are visible in service mode (geservice).

HDTRKneeService and **HDTRK1chService** must be installed for the SNR tool to run.

Check the system coil configuration file for these configurations. If they have not been installed, select them in the appropriate mode, from the list of available coils, and install.

Revision History

Rev	Date	Author	Primary Reason for Change
A	08-04-04	L. Hyler	Initial Release
B	11-10-04	L. Hyler	Changes per Formal Review Log, Kyle Zeller, 11.0/04
C	12-07-04	L. Hyler	QA Section Changes from Julia Klinge
D	12-21-04	L. Hyler	Add External Cable Checks from Asdrubal Mata
E	01-07-05	L. Hyler	Change Name to "1.5T HD Knee Array" Throughout Manual
F	01-13-05	L. Hyler	Implement Changes from J. Klinge e-mail of 01-13-05
01	01-21-05	L. Hyler	Released
02	6-5-08	L. Hyler	Customer Change Request CCN0021 (DCN 3994)
03	10/29/09	M. Vaughn	Updated per Change Request CCN0036 (Unified phantom set, 450/450w addition)