

SIEMENS

MAGNETOM

MR

Software Information for Service

MR A30

**Concerto,
Harmony (Turbo, Ultra, Sprint, Quantum),
Symphony (Turbo, Ultra, Sprint, Quantum),
Harmony Impact/Expert, Symphony Vision,
Sonata, Sonata Vision, Trio, C!, Symphony P,
Novus 1.5 and Allegra**

© Siemens AG 2007

The reproduction, transmission or use of this document or its contents is not permitted without express written authority. Offenders will be liable for damages. All rights, including rights created by patent grant or registration of a utility model or design, are reserved.

Print No.: MR-000.816.27.01.02

Replaces: n.a.

English

Doc. Gen. Date: 03/07

Document revision level

This document corresponds to the version/revision level effective at the time of system delivery. Revisions to hardcopy documentation are not distributed automatically.

Please contact your local Siemens office to order current revision levels.

Disclaimer

The installation and service of equipment described herein is to be performed by qualified personnel who are employed by Siemens or one of its affiliates or who are otherwise authorized by Siemens or one of its affiliates to provide such services.

Assemblers and other persons who are not employed by or otherwise directly affiliated with or authorized by Siemens or one of its affiliates are directed to contact one of the local offices of Siemens or one of its affiliates before attempting installation or service procedures.

Author: CSSDMR Marco Manchinu

	Page
1 Release Summary	1-1
Release summary	1-3
Release Data	1-3
Release enclosure.	1-4
Software label and delivery media.	1-4
Hardware/Software operating system	1-5
Functionality	1-5
Released Cameras	1-7
Test phases	1-8
Service Features for <i>syngo</i> MR A30.	1-8
2 Hints for Service	2-1
Installing software	2-1
Service.	2-1
Imager	2-1
Configuration SRS / DICOM.	2-1
Remote diagnostics	2-1
DICOM.	2-1
Known bugs for service	2-2
3 Changes to previous version	3-1

Page

MAGNETOM

Software Release

***syngo* MR A30
NUMARIS/4 V30A
for the systems:**

**Concerto,
Harmony (Turbo, Ultra, Sprint, Quantum),
Symphony (Turbo, Ultra, Sprint, Quantum),
Harmony Impact/Expert, Symphony Vision,
Sonata, Sonata Vision, Trio, C!,
Symphony P, Novus 1.5 and Allegra**

March 30, 2007

This page intentionally left blank.

Release summary

Release Data

Date of release: March 30, 2007

Products: MAGNETOM systems:
Concerto,
Harmony (Turbo, Ultra, Sprint, Quantum),
Symphony (Turbo, Ultra, Sprint, Quantum),
Harmony Impact/Expert, Symphony Vision,
Sonata, Sonata Vision, Trio, C!, Symphony P,
Novus 1.5 and Allegra

Software version *syngo* MR A30
which consists of
NUMARIS/4 VA30A
including
syngo VE31G
ARGUS VA70B
3D Vesselview VA30B
3D Neuro VA10F
Composer VA71B

The *syngo* MR A30 software has been developed according to valid software engineering and quality regulations. The software has undergone extensive testing through the different phases of internal and external tests.

The contents of the released software are described in chapter 1 as well as in chapter 2 of "Hints for Service."

Release enclosure

Software label and delivery media

Software version: *syngo* MR A30
NUMARIS/4 VA30A
Software baseline: N4_VA30A_LATEST_20070315

The software baseline can be viewed with the function "Help, Info...", the complete naming is:

N4_VA30A_LATEST_20070315
syngo VE31G
N4_VA30A_LATEST_20070315_P1 (only TRIO)

Delivery media:

syngo MR A30 (2 DVDs)
NUMARIS/4 VA30A
MAGNETOM Harmony / Symphony / Sonata / Concerto /
Trio / Novus / C! / Allegra
Software **R640**
with part number 100 92534 K2219

syngo MR A30 (2 DVDs)
NUMARIS/4 VA30A
MAGNETOM Harmony / Symphony / Sonata / Concerto /
Trio / Novus / C! / Allegra
Software **R630**
with part number 102 76973 K2219

syngo MR A30 (1 DVDs)
NUMARIS/4 VA30A
MAGNETOM Trio
VA30A_SP1
with part number 102 77194 K2219

syngo MR A30 (1 DVDs)
NUMARIS/4 VA30A
MAGNETOM Harmony / Symphony / Sonata / Concerto /
Trio / Novus / C! / Allegra
IDEA
with part number 100 92535 K2219

Hardware/Software operating system

Host (MRC or MRSC) / Imager

Operating systems supported by *syngo* MR A30 are:

R630 for host (3 Gbyte) and R630 for imager (1Gbyte)

R640 for host (3 Gbyte) and M450 for imager (1Gbyte)

MPCU:

- MPCU-2 / MPCU-3
- Operating system VxWorks 5.31 (on DVD 1/2)

Functionality

The following systems of the MAGNETOM product line are supported by *syngo* MR A30 NUMARIS/4 VA30A:

MAGNETOM Systems

Concerto,
Harmony (Turbo, Ultra, Sprint, Quantum),
Symphony (Turbo, Ultra, Sprint, Quantum),
Harmony Impact/Expert, Symphony Vision,
Sonata, Sonata Vision, Trio, C!, Symphony P,
Novus 1.5 and Allegra

The main goals for the NUMARIS/4 VA30A software are as follows:

- Future security for Symphony, and C!
Assure customers that they invested in future proof and still well tended system with Syphony and C!
Achieved via features:
 - BLADE for Symphony
 - SWI for Symphony, TRIO
 - DTI Tensor only (12 directions) for Symphony, TRIO
 - VIEWS Breast Protocols for Symphony
 - C! - Applications
 - Single Shot EPI
 - DWI Neuro-Perfusion
 - 3D SPACE
 - DESS
 - Workflow for all systems
 - Same *syngo* version as MR B13 - DVD supportAll new features of Symphony should be available for SymphonyP, but not for NOVUS15.
- Grow market share with C!
Increase the perception of continual improvement and enhancement of the product
- Cover Evolve demands
Use VA30A for software update on non Tim systems to provide Evolve version to customers.

- N4VA30A SW Maintenance (CHARMs)
Please note that problems found in previous versions VA25 and VA26 will be corrected
- Expert-i option. The Expert - i SW allows remote control of *syngo*-speaking MR systems from a standard Windows PC connected to the clinical network.

Details regarding software installation can be found in the "Update instruction" for this software.

Released Cameras

Camera
Agfa Drystar 3000
AGFA Drystar 5500
Agfa LR 5200/MG 3000
AGFA LR5200 / LR DICOM Controller
CODONICS Horizon
FUJI Drypix 3000
FUJI Drypix 7000
FUJI FM DPL
Kodak Dryview 8100
Kodak Dryview 8200
Kodak Dryview 8700 / MIM 200
Kodak Dryview 8900
Kodak Dryview KELI 160
Kodak KELP 1120
Kodak KELP 2180
Kodak MLP 190
Konica DRYPRO 722 + PrintLink
Konica DRYPRO 751/752

Test phases

syngo MR A30 has completed the following test phases:

- Unit test (internally performed at SIEMENS by SW engineering teams)
- Integration test and validation of the *syngo* software
(performed internally at SIEMENS by the SW engineering teams)
- Internal system test
(internally performed at SIEMENS by MR test teams and application specialist teams)
- External system test (Clinical workflow at selected customer sites)

Service Features for *syngo* MR A30

Eventlog:

- different layout with more filter possibilities

Installing software

For a detailed description of how to install the System Software, refer to the Software Installation for NUMARIS/4 VA30A.

Service

For the following service activities: (Configuration, Backup & Restore, TuneUp, Shim, Magnet & Cooling, Install & Maintenance, Quality Assurance, TestTools, DICOM Tests, File & Image Tools, Auto Report, and Utilities) the system provides information, see "Service UI" -> Help -> Online Help.

Host

Make sure to follow the instructions concerning important milestones during the software installation of the MRC / MRSC:

- No user actions are allowed during the automatic installation. Wait until the login window is displayed.

Imager

- The imager loads its license at each system start-up via tftp from the Host (MRC).
When replacing the imager a new Imager license is required.

Configuration SRS / DICOM

Remote diagnostics

Only IE 6.0 has been tested and released as the service browser for remote service on *syngo* MR systems running NUMARIS VA30A.

The "SIEMENS Remote Services TM" will only support router connections for NEWLY delivered systems. See <http://cs.med.siemens.de/cms/en/> "For Service", "SRS Toolkit" for the most recent information regarding the router configuration.

The project manager is responsible for planning a router connection to the Siemens Remote Diagnostic server. In case of questions concerning your pre-configured router from Siemens, please contact mailto:srs_final.med@siemens.com or your remote diagnostics person responsible.

Factory pre-configured Auto reports are available for all newly delivered NUMARIS/4 VA30A *syngo* MR systems. You only need to manually activate the pre-configured auto reports, if your customer agrees to remote service.

DICOM

The "CCC Team" provides a powerful configuration management tool.

Please refer to <http://cs.med.siemens.de/cms/en/> "For Service", "DICOM", in the DICOM Homepage (CCC). Navigate to "Interoperability DB", "Download", "Acquisition Tool".

Known bugs for service

Function	TestTools: Imager Functional Test does not run		
Problem	TestTools: Imager Functional Test 1,2,3,4 does not run		
Workaround	Restart the imager processes via the system control in the main UI. Run the Functional Test again.		
(n.a.)	Version	<i>syngo</i> MR A30	Solution (planned) n.a.

Function	Service UI / Testtools / wrong IOP Status		
Problem	Testtools: IOP status is sporadicly wrongly displayed.		
Workaround	n.a.		
(MR_00356276)	Version	<i>syngo</i> MR A30	Solution (planned) n.a.

Initial version

This page intentionally left blank.