



## PM Report Schedule A

### LightSpeed VCT and Discovery CT750 HD

Customer Name: \_\_\_\_\_ Service Record #: \_\_\_\_\_

System ID/ Serial Number: \_\_\_\_\_ Customer Equipment ID: \_\_\_\_\_

Service Manual Document Number: \_\_\_\_\_ Software Version: \_\_\_\_\_

PM Frequency: \_\_\_\_\_ months

System Type (Select:)

#### PRELIMINARY TASKS

##### Last PM Inspection

Verify last PM Inspection

Last Inspection Date :

Last Schedule Completed:

List any Schedule A tasks that were completed earlier

PM Schedule Sys A

#### Review Error Logs

Check the tasks you completed and list any follow-up work scheduled.

| Task                                     | Completed                | Comments/Follow-up |
|------------------------------------------|--------------------------|--------------------|
| Table Hydraulic Errors                   | <input type="checkbox"/> |                    |
| System Application Stop and Start Errors | <input type="checkbox"/> |                    |
| Linux Reboot and Shutdown Errors         | <input type="checkbox"/> |                    |
| Record Gantry Rotations                  | 0                        |                    |

#### Initial Procedures

##### Check Temperature and Humidity

Record the temperature and humidity. Indicate if these are within spec.

| Room         | Temperature | Humidity | Pass/Fail | Comments/Follow-up |
|--------------|-------------|----------|-----------|--------------------|
| Scan room    |             |          |           |                    |
| Control Room |             |          |           |                    |

#### Console

| Task                                  | Completed                | Comments/Follow-up |
|---------------------------------------|--------------------------|--------------------|
| Console Monitor Cleaning & Inspection | <input type="checkbox"/> |                    |



### General Gantry Inspection

| Task                      | Pass/Fail | Comments/Follow-up |
|---------------------------|-----------|--------------------|
| Listen for Unusual Noises |           |                    |
| Test Display Lights       |           |                    |
| Test E-stops              |           |                    |
| Test Patient Tilt Sensors |           |                    |

### Laser Alignment Inspection

| Task                                      | Pass/Fail | Comments/Follow-up |
|-------------------------------------------|-----------|--------------------|
| Internal Alignment Lights                 |           |                    |
| External Axial to Internal Axial Distance |           |                    |

### System Options

| Task                              | Pass/Fail | Comments/Follow-up |
|-----------------------------------|-----------|--------------------|
| Ceiling Component Bolt Inspection |           |                    |
| Nemoto Injector Inspection        |           |                    |
| Cardiac Option Inspection         |           |                    |

### Gantry Power Off (Covers Off)

#### General Tasks

| Task                           | Pass/Fail | Comments/Follow-up |
|--------------------------------|-----------|--------------------|
| Emergency Off Switch           |           |                    |
| UPS Power Off Inspection       |           |                    |
| Gantry Cover Cables Inspection |           |                    |

### General Console Cleaning & Inspection

| Task                                       | Pass/Fail | Comments/Follow-up |
|--------------------------------------------|-----------|--------------------|
| Clean Component Filters, Fans & Grills     |           |                    |
| Clean, Inspect, Replace Desktop Components |           |                    |
| Visual Fan Inspection                      |           |                    |
| Seismic Anchor Check (if necessary)        |           |                    |

**Slip Ring & Brush Block Inspection and Maintenance**

| Task                          | Pass/Fail | Comments/Follow-up |
|-------------------------------|-----------|--------------------|
| Remove Slip Ring Brush Debris |           |                    |
| Slip Ring Tracks Inspected    |           |                    |
| Clean Brush Blocks            |           |                    |
| Brush Tip Inspection          |           |                    |

**Gantry Filter Cleaning**

| Task                                      | Pass/Fail | Comments/Follow-up |
|-------------------------------------------|-----------|--------------------|
| Clean DAS Detector Plenum Filter          |           |                    |
| Clean Detector Face Plate                 |           |                    |
| Clean Gantry Heater Filter                |           |                    |
| Clean Top Cover Fans                      |           |                    |
| Inspect Rotating Path Components & Cables |           |                    |

**Gantry Power On (Covers Off)****Tube Heat Exchanger & Pump**

| Task                                                             | Select | Comments/Follow-up |
|------------------------------------------------------------------|--------|--------------------|
| Inspect/Clean Tube Heat Exchanger                                |        |                    |
| Inspect/Clean JEDI Inverter Fan                                  |        |                    |
| Tube Heat Exchanger Pump Inspection<br>Mfg date/last replacement |        |                    |

**Gantry Thermistor (Does not apply to Discovery CT750 HD)**

|                              |           |  |      |                    |  |      |
|------------------------------|-----------|--|------|--------------------|--|------|
| Resistance Readings:         | COLD      |  | ohms | HOT                |  | ohms |
|                              | Pass/Fail |  |      | Comments/Follow-up |  |      |
| Gantry Thermistor Inspection |           |  |      |                    |  |      |

**General Tasks**

| Task                                                       | Select | Comments/Follow-up |
|------------------------------------------------------------|--------|--------------------|
| UPS Powered and Operational                                |        |                    |
| Plenum Fans (5) Operational                                |        |                    |
| Grease Gantry Main Bearing<br>No. Revs since last greasing |        |                    |



## Gantry Power On (Covers On)

### General Tasks

| Task                               | Select | Comments/Follow-up |
|------------------------------------|--------|--------------------|
| Scan Window Inspection/Replacement |        |                    |
| Inspect Rotational Clearances      |        |                    |

## Finalization

### System Scanning Test

| Task                          | QA Phantom Serial # | Location  | Gantry Serial # |
|-------------------------------|---------------------|-----------|-----------------|
|                               |                     |           |                 |
| Scans                         | Scan Failures       | Artifacts | Artifact Types  |
| Series 1 (Scout/Auto Voice)   |                     |           |                 |
| Series 2 (Axial/Tilt)         |                     |           |                 |
| Series 3 (Helical/Auto Voice) |                     |           |                 |

### General Tasks

|                                       |                          |  |
|---------------------------------------|--------------------------|--|
| X-ray On Indicators                   |                          |  |
| Scan Control Pushbuttons Test         |                          |  |
|                                       | <b>Completed</b>         |  |
| Site Log Updated                      | <input type="checkbox"/> |  |
| PM Paper work completed for this site | <input type="checkbox"/> |  |

Is follow-up work required?

Service Scheduled Date:

☐ This inspection passed and meets all GE Healthcare PM specifications

Comments:

Characters Remaining:

Tools Used:

| Description | Serial Number | Bar Code | Cal Due Date |                                            |
|-------------|---------------|----------|--------------|--------------------------------------------|
|             |               |          |              | <input type="button" value="Add Tool"/>    |
|             |               |          |              | <input type="button" value="Remove Tool"/> |



**GE Representative:**

FSEs completing this PM

FSE Name: \_\_\_\_\_

SSO #: \_\_\_\_\_

FSE SR #: \_\_\_\_\_

FSE Signature Data...

**Customer Signature:** ☐ Yes ☐ No

☐ Final PM Complete

Date  
Completed: \_\_\_\_\_