



System ID: _____

SR#: _____

Date: _____

PM Report Schedule B

LightSpeed VCT, Discovery CT750 HD, Revolution Discovery CT, Revolution GSI, Revolution HD, and Discovery CT

Customer Name: _____ Service Record #: _____

System ID/ Serial Number: _____ Customer Equipment ID: _____

Service Manual Document Number: _____ Software Version: _____

PM Frequency: _____ months

System Type (Select:)

PRELIMINARY TASKS

Last PM Inspection

Verify last PM Inspection

Last Inspection Date :

Last Schedule Completed:

List any Schedule B tasks that were completed earlier

PM Schedule Sys B

Review Error Logs

Task	Completed	Comments/Follow-up
System Application Stop and Start Errors	☐	
Table Hydraulic Errors	☐	

Initial Procedures

Check Temperature and Humidity

Room	Temperature	Humidity	Pass/Fail	Comments/Follow-up
Scan room Temp & Humidity				

General Gantry Inspection

Task	Pass/Fail	Comments/Follow-up
Test Display Lights		
Test E-stops		
Test Patient Tilt Sensors		



System ID:

SR#:

Date:

Laser Alignment Inspection

Task	Pass/Fail	Comments/Follow-up
Internal Alignment Lights		
External Axial to Internal Axial Distance		

System Options

Task	Pass/Fail	Comments/Follow-up
Nemoto Injector Inspection		

Collimator Inspection

Task	Pass/Fail	Comments/Follow-up
Aperture and Filter Positions and Operation		
Protective Window & Opening, Inspection		

Gantry Power Off (Covers Off)

General Tasks

Task	Pass/Fail	Comments/Follow-up
Gantry Tilt Zero Check		
Emergency Off Switch (on wall)		
Gantry Cover Cables Inspection		
UPS Power Off Inspection		

General Console Cleaning & Inspection

Task	Pass/Fail	Comments/Follow-up
Clean Component Filters, Fans & Grills		
Visual Fan Inspection		
Seismic Anchor Check (if necessary)		

Slip Ring & Brush Block Inspection and Maintenance

Task	Pass/Fail	Comments/Follow-up
Remove Slip Ring Brush Debris		
Slip Ring Tracks Inspected		
Clean Brush Blocks		
Brush Tip Inspection		



System ID:

SR#:

Date:

Gantry Filter Cleaning

Task	Pass/Fail	Comments/Follow-up
Clean DAS Detector Plenum Filter		
Clean Detector Face Plate		
Clean Gantry Heater Filter		
Clean Top Cover Fans		
Inspect Rotating Path Components & Cables		

Hydraulic Tilt Pump (Inspect for Leaks)

Task	Pass/Fail	Comments/Follow-up
Hydraulic Tilt Cylinder & Fittings		
Hydraulic Tilt Pump & Reservoir		

NGPDU

Task	Pass/Fail	Comments/Follow-up
Inspect Electronic Components		
Inspect Seismic Anchor Bolts		
Torque all Connectors		

Gantry Power On (Covers Off)

General Tasks

Task	Select	Comments/Follow-up
Grease Gantry Main Bearing Record Gantry Rotations		
Plenum Fans (5) Operational		
UPS Powered and Operational		

Tube Heat Exchanger & Pump

Task	Select	Comments/Follow-up
Inspect/Clean Tube Heat Exchanger		
Inspect/Clean JEDI Inverter Fan		
Tube Heat Exchanger Pump Inspection		



System ID:

SR#:

Date:

Power On (Covers On)

General Tasks

Task	Select	Comments/Follow-up
Inspect Rotational Clearances		
Scan Window Inspection/Replacement		

Finalization

System Scanning Test

Task	QA Phantom Serial #	Location	Gantry Serial #
Scan	Scan Failures	Artifacts	Artifact Types
Series 1 (Scout/Auto Voice)			
Series 2 (Axial/Tilt)			
Series 3 (Helical/Auto Voice)			

General Tasks

Task	Select	Comments/Follow-up
X-ray On Indicators		
Scan Control Pushbuttons Test		
	Completed	
Update Site Log	☐	
Complete PM Paper Work	☐	
Complete Customer Reports	☐	

Is follow-up work required?

Service Scheduled Date:

☐ This inspection passed and meets all GE Healthcare PM specifications

Comments:

Characters Remaining:



System ID:

SR#:

Date:

Tools Used:

Description	Serial Number	Bar Code	Cal Due Date	
				Add Tool Remove Tool

GE Representative:

FSEs completing this PM

FSE Name: _____

SSO #: _____

FSE SR #: _____

FSE Signature Data...

Customer Signature: ☐ Yes ☐ No

☐ Final PM Complete

Date
Completed: _____