



System ID: _____

SR#: _____

Date: _____

PM Report Schedule C

LightSpeed VCT, Discovery CT750 HD, Revolution Discovery CT, Revolution GSI, Revolution HD, and Discovery CT

Customer Name: _____ Service Record #: _____

System ID/ Serial Number: _____ Customer Equipment ID: _____

Service Manual Document Number: _____ Software Version: _____

PM Frequency: _____ months

System Type (Select:)

PRELIMINARY TASKS

Last PM Inspection

Verify last PM Inspection

Last Inspection Date :

Last Schedule Completed:

List any Schedule C tasks that were completed earlier

PM Schedule Sys C

Review Error Logs

Check the tasks you completed and list any follow-up work scheduled.

Task	Completed	Comments/Follow-up
System Application Stop and Start Errors	☐	

Initial Procedures

General Gantry Inspection.

Task	Pass/Fail	Comments/Follow-up
Test Display Lights		
Test E-stops		
Test Patient Tilt Sensors		

System Options

Task	Pass/Fail	Comments/Follow-up
Cardiac Option Inspection		



System ID:

SR#:

Date:

Gantry Power Off (Covers off)

General Tasks

Task	Pass/Fail	Comments/Follow-up
Emergency Off Switch (on wall)		
Gantry Cover Cables Inspection		
UPS Power Off Inspection		

General Console Cleaning & Inspection

Task	Pass/Fail	Comments/Follow-up
Clean Component Filters, Fans & Grills		
Visual Fan Inspection		
Seismic Anchor Check (if necessary)		

Slip Ring & Brush Block Inspection and Maintenance

Task	Pass/Fail	Comments/Follow-up
Remove Slip Ring Brush Debris		
Slip Ring Tracks Inspected		
Clean Brush Blocks		
Brush Tip Inspection		

CT Belt Tension & Adjustment

Task	Pass/Fail	Comments/Follow-up
Inspect and Adjust Tension Belt		

Gantry Filter Cleaning

Task	Pass/Fail	Comments/Follow-up
Clean DAS Detector Plenum Filter		
Clean Detector Face Plate		
Clean Gantry Heater Filter		
Clean Top Cover Fans		
Inspect Rotating Path Components & Cables		

Gantry Power On (Covers Off)

General Tasks



System ID:

SR#:

Date:

Task	Select	Comments/Follow-up
Grease Gantry Main Bearing		
Record Gantry Rotations		
Plenum Fans (5) Operational		
UPS Powered and Operational		

Tube Heat Exchanger & Pump

Task	Select	Comments/Follow-up
Inspect/Clean Tube Heat Exchanger		
Inspect/Clean JEDI Inverter Fan		
Tube Heat Exchanger Pump Inspection		

Gantry Verifications (If required)

Task	Select	Comments/Follow-up
HV Tank Resistor Verification		
Meter Verification		
HHS Scans		
Task (CT750 HD Only)	Select	Comments/Follow-up
GSI Meter Verification		

Power On (Covers On)

General Tasks

Task	Select	Comments/Follow-up
Inspect Rotational Clearances		
Scan Window Inspection/Replacement		
Save State Completed	☐	

General Table Cleaning & Inspection

Table and Tilt Functional Checks: Table Height Min Max Gantry Tilt S30 and Gantry Tilt I30		
--	--	--

Finalization

System Scanning Test



System ID:

SR#:

Date:

Scans	Scan Failures	Artifacts	Artifact Types
Task	QA Phantom Serial #	Location	Gantry Serial #
Scans	Scan Failures	Artifacts	Artifact Types
Series 1 (Scout/Auto Voice)			
Series 2 (Axial/Tilt)			
Series 3 (Helical/Auto Voice)			

General Tasks

Task	Select	Comments/Follow-up
X-ray On Indicators		
Scan Control Pushbuttons Test		
Quality Assurance Test		
GE Phantom Serial #		
	Completed	
Complete PM Paperwork	☐	
Update Site Log	☐	
Complete Customer Reports	☐	



System ID:

SR#:

Date:

Appendix A: Data Forms

Meter Verification

	DMM Reading	A/D	Delta	Limit
mA Low Scale				± 3.0 mA
mA High Scale				± 3.0 mA
kV				Not Applicable

HV Tank Feedback Resistor Verification

Selected kV	On Screen kV	Measured kV
80		
100		
120		
140		
Scan Timer Value		



System ID:

SR#:

Date:

HHS Data - Tube, HV Tank or Collimator Replacement

CT System

Small Focal Spot		On Screen Data			
kV	mA	kV	mA	5ms mA	Time
80					
100					
120					
140					

Large Focal Spot		On Screen Data			
kV	mA	kV	mA	5ms mA	Time
80					
100					
120					
140					

GSI Meter Verification - Complete for CT750 HD Only

Selected Protocol	Measured kV	Measured mA
GSI Protocol 1		
GSI Protocol 2		
GSI Protocol 3		

Is follow-up work required?

Service Scheduled Date:

☐ This inspection passed and meets all GE Healthcare PM specifications



System ID:

SR#:

Date:

Comments:

Characters Remaining:

Tools Used:

Description	Serial Number	Bar Code	Cal Due Date	
				Add Tool Remove Tool

GE Representative:

FSEs completing this PM

1

FSE Name:

SSO #:

FSE SR #:

FSE Signature Data...

Customer Signature: ☐ Yes ☐ No

☐ Final PM Complete

Date

Completed: