

Services Knowledge Base

ID: SKB0098949

Important note: The originating document might have changed - please check the online document prior to taking any action <https://skb.healthcare.siemens.com/Pages/SKBAdvancedViewKBDocument.aspx?docid=SKB0098949>

Exchange of MSUP (part-no. 10523587 or 10892850) after a quench

Document Type: Information

Products: Aera, Skyra, Verio, Avanto, Espree, Aera24, Amira, Avanto fit, Avanto Dot, Verio Dot, Prisma, Prisma fit, Sempira, Skyra fit, Sola, Spectra, Terra, Vida

Software Version:

Description

Every now and then returned MSUPs are received with the fault diagnosis “exchange after quench”. The manufacturing test typically results in a “no fault found”.

We recommend not to swap the MSUP, the alarm box, the ERDU battery or the DAFID after a quench.

Typically these components are not responsible for a quench.

To actively quench a magnet (due to energy supplied by via the MSUP) you have to press one of the magnet stop buttons on the alarm box or the one in the RF cabin next to the door. Then the quench and switch heaters in the magnet are directly supplied with energy, either from mains power (which is normally available) or in very few cases from the ERDU battery.

On the above mentioned MSUP types there is no active element involved in the quench circuit, so there is nothing that can go defective.

The ERDU battery (part-no. 10128924/10128925) is a passive device which is a backup for the case that the mains power is off. The battery cannot cause a magnet quench without activating one of the magnet stop buttons. As long as there is mains power, this source will be used to trigger the quench, not the battery.

Please only exchange the ERDU battery if Magnet & Cooling Magnet Status shows a battery voltage which is below the warning limit.

Root Cause

Resolution

Therefore only change the MSUP, ERDU battery, Alarm Box or DAFID under instruction from the HSC.

When quenches are reported to the HSC they review the logs and other information obtained to ensure only defective parts are changed.

Information For Customer

Applies to

All systems with MSUP K2306 pn 10523587 or 10892850 installed.

Document links

Copyright

"© Siemens Healthcare GmbH, 2016" refers to the copyright of a Siemens entity such as Siemens Aktiengesellschaft - Germany, Siemens Shenzhen Magnetic Resonance Ltd. - China, Siemens Shanghai Medical Equipment Ltd. - China, Siemens Medical Solutions USA Inc. - USA, Siemens Healthcare Diagnostics Inc. - USA and/or Siemens Healthcare Diagnostics Products GmbH - Germany.

Disclaimer

October 2015

This database contains copyrighted Siemens Healthcare information, trade secrets and other confidential information. The Services Knowledge Base is intended as a repository for frequently asked support questions, application fixes, workarounds, other issues, and for no other purpose. Services Knowledge Base solutions can be used as assistance for error diagnostics and error solving. Each activity/action by the user (e.g. Service Engineers, Administrators, Clinical Educational Specialists) as part of maintenance/repair must to be verified against Siemens' technical documentation.

All information contained herein is proprietary to Siemens and may not be used for any purpose other than as identified above, nor may it be disclosed to any unauthorized user without Siemens' advance written consent. Please note that while every reasonable attempt is made to assure the correctness of information provided herein, Siemens cannot guarantee that any particular workaround or solution will solve a facility's particular issue. Since each Siemens customer has different systems installed and working together, no one solution can be offered for each problem.

Attention: Copies are uncontrolled documents!