

Service Request Number _____

Site Name: _____

Planned Maintenance Checklist

3200 NW 27th Ave Ste. 100
Pompano Beach, FL

FE Name: _____

Month: _____

33069
Office: 888-400-3979
aisserv.com

Procedure

OK
Cleaned
Adjusted
Needs
Correction

PM Checklist

Talk With Technologist About System Operation..... ☐ ☐ ☐ ☐

Look at Image Quality of Previous Exams..... ☐ ☐ ☐ ☐

Check and Set System Date and Time..... ☐ ☐ ☐ ☐

SAFETY AND REGULATORY

SYSTEM - Check Gating Cable/ Respiratory Bellows / PG Probe..... ☐ ☐ ☐ ☐

SYSTEM - Check Pneumatic Patient Alert System..... ☐ ☐ ☐ ☐

SYSTEM - Check Emergency Release of Cradle and Patient Transport..... ☐ ☐ ☐ ☐

SYSTEM - Check Patient Transport Casters and Arm board Set Screws..... ☐ ☐ ☐ ☐

MAGNET - Check Magnet Bore for Metal Objects and Remove..... ☐ ☐ ☐ ☐

MAGNET - Test Magnet Emergency Rundown (Test Battery, Htr LED)..... ☐ ☐ ☐ ☐

MAGNET - Check for Condensation or Leaks on Top and Around Magnet..... ☐ ☐ ☐ ☐

MAGNET - Record Cryostat Pressure and Flow Rates..... ☐ ☐ ☐ ☐

He Level 1 _____ Pressure _____ Flow One _____ Flow Two _____ (HFO ONLY) He Level 2 _____

MAGNET - Compressor Run-Time / Pressure (Dynamic - Static)..... ☐ ☐ ☐ ☐

Compressor One -Hrs _____ DY _____ -ST _____ = _____ Compressor Two - Hrs _____ DY _____ -ST _____ = _____

MAGNET - Coldhead and Magnet Temps..... ☐ ☐ ☐ ☐

Coldhead_Ru0 _____ Recon_Si410 _____ Recon_Ru0 _____ Shield_Si410 _____ Water Temp _____ Water flow _____

Clean, Test and Check

Check Correlated Noise (RF Integrity / RF Door Seal, top hatch, clean, Check)..... ☐ ☐ ☐ ☐

Clean Light-Weight Cradle Wheels..... ☐ ☐ ☐ ☐

Check Table Functionality (Clutch, Alignment, Undock)..... ☐ ☐ ☐ ☐

Check Body Chiller Temp and Fluid..... ☐ ☐ ☐ ☐

Water Temp and flow _____

Clean Patient Blower Filter..... ☐ ☐ ☐ ☐

Clean All Cabinet Filters..... ☐ ☐ ☐ ☐

Clean Console area (Vacuum PC, Octane, or GOC)..... ☐ ☐ ☐ ☐

Semi-Annual or First PM Safety Checks

Month: _____

Test Emergency off Buttons..... ☐ ☐ ☐ ☐

Test Emergency Stop Buttons..... ☐ ☐ ☐ ☐

Inspect all System Cables and Connections for Tightness and Security..... ☐ ☐ ☐ ☐

Check Condition and Placement of Warning Labels..... ☐ ☐ ☐ ☐

Check Condition and Placement of 0.5 mT Exclusion Warnings..... ☐ ☐ ☐ ☐

Check Magnet Exhaust Vent and Vent Glass..... ☐ ☐ ☐ ☐

Clean Water Chiller Filter - If filter is really dirty then need to clean on next PM and note status below..... ☐ ☐ ☐ ☐

Tip Seal and Solenoid Valve Replacement (Twinspeed)..... ☐ ☐ ☐ ☐

Change Out All Patient Cradle Wheels..... ☐ ☐ ☐ ☐

Inspect and Clean Erbttec Interior..... ☐ ☐ ☐ ☐

Check PDU power and cables..... ☐ ☐ ☐ ☐

Check gradient cables..... ☐ ☐ ☐ ☐

Check gradient filters and PHPS..... ☐ ☐ ☐ ☐

Service Request Number _____

Site Name: _____

FE Name: _____

Planned Maintenance Checklist

Month: _____

3200 NW 27th Ave Ste. 100
Pompano Beach, FL
33069
Office: 888-400-3979
aisserv.com

Procedure

IMAGE QUALITY

	OK	Cleaned	Adjusted	Needs Correction
Check ISO-Center (Check Alignment Lights).....	☐	☐	☐	☐
TG _____ AX _____ iso Z _____				
Check Gradient Calibration: X,Y,Z-axis.....	☐	☐	☐	☐
Cfxfull - Old _____ New _____ Cfyfull - Old _____ New _____ Cfzfull - Old _____ New _____				
Check Gradient Calibration zoom mode (TRM ONLY): X,Y,Z-axis.....	☐	☐	☐	☐
Cfxfull - Old _____ New _____ Cfyfull - Old _____ New _____ Cfzfull - Old _____ New _____				
Check SNR and RF output- Head.....	☐	☐	☐	☐
TG _____ Signal Noise SNR				
Check SNR and RF output- Body.....	☐	☐	☐	☐
TG _____ Signal Noise SNR				
Check (CTLTOP)..... Display CV (saveinter) check for grainy images	☐	☐	☐	☐
TG _____				
Check (CTLBOT)..... Display CV (saveinter) check for grainy images	☐	☐	☐	☐
TG _____				
Check LVShim Localizer (DC off-sets before centering phatom).....	☐	☐	☐	☐
TG _____				

LVshim Calibrations

Shim Values After

Shim Values Before

X _____ Y _____ Z _____

	X	Y	Z
Inhomogeneity	Z1(1,0)	X(1,1)	Z2X(3,1)
Spine	Z2(2,0)	Y(1,-1)	Z2Y(3,-1)
Shoulder	Z3(3,0)	ZX(2,1)	Zxy2(3,2)
	Z4(4,0)	ZY(2,-1)	ZXY(3,-2)
Average B field	Z5(5,0)	X2-Y2(2,2)	X3(3,3)
Std B field	Z6(6,0)	XY(2,-2)	Y3(3,-3)

**T
R
M

O
N
L
Y**

Zoom Mode (After)

	X	Y	Z	Z2
Inhomogeneity	Z1(1,0)	X(1,1)	Z2X(3,1)	
Spine	Z2(2,0)	Y(1,-1)	Z2Y(3,-1)	
Shoulder	Z3(3,0)	ZX(2,1)	Zxy2(3,2)	
	Z4(4,0)	ZY(2,-1)	ZXY(3,-2)	
Average B field	Z5(5,0)	X2-Y2(2,2)	X3(3,3)	
Std B field	Z6(6,0)	XY(2,-2)	Y3(3,-3)	

Supercon Shim Currents

	Ax1	Ax2	Ax3	Ax4	Ax5	Ax6
	T1-1	T1-2	T1-3	T1-4	T1-5	T1-6
Center Frequency	T2-1	T2-2	T2-3	T2-4	T2-5	T2-6

FINAL CHECKS

Check and Delete Error Message Log.....	☐	☐	☐	☐
Do Save INFO, Save GUI and Save Protocol.....	☐	☐	☐	☐
RESTORATION - DQA Phantom scan, Remove Test Scans, and Check Cabinet Doors / Covers.....	☐	☐	☐	☐
RESTORATION - Shut down system pre site procedure.....	☐	☐	☐	☐
RESTORATION - Update site log book, E-mail update with any follow up's and status.....	☐	☐	☐	☐

N
o
t
e
s