



CONTROL

8

7

2

65

64

62

COUNTER

5 G*

63

11

1 G*

3' G*

10 G*

9

6

5

4

66

12

3

1

90

0.5 G*

TABLE

MIN. 9000 SIZE

10

60

MAGNET ROOM

8'-0"

8'-0"

9'-2"

14'-1"

EQUIPMENT

MOV. CAP. AND

NOT MOV. MET.

MOV. BUS. TRA.

SEE APPENDIX "A" FOR POWER SPECIFICATIONS

ROOM PLAN LOCATED ON PAGE 8-68

ITEM	QUANTITY	ITEM DESCRIPTION	WEIGHT	HEAT OUTPUT
①	1	REGULATED STEP DOWN TRANSFORMER	861 lbs	7003 btu
②	1	UNINTERRUPTIBLE POWER SUPPLY	44 lbs	61 btu
③	1	GRADIENT POWER SUPPLY	881 lbs	15201 btu
④	1	SHIELD COOLER CABINET	313 lbs	28344 btu
⑤	1	SCAN ROOM UNIT	264 lbs	3412 btu
⑥	2	SCAN ROOM UNIT COVERS		
⑦	1	OPERATOR'S CONSOLE	555 lbs	2730 btu
⑧	1	OPERATOR'S CHAIR		
⑨	1	MAGNET RUNDOWN UNIT	8 lbs	
⑩	1	0.5T VMX-II MAGNET	9045 lbs	12389 btu
⑪	1	PATIENT TABLE	793 lbs	
⑫	1	INTEGRATED POWER SYSTEM	661 lbs	3412 btu
***** CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS *****				
60		MINIMUM 7 FT. -0 IN. [2133 MM] X 8 FT. -0 IN. [2438 MM] REMOVABLE WALL SECTION FOR MAGNET DELIVERY.		
61		MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 43 IN. W x 82 IN. H [1092mm x 2083mm], CONTINGENT ON A 96 IN. [2438mm] CORRIDOR WIDTH		
62		BASE CABINET FOR STORAGE OF SURFACE COILS, PATIENT POSITIONING PADS, ETC.		
63		MAGNET ROOM EXHAUST FAN		
64		RF SCREEN, INCLUSIVE OF WALLS, FLOORS, DOOR, ETC. GROUND IMPEDANCE GREATER THAN 1000 OHMS. ATTENUATION 100 DB AT 10-100MHZ PLANEWAVE.		
65		METAL DETECTOR (HAND HELD)		
66		ACCESS FLOORING		
90		MAIN DISCONNECT CONTROL - 94 LBS. <43 KG.>, 900 BTU/HR (264W) CAT NO. R4503AT FOR 480-3 WYE.		

ENVIRONMENT

Ambient Operation Temperature: 59-86° F [15-30° C]
 Allowable Temperature Change: 5° F [3° C] per hour
 Humidity: 30-70% non condensing
 Allowable Humidity Change: 5% per hour
 Max - Altitude: 8,000 ft. [2,438m] above sea level