

Service Request Number \_\_\_\_\_

# PM Checklist

Site Name: \_\_\_\_\_

SwMR CT Planned Maintenance Checklist

FE Name: \_\_\_\_\_

Month: \_\_\_\_\_

| Procedure                                                                                                                                              | 250k   | 500k | 750k | 1000k |
|--------------------------------------------------------------------------------------------------------------------------------------------------------|--------|------|------|-------|
| <b>GANTRY</b>                                                                                                                                          |        |      |      |       |
| Check number of gantry revolutions                                                                                                                     |        |      |      |       |
| Check scan and laser windows                                                                                                                           |        |      |      |       |
| Clean ventilation holes and check fans                                                                                                                 |        |      |      |       |
| Handling and Removal of Slip Ring Brush Debris                                                                                                         |        |      |      |       |
| Slip Ring Vacuuming                                                                                                                                    |        |      |      |       |
| Clean or Replace tube fan filter                                                                                                                       |        |      |      |       |
| SDAS DAS/Duct Cleaning (Only for assemblies with filter)                                                                                               |        |      |      |       |
| Inspect the axial motor brushes                                                                                                                        |        |      |      |       |
| Grease main bearing                                                                                                                                    |        |      |      |       |
| Check drive belt for wear                                                                                                                              |        |      |      |       |
|                                                                                                                                                        |        |      |      |       |
|                                                                                                                                                        |        |      |      |       |
| <b>PDU/GPDU/CRPDU</b>                                                                                                                                  |        |      |      |       |
| Inspect / Clean                                                                                                                                        |        |      |      |       |
| Inspect fan operation                                                                                                                                  |        |      |      |       |
| Inspect 7 green indicator lamps on GPDU                                                                                                                |        |      |      |       |
| Inspect 4 green lamps (L1, L2, L3, DCRGS status) PDU                                                                                                   |        |      |      |       |
| Inspect 6 green lamps on transformer assembly PDU                                                                                                      |        |      |      |       |
| Check tightness of wire lugs                                                                                                                           |        |      |      |       |
| CRPDU PM checks                                                                                                                                        |        |      |      |       |
| <b>SYSTEM</b>                                                                                                                                          |        |      |      |       |
| Slice/mA Count                                                                                                                                         | Enter: |      |      |       |
| Verify emergency off buttons (table/console)                                                                                                           |        |      |      |       |
| Check X-Ray on lights (HHS)                                                                                                                            |        |      |      |       |
| Check operation of final scan abort push buttons                                                                                                       |        |      |      |       |
| Monitor Sys Logs                                                                                                                                       |        |      |      |       |
| General Cleaning                                                                                                                                       |        |      |      |       |
| Image Performance QA                                                                                                                                   |        |      |      |       |
| <b>HIGH VOLTAGE</b>                                                                                                                                    |        |      |      |       |
| Check Oil in Tube and Tank Wells                                                                                                                       |        |      |      |       |
| Auto mA cal                                                                                                                                            |        |      |      |       |
|                                                                                                                                                        |        |      |      |       |
|                                                                                                                                                        |        |      |      |       |
| <b>CONSOLE / TABLE</b>                                                                                                                                 |        |      |      |       |
| Check for smooth rolling trackball / mouse                                                                                                             |        |      |      |       |
| Upper Assembly PM                                                                                                                                      |        |      |      |       |
| Lower Assembly PM (CT/I only)                                                                                                                          |        |      |      |       |
| Check for wear or damage on table extensions                                                                                                           |        |      |      |       |
| <b>FINAL</b>                                                                                                                                           |        |      |      |       |
| Perform Site Restoration - Check DQA Phantom, Check Alignment Lights, Remove Test Scans, and Check Cabinet Doors and Covers, Update vBLOG and site log |        |      |      |       |