GE Medical Systems

gehealthcare.com

Technical Publication

Direction 5113145-100

Revision 03

GE Medical Systems LightSpeed 5.X Functional Interconnects

ADVANCED SERVICE DOCUMENTATION
PROPERTY OF GE
FOR GE SERVICE PERSONNEL ONLY
NO RIGHTS LICENSED
DO NOT USE OR COPY
DISCLOSURE TO THIRD PARTIES PROHIBITED

Copyright © 2004, 2005 by General Electric Company.
All rights reserved.

LEGAL NOTES

PROPRIETARY TO GE MEDICAL SYSTEMS

Permission to use this Advanced Service Software and related documentation (herein called the material) by persons other than GE Medical Systems employees is provided only under an Advanced Service Package License relating specifically to this Proprietary Material. This is a different agreement from the one under which operating and basic service software is lessened. A license to use operating or basic service software does not extend to or cover this software or related documentation.

If you are a GE Medical Systems employee or a customer who has entered into such a license agreement with GE Medical Systems to use this proprietary software, you are authorized to use this Material according to the conditions stated in your license agreement. However, you do not have the permission of GE Medical Systems to alter, decompose or reverse-assemble the software, and unless you are a GE employee, you may not copy the Material. The Material is protected by Copyright and Trade Secret laws; the violation of which can result in civil damages and criminal prosecution.

If you are not party to such a license agreement or a GE Medical Systems Employee, you must exit this Material now.

COPYRIGHTS

All Material, Copyright © 2004, 2005 by General Electric Company. All rights reserved.

IMPORTANT PRECAUTIONS

LANGUAGE

WARNING

- THIS SERVICE MANUAL IS AVAILABLE IN ENGLISH ONLY.
- IF A CUSTOMER’S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER’S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR OR PATIENT FROM ELECTRIC SHOCK, MECHANICAL OR OTHER HAZARDS.

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DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation “damage in shipment” written on all copies of the freight or express bill before delivery is accepted or “signed for” by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report:

- Call 1-800-548-3366 and use option 8.
- Fill out a report on <http://us44hdd21/sctq/InstallFulfill/InstalFulfillment.htm>
- Contact your local service coordinator for more information on this process.

Rev. Jan. 5, 2005

CERTIFIED ELECTRICAL CONTRACTOR STATEMENT

All electrical Installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE Medical personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE’s electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE’s field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

IMPORTANT...X-RAY PROTECTION

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The General Electric Company, Medical Systems Group, will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the General Electric Company, Medical Systems Group, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective materials and devices are available. It is urged that such materials or devices be used.

CAUTION
Risk of
Explosion

ATTENTION
Danger
d’Explosion

LITHIUM BATTERY CAUTIONARY STATEMENTS

Danger of explosion if battery is incorrectly replaced. Replace only with the same or equivalent type recommended by the manufacturer. Discard used batteries according to the manufacturer’s instructions.

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OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the GEMS PQR Process to report all omissions, errors, and defects in this publication.

Revision History

Rev	Date	Reason for change
1	06/07/04	Initial Release
2	04/21/05	PQR 13034160; PQR 13034162
3	07/21/05	PQR 13036875

Table of Contents

Chapter 1	
System Communications.....	7
Section 1.0 — Slip Ring Communications Interface - 2399621	8
Section 2.0 — Final Scan Control, Console Ethernet - 2399622	9
Section 3.0 — Rotating Controller Interface Bus (RCIB) - 2399623.....	10
Chapter 2	
X-Ray Generation	11
Section 1.0 — Tube Fans / Pump / SmartID Interface - 2399624.....	12
Section 2.0 — HP60 JEDI Tube Fans/Pump/Smart ID Interface - 5141517..	13
Section 3.0 — JEDI Interface - kV Loop - 2399625	14
Section 4.0 — HP60 JEDI kV Loop - 5141524	15
Section 5.0 — JEDI Interface - mA Loop - 2399628.....	16
Section 6.0 — HP60 JEDI mA Loop - 5141527	17
Section 7.0 — X-Ray Collimation & Filtration - 2399626.....	18
Section 8.0 — Exposure Enable and Final Command - 2399629	19
Section 9.0 — Tube Rotor Control - 2399630.....	20
Section 10.0 — HP60 JEDI Tube Rotor Control - 5141528.....	21
Section 11.0 — HP X-Ray Tube Coolant Flow - 2399628	22
Chapter 3	
Data Acquisition	23
Section 1.0 — Trigger Generation - 2399724	24
Section 2.0 — Detector Heater (DAS Item) - 2399725	25
Chapter 4	
Patient Positioning	27
Section 1.0 — Alignment Lights - 2399726	28
Section 2.0 — Gantry Tilt, from Gantry, Console - 2399727	29
Section 3.0 — Table Cradle Sync - 2399728	30
Section 4.0 — Table Elevation - 2399729	31
Section 5.0 — Table Cradle Drive - 2399730	32

Chapter 5	
Axial Control	33
Section 1.0 — Axial Servo Control / Dynamic Brake - 2399731	34
Section 2.0 — Axial Power Distribution - 2399732	35
Chapter 6	
Operator I/O.....	37
Section 1.0 — Gantry Operator Interface and Displays - 2399733.....	38
Section 2.0 — Gantry Reset from Console - 2399734	39
Section 3.0 — Cardiac Options - 2399735	40
Section 4.0 — Pulmonary Options - 2399736.....	41
Section 5.0 — Scan Footswitch and IHHC Controller - 2399737	42
Section 6.0 — Site X-Ray Light, System Light - 2399738.....	43
Section 7.0 — Autovoice / Intercom - 2399739	44
Chapter 7	
System Monitoring	45
Section 1.0 — Gantry Temperature and Balance Sensors - 2399740	46
Section 2.0 — Gantry Cover Fans - 2399741	47
Chapter 8	
System Power Control	49
Section 1.0 — System E-Stop and Power Interlocks - 2399742	50
Section 2.0 — PDU - Gantry 120 VAC Control - 2399743	51
Section 3.0 — HV_DC_Control - 2399744.....	52
Section 4.0 — HP MDAS/GDAS Gantry Power Distribution - 2399627	53

THIS PAGE INTENTIONALLY BLANK.

Chapter 1

System Communications

Section 1.0 — Slip Ring Communications Interface - 2399621


Figure 1-1 2399621 Sheet 1, Rev. 1

Section 2.0 — Final Scan Control, Console Ethernet - 2399622


Figure 1-2 2399622 Sheet 1, Rev 1

Section 3.0 — Rotating Controller Interface Bus (RCIB) - 2399623


Figure 1-3 2399623 Sheet 1, Rev. 1

Chapter 2

X-Ray Generation

Section 2.0 — HP60 JEDI Tube Fans/Pump/Smart ID Interface - 5141517


Figure 2-2 5141517 Sheet 1, Rev 1

Section 3.0 — JEDI Interface - kV Loop - 2399625


Figure 2-3 2399625 Sheet 1, Rev. 1

Section 4.0 — HP60 JEDI kV Loop - 5141524


Figure 2-4 5141524 Sheet 1, Rev 1

Section 5.0 — JEDI Interface - mA Loop - 2399628


Figure 2-5 2399628 Sheet 1, Rev. 1

Section 6.0 — HP60 JEDI mA Loop - 5141527



Figure 2-6 5141527 Sheet 1, Rev 1

Section 7.0 — X-Ray Collimation & Filtration - 2399626

Section 8.0 — Exposure Enable and Final Command - 2399629


Figure 2-8 2399629 Sheet 1, Rev. 1

REV	---	PCN	-----	DESCRIPTION	-----						MADE BY:	MICHAEL MILLER	 GE Medical Systems Milwaukee, WI USA	DESIGN TITLE		EXPOSURE ENABLE + FINAL COMMAND FUNCTIONAL						
MADE BY	-----			DATE	-----	APP BY	-----	DATE	-----		DATE:	20 JUNE 2001		FIRST MADE FOR		H-POWER FUNC						
REV	---	PCN	-----	DESCRIPTION	-----									 GE Medical Systems Milwaukee, WI USA	2399629SCH		REV	1	SIZE	B	SHEET	1 OF 1
MADE BY	-----			DATE	-----	APP BY	-----	DATE	-----													

Section 9.0 — Tube Rotor Control - 2399630



Figure 2-9 2399630 Sheet 1, Rev. 1

Section 10.0 — HP60 JEDI Tube Rotor Control - 5141528


REV --- PCN --- DESCRIPTION --- MADE BY --- DATE --- APP BY --- DATE ---	MADE BY: KENICHI ODAGIRI DATE: MAY/10/05	 GE Medical Systems Milwaukee, WI USA	DESIGN TITLE HP60_JEDI_TUBE_ROTOR_CONTROL FIRST MADE FOR H-POWER 60 FUNC			
REV --- PCN --- DESCRIPTION --- MADE BY --- DATE --- APP BY --- DATE ---	CHECK BY: NAOKI ARIYAMA DATE: MAY/10/05		5141528SCH	REV 1	SIZE B	SHEET 1 OF 1

Figure 2-10 5141528 Sheet 1, Rev 1

Section 11.0 — HP X-Ray Tube Coolant Flow - 2399628



Figure 2-11 2399628 Sheet 1, Rev. 1

Chapter 3

Data Acquisition

Section 1.0 — Trigger Generation - 2399724

REV --- PCN ----- DESCRIPTION -----	MADE BY: NATHANAEL HUFFMAN	DESIGN TITLE HP TRIGGER GENERATION
MADE BY ----- DATE ----- APP BY ----- DATE -----	DATE: 12-MAR-2004	FIRST MADE FOR H-POWER FUNC
REV --- PCN ----- DESCRIPTION -----	GE Medical Systems Milwaukee, WI USA	2399724SCH
MADE BY ----- DATE ----- APP BY ----- DATE -----		REV 1 SIZE B SHEET 1 OF 1

Figure 3-1 2399724 Sheet 1, Rev. 1

Section 2.0 — Detector Heater (DAS Item) - 2399725



Figure 3-2 2399725 Sheet 1, Rev. 1

THIS PAGE INTENTIONALLY BLANK.

Chapter 4

Patient Positioning

Section 1.0 — Alignment Lights - 2399726

Section 2.0 — Gantry Tilt, from Gantry, Console - 2399727


Figure 4-2 2399727 Sheet 1, Rev. 1

Section 3.0 — Table Cradle Sync - 2399728



Figure 4-3 2399728 Sheet 1, Rev. 1

Section 4.0 — Table Elevation - 2399729

Section 5.0 — Table Cradle Drive - 2399730


Figure 4-5 2399730 Sheet 1, Rev. 1

Chapter 5

Axial Control

Section 1.0 — Axial Servo Control / Dynamic Brake - 2399731


Figure 5-1 2399731 Sheet 1 Rev. 1

THIS PAGE INTENTIONALLY BLANK.

Chapter 6

Operator I/O

Section 1.0 — Gantry Operator Interface and Displays - 2399733


Figure 6-1 2399733 Sheet 1, Rev. 1

Section 2.0 — Gantry Reset from Console - 2399734



Figure 6-2 2399734 Sheet 1, Rev. 1

Section 3.0 — Cardiac Options - 2399735

Section 4.0 — Pulmonary Options - 2399736


Figure 6-4 2399736 Sheet 1, Rev. 1

Section 5.0 — Scan Footswitch and IHHC Controller - 2399737


Figure 6-5 2399737 Sheet 1, Rev. 1

Section 6.0 — Site X-Ray Light, System Light - 2399738


Figure 6-6 2399738 Sheet 1, Rev. 1

Section 7.0 — Autovoice / Intercom - 2399739


Figure 6-7 2399739 Sheet 1, Rev. 2

Chapter 7

System Monitoring

Section 1.0 — Gantry Temperature and Balance Sensors - 2399740


Figure 7-1 2399740 Sheet 1, Rev. 1

Section 2.0 — Gantry Cover Fans - 2399741


Figure 7-2 2399741 Sheet 1, Rev. 1

THIS PAGE INTENTIONALLY BLANK.

Chapter 8

System Power Control

Section 1.0 — System E-Stop and Power Interlocks - 2399742


Figure 8-1 2399742 Sheet 1, Rev. 1

Section 2.0 — PDU - Gantry 120 VAC Control - 2399743



Figure 8-2 2399743 Sheet 1, Rev. 3

Section 3.0 — HV_DC_Control - 2399744

THIS PAGE INTENTIONALLY BLANK.

GE MEDICAL SYSTEMS

GE MEDICAL SYSTEMS-AMERICAS: FAX 262.312.7434
3000 N. GRANDVIEW BLVD., WAUKESHA, WI 53188 U.S.A.
GE MEDICAL SYSTEMS-EUROPE: FAX 33.1.40.93.33.33
PARIS, FRANCE
GE MEDICAL SYSTEMS-ASIA: FAX 65.291.7006
SINGAPORE