GE Healthcare

Accessory Interface Hardware Installation Manual

5321858-1EN

Rev 1

© 2008, General Electric Company,
All Rights Reserved.

LEGAL NOTES

TRADEMARKS

Adobe, the Adobe logo, Acrobat, the Acrobat logo, Exchange, and PostScript are trademarks of Adobe Systems Incorporated or its subsidiaries and may be registered in certain jurisdictions.

Microsoft is a registered trademark and Windows is a trademark of Microsoft Corporation.

All other products and their name brands are trademarks of their respective holders.

COPYRIGHTS

All Material, Copyright © 2008 by General Electric Company, Inc. All rights reserved.

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DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "damage in shipment" written on all copies of the freight or express bill before delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report:

- Call 1-800-548-3366 and use option 8.
- Fill out a report on <http://egems.med.ge.com/edq/home.jsp>
- Contact your local service coordinator for more information on this process.

Rev. June 13, 2006

CERTIFIED ELECTRICAL CONTRACTOR STATEMENT

All electrical Installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE Medical personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

IMPORTANT...X-RAY PROTECTION

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The General Electric Company, Medical Systems Group, will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the General Electric Company, Medical Systems Group, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective materials and devices are available. It is urged that such materials or devices be used.

LITHIUM BATTERY CAUTIONARY STATEMENTS

CAUTION
Risk of
Explosion

Danger of explosion if battery is incorrectly replaced. Replace only with the same or equivalent type recommended by the manufacturer. Discard used batteries according to the manufacturer's instructions.

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OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the GEMS CQA Process to report all omissions, errors, and defects in this publication.

Revision History

Revision	Date	Reason for change
1	11/13/08	Initial release.

Table of Contents

Chapter 1

Accessory Interface Hardware Installation 17

Section 1.0

Introduction 17

1.1 Overview 17

1.2 Prerequisites 17

Section 2.0

Material List 18

Section 3.0

Installation Procedure..... 20

3.1 Base Cover Removal..... 20

3.2 Install LAN Switch..... 21

3.2.1 LAN Switch (5115334) Assembly 21

3.2.2 LAN Switch Adapter (5136650) Assembly..... 21

3.3 Install IPC and Accessory Panel..... 22

3.4 Install Cable Harness..... 25

Section 4.0

Finish Installation..... 32

4.1 LED Inspection of IPC Board..... 32

4.2 LED Inspection of LAN Switch..... 33

4.3 Finalization..... 33

Chapter 1

Accessory Interface Hardware Installation

This option is only applied to LightSpeed RT16 system.

Section 1.0

Introduction

1.1 Overview

The Accessory Interface Hardware option provides interface for the following options.

- Prospective Respiratory Gating
- Cardiac Gating
- Advantage 4D Respiratory Gating
- Smart Step
- Xtream Injector

The Accessory Interface Hardware option is shipped ready for immediate installation. The following installation procedure describes how to install that hardware on the Gantry.


Refer to [Table 1-1](#) for a listing of the parts in the kit.

1.2 Prerequisites

LightSpeed RT16 Gantry with New Power Pan (5269061)

Section 2.0

Material List

MARK	PART NUMBER	QUANTITY	DESCRIPTION
1	5136002	1.0	Integrated Peripheral Controller PWA
2	5115334	1.0	Switch, Ethernet 8-port GBe
3	5136650	1.0	ups-120vac to 12vdc,
4	46-136246P18	2.0	1~ WIDE, WHT, ADH BACK, VELCRO HOOK 75 FT PER RL
5	2315482	1.0	VELCRO TAPE LENGTH IS 320
6	2315482-3	1.0	VELCRO TAPE LENGTH IS 830
7	5316347	1.0	IPC option kit cable harness, RT Vision
8	5306106	1.0	Accessory Panel LH ASM, CT HD
9	5306107	1.0	Accessory Panel RH ASM, CT HD
10	5306201	1.0	Cable, Cardiac Pulmonary, CT HD
11	5309358	1.0	BSD V3 IPC MNT brkt cushion
12	5309359	1.0	BSD V3 IPC MNT brkt
13	5309361	1.0	BSD V3 connector MNT brkt

Table 1-1 5323070 Material List

MARK	PART NUMBER	QUANTITY	DESCRIPTION
15	5304664	1.0	Cover IPC Board, Kitty Hawk
16	46-136246P16	2.0	1~ WIDE, WHT, ADH BACK, VELCRO LOOP 75 FT PER RL
104	46-328417P3	18.0	SCREW HEXAGON SOCKET 4 MM 10 MM ZINC PLATED, STEEL
105	46-328432P2	25.0	WASHER LOCK - SPRING 4.1 MM 7.6 MM ZINC PLATED, SPRING STEEL
106	46-328430P2	25.0	WASHER PLAIN - NORMAL 4.3 MM 9 MM ZINC PLATED, STEEL
107	U0099AG	3.0	STANDOFF/SPACER MALE/FEMALE 4 MM 15 FF15 STUD
108	46-328425P1	4.0	NUT HEXAGON 4 MM
	5321858-1EN	1.0	Option Installation Manual (This Manual)

Table 1-1 5323070 Material List

Section 3.0 Installation Procedure

DANGER RISK OF ELECTRIC SHOCK. USE PROPER LOCKOUT/TAGOUT PROCEDURES.

3.1 Base Cover Removal

- 1.) Remove the two gantry side covers. Refer to your System Service Manual for any needed instruction.
- 2.) Turn **OFF** power of PDU and all 3 switches (Axial Drive, HVDC, 120 VAC) on the service switch panel on Gantry.
- 3.) Remove all Gantry base covers.


Figure 1-1 Gantry Base Cover

3.2 Install LAN Switch

3.2.1 LAN Switch (5115334) Assembly

- 1.) Clear the contact surface of 5316196 and 5115334 using alcohol before adhering 46-136246P18 and 46-136246P16.
- 2.) Adhere 46-136246P18 (two slices) to Bracket 5316196 which is at gantry leg of right side.
- 3.) Adhere 46-136246P16 (two slices) to LAN switch (5115334)
- 4.) Press LAN switch on 5316196 and use Velcro tape (2315482-3) to fix it. Please refer to [Figure 1-2](#)


Figure 1-2 LAN Switch Assembly

3.2.2 LAN Switch Adapter (5136650) Assembly

- 1.) Fix the adapter (5136650) using Velcro tape (2315482) to the bracket.

3.3 Install IPC and Accessory Panel

- 1.) Assemble the IPC and the left accessory panel.


Figure 1-3 IPC Board and Left Accessory Panel

- 2.) Install the IPC assembly and left panel to the Gantry rear base.


Figure 1-4 IPC and Left Accessory Panel Installation

- 3.) Connect 5306201 cable (Connector Bracket <-> IPC, J11 <-> Left Accessory Panel).


Figure 1-5 Cable (5306210) Connection

- 4.) Assemble the right accessory panel, and install it to the Gantry rear base.


Figure 1-6 Right Accessory Panel Installation

3.4 Install Cable Harness

- 1.) Carefully route the cable harness (5316347) as shown in illustration [Figure 1-7](#).



Figure 1-7 Cable Routing

- 2.) Connect all cable connectors of cable harness as shown in [Figure 1-8](#).
(Refer to [Figure 1-9](#) through [Figure 1-14](#) for further information)

Cable Harness Assembly (5316347) Connection

Item	Part No.	Description	Cable Connection
1	5189456	Cable Asm - IPC power	IPC, J6 <-> TGP, J1 <-> Power Connector
2	5193967	Cable - TGPU J5 to LAN Switch 1	TGPU, J5 <-> LAN Switch 1
3	5169215	cable, cardiac-pulmonary_PathFinder	IPC, J13 <-> TGPU, J22
4	2286150-2	Smart step cable for BSD Delight Lite.	R-panel, J2, J3 <-> TGPU, J27 <-> L-panel, J2, J3
5	5306481	Power Cable, IPC Accessory Panel, VCT HD	R-panel <-> Power Pan <-> L-panel
6	5193969-10	Cable -LAN	IPC, J5 <-> LAN Switch 2
7	5265804	Cable Asm - Network Switch Power, Kitty Hawk	Service Outlet <->Harness Outlet <->LAN Switch
8	5193969-2	Cable - LAN	Power Pan <-> LAN Switch 3

Figure 1-8 Cable Connection

- 3.) Connect cables 5189456, 5193967, 5169215 and 2286150-2 to DC Power, J5, J22 and J27 of TGPU BD according to their labels. as shown in [Figure 1-9](#).


Figure 1-9 TGPU Connection

4.) Connect cables 2286150-2, 5306481 and 5306201 to Left Panel, as shown in [Figure 1-10](#).

Note: Brown wire of 5306481 should be connected to L pin of the power outlet, blue wire of 5306481 should be connected to N pin of the power outlet.


Figure 1-10 Left Panel Connection

- 5.) Connect cables 5189456, 5169215 and 5193969-10 to J6, J13 and J5 of IPC BD according to their labels. Connect cable 5306201 to J11 of IPC BD and IPC bracket according to its labels, as shown in [Figure 1-11](#).


Figure 1-11 IPC BD Connection

- 6.) Connect cables 5193967, 5193969-10 and 5193969-2 to Port 1, 2 and 3 of Lan Switch according to their labels. Connect cable 5265804 to Service Outlet Power and Lan Switch adapter according to its labels as shown in [Figure 1-12](#).


Figure 1-12 LAN Switch Connection

7.) Connect cables 5306481 and 2286150-2 to Right Panel, as shown in [Figure 1-13](#).

Note: Brown wire of 5306481 should be connected to L pin of the power outlet, blue wire of 5306481 should be connected to N pin of the power outlet.


Figure 1-13 Right Panel Connection

8.) Connect cable 5193969-2 to the bulkhead of Power Pan according to its label, as shown in [Figure 1-14](#).

Figure 1-14 Network Cable Connection

- 9.) Connect cable 5306481 to Power Pan, Gantry should be tilt 25° forward to access the connector, as shown in [Figure 1-15](#).


Figure 1-15 Cable to Power Pan

- 10.) Secure the cables with tie wraps.
- 11.) Restore power to the system.

Section 4.0

Finish Installation

4.1 LED Inspection of IPC Board

- 1.) To make sure the IPC board is functioning, an LED inspection should be made on the board. The LED configuration on the board should be as follows with power on but no options connected to the options interface panels:

LED	Color	Label	Status (without options connected)
DS1	Green	24V_PWR	On Solid
DS2	Green	PHY_FW_LED	Blinking with network activity
DS3	Green	HRTBT	Blinking
DS4	Green	FW_LED_0	Off
DS5	Green	FW_LED_1	Off
DS6	Green	FW_ERROR	Off
DS7	Green	24V_FUSE	On Solid
DS8	Green	RESET	On Solid
DS9	Yellow	RESET	Off
DS10	Green	R_PULSE_FILTERED	Off
DS11	Green	RESP_PULSE_FILTERED	Off
DS12	Green	R_PULSE_UNFILTERED	Off
DS13	Green	RESP_PULSE_UNFILTERED	Off
DS14	Yellow	FPGA_LOAD	Off
DS15	Green	FPGA_LOAD	On Solid
DS16	Yellow	B/A_CODE	Off
DS17	Green	B/A_CODE	On Solid
DS18	Green	+1.2V_OK	On Solid
Ethernet	Orange/Green	Ethernet-Hub	Blinking with network activity

Table 1-2 IPC Board LED Status


Figure 1-16 IPC Board LED locations

4.2 LED Inspection of LAN Switch

- 1.) Visually verify the front panel Power LED of the switch is illuminated.
- 2.) Visually verify that all appropriate Port LEDs show active links.

4.3 Finalization

- 1.) Perform an application scan to ensure proper system operation before turning the system over to the customer.
- 2.) Turn OFF all 3 switches (Axial Drive, HVDC, 120 VAC) on the service switch panel.
- 3.) Reinstall both rear base covers (see [Figure 1-1](#)).
- 4.) Reinstall the gantry side covers. Refer to your System Service Manual for any needed instruction.

www.gehealthcare.com

Imagination at work