

**ADMINISTRATIVE INFORMATION****FMI No: 25366****FMI Title: VCT M3 Patch 3**

| FMI TYPE                                      | CHARGE CODES | IMPLEMENTATION DATES          |
|-----------------------------------------------|--------------|-------------------------------|
| <input type="checkbox"/> Mandatory Safety     | GEMS - AM:   | Issue Date: May 20, 2005      |
| <input checked="" type="checkbox"/> Mandatory | GEMS - E:    |                               |
| <input type="checkbox"/> Mandatory on Request | GEMS - A:    | Close Date: September 2, 2005 |

- ☒ **Staged FMI** **Customer List** provided by FMI Administration. If you identify additional units within effectivity range, but they are not shown on the customer list, contact FMI Administration in your Pole.
- ☒ **Parts** and any special tools are furnished automatically. To address shortages and ordering exceptions, contact FMI Administration in your Pole.
- ☐ **Parts Ordering:** GEMS - Americas, Contact your Region logistics coordinator  
GEMS - Europe, Order FMI kits from GPO/EPL Buc  
GEMS - Asia, Contact your region logistics coordinator
- ☐ **No Parts** are required to complete this FMI.
- ☐ **Batch with** \_\_\_\_\_
- To minimize costs, combine implementation of this FMI with the FMI's listed for all customers where the effectivity is the same (Please see your customer lists).
- ☐ **Non-Staged** **No Customer List is provided.** Each service location must identify affected units, and place part orders as indicated in the Special Instructions.

**Technical Support**

| GEMS - Americas                                                                                   | GEMS - Europe                                                                                                    | GEMS - Asia                                                         |
|---------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------|
| GE Medical Systems<br>Insite Center<br>(414) 524-5300<br>800-321-7973<br>Milwaukee, WI 53201, USA | GEMS - E / ESC<br>BP34<br>F 78533 BUC CEDEX, FRANCE<br>Tel: (33 1) 39 20 00 07 or ESC<br>Fax: (33 1) 30 70 99 70 | Asia Support Center<br>Phone: 81-426-56-0033<br>Fax: 81-426-48-2905 |

**Special Instructions:****Unit Tracked:****FMI Instructions Part No.: 5141223-100****Do NOT leave unsecured on site**

# **FMI** Completion Certification Report

## Mail To

|                                                                                       |                                                                     |                                                         |
|---------------------------------------------------------------------------------------|---------------------------------------------------------------------|---------------------------------------------------------|
| <b>GEMS - Americas</b>                                                                | <b>GEMS - Europe</b>                                                |                                                         |
| G.E. Medical Systems<br>Attn: FMI Admin. W-523<br>P.O. Box 414<br>Milwaukee, WI 53201 | GEMS - Europe<br>FMI Admin. 322<br>BP34<br>F78533 Buc Cedex, France | Send to your<br>Country FMI or Service<br>Administrator |

FMI No.: 25366

FMI Title: VCT M3 Patch

Customer Full Name:

Number and Street Address:

City and Country:

Country Code

Site Number

System I.D.

### Unit Tracked:

| Model Number | Serial Number | Compl. Date<br>DD/MM/YY | FMI Comp.<br>Code | Service Hours |        |
|--------------|---------------|-------------------------|-------------------|---------------|--------|
|              |               |                         |                   | Labor         | Travel |
|              |               |                         |                   |               |        |
|              |               |                         |                   |               |        |
|              |               |                         |                   |               |        |
|              |               |                         |                   |               |        |

### FMI Completion Code

1. Modification Complete or previously modified.
2. FMI not required per FMI instruction.
3. Not modified. Unit in storage. FMI kit returned to stock.
4. Modification refused or equipment altered. Customer signature and title required.
5. Unit location unknown or unit not in area. Product Locator notified.

### Customer Refusal or Equipment Altered - Modification Pending.

Installer

Comments

## Effectivity

The following CT Equipment is affected:

| Pole   | Customer / Location            | Primary FE                          | Shipped as | System ID    |
|--------|--------------------------------|-------------------------------------|------------|--------------|
| AM     | Froedtert                      | Lee Simplaski                       | ME         | 414805CT4    |
| AM     | Education Center               | Gary Juedes /                       | ME         | 414541VCT1   |
| Europe | CCN                            | Francis Bonin /<br>Christian Garcia | ME         | M2322223     |
| AM     | JFK                            | Joe Ward                            | ME         | 732321VCT    |
| Japan  | Keio University Hospital       | Jiro Sato                           | ME         | XG0009       |
| China  | Hong Kong Sanatorium           | Benny Liu                           | ME         | 0828210048   |
| AM     | Univ. of Washington            | Bob Flye                            | ME         | 206598VCT2   |
| Europe | Charite                        | Andreas Kuring                      | ME         | CC0618CT07   |
| China  | AsiaMedic                      | Newton Lee                          | ME         | 0847210010   |
| Japan  | Toki Municipal Hospital        | Shuji Suzuki                        | ME         | QS0001       |
| Japan  | Handai Fuzoku (Osaka)          | Jiro Sato/Shigeya<br>Yoshimura      | ME         | QX0008       |
| Europe | Dortmund                       | Juergen Herstell                    | ME         | HC0074CT04   |
| AM     | Huntsville Medical Center, AL  | James Stewart                       | M3         | 256MAINVCT   |
| AM     | Cardiology Associates, AL      | Ned (Butch) Cobb                    | M3         | 251340VCT    |
| AM     | West Suburban, IL              | Luke Harvey/Anthony<br>Lavareto     | M3         | 630972VCT    |
| Europe | CHUV Lausanne                  | Jean-Louis Loconte                  | ME         | 260010CT15   |
| Europe | RIKSHOSPITALET HF              | Anders Hersvik                      | M3         | N01016CT06   |
| AM     | University of Michigan         | Bill Marlette                       | M3         | 734936VCT    |
| Europe | DI ROMA TOR VERGATA, Italy     | Fabrizio Franceschini               | M3         | A5333029     |
| China  | BEIJING PLA NO.301 HOSPITAL    | Zhong Ren/Zhang<br>Guojon           | M3         | 082421040148 |
| AM     | AMI                            | Nick Iannotti                       | M3         | 609652VCT    |
| China  | Beijing Xuanwu Hospital        | Zhong Ren/Song<br>Minghe            | M3         | 082421040151 |
| AM     | Miami Baptist, FL              | Tony Mancini                        | M3         | 305596VCT    |
| AM     | Boston Med Ctr, MA             | David Goodline                      | M3         | 617414VCT1   |
| AM     | North Shore Univ., NY          | Michael Bozeat                      | M3         | 516562VCT    |
| AM     | Bronson Methodist Hospital, MI | David Newberry                      | M3         | 269341VCT    |

| Pole   | Customer / Location                         | Primary FE         | Shipped as | System ID  |
|--------|---------------------------------------------|--------------------|------------|------------|
| AM     | Mobile Infirmary                            | Ned (Butch) Cobb   | M3         | 251435CT3  |
| AM     | Morristown, NJ / Mountainside               | Herb Tunison       | M3         | 973429VCT  |
| AM     | East River Medical imaging, NY              | Eric Wesen         | M3         | 212288CT64 |
| AM     | St. Mary's Hospital, MI                     | Paul Vanslambrouck | M3         |            |
| AM     | Methodist Hospital, TX                      | Don Swisher        | M3         | 713441MCT9 |
| Europe | HOPITAL BEAUJON                             |                    | M3         | M4016557   |
| China  | TZU CHI BUDDHIST GENERAL HOSPITAL HSIN DIAN |                    | M3         | 0853210161 |
| AM     | LHSC - UNIVERSITY HOSPITAL                  | Canada             | M3         | 519685VCT  |
| Japan  | IRYOU HOUJIN SHADAN HOKUTO HOSPITAL         | Sasaki Takeshi     | M3         | QV0030     |
| Europe | HOPITAL NECKER ENFANTS MALADES              | Michel Joseoh      | M3         | M4033088   |
| Japan  | RENGOU HIROSHIMA SOUGOU HOSPITA             | Tanigawa Kosuke    | M3         | QS0002     |
| Europe | HOPITAL EUROPEEN GEORGES POMPIDOU           | Gilbert Huet       | M3         | M5460685   |
| Japan  | KENKOU HOKEN ISAHAYA SOUGOU HOSPITAL        | Imogawa Keiji      | M3         | QS0005     |
| AM     | UNIV OF WASHINGTON MED CTR                  | Bob Flye           | M3         | 206598VCT1 |
| Europe | MC Morice (Massy)                           |                    | M3         | M9227209   |
| AGM    | Metropolitan Health Service                 |                    | M3         |            |
| Europe | Huddinge (Karolinska)                       |                    | M3         | A007CT01   |
| Japan  | KOUBE SHIRITSU CHUO SHIMIN HOSPITAL         | Sasaki Masaki      | M3         | QS0007     |
| China  | SHOW-CHWAN MEMORIAL HOSPITAL                |                    | M3         | 0853210162 |
| AM     | INOVA MT VERNON HOSPITAL                    |                    | M3         | 703664LS64 |
| AM     | EVANSTON NORTHWESTERN HEALTHCARE CORP       |                    | M3         | 847570VCT  |
| AM     | PARK RIDGE HOSPITAL                         |                    | M3         | 828681CT   |

## ***Purpose***

The purpose of this FMI is to update the current VCT installed base with the newest M3 software patch.

## ***Furnished Materials***

FMI 25366 Kit (P/N 5141836) contains the following materials:

| Qty | Description                        | Part #      |
|-----|------------------------------------|-------------|
| 1   | FMI Document                       | 5141223-100 |
| 1   | LightSpeed 7.X operating software. | 5139419     |

## ***Tool Required***

| Qty | Description  | SW Version | Part #      |
|-----|--------------|------------|-------------|
| 1   | FMI Document | 05MW10.2   | 5141223-100 |

## ***Personnel Requirements***

1 hour labor and 1 hour of travel

## Procedure

### Installation Instructions for the 05MW10.2 VCT M3 Patch 3

1. Remove Auto mA option before starting patch installation.
2. If applicable, Shutdown Applications from the Netscape: Service Desktop:
  - Click on the Service Desktop Icon
  - 1, 2, 3 as required
  - Click on Utilities/Tools, Utilities
  - Click on Application Shutdown from the menu
3. Insert the 05MW10.2VCT M3 Patch 3 (5138492) CD-ROM into the DVD-ROM drive on the HP Host Computer.
4. When Application Software is down open a Unix Shell.

Type the following:

- `{ctuser@hostname} su - <Enter>`
- Password: `#bigguy <Enter>`
- `hostname# mount /mnt/cdrom <Enter>`
- `hostname# cd /mnt/cdrom <Enter>`
- `hostname# ./install_patch -p 05MW10.2VCTM3PATCH.3-1-1.i386.rpm <Enter>`

[Press the Escape key after 05MW10.2 instead of typing the whole string]

Installing 05MW10.2VCTM3PATCH.3 -1-1

*<messages indicating the progress>*

- `hostname# cd (space) / <Enter>`
- `hostname# umount /mnt/cdrom <Enter>`

5. Manually eject the CD.
6. Remove the 05MW10.2VCT M3 Patch 3 CDROM from the DVD ROM drive on the HP Host Computer.
7. Reboot before starting Application Software:
  - `hostname# sync <Enter>`
  - `hostname# sync <Enter>`
  - `hostname# reboot <Enter>`

Message appears:

The system is going down for reboot NOW!

8. Complete a flash download when the system comes back up. On the first attempt, the flash load may fail and cause the system to reboot. Try the flash download again after the reboot.
9. After completing a successful flash download, install the Auto mA option.
10. Complete a fastcal.

Installation of 05MW10.2VCT M3 Patch 3 process is complete.

## ***Customer communication***

**Please communicate the following to your customers:**

1. The maximum mA allowed for Scout scan is 80 mA. Use of higher mA values may cause recon to fail.
2. Auto mA for Plus mode has been corrected in this release of software.
3. To enable retro recon for 0.625mm, the Recon mode must be set to "Plus" in the Recon Options menu.
4. If the left head LCD goes blank, loses the view edit or Dynaplan information; select Retro Recon then Cancel when the Retro Recon interface is displayed. This will force the screen to reorder and the viewedit or Dynaplan screen should then be visible.