



Copyright © 2006 by General Electric Company, Inc.

FMI Completion Certification Report

Mail To

GEMS - Americas

GEMS - Europe

G.E. Medical Systems
Attn: FMI Admin. W-523
P.O. Box 414
Milwaukee, WI 53201

GEMS - Europe
FMI Admin. 322
BP34
F78533 Buc Cedex, France

Send to your
Country FMI or Service
Administrator

FMI No.: FMI 25370A

FMI Title: VCT EMI Gasket Upgrade

Customer Full Name:

Number and Street Address:

City and Country:

Country Code

Site Number

System I.D.

Unit Tracked:

Model Number	Serial Number	Compl. Date DD/MM/YY	FMI Comp. Code	Service Hours	
				Labor	Travel

FMI Completion Code

1. Modification Complete or previously modified.
2. FMI not required per FMI instruction.
3. Not modified. Unit in storage. FMI kit returned to stock.
4. Modification refused or equipment altered. Customer signature and title required.
5. Unit location unknown or unit not in area. Product Locator notified.

Customer Refusal or Equipment Altered - Modification Pending.

Installer
Comments

Effectivity

VCT systems installed before July 2006 that were missed on FMI 25370.

Tracked Model numbers: 5124069 VCT and 5123200 VCT Select 32 Slice.

Purpose

The purpose of this FMI:

- Install EMI gasket on the Slipring receivers.

NOTE: FMI 25370 – VCT Intake Filter and EMI Gasket Upgrade was distributed regarding this fix. It was realized in Engineering that a select set of systems was missed in the effectivity list. This FMI addresses those systems that were missed. If you already performed FMI 25370, **DO NOT** execute this FMI.

Close your dispatch as code 2 and return you FMI kit to:

HM Product Solutions

3225 Gateway Road

Brookfield, WI 53045

Attn: David Schultz

Prerequisite FMIs

None

Related FMIs

FMI 25370 – VCT Intake Filter and EMI Gasket Upgrade

Furnished Materials

FMI kit (5194285)

Qty	Description	Part #
1	FMI Document	5194285-100
1	Gasket Kit for VCT Slipring Receivers with installation instruction.	5144787

Tool Required

Qty	Description
1	Standard field service kit

Personnel Requirements

"2" hours labor and "1" hours of travel

Procedure

1. Install EMI gasket on slipring receivers

Install EMI gasket using instruction included in kit 5144787. Kit 5144787 contains two gaskets, two rubber spacers and the **instruction sheet**.

2. Finalization

1. Perform the **Gantry Rotation Safety Check**, found under Functional Checks on the Common Service Desktop.
2. Perform a test scan using the 20cm QA phantom to verify system operation.
3. Close the FMI.