

**ADMINISTRATIVE INFORMATION****FMI No: 25391****FMI Title:** VCT and Helios Collimator Inspection Procedure

| FMI TYPE                                                    | CHARGE CODES      | IMPLEMENTATION DATES |
|-------------------------------------------------------------|-------------------|----------------------|
| <input checked="" type="checkbox"/> <b>Mandatory Safety</b> | <b>GEMS - AM:</b> | <b>Issue Date:</b>   |
| <input type="checkbox"/> <b>Mandatory</b>                   | <b>GEMS - E:</b>  |                      |
| <input type="checkbox"/> <b>Mandatory on Request</b>        | <b>GEMS - A:</b>  | <b>Close Date:</b>   |

|                                                       |                                                                                                                                                                                                                                                                                                                                                |
|-------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input checked="" type="checkbox"/> <b>Staged FMI</b> | <b>Customer List</b> provided by FMI Administration. If you identify additional units within effectivity range, but they are not shown on the customer list, contact FMI Administration in your Pole.                                                                                                                                          |
| <input type="checkbox"/>                              | <b>Parts</b> and any special tools are furnished automatically. To address shortages and ordering exceptions, contact FMI Administration in your Pole.                                                                                                                                                                                         |
| <input type="checkbox"/>                              | <b>Parts Ordering:</b> GEMS - Americas, Contact your Region logistics coordinator<br>GEMS - Europe, Order FMI kits from GPO/EPL Buc<br>GEMS - Asia, Contact your region logistics coordinator                                                                                                                                                  |
| <input checked="" type="checkbox"/>                   | <b>No Parts</b> are required to complete this FMI.                                                                                                                                                                                                                                                                                             |
| <input type="checkbox"/> <b>Batch with</b>            | <hr/> To minimize costs, combine implementation of this FMI with the FMIs listed for all customers where the effectivity is the same (Please see your customer lists).<br><br>To minimize costs, combine implementation of this FMI with the FMIs listed for all customers where the effectivity is the same (Please see your customer lists). |
| <input type="checkbox"/> <b>Non-Staged</b>            | <b>No Customer List is provided.</b> Each service location must identify affected units, and place part orders as indicated in the Special Instructions.                                                                                                                                                                                       |

**Technical Support**

| GEMS - Americas                                                                 | GEMS - Europe                                                                                                    | GEMS - Asia                                                         |
|---------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------|
| GE Medical Systems<br>Insite Center<br>800-321-7937<br>Milwaukee, WI 53201, USA | GEMS - E / ESC<br>BP34<br>F 78533 BUC CEDEX, FRANCE<br>Tel: (33 1) 39 20 00 07 or ESC<br>Fax: (33 1) 30 70 99 70 | Asia Support Center<br>Phone: 81-426-56-0033<br>Fax: 81-426-48-2905 |

**Special Instructions:**

If during the course of performing this FMI a defective collimator is discovered, follow these steps.

1. Inform the customer the site is down, until a new collimator is installed.
2. Submit a PSR in iTrak to report site.
3. Hold onto defective collimator until give instructions as to where to return it for evaluation.
4. Order replacement collimator per instruction include in FMI Inspection Procedure.

**Unit Tracked:** Refer to Effectivity list on page 2.**FMI Instructions Part No.:** 5193689-100, Rev 3**Do NOT leave unsecured on site**

# **FMI** Completion Certification Report

## Mail To

|                                                                                       |                                                                     |                                                         |
|---------------------------------------------------------------------------------------|---------------------------------------------------------------------|---------------------------------------------------------|
| <b>GEMS - Americas</b>                                                                | <b>GEMS - Europe</b>                                                |                                                         |
| G.E. Medical Systems<br>Attn: FMI Admin. W-523<br>P.O. Box 414<br>Milwaukee, WI 53201 | GEMS - Europe<br>FMI Admin. 322<br>BP34<br>F78533 Buc Cedex, France | Send to your<br>Country FMI or Service<br>Administrator |

FMI No.: 25391

FMI Title: VCT and Helios Collimator Inspection Procedure

Customer Full Name:

Number and Street Address:

City and Country:

Country Code

Site Number

System I.D.

### Unit Tracked:

| Model Number | Serial Number | Compl. Date<br>DD/MM/YY | FMI Comp.<br>Code | Service Hours |        |
|--------------|---------------|-------------------------|-------------------|---------------|--------|
|              |               |                         |                   | Labor         | Travel |
|              |               |                         |                   |               |        |
|              |               |                         |                   |               |        |
|              |               |                         |                   |               |        |
|              |               |                         |                   |               |        |

### FMI Completion Code

1. Modification Complete or previously modified.
2. FMI not required per FMI instruction.
3. Not modified. Unit in storage. FMI kit returned to stock.
4. Modification refused or equipment altered. Customer signature and title required.
5. Unit location unknown or unit not in area. Product Locator notified

### Customer Refusal or Equipment Altered - Modification Pending.

Installer

Comments

## Table of Contents

|                                                                       |    |
|-----------------------------------------------------------------------|----|
| Table of Contents .....                                               | 1  |
| Revision History .....                                                | 2  |
| Effectivity .....                                                     | 2  |
| Purpose .....                                                         | 3  |
| Prerequisite FMIs.....                                                | 3  |
| Related FMIs .....                                                    | 3  |
| Furnished Materials .....                                             | 4  |
| Tool Requirements .....                                               | 4  |
| Personnel Requirements .....                                          | 4  |
| Before You Begin.....                                                 | 4  |
| Safety.....                                                           | 4  |
| INSPECTION PROCESS .....                                              | 5  |
| I.    HELIOS COLLIMATOR (LIGHTSPEED 1X-5X) INSPECTION PROCEDURE ..... | 6  |
| II.   VCT (LIGHTSPEED 7X) COLLIMATOR INSPECTION PROCEDURE.....        | 10 |
| Finalization.....                                                     | 12 |

## Revision History

| Revision | Date              | Reason for Change                |
|----------|-------------------|----------------------------------|
| 1        | December 12, 2006 | Initial release                  |
| 2        | January 30, 2007  | Updated with Pilot site feedback |
| 3        | February 23, 2007 | Updated Purpose section          |

## Effectivity

The effectivity for this FMI is a select list of installed base system, tracked by either collimator or gantry serial number.

This FMI applies to all CT and PET systems that utilize the VCT (LightSpeed 7X) or Helios (LightSpeed 1X-5X) collimators that have ordered replacement FRUs (Spare Parts) between Jan. 1<sup>st</sup>, 2006 to Sept. 1<sup>st</sup>, 2006.

FRUs include complete collimator assemblies and filter sub-assemblies with the following model numbers:

| Collimators           | Description             |
|-----------------------|-------------------------|
| 5130001 (2383390)     | VCT Failsafe Collimator |
| 5140001 (2214768)     | Helios Collimator       |
| 5140001-2 (2214768-2) | Helios 2 Collimator     |
| 5140001-3 (2214768-3) | Helios 3 Collimator     |
| 5140001-4 (2214768-4) | Helios 4 Collimator     |

| Collimator Filter Assemblies (sub-FRUs) | Description                    |
|-----------------------------------------|--------------------------------|
| 5148339                                 | VCT Filter Assembly            |
| 2126674                                 | Helios & -2 Filter Assembly    |
| 2291958                                 | Helios -3 & -4 Filter Assembly |

## Purpose

GE Healthcare is in the process of releasing a field inspection FMI for VCT and Helios collimators used on GE CT and PET/CT scanners. This inspection is being conducted as a result of an investigation into two isolated reports out of thousands of units, from GE manufacturing, regarding collimator internal shielding.

GE Healthcare inspection and test methods for collimators include signed-off collimator manufacturer control steps, Spherical Radiation Monitoring (SRM) sample testing of collimators, and system bay staging radiation leakage testing. These redundant processes effectively identified the issue prior to product release from GE CT manufacturing.

GE Healthcare analysis of these events determined that there is not a systemic issue with collimator design or manufacture.

GE Healthcare has not received any field reports related to anomalies in collimator internal shielding from customers, service engineers, or repair centers. Collimators used for service spares are subject to the same manufacturer control sign-off as well as SRM sample testing. The vast majority of collimators are processed through manufacturing with its additional testing. GE Healthcare does not believe any units in the field are affected.

GE is proactively performing this inspection FMI on a sample of installed service spare collimators as an additional step to reconfirm and document that no shielding anomalies are in this population.

## Prerequisite FMIs

| FMI No. | Title                                                            |
|---------|------------------------------------------------------------------|
| 25360   | Helios 221478-3 & 221478-4 Collimator X-ray Scatter Filter Bar   |
| 25367   | Helios 2214768 & 2214768 - 2 Collimator X-ray Scatter Filler Bar |

## Related FMIs

See *Prerequisite FMIs* section above.

## ***Furnished Materials***

None.

## ***Tool Requirements***

| Quantity | Description             |
|----------|-------------------------|
| 1        | 3 mm Hex key socket     |
| 1        | Phillips #2 screwdriver |

## ***Personnel Requirements***

|          |       |
|----------|-------|
| Resource | 1 FE  |
| Labor    | 2 hrs |
| Travel   | 1 hr  |

## ***Before You Begin***

Before beginning this procedure, please read the gantry safety information.

## ***Safety***



### **⚠ WARNING**

POTENTIAL FOR SHOCK.

VOLTAGE MAY BE PRESENT.

USE PROPER LOCKOUT/TAGOUT PROCEDURES BEFORE REPLACING COMPONENTS.

## INSPECTION PROCESS

This FMI is for inspection of CT and PET VCT and Helios Collimators with a set effectivity. Follow all applicable instruction and return the system to normal operation.

Upon discovery of missing lead or tungsten shields in the following instructions, follow the steps below.

|   |                                                                                                                                              |
|---|----------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | Contact:<br>Michael L. Bauer<br>MI&CT Service Engineering<br>262-312-7004<br>michael.l.bauer@med.ge.com                                      |
| 2 | Submit a PSR in iTrak to report site.                                                                                                        |
| 3 | Inform customer that the system is down until a replacement collimator is installed.                                                         |
| 4 | Hold defective collimator until given instructions as to where to send the collimator for evaluation.                                        |
| 5 | Order replacement collimator:<br><br>For VCT Collimators - Order FMI 25374<br>For Helios Collimators – Order through FE Tools (Spares Stock) |

## Procedures


It was in a VCT collimator that the missing pieces of lead and tungsten were first discovered. The VCT collimator is similar to the Helios collimator, which is made by same supplier. While there have been no similar problems reported with the Helios collimator, actions are being taken for the Helios collimators to be inspected as well to avoid any potential issues.

*Due to certain differences in configuration, separate inspection procedures follow for the Helios and VCT collimators. If the system of interest is a VCT, proceed to the VCT section of this instruction on Page 10. Otherwise perform the following steps.*

### I. Helios Collimator (LightSpeed 1X-5X) Inspection Procedure

This procedure describes the process of inspecting the Helios collimator assembly.

1. Shutdown system console / software.
2. Turn off system power and follow appropriate LOTO procedures.
3. Remove the gantry side, top, and front cover.
4. Position gantry so that the collimator is in the 6 o'clock position.
5. Engage gantry Rotation Lock mechanism.
6. Position a drop cloth below the collimator to collect any hardware that might be dropped.
7. Inspect for the presence of the Tungsten Filler Bar. (If found missing, return to prerequisite section and install FMI appropriate for this system.)







8. Follow the detailed procedure in the service manual for removing the Collimator Control Board (CCB) chassis and Collimator Filter assembly. The basic steps are these:
  - a. Remove Filler Bar from collimator frame to allow removal of CCB chassis cover
  - b. Remove CCB cover by removing the five pan head M4 mounting screws and set aside.
  - c. Remove the CAM A and B encoder cables (J1, J2& J3) including grounding ring terminals and washers from CCB and frame mount. Return the screws and washers to their original location after removing the ring terminals.
  - d. Slide the grommet on the CAM A and B encoder cables out of the frame mount and cut tie wraps to fully disconnect the encoder cables from the frame.
  - e. Remove the cam motor drive covers by removing the two pan head M4 mounting screws and set aside.
  - f. Disconnect the Cam A and B motor drive cables from motor drivers.
  - g. Remove the two hex socket M4 frame mounting screws and carefully slide the CCB chassis and Collimator Filter assembly from the collimator base.

**NOTE:** Do not damage the fragile aluminum, graphite, and copper filter assembly that extends into the collimator base during the filter assembly removal or installation.

9. With the filter assembly removed from the collimator base, check for the presence of the lead shield in the back of the cavity in the collimator base.



10. Inspect the filter assembly for the presence of lead shields, as shown in the following illustrations.



11. Follow the detailed procedure in the service manual for installing the Collimator Control Board (CCB) chassis and Collimator Filter assembly. The basic steps are these:

- a. Reinsert the CCB chassis and Collimator Filter assembly in the collimator base assembly. On the filter assembly, insert the two hex socket M4 frame mounting screws and torque to  $3 \pm 0.3$  N-m.

**NOTE:** Conversion Factor: 1 N-m = 1.356 ft-lb

- b. Connect the Cam A and B motor drive cables to motor drivers.
- c. Replace the motor driver covers on the motor mount assembly; insert the two pan head M4 mounting screws and torque to  $3 \pm 0.3$  N-m. Repeat for other motor drive cover.
- d. Connecting the Cam A and B encoder cables (J1, J2 & J3) to the CCB.
- e. Remove the screw and washer for the grounding ring terminals from the frame mount and use them secure the CAM A and B encoder grounding ring terminals and tighten.
- f. Replace the top cover of the CCB chassis using the five pan head mounting screws with spring washers and torque to  $3 \pm 0.3$  N-m.
- g. Install the filler bar using the two M4 frame mounting screws and torque to  $3 \pm 0.3$  N-m

12. Remove the drop cloth and follow the procedures for closing the gantry.

13. Proceed to Finalization section of this instruction on Page 12.


## II. VCT (LightSpeed 7X) Collimator Inspection Procedure

This procedure describes the process of inspecting the VCT collimator assembly.


1. Shutdown system console / software.
2. Turn off system power and follow appropriate LOTO procedures.
3. Remove the gantry side, top, and front cover.
4. Position gantry so that the collimator is in the 6 o'clock position.
5. Engage gantry Rotation Lock mechanism.
6. Position a drop cloth below the collimator to collect any hardware that might be dropped.
7. Follow the detailed procedure in service manual for removing the Collimator Control Board (CCB) chassis and Collimator Filter assembly. The basic steps are these:
  - a. Remove the CCB cover by removing the six pan head M4 mounting screws and set aside.
  - b. Remove the CAM A and B encoder and motor drive cables from CCB chassis.
  - c. Remove the four hex socket M4 frame mounting screws and carefully slide the CCB chassis and Collimator Filter assembly from the collimator base.

**NOTE:** Do not damage the fragile aluminum, graphite and copper filter assembly that extends into the collimator base during the filter assembly removal or installation.

8. With the filter assembly removed from the collimator base, check for the presence of the lead shields indicated in the following illustrations.



9. Inspect the filter assembly for the presence of lead and tungsten shields, as shown in the following illustrations.



10. Follow the detailed procedure in service manual for installing the Collimator Control Board (CCB) chassis and Collimator Filter assembly. The basic steps are these:
  - a. Reinsert the CCB chassis and Collimator Filter assembly in the collimator base assembly. On the filter assembly, insert the four hex socket M4 frame mounting screws and torque to  $3 \pm 0.3$  N-m.  
**NOTE:** Conversion Factor: 1 N-m = 1.356 ft-lb
  - b. Connect the CAM A and B encoder and motor drive cables to CCB chassis.
  - c. Replace the top cover of the CCB chassis using the six pan head mounting screws with spring washers and torque to  $3 \pm 0.3$  N-m.
11. Remove the drop cloth and follow the procedures for closing the gantry.
12. Proceed to Finalization section of this instruction

## ***Finalization***

1. Inspect gantry for loose hardware or tools left behind.
2. Disengage Rotation Lock mechanism.
3. Replace gantry covers.
4. Remove LOTO per detailed instructions in the Service Manual and power up system.
5. With covers in place, rotate the gantry to ensure that there is no mechanical interference between the collimator and gantry cover components. Adjust covers, if required (no cover adjustment expected).
6. Perform test scan to confirm system performance before returning system for patient scanning.
7. Remove document from customer site and close FMI RFS per regional procedures.