

**Restricted to General Electric Medical Systems
(GEMS) personnel only.**

Permission to use this Restricted Service Software and related documentation (hereinafter called the "Material") by persons other than GE Medical Systems employees is prohibited. This Material is the property of GEMS and contains information which is restricted to use by GEMS Manufacturing, Engineering and Service personnel. Any re-production in whole or in part is strictly prohibited. Neither this material nor any information contained herein shall be disclosed to anyone outside of GEMS Manufacturing, Engineering and Service personnel.

Effectivity:

LightSpeed RT

Purpose:

The purpose of this service note is to remind FE the existence of the extending pin on the cover. And inform them to use this extending pin (PN:2383599) during the assembly and the disassembly of the front gantry cover

Solution/Details:

There is a design change for the LightSpeed RT tilt dollies pins. An extending pin on each side of the gantry is placed on the Gantry cover of each side as indicated in Fig1 during shipment. When disassembling and/or assembling the front cover, this pin should be connected to the lower one of the dolly pin (refer to Fig2). Without this extending pin, the tilt dolly will not be able to close enough and is not useable. Related manual will be updated to indicate this change.



Fig 1



Fig 2

LI GANG

GEHW, CT Engineering