

FCT SERVICE NOTE

SN05102815a

Class A Service Information Enclosed.

Date: October 28, 2005

Subject: Error Code 290012000 is misleading

Effectivity: LightSpeed VCT, All Versions

Details: The following gesyslog error code entry has been causing confusion and driving unnecessary tube replacements.

```
Thu Oct 13 13:47:59 2005      Error: 290012000
Host: ig2                    Process: image_generation
File: MissedView.cpp         Line: 102
View corrections were applied in view 370.
exam/series/scan            = 3/2/3
```

Error code 290012000 is a status notification that Recon has made a correction to maintain image quality. This error is transparent to the Customer and has no effect on system performance.

This condition is generated when the IG reads in the scan data and checks the KV signal for each view. If the KV signal drops below a certain threshold, a correction is made by design. The low KV signal can be defined as a "High Voltage Disturbance".

Do not assume the tube is the source of the problem. A healthy tube will typically generate 2000 to 3000 of these events during its normal life.

Other potential sources for high voltage disturbances are the Jedi Generator, HV Tank, and cabling to the DAS Backplane.

- Service action might be required if Scan Aborts accompany this error.
- Diagnose the source of the Scan Aborts, not this error.
- Examine the log for Spit Max (SME=260118213) or Spit Ratio (SRE=260118214) errors.
If no SME or SRE errors exist in the gesyslog but the Spit Rate for the entire day (not individual exams) exceeds 1.0 spits per 10 scan seconds, a tube heatsoak can help.

Erwin Sulma
FCT Service Engineering