

## **FCT SERVICE NOTE**

### **06062207b**

**Class C information is Property of GE.**

**Date:** June 23, 2006

**Subject:** TCB Flash Download Patch for 06MW03.4 VCT M4

**Effectivity:** VCT systems (7.x systems) with 06MW03.4

**Problem:**

Some (approx. 5%) VCT M4 software upgrade sites have reported GTCB failure during the Flash Download process even though the workaround described in 7X Software Load Procedure rev 8 is correctly followed.

During the TCB flash download, Two Binary files (tcb.bin and tcbFpga.bin) are flashed to the GTCB. The GTCB failure occurs when the TCB flash download is interrupted for some reason, this depends on the timing, a mismatch of the versions of the Two Binary files. If the Two Binary files do not match then the GTCB will never boot up and is rendered inoperable.

The same failure could occur if GTCB was replaced and flash download was executed. Once GTCB is completely flashed to M4 version, this issue will never happen.

**Resolution:**

A software patch was prepared to resolve this TCB flash download issue.

**Please install the patch CD between the following Steps:**

**Direction 5116412-100 rev8 ( 7X Software Load Procedure )**

**1.3.2.14 Changing DARC Subnet Using the "reconfig" Command**

-----> **INSTALL PATCH CD HERE**

**1.3.3 CT Applications Start Up**

**Or when using the Service Engineering Web page Direction 5116424-800 rev2 (LS 7.X LFC Procedure - 06MW03.4 SW), install the patch CD between the following steps:**

**3.2.14 Changing DARC Subnet Using the "reconfig" Command**

-----> **INSTALL PATCH CD HERE**

**3.3 CT Applications Start Up**

After installing the patch, you can run the Flash Download normally. (No need to power off the GTCB to avoid flash download failure.)

- **Patch Install/Uninstall Procedure**  
Refer to the Appendix in this service note.
- **How to get the patch**  
The patch ISO file can be obtained from CT On Line SupportCentral.

GE Healthcare Functional CT (FCT) Class C eService Notes are copyright (c) to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.

Duplication and distribution of any class of FCT eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

[http://supportcentral.ge.com/bookupload/sup\\_book\\_default.asp?prod\\_id=15768&trng\\_doc\\_id=281360&vendor\\_flag=11&doc\\_cat\\_id=115255](http://supportcentral.ge.com/bookupload/sup_book_default.asp?prod_id=15768&trng_doc_id=281360&vendor_flag=11&doc_cat_id=115255)

Go to **LightSpeed 7.x (VCT64, VCT32) > LightSpeed VCT64**

Please follow the *Burning ISO File.doc* to download the ISO file and burn the CD.

## Appendix Patch Install/Un-install Procedure

### 1. Patch Installation

- 1) Insert the patch CD into the DVD Drive of the SCSI tower.
- 2) Shutdown application from common service desktop.
- 3) Open a Unix Shell.
- 4) Type the following:

```
su - <ENTER>
#bigguy <ENTER>
mount /mnt/cdrom1 <ENTER>
cd /mnt/cdrom1 <ENTER>
./install_patch -p VCT_M4_Flash_Download_patch-1-1.i386.rpm <ENTER>
```

(The following installation message appears.)

```
*****PRE-INSTALL*****
```

```
No pre-install activity required
```

```
*****PRE-INSTALL*****
```

```
*****INSTALL*****
```

```
Installing VCT_M4_Flash_Download_patch-1-1
```

```
##### [100%]
##### [100%]
```

```
*****INSTALL*****
```

```
*****POST-INSTALL*****
```

```
No post-install activity required
```

```
*****POST-INSTALL*****
```

```
cd <ENTER>
umount /mnt/cdrom1 <ENTER>
```

- 4) Eject and remove the patch CD from the DVD drive of SCSI tower.
- 5) Type the following and verify the Flash Download patch is installed:  
**rpm -qa|grep patch <ENTER>**

```
VCT_M4_Flash_Download_patch-1-1
patch-2.5.4-19
```

- 6) Type the following to exit from Unix shell:

```
exit <ENTER>
exit <ENTER>
```

- 7) Return to the LFC procedure.

Duplication and distribution of any class of FCT eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

## 2. Patch Un-Installation

**Note: This is just for your reference. You don't have to uninstall the patch unless any problem occurs on TCB Flash download.**

- 1) Insert the patch CD into the DVD drive of SCSI tower.
- 2) Shutdown application from common service desktop.
- 3) Open a Unix Shell, and type the following:

```
su - <ENTER>
#bigguy <ENTER>
mount /mnt/cdrom1 <ENTER>
cd /mnt/cdrom1 <ENTER>
./uninstall_patch <ENTER>
```

Following messages appear on the shell, select 1 to un-install the patch.

```
The installed patch are listed below.
Please selected the number associated with the patch to be un-
installed.
```

```
0 - exit
1 - VCT_M4_Flash_Download_patch-1-1
Enter a number: 1 <ENTER>
```

```
*****PRE-UNINSTALL*****
No pre-uninstall activity required
```

```
*****PRE-UNINSTALL*****
```

```
*****UNINSTALL*****
Uninstall VCT_M4_Flash_Download_patch-1-1
```

```
Replacing files to the state before the patch
```

```
Remove the /usr/g/patch/VCT_M4_Flash_Download_patch-1-1 and backup
files
```

```
*****UNINSTALL*****
```

```
*****POST-UNINSTALL*****
No post-uninstall activity required
```

```
*****POST-UNINSTALL*****
```

```
cd <ENTER>
umount /mnt/cdrom1 <ENTER>
```

- 4) Eject and remove the patch CD from the DVD drive of SCSI tower.
- 5) Type the following and verify the Flash Download patch is uninstalled:  
**rpm -qa|grep patch <ENTER>**

```
patch-2.5.4-19
```

GE Healthcare Functional CT (FCT) Class C eService Notes are copyright (c) to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.

Duplication and distribution of any class of FCT eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

- 6) Type the following to exit from Unix shell:  
**exit <ENTER>**  
**exit <ENTER>**
- 7) Open a new Unix Shell and type the following to start up application:  
**st <ENTER>**
- 8) Update firmware by selecting Service desktop -> Utilities -> Flash Download. Select Update in Flash Download GUI - tcb.bin and tcbFpga.bin files should be downloaded.

Takashi Uehara,  
Lead Service Integrator Global Table