

FCT SERVICE NOTE SN06101011A

Class A Service information enclosed.

Date: Oct 20, 2006

Subject: Limited Access / Small Room CT Installations & Service Changes

Effectivity: BrightSpeed 4/8/16, LightSpeed 5X HP100/16, and 7X VCT 32 and 64-Slice Systems. Available Q4 2006 for ship & install for the above products.

Introduction: CT systems may be installed into “Limited Access,” or Small Room environments to satisfy USA OSHA (Occupational Safety and Health Administration) All installation must meet the regulatory and service minimums to be installed. No Energized Left Side Service.

Minimum Room Size (Limited Access)

Please note **all** systems installed in the United States must comply with all Federal and local regulations. This includes room moves, room upgrades, and Gold Seal products. For installations outside the United States, **all** country-specific or local regulatory clearance requirements must be met. The requirements listed are the GE Healthcare minimums for servicing these products. Limited access is a system design for production products and **is not** a field **add-on** for existing products.

Regulatory Caution

Site prints are required for all system installations including upgrades, relocation and moves and each shall meet all regulatory requirements. System room layouts as shown on your site print shall meet all regulatory requirements as described in the pre-installation manual. This also includes placing equipment as shown on the site print. Customer equipment such as additional cabinets, not shown on the site print may void the caution statement, making the room non-compliant. Room measurements will be taken before the installation for to determine actual room size and regulatory compliance. The project manger of installation is reasonable for confirming that all requirements are met and the site is ready prior to the installation.

For minimum size rooms, where the gantry distance between the wall and side cover is between 14 in. – 27 in. (256mm x 686mm) the left side cannot be considered or used for egress. In this configuration, egress around the left side of the gantry will be restricted. Unobstructed egress route shall be provided around the table foot end or around the back of the gantry. If millwork is required or will be installed later, consider room placement locations that allows for regulatory compliance and clearances. Adding additional millwork may void regulatory and service compliance. One-man service is required for ALL systems installed. All rooms including minimum size rooms shall provide a unrestricted route to move service components to the gantry right side and provide space for the front and back cover removal. The project manger of installation is reasonable for confirming that these conditions are met.

Operational Caution

In this minimum room layout 14 in. to 27 in. (256mm x 686mm) the customer should consider workflow, access for patient care and critical care operations space requirements. Additionally, there may be limited equipment access on the gantry left side when loading patients, or when positioning patient equipment in the room between the gantry and the wall. Customer installation tasks are detailed in the products PIM chapters 1-4.

George Farrington SMD, Molecular Imaging & Computed Tomography Service Engineering

Suggested Room Size

This room configuration offers the most flexibility for future upgrades, has sufficient workspace, with space to add millwork and meet all regulatory and service requirements. This room would be compatible with most two-step future installations.

Recommended Room Size

This room configuration allows for some future upgrades, sufficient workspace, limited space to add millwork and meet all regulatory and service requirements. This room may be compatible with some two-step future installations.

Minimum Room Sizes A 35" – 28" (889mm x 711mm) and B 27" - 14" (256mm x 686mm)

This room configuration allows for no future upgrades, has limited workspace, no in room millwork, but meet all regulatory and service requirements. This room size is not compatible with two-step future installations.

Purpose:

The CT Gantry Left Side Limited Access Initiative provides the capability to improve sales opportunities by reducing the minimum room size for CT Systems. Limited Access adds left side siting flexibility while meeting all service and regulatory requirements. Doing so allowing a CT system to be sited in rooms with widths 22 in. (559mm) smaller than the current minimum room widths.

In the new minimum room configuration, left-side access and egress will be restricted for sites with less than 28 in (711mm) of left side clearance. Refer to your site's installation print for your room's detail.

Limited Access Changes

Pre-Installation

- All Installations Requires a Site Prints
- Limited Access B-Cat added to Print and order by the design center
B-Cats
- **B7800SA** Identifier Only See room "A" 35" to 28" (889mm x 711mm)
- **B7800SB** Identifier with parts Includes:
 - Installation Positioning Dolly
 - * Front Cover dollies
 - * Q4 Top Cover Tools
- US Systems Shipped Between Wk 42 – Wk 48 will be engineering supported.

Note: Limited Access is Not field Upgradeable as an add-on to existing systems!

Limited Access Service Requirements

- 06MW29.7 System Software Loaded
- Exceed & Hummingbird software on your service laptop
- Limited Access System Troubleshooting patch CD ROM
- Laptop service cable

Training:

Limited Access Changes Will Be Integrated Into Existing Product Training Classes

- ERT – T-Cons weekly during October 2006
- ERT Documents
 - **Regulatory Clearance Quick Reference Guide**
 - **Service Limited Access ZSE**
 - **Small Room Deliverables Doc**

New Service Tools:

Laptop service cable
Front Cover Service Dollies
Top Cover Removal Tool

George Farrington SMD, Molecular Imaging & Computed Tomography Service Engineering

Minimum room installation dolly

Revised Service Docs:

- IPL
- Pre-Installation
- Installation Manual
- System service CD (new part number)
- Patch CD ROM

Tips and Work-Around:

When loading the Limited Access Patch be sure to load Class "C" software before loading the Limited Access Patch. Loading the Limited Access Patch first will result in the service tools not being displayed on the common service desktop. Limited Access Service Tools will work in both class "A" and class "C".

"No Energized" Service Cut-In Dates

Limited Access Service Changes Applies To All New Manufactured Products*

[You can identify "Limited Access" systems by the yellow sticker on the TGPu cover.](#)

*VCT Family	After Wk 45 -06
*Pro 100 16	After Wk 45 -06
*Brightspeed Family	After Wk 48 -06
*RT 16 & Xtra 16	After Wk 07 -07

Products Not included

RT 4 Systems - **Must meet ALL of the current minimum requirements**

HP60 Systems - **Must meet ALL of the current minimum requirements**

Pro 80 Systems - **Must meet ALL of the current minimum requirements**

All Goldseal Systems - **Must meet ALL of the current minimum requirements**

Other Global Installation Requirements

- All systems installed will meet **ALL** minimum regulatory requirements.
- Limited Access changes applies to **All** service activities.
- **No Energized** Left Side Service applies to all production systems listed above and is not room size dependant.
- **All** CT Systems Shall be Installed with A LOTO compatible A1 disconnect.
- A1 Disconnect, PDU & UPS cannot be remote located outside of the CT scan suite. Adjacent equipment room on the same floor but next to the scan room or control room is acceptable.

New FRUs:

P/N Limited Access FRU	Part Description	Qty / System	Replaces P/N
5182284	Encoder, Axial Drive Assembly	1	2369747-2
2373436-7	Cable, Network Cable, 23000mm long, FE Laptop to Host Console	1	N / A
5176413	Clamping Lever, Cover Dolly	4	N / A
5179183	Cover Dolly Right, Limited Access	1	5101127
5179183-2	Cover Dolly Left, Limited Access	1	5101127-2

Planned Maintenance:

- No Changes

George Farrington SMD, Molecular Imaging & Computed Tomography Service Engineering