

## MI&CT SERVICE NOTE SN07012210a

Class C Service information enclosed.

**Date:** January 22, 2007

**Subject:** LightSpeed QXi Network Service Tool

**Effectivity:** VCT and GRE CT Systems using software 06MW3.4, 06MW3.5 or 06MW29.7

**Issue:**

The issue occurs when networks contain systems running newer software 06MWxx.x AND LightSpeed QXi systems running previous software versions. When systems running 06MW03.4, 06MW03.5 or 06MW29.7 software boot up, the LightSpeed QXi systems running the previous software can lock-up, QXi system being on the same hospital network.

**Cause:**

This issue was originally fixed by the script fix\_irix\_sd, however, with the release of the 06MWxx.x software, a sort command was introduced which made the fix\_irix\_sd script ineffective.

**Solution:**

A software tool has been developed to fix the issue. The tool modifies the fix\_irix\_sd script by reversing the order of the startup and shutdown of network interfaces and eliminating the arpings sent out when starting the network.

The part number for the software tool is **5194766** and can be ordered from Coakley Tech.

The service tool is also available via CT Online Support Central and can be downloaded for installation in the system. The following is the link to the iso file:

[http://supportcentral.ge.com/bookupload/sup\\_book\\_access.asp?book\\_id1=3400505&chap\\_id=3019134&book\\_id2=347853&doc\\_id1=281360&prod\\_id1=15768](http://supportcentral.ge.com/bookupload/sup_book_access.asp?book_id1=3400505&chap_id=3019134&book_id2=347853&doc_id1=281360&prod_id1=15768)

Link to the read me file is:

[http://data.supportcentral.ge.com/books/15768/B19097/doc\\_3400509.pdf](http://data.supportcentral.ge.com/books/15768/B19097/doc_3400509.pdf)

N. Srinivasan  
MI & CT Service Engineering

**Manager's Signature:**

