

GE Healthcare Functional CT (FCT) Class C eService Notes are copyright (c) to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.

Duplication and distribution of any class of FCT eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

FCT SERVICE NOTE SN07012514a

Class C information is Property of GE.

Date: January 25, 2007

Subject: VCT DIFB slot 1 issue

Effectivity: All VCT systems

Details: There is a firmware issue that interferes with the systems ability to update the flash code in the DIFB in slot 1. When a new DIFB is installed in slot 1 that doesn't have the latest firmware code, flash download fails. This issue is known to occur in 06W03.4 (VCT M4) and 06MW29.7 software.

Resolution:

For any DIFB issue with slot 1, move DIFB from slot 3 to slot 1 and install the new DIFB in slot 3. This bypasses the flash download issue for slot 1 and flash download will work properly. The firmware issue is resolved in all software releases after 06MW29.7.

Steve Merritt
FCT Service Engineering