

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.  
Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

**MI&CT ESERVICE NOTE  
SN SN07092012a**

**Class C information is Property of GE.**

This short-term informational note is invalid 6 months after the issue date.

**Date:** September 20, 2007

**Subject:** Service Process for IVY Cardiac Trigger Monitors

**Effectivity:** Ivy 3150-A Cardiac Trigger Monitor E8007RD and all MICT Accessories

**Details:** Field Service is unable to find information on how to handle the service activities for the Ivy 3150-A Cardiac Trigger Monitor and associated MICT Accessories.

**Resolution:** See the information enclosed for details concerning the IVY Trigger Monitor. This information will be added to the Accessories and Supplies Catalog upon its next release.

**Affected Service Publication:** None

**PQR/SPR #:** PQR 13077008

Please See Next Page for Details

## **Accessories Service Process:**

### **1. How does an FE find out if an Accessory item is covered under a service contract?**

If an issue with an accessory comes up and the FE needs to determine if it is under contract or not, the FE is responsible for contacting the Local Sales and Marketing Team and either the Local Sales and Marketing Team take over the conversation with the customer or the FE works with the Local Sales Team to determine.

1. Contract Status of the Accessory
2. Warranty information defined by the contract and by the Accessory Supplier. The information regarding the Warranty from the Accessory Supplier can be found in the product manual that the Customer would have or in the Product Solution Database.

Based on #1 and #2 above, the FE would communicate this information to the customer and based on the customers response proceed with a repair action if requested.

### **2. For further Service Information, the Accessories Service Plan is defined as follows:**

- a. Refer to the "Troubleshooting Section" of the product manual for further Service information. This manual is the property of the customer. You may need to consult with the Customer or locate the product by some other means, such as Web Search, PSDB, or contacting Accessories at 1- 877-248-5355.
- b. If unable to determine the problem, contact the Accessory Supplier(Contact information is in the Product Manual or by working with Accessory Inside Sales Team at 1-877-248-5355) for further assistance in identifying the problem.
- c. Consult the Problem Solution Database for further information

### **3. Escalation Process if you have not resolved the issue with the Accessory Product:**

- a. **LEVEL 1:** After 4 hours elevate to OLC Support. If the problem is not resolved, provide status and your action plan to your customer and escalate to Level 2.
- b. **LEVEL 2:** After 8 hours elevate to Accessories calling 1-877-248-5355. If the problem is not resolved after requesting Accessory Product Manager involvement, provide status and your action plan to your customer and local sales leadership, and escalate to Level 3.
- c. **LEVEL 3:** If problem is not resolved, generate a CSO and provide status and your Action plan to your customer. Work with CSO team to resolve the issue.

### **4. Specific Information for the Ivy 3150 or 3150-A Cardiac Trigger monitor**

- a. Replacement Monitors can be ordered through FEMC (FE TOOLS)
- b. Replacement parts are available through FEMC or as Accessories (see next page for non-contract customers)
- c. Manufacture Warranty Coverage is 1 year from the date of purchase, which is covered in the operator's manual on page iii.

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.  
Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.



## 3150-A FRU 5235192

- Cardiac Trigger monitor with P Lock algorithm to help prevent false triggers
- 3150-A monitor is mandatory for Snap Shot Pulse option



## Patient Cable E8007RE

- Low noise IVY ECG cable
- Fully shielded, AAMI 6 pin connector, Length: 10 ft(3.3m)



## Set of 3 lead wires E8007RF or RH

- Radio translucent and Fully shielded
- E8007RF- 30 inches long E8007RH – 24 inches long



## Radio translucent Electrodes E8007RG

- Radio translucent, latex free 600/case, bags of 30 (10 pouches of 3)
- Gel: 10% KCL wet gel
- Snap and Sensor: Ag/AgCl plated plastic



## ECG Chart Paper E8007RJ

- Case of 100 rolls (10 packs of 10 rolls)
- Thermal paper, No Grid, 50 mm x 80 feet



## Ivy Monitor Roll Stand E8007RK

- Load Capacity up to 12 lbs
- Height: 44 inches, 21 inch diameter aluminum base
- Mounting plate, basket, and handle with padded grip

**GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.**

**Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.**

Rob Schaefer

Accessories Global Product Manager

MI&CT Engineering Manager Reviewer: \_\_\_\_\_ Mark Nyman \_\_\_\_\_

| Service Note Number | Service Note Title                          | Unit # or Department | Signature    | Release Date |
|---------------------|---------------------------------------------|----------------------|--------------|--------------|
| SN07092012a         | IVY 3150-A FRU and Accessories Service Plan | Accessories          | Rob Schaefer | 09/20/07     |