



GE Medical Systems

Technical Publication

Direction 2333263-100

Revision 0

Octane I Computer Cleaning Service Note T-1445

THIS MATERIAL IS THE PROPERTY OF GEMS AND CONTAINS INFORMATION, WHICH IS RESTRICTED TO USE BY GEMS MANUFACTURING, ENGINEERING AND SERVICE PERSONNEL. ANY REPRODUCTION IN WHOLE OR IN PART IS STRICTLY PROHIBITED. NEITHER THIS MATERIAL NOR ANY INFORMATION CONTAINED HEREIN, SHALL BE DISCLOSED TO ANYONE OUTSIDE OF GEMS MANUFACTURING, ENGINEERING AND SERVICE PERSONNEL.

*Copyright © 2002 by General Electric Company, Inc
All rights reserved*



Page Left Intentionally Blank

LEGAL NOTES

GE EXCLUSIVE USE ONLY

Restricted To General Electric Medical Systems (GEMS)

This Material is the property of GEMS and contains information, which is restricted to use by GEMS Manufacturing, Engineering and Service personnel. Any reproduction in whole or in part is strictly prohibited. Neither this material nor any information contained herein shall be disclosed to anyone outside of GWMS Manufacturing, Engineering and Service personnel.

If you are not a GE Medical Systems Employee, you must exit this material now.

TRADEMARKS

All products and their name brands are trademarks of their respective holders.

OMISSIONS & ERRORS

GE personnel, please use the GEMS CQA Process to report all omissions, errors, and defects in this publication.

COPYRIGHTS

All Material Copyrighted(c) 2002 by the General Electric Company Inc. All rights reserved.

REVISION HISTORY

REVISION HISTORY

Revision	Date	Reason for change
0	02/07/02	Initial Release, CQA10110931

TABLE OF CONTENTS

Legal Notes	3
GE EXCLUSIVE USE ONLY	3
TRADEMARKS.....	3
OMISSIONS & ERRORS.....	3
COPYRIGHTS.....	3
REVISION HISTORY	3
Table of Contents.....	4
Octane I Computer Cleaning.....	5
PURPOSE OF SERVICE NOTE.....	5
BEFORE YOU BEGIN	5
REQUIRED TOOLS	5
CLEANING PROCEDURE	5

OCTANE I COMPUTER CLEANING

PURPOSE OF SERVICE NOTE

This Service Note describes precautions and ordered steps for removing dust/dirt from an Octane I computer. It has been reported; dirt and dust accumulation can restrict airflow in an Octane I. To remove these obstructions, special attention must be paid to the handling of computer parts.

BEFORE YOU BEGIN

1. The Octane computer must be removed from the console and cleaned in a location suitable for the tasks described below. A suitable place provides:
 - ~~///~~ A level surface for disassembly
 - ~~///~~ A static free work environment
 - ~~///~~ Adequate ventilation
 - ~~///~~ Good lighting
2. Arrangements should be made for the disposal of all cleaning materials, in accordance with GE policy.

REQUIRED TOOLS

- ~~///~~ Safety Glasses
- ~~///~~ Anti-Static Bags (Large enough to accommodate largest parts)
- ~~///~~ Grounded Static-free Wrist Strap
- ~~///~~ Grounded Static-free Work Mat
- ~~///~~ Compressed Air (Aero Duster Part # MS-222N)

CLEANING PROCEDURE

The following steps outline the tasks required to clean an Octane I computer. Consult the your System Service Manual, if you have any difficulty performing the tasks described below.

1. Shutdown CT applications and Operating System software
2. Switch off power to the console. We suggested that you LOTO power.
3. Remove front cover
4. Rotate computer and disconnect all cables.
5. Remove computer and take to work area.

**OBSERVE THE FOLLOWING PRECAUTIONS BEFORE BEGINING:**

 **Wear safety glasses**



 **Wear and use ESD protective equipment. Always wear a grounded ESD strap.**

 **Do not touch the pads of the compression connector. Whenever a module or board is not in the chassis, place the protective cap on the compression connector. A compression connector should never need to be cleaned if you keep the protective cover on whenever the module or board is not in the chassis.**

 **Before laying a part on a surface, make sure that the surface is free of dust, lint, powder, metal filings, oil, and water.**

 **Do not blow dust, dirt, or powder anywhere near a board when it is not inside an antistatic bag.**

 **Do not clean any contacts with aerosol or liquid cleaners. These substances can cause irreparable damage to the connector's surface.**



 **Allow sufficient time for parts to cool before handling. Wait five minutes after power is off before you begin.**

6. Disassemble Computer
 - a. Remove Front Bezel.
 - b. Remove Disk Drives
 - c. Remove System Module
 - d. Remove XIO module.
 - e. Remove PCI Module.
 - f. Remove Power Supply.
 - g. Remove Front plane.
 - h. Remove Power Supply Fan
7. Use compressed air to blow dirt and dust from computer components.

 **Do not touch circuit board contacts. Potential for damage is too great.**
8. Assemble Computer in reverse of order shown in step 6.
9. Place computer back in console
10. Reconnect cables
11. Replace covers
12. Turn on power
13. Start software.

END DOCUMENT

g

GE Medical Systems

P.O. BOX 414

MILWAUKEE WISCONSIN 53201-0414 U.S.A.