



GE Medical Systems

gemedical.com

Technical Publication

Direction 2341063-100

Revision 0

GE Medical Systems

LightSpeed MDAS DCB Corrupted Firmware

Service Note

T-1450

Restricted To General Electric Medical Systems (GEMS)

This Material is the property of GEMS and contains information that is restricted to use by GEMS Manufacturing, Engineering and Service personnel. Any reproduction in whole or in part is strictly prohibited. Neither this material nor any information contained herein shall be disclosed to anyone outside of GEMS Manufacturing, Engineering and Service personnel.



LEGAL NOTES

GE EXCLUSIVE USE ONLY

Restricted To General Electric Medical Systems (GEMS)

This Material is the property of GEMS and contains information that is restricted to use by GEMS Manufacturing, Engineering and Service personnel. Any reproduction in whole or in part is strictly prohibited. Neither this material nor any information contained herein shall be disclosed to anyone outside of GEMS Manufacturing, Engineering and Service personnel.

If you are not a GE Medical Systems Employee, you must exit this material now.

TRADEMARKS

Adobe, the Adobe logo, Acrobat, the Acrobat logo, Exchange, and PostScript are trademarks of Adobe Systems Incorporated or its subsidiaries and may be registered in certain jurisdictions.

Microsoft is a registered trademark and Windows is a trademark of Microsoft Corporation.

All other products and their name brands are trademarks of their respective holders.

COPYRIGHTS

All Material, Copyright © 2002 by General Electric Company, Inc. All rights reserved.

REVISION HISTORY

Rev.	Date	Reason for change
0	02/25/02	CQA 1021530

LightSpeed MDAS DCB Corrupted Firmware

Effectivity:

LightSpeed Systems with MDAS Hardware and 2002 LightSpeed Apps ver 303.1_H3M4 software.

Purpose:

To identify root cause and aid the Field Engineer in quick recovery of a down system situation.

Problem Statement:

- LightSpeed Plus 2.1 software was inadvertently loaded on a LightSpeed Ultra system with MDAS hardware.
- The net result is the MDAS DCB board was flashed with SDAS DCB code.
- This results in the MDAS DCB board failing its power up self test. No further communication is possible with the MDAS and hardware reset will fail.
- No useful error information is presented.
- Replacement DCB boards will Flash without error and the same condition exists.
- After reloading the correct software the corrupted DCB will not flash and the same condition exists.

Solution:

There are 2 possible actions the Field Engineer can take. Both require an LFC with the correct software.

- A.) Have the OLC download a script to flash the corrupted MDAS DCB.
 - i.) Requires InSite Connectivity.
 - ii.) Recovery with minimum down time and no material cost.
- B.) Replace MDAS DCB and follow the specific instruction in procedure "B" on the next page.
 - i.) Cost of DCB incurred against your System.
 - ii.) Additional down time while waiting for parts.

Procedure: A

Comment: Requires InSite Connectivity for script download.

- 1.) Load correct version of software.
Helios_H2_O2LESS_OC_Product 02/15/2002 LightSpeed Apps ver 303.1_H3M4
- 2.) Ensure the correct hardware configuration is selected in Reconfig screens.
- 3.) With Applications up, run Flash Download. Wait for completion... it will fail.
- 4.) Execute *dasInstall* from a shell window.
- 5.) This was downloaded by the OLC.
- 6.) CTUSER prompt: *cd /tmp*
- 7.) CTUSER prompt: *dasInstall*
- 8.) Script output will indicate "Loading ApplicationOK"
"Waiting for node resetOK"
Script completes after the node reset OK.
- 9.) Perform Flash Download a second time.
- 10.) Perform Scan Hardware Reset. Verify successful completion.
- 11.) Perform test scan in New Patient.

Procedure: B

Comment: Requires new Replacement MDAS DCB Board Part # 2253562



NOTICE DO NOT PERFORM FLASH DOWNLOAD with the wrong software loaded. This is how you got here in the first place!

- 1.) Load correct version of software.
Helios_H2_O2LESS_OC_Product 02/15/2002 LightSpeed Apps ver 303.1_H3M4
- 2.) Ensure the correct hardware configuration is selected in Reconfig screens.
- 3.) Bring system to Applications level.
- 4.) Install new uncorrupted MDAS DCB board.
- 5.) Perform Flash Download.
- 6.) Perform Scan Hardware Reset. Verify successful completion.
- 7.) Perform test scan in New Patient.

Note: You can perform an LFC while you are awaiting the arrival of the new MDAS DCB board.



GE MEDICAL SYSTEMS

***GE MEDICAL SYSTEMS-AMERICAS: FAX 414.544.3384
P.O. BOX 414; MILWAUKEE, WISCONSIN 53201-0414, U.S.A.***

***GE MEDICAL SYSTEMS-EUROPE: FAX 33.1.40.93.33.33
PARIS, FRANCE***

GE MEDICAL SYSTEMS-ASIA: FAX 65.291.7006