GE Healthcare

Technical Publications

Direction 5198939-100
Revision 4

BrainwavePA 1.6 for AW Installation/Service Manual

Copyright© 2007 by General Electric Company Inc.
All Rights Reserved

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation *Damage in shipment* written on **all** copies of the freight or express bill **before** delivery is accepted or signed for by a GE Healthcare representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier immediately upon discovery or within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14-day period.

To file a report:

- Call 1-800-548-3366, and use option 8.
- Fill out a report at <http://egems.med.ge.com/edq/home.jsp>.

Contact your local service coordinator for more information on this process.

Rev. 09/17/2006

OMMISSIONS & ERRORS

Customers, contact your GE Healthcare Sales or Service representatives. GE personnel, use the GEMS iTrak Process to report all omissions, errors, and defects in this publication.

WARNING

(EN)

- THIS SERVICE MANUAL IS AVAILABLE IN ENGLISH ONLY.
- IF A CUSTOMER'S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER'S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR, OR PATIENT FROM ELECTRIC SHOCK, OR FROM MECHANICAL OR OTHER HAZARDS.

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REVISION HISTORY

REV	DATE	PRIMARY REASON FOR CHANGE
1	January 15, 2007	Initial Version
2	June 20, 2007	Added Korean warning, changed cd/ to cd / in section 2-2-2, updated FRU number
3	Nov 07,2007	Updated rev # and FRU
4	April 30, 2008	Changed title and revision number to BrainWavePA 1.6 for AW Installation/ Service Manual; Revision 4. Updated Damage in Transportation and Warnings sections. Modified Introduction (AW models) and License Key Number section (URL)

TABLE OF CONTENTS

DAMAGE IN TRANSPORTATION	2
OMMISSIONS & ERRORS.....	2
1- INTRODUCTION	12
1-1 BrainWavePA Advanced Application “Clinical” Summary	12
2- PRELIMINARY REQUIREMENTS	12
2-1 License Key Number	12
2-2 Save Configuration	12
2-2-1 Reading Config Media → PC.....	13
2-2-2 Reading Config Media → AW.....	13
3- INSTALLATION PROCEDURE	13
4- AW AND BRAINWAVEPA ERRORS	13
5- ADVANCED SOFTWARE APPLICATION OPERATION AND DEMO EXAM INFORMATION	14
6- SPECIAL CONSIDERATIONS.....	14
7- FRU LIST.....	14

1- INTRODUCTION

The following manual gives a brief overview of the installation process for the BrainWavePA 1.6 advanced application software option. The BrainWavePA 1.6 application is designed to run on the AW4.3 and forward AW platforms only. Note that the previous version of application BrainWavePA 1.5 provides the same functionality but is designed to run on AW4.2 and AW4.2P platforms only. This guide makes references to all the appropriate standard Advantage Workstation service manuals for more detailed information about installing an advanced application software option on the AW. The current revision of the appropriate AW Service Manual should be used when installing this software. Unique installation processes for BrainWavePA and NOT found in the AW service manual can be found in this document under the section entitled SPECIAL CONSIDERATIONS.

1-1 BrainWavePA Advanced Application “Clinical” Summary

BrainWavePA is the post processing and visualization software for functional Magnetic Resonance Imaging (fMRI) data collected with Brainwave Real Time (RT). BrainWavePA allows the user to interactively view and edit fMRI data.

This tool allows the transformation of the three-dimensional (3D) volume fMRI data into an image. The user can then use volume rendered data to interactively view the fMRI images. The key advantages of BrainWavePA include:

- ❑ Elaborate processing of fMRI data sets (motion correction and alignment with segmented data).
- ❑ 3D activation display and display of fMRI data from multiple paradigms against the rendered high-resolution 3D volume.

2- PRELIMINARY REQUIREMENTS

2-1 License Key Number

The BrainWavePA option requires a license key # based upon the workstation's Host ID. This key number may be provided by:

- ❑ Paperwork shipped with the option that will include one or more 8 to 10 digit license key number(s) – created based on the workstation's/application's FDO and the workstation's Host ID/License ID.
- ❑ The software package may include a portable media (i.e. floppy disk) with the license key information on it.
- ❑ The license key information will also be available online at the eLicense website link by searching under the order number or the workstation's Host ID/License ID: <http://elicense.gehealthcare.com>.

All of these elements may exist (hard copy license #(s), media, and eLicense data), and at least one of them is necessary to load the software.

2-2 Save Configuration

Create an AW 'CONFIG' media before and after the application software installation for use as needed throughout this and other AW service activities.

2-2-1 Reading Config Media → PC

Configuration information can be read on a PC by viewing the config.txt file on the media using text readers like WordPad, Internet Explorer or MS Word (MS NotePad may not work).

2-2-2 Reading Config Media → AW

To read a Save Configuration using a floppy disk insert the floppy and type the following in a c shell:

- **su**
- Password '**operator**'
- **mount /mnt/floppy**
- **cd /mnt/floppy**
- **more config.txt**
- **cd /**
- **umount /mnt/floppy**

To read a Save Configuration using a CDROM disk insert the CD and type the following in a c shell:

- **su**
- Password '**operator**'
- **mount /mnt/cdrom**
- **cd /mnt/cdrom**
- **more ("backup directory name")**
- **cd /**
- **umount /mnt/cdrom**

3- INSTALLATION PROCEDURE

The 'Standard' AW Software Options Installation – Easy Install process applies. For ongoing process revision viability, the details will not be printed here, but they can be found in the AW Service Manual.

The current revisions for this manual can be viewed online at: <http://aw-ib.euro.med.ge.com>.

1. Open up the appropriate Advantage Workstation (AW) Service Manual.
2. Proceed to Job Card entitled – Software Options Installation.
3. Follow the procedure using the "Easy Install" method.
4. When complete, perform a save configuration.

4- AW AND BRAINWAVEPA ERRORS

BrainWavePA and AW System errors can be examined and investigated in either the usual AW error logs – such as /var/adm/messages, and/or ~sdc/logfiles/easyInstallLog – and also in the application's dedicated logfile:

/export/home/sdc/logfiles/NeuroNav.log.

5- ADVANCED SOFTWARE APPLICATION OPERATION AND DEMO EXAM INFORMATION

To test the functionality of the software, first you must load the DEMO exam CD shipped with the product. The 'Standard' AW Software Demo Exam Installation process applies. For ongoing process revision viability, the details will not be printed here, but they can be found in the AW Service Manual.

The current revisions for this manual can be viewed online at: <http://aw-ib.euro.med.ge.com>.

1. Open up the appropriate Advantage Workstation (AW) Service Manual.
2. Proceed to Job Card entitled – Advantage Workstation Demo Exam Reinstallation.
3. Follow the procedure to load the Demo Exams from the CD provided.

To test this application, after the application is installed, logout and log back into the AW (as suggested by the messages after the installation). Then select any data set and click on BrainWavePA in the "Applications" menu. The BrainWave window will open (to the left side of the screen). This ensures that Brainwave is installed correctly. Click on the Quit button to exit BrainWavePA.

6- SPECIAL CONSIDERATIONS

There are no special considerations to deal with when installing this advanced application option.

7- FRU LIST

The replacement parts for this option are as follows:

Part # 5213094-4 — BrainWavePA 1.6 Software and Documentation CD-ROM