SERVICE MANUAL

Invivo Corporation
3545 SW 47th Avenue
Gainesville, FL 32608

Document TR0528S

Revision 03

GE 1.5T 32-Channel Cardiac Coil

GE Part Number: M7000LB

Invivo Part Number: 800528

Manufactured By:

© Invivo Corporation, 2009

Distributed By:

GE Healthcare

Damage in Transportation

All packages should be closely examined at time of delivery. If damage is apparent, have notation “**damage in shipment**” written on **all** copies of the freight or express bill **before** delivery is accepted or “signed for” by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier **immediately** upon discovery, or in any event, within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this **14** day period.

File a report with

- Call 1–800–548–3366 and use option 8.
- Fill out a report on <http://egems.med.ge.com/edq/home.jsp>
- Contact your local service coordinator for more information on this process.

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TABLE OF CONTENTS

SECTION 1 – INTRODUCTION.....	9
1.1 Product Identification and Shipping List	9
1.2 Shipping Contents.....	10
1.3 Compatibility	10
1.4 Related Documentation.....	10
1.5 Environmental Requirements	10
SECTION 2 – INSTALLATION PROCEDURE	11
2.1 Personnel Requirements.....	11
2.2 Tools and Test Equipment.....	11
2.3 Safety – Not applicable.....	11
2.4 Required Conditions.....	11
2.5 Procedure.....	11
2.6 Finalization	15
SECTION 3: 32CH CARDIAC COIL MCQA TEST	16
3.1 Tools Required	16
3.2 MCQA Setup	16
3.3 MCQA Test	18
3.4 Periodic Quality Assurance Check	19
SECTION 4: TROUBLESHOOTING	19
4-1 Mechanical Hardware Check	19
4-2 Electrical Hardware check.....	19
4-3 Coil ID Troubleshooting Tips	21
SECTION 5 – MAINTENANCE	22
5-1 Coil Care.....	22
SECTION 6 – RENEWAL PARTS	22
6.1 Field Replaceable Units	22
6.2 Other Replaceable Accessories	23

ATTENTION!

A GE field engineer **MUST** perform installation of this coil. Please contact your field engineer to arrange for installation.

Attempts to perform installation of this coil without proper system configuration will result in system errors.

SECTION 1 – INTRODUCTION

1.1 Product Identification and Shipping List

This is a service manual for the GE 1.5T 32-Channel Cardiac Coil. The 32-channel Cardiac Coil consists of two 16-channel coil halves (anterior & posterior). Each half is uniquely identified by using the labels below:


Illustration 1: Cardiac Coil Labels

Illustration 2 shows images of each half of the 1.5T 32 Channel Cardiac Coil


Illustration 2

1.2 Shipping Contents

SHIPPING LIST BOX 1– TABLE 1-2-1

Description	Invivo Part #	GE Part #	Qty
1.5T 32Ch Cardiac Coil - Anterior	4535-302-48071	5366664-2	1
1.5T 32Ch Cardiac Coil - Posterior	4535-302-48081	5366664-3	1
32Ch Cardiac Coil - Phantom Positioner (2 pieces)	4535-302-66001	5366627-4	1
32Ch Cardiac Coil - Coil Definition Disk	N/A	5366627-5	1
32Ch Cardiac Coil - Posterior Pad, Bottom	4535-302-58891	N/A	1
32Ch Cardiac Coil - Posterior Pad, Top	4535-302-58901		1
32Ch Cardiac Coil - Head Pad	4535-302-58911	N/A	1
32Ch Cardiac Coil - Foot Pad	4535-302-58921	N/A	1
32Ch Cardiac Coil - Anterior Pad	4535-302-60711	N/A	1
32Ch Cardiac Coil - Operator Manual Insert	4535-302-85131	N/A	1
32Ch Cardiac Coil - Service Manual	TR0528S	N/A	1

1.3 Compatibility

This coil is compatible with the following hardware configurations:

Discovery MR450 and Optima MR450w 1.5T Systems

1.4 Related Documentation

- 5500110 Discovery MR450 1.5T Service Methods DVD
- 5670000 Optima MR450w 1.5T Service Methods DVD
- 1.5T & 3.0T 32-Channel Cardiac Coil Operator's Manual, 4535-302-58991

1.5 Environmental Requirements

Storage

Store the coil in an air-conditioned scan room or equipment room. To store the coils, a storage space of 48.3 cm/19.0 in (width) x 16.5 cm/6.5 in (depth) x 101.6 cm/40.0 (length) is required.

Dimensions

Coil (with posterior pad):	48.3 cm x 16.5 cm x 101.6 cm	(19.0 in x 6.5 in x 40.0 in)
Phantom Positioner:	50.6 cm x 50.6 cm x 17.7 cm	(20 in x 20 in x 7 in)

Weight

Coil half:	2.95 Kg	(6.5 lb)
Phantom Positioner:	0.68 Kg	(1.5 lb)

SECTION 2 – INSTALLATION PROCEDURE

2.1 Personnel Requirements

Required Persons	Procedure	Finalization
1	15 mins	45 minutes

2.2 Tools and Test Equipment

Item	Qty	Part#
32Ch Cardiac Coil – Coil Definition Disk	1	5366627-5
Phantoms and phantom positioners are needed for the MCQA test.	N/a	See section 3

2.3 Safety – Not applicable

2.4 Required Conditions

Condition
For 450: The minimum software build to install this coil is 20.1 M3 or 20.1_IB1
For 450w: The minimum software build to install this coil is 21.0 M3

2.5 Procedure

NOTE: Do not plug in the coil and use AutoCoil to install the coil. This method will not install the coil protocols.

- 2.5.1 Make sure that the latest software build is loaded on the system prior to installing the Cardiac coil. The 32Ch Cardiac protocols will only install if the required conditions are met.
- 2.5.2 Check to confirm if the 32chCardiac coil is already installed on the system software by opening the Service Browser under the Service Desktop Manager. In the Service Browser click on the Configuration tab, then the Coil Database Explorer under the configuration folder on the left.
- 2.5.3 Select “Click here to start this tool”. The Coil DB Explorer window will open and under the "Coil Name" tab there is a "Clinical Coil Name" control. Click on this and there will be a dropdown list of all the logical coil names. There should be three names that start with "C-GE_32Ch Cardiac". This indicates the coil is present in the Coil DB. Reference Illustration 3. If the coil is not present on the list proceed to step 2.5.6.


Illustration 3: Confirmation of coil configurations installed

- 2.5.4 The 32Ch Cardiac coil also has two coil protocols that must be on the system. Open the Protocol Management tab. Confirm the coil protocols are installed by viewing the “Chest” protocols under the GE Protocol Library. Two protocols that end in “- 32ch Cardiac” must be present as shown in Illustration 4.

If protocols are not present and the coil CD was just installed one may have to exit the “Chest” protocols by selecting another area such as “Abdomen” and then reselecting “Chest” for the two protocols to appear.


Illustration 4: Confirmation of 32ch Cardiac coil protocols installed

- 2.5.5 If the cardiac coil configurations and protocols are on the system skip to the finalization steps in section 2-6.
- 2.5.6 To install the cardiac coil configurations and protocols, insert the Cardiac Coil Definitions Disc into the default (typically upper) host DVD drive.
- 2.5.7 Open a C-Shell from the Service Desktop manager and type the following commands:

su

```
<root password>  
mount /mnt/cdrom  
/mnt/cdrom/autorun
```

- 2.5.8 A Coil Data Disk Install window will open. Type “i” to install the CD contents.
- 2.5.9 When the contents of the CD installation is complete, the cursor will appear and you will see the message as shown in Illustration 5. Press enter to exit the Coil CD installation.


Illustration 5: CD Contents Installation Complete

- 2.5.10 Type <eject> in the original c shell. Remove the CD from the host DVD drive.
- 2.5.11 To confirm the coil has been properly installed from the coil CD follow steps 2.5.2 – 2.5.4.

2.6 Finalization

- 2.6.1 Perform a SaveInfo under Guided Install to make the changes permanent.
- 2.6.2 Perform the Multi-Coil Quality Assurance (MCQA) test with this coil. Refer to Section 3.
- 2.6.3 Complete Product Locator information for serialized components. Failure to do this may result in failure of your site to receive future FMI's. Instructions are below.
 - Go to GIB: http://gib.gehealthcare.com/gib/gib_MainMenu.jsp
 - After the log in screen, click on [Global Installed Base entries]
 - Under "Serial Changes to Installed Base", enter the system ID of the system on which this coil is being installed. Click [Submit]
 - Click on [New Entries]
 - For both the Anterior and Posterior portions of the coil, enter the Model (Part number), Serial number, Installation Date and "A", then click [Submit].

SECTION 3: 32CH CARDIAC COIL MCQA TEST

Follow this procedure to perform the MCQA test. There are two phantom options that can be used for MCQA on the MR450 system. Only the unified phantom can be used on MR450w.

3.1 Tools Required

TOOLS REQUIRED - TABLE 3-1-1

Description	GE Part #	Qty
Legacy TL Phantom Brick (Not an option on MR450w) Or Unified TL Phantom Brick	2264744-8 Or 5343347	1
32Ch Cardiac Coil - Phantom Positioner (2 pieces)	5366627-4	1

3.2 MCQA Setup

The MCQA tool will be run separately for the posterior and the anterior coils. The coil positioning and the phantom setup for the MCQA is the same for both the coil components.

TABLE 3-2-1 Steps for Coil and phantom setup (Anterior or Posterior coil MCQA test)

Visual	Description
	1. Remove any other surface coils (if present) from the cradle.
 Illustration 6	2. Remove any pads present on the coil. The cable must point towards the magnet and plug into ports P1 or P2. Place lower half of phantom positioner on table. Ensure that the positioner is centered in the R/L direction. Insert coil half (anterior or posterior) and place top part of the positioner on the coil as shown in Illustration 6 and Illustration 7.

Illustration 7


Illustration 8


Illustration 9: MCQA SETUP WITH UNIFIED TL BRICK

3. Turn alignment button on

Move table into bore until landmark lasers intersect with crosshairs on the coil as shown in Illustration 8.

Press 'Landmark' button

Place phantom (legacy or unified TL brick) on the coil. The legacy phantom option is not available on MR450w systems. Ensure that the phantom is contained within the phantom positioner as shown in Illustration 9 and Illustration 10.

Press 'Advance to Scan' button.

Illustration 10: MCQA setup with legacy TL Brick

3.3 MCQA Test

3.3.1 At the operator workstation, end the previous exam (if required).

3.3.2 From the Common Service Desktop, select the [Image Quality] tab, choose “Multi-Coil QA Tool”, and then “Click here to start this tool” as shown in Illustration 11.



Illustration 11: MCQA from the Service Desktop

3.3.2 Further instructions on using the MCQA tool can be found on the MR450 or MR450w Service Methods CD or the web at: [Troubleshooting -> System -> Multi-Coil Quality Assurance Tool](#).


For troubleshooting: The relationship between MCQA image order, coil Receive Pin and coil element are the same (IMG 1 = Rx Pin 1 = Element 1). *Illustration 12* below shows the element layout for the anterior and posterior sections.

Posterior Section Viewed From Top (In positioning pad!)

Foot (cable side, also exits along patient right)
Hexagon diameter = 3.65" measured from center of trace

Anterior Section Viewed From Top

Foot (cable side also exits along patient right)

Illustration 12: Element layout for 32Ch Cardiac Coil

3.4 Periodic Quality Assurance Check

QA checks should be performed during installation and during planned maintenance. Perform the quality assurance checks as outlined below to ensure the coil is operating properly.

3.4.1 Check external cable for cracks or cuts.

3.4.2 Perform Section 3: 1.5T 32-CHANNEL CARDIAC COIL MCQA TEST

SECTION 4: TROUBLESHOOTING

Note: The Invivo 32-channel cardiac coil is not to be opened for any reason!

4-1 Mechanical Hardware Check

- Verify that there are no loops in the cables.
- Verify that there are no metal or ferromagnetic objects close to the coil, patient or magnet (i.e., safety pin, hair pin).
- Verify that the coil is properly positioned.


4-2 Electrical Hardware check

Visually inspect all coaxial connectors and pins on the system plug. If there are any broken, deformed or recessed pins or coaxial connectors, replace the coil component.

Coil Connector Illustration 13

Connector Pin Out Illustration 14

Connector pin out shown from the pin side of connector.

NOTE: The center pins of the System Connector *receive coaxial connectors* are extremely fragile [Connector Pin Out Illustration 14, R1-R8]. **Do not attempt to make any measurements at these center pins.**

For GND measurements, use any of the receive coaxial pins' outer shields (R1-R8).

RESISTANCE MEASUREMENTS – TABLE 4-2-1

System Connector		Reading
Positive Lead	Negative Lead	
1	2	< 8Ω
3	4	< 8Ω
5	GND	Open
6	GND	Open
20	GND	140 ±40Ω
21	GND	< 5Ω
22	GND	140 ±40Ω
23	GND	< 8Ω

If one or more of the readings are out of specification, replace the coil. If all of the above readings are correct, redo *Section 3: 32CH CARDIAC COIL MCQA TEST*. If the problem still exists, replace the 1.5T Cardiac Coil.

4-3 Coil ID Troubleshooting Tips

If the following error is displayed, there may be a problem with the Coil ID coil or system hardware: *The prescribed RF coil and detected RF coil do not match.*

Troubleshooting instructions are available through the Service Methods DVD or online via the path: [Troubleshooting](#) -> [Diagnostics](#) -> [Hardware Location](#) -> [Magnet Room](#) -> [Coil ID Functional Diagnostic](#)

The diagnostic will display the contents of the coil ID chip whenever the connector is inserted. Follow the troubleshooting instructions if the correct ID is not displayed. When the troubleshooting is complete, and if the correct ID is still not displayed, replace the coil.

SECTION 5 – MAINTENANCE

5-1 Coil Care

WARNING!

Detach coil connector from scanner before attempting to clean. Do not reattach after cleaning until coil has dried completely. Having the coil attached to the system during cleaning or when it is wet may result in electrical shock.

CAUTION

Do not spray or pour cleaning solution directly on coil. Do not submerge coil in solution. The coil contains sensitive electronics components that could be damaged by the solution.

The 1.5T 32-Channel Cardiac Coil, posterior pads must be cleaned and stored using the following procedure:

- Wipe with a cloth that has been dampened in a solution of 10% bleach and 90% tap water, or 30% isopropyl alcohol and 70% tap water.
- **Do not pour any cleaning solution directly on the coil!**
- Let the coil housing and pads dry before use.
- **Under no circumstances should the coil be placed into any type of sterilizer.**

SECTION 6 – RENEWAL PARTS

6.1 Field Replaceable Units

FIELD REPLACEABLE UNITS LIST – TABLE 6-1

Description	GE Part #	Invivo #
1.5T 32Ch Cardiac Coil - Anterior	5366664-2	4535-302-86221
1.5T 32Ch Cardiac Coil - Posterior	5366664-3	4535-302-86231
32Ch Cardiac Coil Phantom Positioner	5366627-4	4535-302-86261

6.2 Other Replaceable Accessories

The following table lists replaceable accessories that may be purchased from Invivo Corporation. The operator manual contains the contact information for Invivo.

OTHER REPLACEABLE ACCESSORIES LIST – TABLE 6-2

Description	Invivo #
32Ch Cardiac Coil Head Pad	4535-302-86281
32Ch Cardiac Coil Foot Pad	4535-302-86291
32Ch Cardiac Coil Anterior Pad	4535-302-86301
32Ch Cardiac Coil – Posterior Pad Includes: top posterior pad (Invivo part# 4535-302-58901), bottom posterior pad (Invivo part# 4535-302-58891), and straps (Invivo part# 4535-302-59791)	4535-302-86271
Straps Set (Includes both Anterior and Posterior Velcro straps)	4535-302-86311

Revision History

Rev	Date	Author	Primary Reason for Change	DCN Number
A	10-02-09	J. Ciesla	Initial Release	
B	10-12-09	Hohmuth	Updates	
C	10-28-09	Hohmuth	Update pics, index table	
D	12-3-09	J. Ciesla	Added install and MCQA sections	
E	12-10-09	J. Ciesla	Added element layout, edited required conditions for install	
F	02/05/10	J. Ciesla	Corrected part numbers	6809
01	02/05/10	J. Ciesla	Issued	6809
02	02/05/10	J. Ciesla	Corrected FRU numbers & label pictures	6809
03	04/01/10	J.Ciesla	Corrected Resistance Measurement in Table 4-2-1	7520