GE Signa[®] 1.5T 8-Channel Cardiac Array

Service Manual

Revision History

Date	Rev	Change Description
10/03	1	First Issue
10/07	2	Added, "Please see the locations of the preamps shown in Figure 9A,B below. The wrist strap MUST be worn all times when handling Preamps" on page 16. Added figure 9A & B to page 17. Re-numbered figures.
11/07	3	Applied New Template, Added preamplifier to table 6.1 and revised Manual, Replaced the word "coil" to "image" on section 1-5 /line 2, Added the line "Note: The wrist strap MUST be worn all times when handling Coils" to sections 3-3, 5-1 and 5-2, Removed 2 channels from the table on section 3-5

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "damage in shipment" written on all copies of the freight or express bill before delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

File a report with

- Call 1-800-548-3366 and use option 8.
- Contact your local service coordinator for more information on this process.

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WARNING

(EN)

- IF A CUSTOMER'S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER'S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR, OR PATIENT FROM ELECTRIC SHOCK, OR FROM MECHANICAL OR OTHER HAZARDS.

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TABLE OF CONTENTS

SECTION 1 – INTRODUCTION	14
1-1 PRODUCT IDENTIFICATION AND SHIPPING LIST	14
1-2 COMPATIBILITY	15
1-3 RELATED DOCUMENTATION	15
1-4 ENVIRONMENTAL REQUIREMENTS	15
1-5 THEORY OF OPERATION	16
SECTION 2 – SETUP AND CALIBRATION	19
2-1 COIL INSTALLATION.....	19
2-1-1 SPECIAL INSTALL NOTES	19
2-1-2 CONFIGURATION FOR 11X SOFTWARE.....	19
2-1-3 CONFIGURATION FOR 10X SOFTWARE.....	19
2-2 INSTALLATION FUNCTIONAL CHECKS.....	19
2-3 PERIODIC QUALITY ASSURANCE CHECK	20
SECTION 3 – FUNCTIONAL CHECKS	21
3-1 SCANNER VERIFICATION	21
3-2 COIL IMAGING PERFORMANCE VERIFICATION.....	21
3-2-1 TOOLS REQUIRED	21
3-2-2 EXPLANATION OF PROCEDURE	21
3-2-3 SNR SETUP AND TEST EXECUTION FOR 11.0 M3A AND GREATER SOFTWARE VERSIONS.....	21
3-3 EXTERNAL CABLE CHECK	23
3-4 PIN DIODES CHECK	24
SECTION 4 – MAINTENANCE	29
4-1 COIL CARE	29
4-2 SPECIAL CARE REQUIREMENTS.....	29
SECTION 5 –REPLACEMENT	30
5-1 ANTERIOR CABLE REPLACEMENT	30
5-2 POSTERIOR CABLE REPLACEMENT	31
SECTION 6 – RENEWAL PARTS	34
6-1 FIELD REPLACEABLE UNITS	34
6-2 OTHER REPLACEABLE ACCESSORIES	34
SECTION 7 – APPENDIX	35
7-1 SNR DATA SHEET	35
7-2 SCHEMATIC	36

SECTION 1 – INTRODUCTION

1-1 Product Identification and Shipping List

To identify the 8-Channel Cardiac Array, refer to the coil label (see *Figure 1*).


Figure 1: Coil label.

Figure 2 shows a picture of the 8-Channel Cardiac Array.

Figure 2: The 8-Channel Cardiac Array.

SHIPPING LIST – TABLE 1-1

Box #	Part Name	GE Part #	USAI Part #	Qty
1	Coil	2304255-2	100391	1
2	Phantom Kit	2304255-3	150511	1
3	Patient Comfort Pad	2304255-6	150478	1
3	Patient Headrest Pad	N/A	150497	1
3	Patient Lower Body Pad	N/A	150477	1
1	Operator Guide	2304255-4	770012	1
1	Service Manual	2304255-5	780032	1

1-2 Compatibility

The 8-Channel Cardiac Array is compatible with the Excite 1.5T system and requires 10.0 M4 or greater software. 11.0 M4 software is recommended.

The automated QA Tool is only available for software 11.0 M3 and greater.

1-3 Related Documentation

Operator's Guide, GE Part Number 2304255-4.

Signa Excite Service Methods CD, GE Part Number 2333500-3.

1-4 Environmental Requirements**Storage Requirements**

Coil should be stored in the scanner room.

Dimensions

19.4" L x 16.5" W x 4.7" H (49.4cm x 42.0cm x 11.9cm)

Weight

12 lbs (5.4 kg)

1-5 Theory of Operation

The block diagram of 8-Channel Cardiac Array is shown in *Figure 3*. The coil is designed to give a high-resolution image for the cardiovascular system. It works at 1.5T field strength and is capable of doing Array Spatial Sensitivity Encoding Technique (ASSET). This design reduces the breath-hold time, which will aid the patient to successfully hold their breath for the designated time while scanning. The coil is composed of two separate pieces, anterior and posterior. The posterior piece is rigid. The anterior piece is flexible. The flexible design of the coil can accommodate different patient habits. Thus the optimum Signal to Noise Ratio (SNR) can be achieved.

The 8-Channel Cardiac Array is a receive-only coil, as indicated in the schematic shown in the Appendix. The coil is actively decoupled from the RF transmit coil during transmit by means of a RF choking circuit. The choking circuit elements are switched on and off by pin diodes or by small-signal diodes. During body coil transmission, the pin diodes are turned on, as shorted circuits, with a forward bias DC source provided by the MR system (active decoupling). The small-signal diodes are turned on by the induced RF voltage in the chokes, coupled from the transmit field (passive decoupling). When either the pin diodes or the small-signal diodes are turned on, the RF choking elements become very high in impedance (typically above 2 kilo-ohm), compared with the other circuit elements (less than 70 ohm). These high impedance elements in the coil segregate the coil circuitry into several isolated electrical segments, preventing any current flow in the coil circuit.


Figure 3: Block diagram of the coil.


Figure 4: Block Diagram of the coil to connector interface.


Figure 5: Block diagram of the circuit layout.

SECTION 2 – SETUP AND CALIBRATION

2-1 Coil Installation

2-1-1 Special Install Notes

None.

2-1-2 Configuration for 11x Software

Use the Configuration File Manager to install the coil configurations listed in Table 2-1-3.

2-1-3 Configuration for 10x Software

The 8-Channel Cardiac Array requires 10.0 M4 or greater software. 11.0 M4 software is recommended.

Updated coil configurations in the system are required for software versions prior to 11.0 M4. Updated coil configurations are recommended for 11.0 M4 and greater software revisions.

To update the coil configuration files:

1. Using the configuration manager, remove all currently installed (if any) configurations for the 8-Channel Cardiac Array.
2. Update the coil configurations to the latest version (Refer to Section 9 of “Host Software Install Release 10.X”, direction OW2SCA8G. as found on the Service Methods CD, GE Part Number 2333500-3).
3. Install the new configurations listed in table 2-1-3.

CONFIGURATION MODES – TABLE 2-1-3

8CARDIAC	8CARDIACANT	8CARDIACPOST
----------	-------------	--------------

2-2 Installation Functional Checks

1. From the Scan Desktop, start new scan by selecting [**New Pt**]; set **Patient ID** to “geservice” and **Patient Weight** to “111” pounds. Click [**Patient Position**] to open protocols window.

2. Plug the coil into the patient table.
3. At the console, verify the Coil has been properly identified by the system: correct name in **Coil** field and correct picture on the screen.
4. Perform Section 3-2 Coil Imaging Performance Verification.

2-3 Periodic Quality Assurance Check

On a periodic basis, such as during planned maintenance, perform the quality assurance checks as outlined below to ensure the coil is operating properly.

1. Check external cable for cracks or cuts.
2. Perform Section 3-2 Coil Imaging Performance Verification and record data values in Data Sheet (See section 7).

SECTION 3 – FUNCTIONAL CHECKS

3-1 Scanner Verification

Perform system level Signal to Noise Check. Refer to Service Methods CD; System Level Procedures; Functional Checks; Signal to Noise Check.

3-2 Coil Imaging Performance Verification

3-2-1 Tools Required

TOOLS REQUIRED – TABLE 3-2-1

Description	GE Part #	USAI Part #	Qty
Phantom Kit	2304255-3	150511	1

3-2-2 Explanation of Procedure

The image quality check uses two consecutive acquisitions for the same slice. The signal is measured at the center of the phantom with a ROI of about 60% of the phantom size. The noise is measured on the subtracted image of the two consecutive signal images using the same ROI. SNR is calculated by dividing the Signal mean by the noise standard deviation. All measurements are performed on composite images. Interim images are acquired for reference only.

3-2-3 SNR Setup and Test Execution for 11.0 M3A and Greater Software Versions

1. Remove any other surface coils (if present) from the cradle.
2. Place the coil on the patient table. Position the phantom set on the coil without the pad, lining up the pegs on the phantom to the holes in the posterior section of the coil (see *Figure 6*).


Figure 6: Positioning the phantom on the coil.

3. Using the raised notches on top of the phantom, align the anterior section on the phantom and fasten the Velcro straps.
4. Landmark the coil at the indicated line and advance to scan (see *Figure 7*).


Figure 7: Landmarking the coil.

5. From the Common Service Desktop, select **Image Quality** and **Coil SNR**, then **Click here to start this tool.**
 6. Select the coil type you want to run, enter the serial number and select **Start Scan**.
- Note: Additional documentation can be accessed by selecting **File, Help, [coil type]**.

3-3 External Cable Check

Diode Check from the Connector:

Note: The wrist strap **MUST** be worn all times when handling Coils.

Step 1 – Using a multi-meter set to the diode check function, check the coil diodes in the following manner.

Step 2 – Orient the connector as shown in *Figure 8*.

Step 3 – Follow Table 3-3 for the correct pins and lead locations to verify the function of the diodes.

EXTERNAL CABLE EXPECTED READINGS – TABLE 3-3

PIN Number	Positive Lead	Negative Lead	DMM Reading
1	B12	A12	0.5 +/- 0.2V
2	B14	A14	0.5 +/- 0.2V
3	B16	A16	0.5 +/- 0.2V
4	B18	A18	0.5 +/- 0.2V
5	B20	A20	0.5 +/- 0.2V
6	B22	A22	0.5 +/- 0.2V
7	B24	A24	0.5 +/- 0.2V
8	B26	A26	0.5 +/- 0.2V


Figure 8: 60 pin Bendix connector.

Note: This connector is system side. Coil side has rows A and B reversed.

3-4 PIN Diode Check

Same as External Cable Check.

The posterior section of the coil must be plugged into the system connector for the pin diode check.

3-5 Troubleshooting Tips

Symptom #1: The system reports a coil fault during prescan or does not recognize the coil connection to the system when selected in the software.

Probable Cause	Suggested Actions	Resolution
The coil connector has become disconnected from the system interface.	Check to make sure the coil connector is fully engaged.	Engage connector and try the scan again. Make sure the indicator light on the patient table turns green with the proper coil selected on the console.
One or more PIN diodes have failed.	From inside the coil, check the PIN diodes using a DVM. Use diode mode for forward-bias measurements and resistance mode for reverse-bias measurements. Note that there should be no reverse-bias leakage. Bad PIN diodes are usually shorted or exhibit reverse-bias leakage current.	Replace the coil.
There is a DC short or open somewhere within the DC path inside the coil.	Check the coil and associated DC paths with a DVM, in both diode check and resistance modes.	If a short or open is observed, check the output cable by itself (see below). If the cable and PIN diodes are okay, then replace the coil.
The output cable has a short or an open.	Disconnect the cable at the coil and check each coaxial line with a DVM (resistance). Try moving and twisting the cable as you watch the meter to find any intermittent connections. Inspect the center pins of the SMB and Bendix connections for wear and proper engagement into the mating connection.	Replace the coil.
There is a problem with the DC bias from the system interface or the coil ID circuitry.	Check to ensure the correct DC bias is being provided to the coil during transmit and receive per GEMS System Service Manual. Try a similar coil to see if the same problem occurs.	Correct MRI system problem.
There is a problem with the coil ID; <i>"The prescribed RF coil and detected RF coil do not match."</i>	Check online instructions by accessing the service browser at Diagnostics; Peripherals; SRI Functional Tests; Coil ID	Start the diagnostic by selecting SRI Coil ID and Run . The diagnostic will display the contents of the coil ID chip whenever the connector is inserted. Follow the troubleshooting instructions if the correct ID is not displayed. When the troubleshooting is complete and the ID is still not displayed, replace the cable.

Symptom #2: The coil does not pass SNR tests or exhibits poor image quality on patient scans.

Probable Cause	Suggested Actions	Resolution
One or more of the coil channels has a high noise level.	Look at the uncombined images to determine which channel has the problem. Compare the noise standard deviation measurements between channels and system performance logs; they should be approximately the same. Try swapping preamplifiers between channels to see if the problem follows the preamp. If it does, then the preamplifier is bad. Please see the locations of the preamps shown in Figure 9A, B below. The wrist strap MUST be worn all times when handling Preamps.	Replace the coil.
One or more of the coil channels has low signal.	Look at the uncombined images to determine which channel has the problem. Compare signal mean measurements between channels and system performance logs.	Replace the coil.
There is a problem with the DC bias from the system interface or the coil ID circuitry.	Check to ensure the correct DC bias is being provided to the coils during transmit and receive per GEMS System Service Manual.	Correct MRI system DC bias problem.

Preamplifier
Locations

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
Figure 9B: Posterior Section

3-6 Coil Element Information

Mode	Rcv 1	Rcv 2	Rcv 3	Rcv 4	Rcv 5	Rcv 6	Rcv 7	Rcv 8	PIN 1	PIN 2	PIN 3	PIN 4	PIN 5	PIN 6	PIN 7	PIN 8
8CARDIAC	ON	ON	ON	ON	ON	ON	ON	ON	-5	-5	-5	-5	-5	-5	-5	-5
8CARDIAC POST					ON	ON	ON	ON	+7	+7	+7	+7	-5	-5	-5	-5
8CARDIAC ANT	ON	ON	ON	ON					-5	-5	-5	-5	N/A	N/A	N/A	N/A

SECTION 4 – MAINTENANCE

4-1 Coil Care

Detach the coil connector from the scanner before attempting to clean the coil. Do not touch the connectors with bare fingers. Never press a sharp object against the surface of the connector. Do not reattach the connector after cleaning the coil until the coil has dried completely. Electrical shock may result if the coil is attached to the system during cleaning or when it is wet.

Do not spray or pour cleaning solution directly on coil. Do not submerge coil in solution. The coil contains sensitive electronic components that could be damaged by the solution. The coil cannot be sterilized and should be cleaned only according to the procedure outlined in this section.

The following solutions are recommended for the coil and pad surfaces: (1) a ten percent bleach solution (some discoloration may occur), (2) one ounce commercial dishwashing liquid mixed with one gallon of water or (3) warm water. Apply cleaning solution to a soft cotton cloth and proceed to clean. To prevent soiling of the coil, the user should place a cotton sheet over the coil before positioning the patient. If the coil is soiled, clean the coil as described above.

4-2 Special Care Requirements

Prior to returning a coil for service, use a ten percent bleach solution (as described above) to eliminate risk of exposure to potentially infectious materials.

SECTION 5 –REPLACEMENT

Simple removals that are clearly obvious are not described here.

Unless otherwise noted, the steps for re-assembly are simply the reverse order of the steps described for disassembly.

5-1 Anterior Cable Replacement

Note: The wrist strap **MUST** be worn all times when handling Coils.

Step 1 – Remove coil cover.

Step 2 – Remove the two screws that hold down the cable assembly.

Step 3 – Disconnect SMB plug and the DC Plug on the cable end from the interface board.

Step 4 – Remove the old cable assembly.

Step 5 – Fix a new cable assembly on the anterior section of the coil and reattach the cable lock with two screws. Connect all four SMB plugs and the DC Plug to the corresponding receptacles on the interface board (see *Figure 12*).

Step 6 – Close the coil cover.


Figure 10: Anterior cable assembly.


Figure 11: Anterior DC plug assembly.


Figure 12: Anterior SMB connector assembly.

5-2 Posterior Cable Replacement

Note: The wrist strap **MUST** be worn all times when handling Coils.

Step 1 – Remove coil cover.

Step 2 – Remove the two screws that hold down the cable assembly.

Step 3 – Disconnect SMB plug and the DC Plug on the cable end from the interface board.

Step 4 – Remove the old cable assembly.

Step 5 – Fix a new cable assembly on the posterior section of the coil and reattach the cable lock with two screws. Connect all four SMB plugs and the DC Plug to the corresponding receptacles on the interface board (see *Figure 15*).

Step 6 – Close the coil cover.


Figure 13: Posterior cable assembly.


Figure 14: Posterior DC plug assembly.


Figure 15: Posterior SMB connector assembly.

SECTION 6 – RENEWAL PARTS

6-1 Field Replaceable Units

FIELD REPLACEABLE UNITS – TABLE 6-1

Part Name	GE Part #	USAI Part #
Coil	2304255-2	165148
Cable Assembly	2304255-7	155156
Phantom Kit (Phantom, Rings and Positioners)	2304255-3	155158
Preamplifier	2377426-10	N/A

6-2 Other Replaceable Accessories

OTHER REPLACEABLE ACCESSORIES LIST – TABLE 6-2

Part Name	GE Part #	USAI Part #
Patient Comfort Pad	N/A	155164
Patient Headrest Pad	N/A	155165
Patient Liner Pad	N/A	155166
Velcro Straps (Set of 2)	N/A	155198


Figure 16: Parts of the coil.

SECTION 7 – APPENDIX

7-1 SNR Data Sheet

Use the table provided below to record the calculated signal to noise ratio (SNR) data obtained from the Functional Checks (see section 3).

Date	Configuration	Slice Location	Corresponding Images in Series	Signal Mean	Noise Std. Dev.	SNR

Comments	

Date	Configuration	Slice Location	Corresponding Images in Series	Signal Mean	Noise Std. Dev.	SNR

Comments	

GE Healthcare Coils
USA INSTRUMENTS

1515 Danner Drive
Aurora, Ohio 44202
USA

www.gehealthcare.com

Imagination at work