GE Healthcare

Signa 16 Channel RF Electronics Upgrade

Technical Publication

5177002
Rev 3

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation “**damage in shipment**” written on **all** copies of the freight or express bill **before** delivery is accepted or “signed for” by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier **immediately** upon discovery, or in any event, within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this **14** day period.

File a report with

- Call 1–800–548–3366 and use option 8.
- Contact your local service coordinator for more information on this process.

Rev. 08/15/2003

© 2006, 2009 General Electric Company
GE Medical Systems, a General Electric Company, going to market as GE Healthcare
3002 N. Grandview Boulevard
Waukesha, Wisconsin 53149
USA

www.gehealthcare.com

Table of Contents

DESCRIPTION	PAGE
Table Of Contents	1
Language Policy For Service Documentation	3
Revision History	5
SECTION 1 – Overview	7
SECTION 2 – Upgrade Procedure	9

This page left blank intentionally.

WARNING

- THIS SERVICE MANUAL IS AVAILABLE IN ENGLISH ONLY.
- IF A CUSTOMER'S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER'S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR OR PATIENT FROM ELECTRIC SHOCK, MECHANICAL OR OTHER HAZARDS.

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Revision History

Date	Rev	Change Description
14-Feb-2006	1	Initial Release
14-Jul-2006	2	Section 1: Removed information that a second MUX board was to be installed. Section 2: Removed instructions in for removing and installing MUX board in slots 13 and 14. (MRIhc16250) Added step for connecting cable to J5 on MUX board.
10-Jun-2009	3	Section 2: Added Mux board 5250106 to step 2 in Section 2-2 (PQR 13238441)

This page left blank intentionally.

SECTION 1 – OVERVIEW

1-1 SIGNA HORIZON SYSTEM CONFIGURATION

This manual provides instructions to upgrade a Signa EXCITE HD or HDx system to a 16 Channel system.

This Direction covers the following catalogs:

M3000PE – 16 Channel RF Electronics Upgrade

This upgrade adds:

- RF-DIF, MUX, UTNS3, and RCVR II boards.
- Cables

1-2 SITE READY CHECK

Pre-installation work must be completed before equipment is delivered to avoid delays and confusion.

1-3 DISCARDED MATERIAL RETURN POLICY

The equipment or material being removed by the upgrade have a potential value for service of the installed base. The disposition of such must be controlled by GE. Disposition of the removed material is per GE Medical Systems Policy and Procedures.

1-4 PRODUCT LOCATOR

The Global Install Base Database (also known as Product Locator System) tracks shipment, trans-shipment, and field location of the serialized models. There are now two methods for submitting Product Locator information.

1. At this time, for **U.S. ONLY**, the preferred method of submitting information is the *FE Site Verification Web Site* at:
<http://gein2.med.ge.com/gib>

The FE Site Verification consists of three components that are available on the web from the main menu. They are:

- Install/deinstall product locator model and serial numbers
- Add/modify ship to address information
- Update CARES FE data for primary/secondary FEs

To obtain a copy of the FE training tutorial for using this website, a downloadable copy is available at the following:
Product Locator Support Central Page – <http://supportcentral.ge.com/15563>

2. One "Shipping Card" is filled out and submitted when shipped (extra cards are supplied for trans-shipments between storage and distribution points), and the "Installation Card" and extra shipment cards are attached.

Verify that serial and model number on each rating plate matches installation card numbers before removing installation card. Note that there may be one or more shipment cards and bar code labels with the installation card. These shipment cards are used to trace the transfer of serialized units between various inventory storage and distribution points until the product reaches its final installation destination. Process just the installation card and discard any extra shipment cards and labels.

1-5 SOFTWARE LOAD FOR UPGRADE

Install 16 Channel option key.

1-6 FINAL COORDINATION

1. Replace all covers that have not been previously replaced.
2. Record and enter applicable data into applicable site configuration files and records.
3. Process product locator installation cards for all serialized components added to the system to include product locator card supplied with the upgrade serial number rating plate. Return to:

**Product Locator File
P.O. Box 414, W-523
Milwaukee, WI 53201-0414**

Note

Failure to fill out and return Product Locator Cards may result in failure of your site to receive future FMI's.


SECTION 2 – UPGRADE PROCEDURE

2-1 REMOVE EXISTING COMPONENTS FROM RRF CHASSIS

Note: Module is a static sensitive component and needs to be handled as such (wrist strap).


- ① Remove blank panel from slot 1 at front of chassis.

FRONT VIEW

Note: It is recommended that the existing MUX board be at the same part number and revision level as the new MUX board being installed in slots 4 & 5.

- ② Remove blank panels from slots 1 to 5 at rear of chassis.


REAR VIEW

2-2 INSTALL MODULES IN RRF CHASSIS


Note: Module is a static sensitive component and needs to be handled as such (wrist strap).


FRONT VIEW

- ① Install RF-DIF Board (5250028) in slot 1.

REAR VIEW


- ② Install MUX Board (2274188-3 or 5250106) in slots 4 & 5 .
- ③ Install UTNS3 Board (2366218-2) in slot3.
- ④ Install RCVR II Board (2349808-2) in slots 1 and 2.

2-3 INSTALL CABLES AT REAR OF RRF MODULE

- ① Connect cable 5123706 as marked.
- ② Connect cable 5113193-2 as marked.
- ③ Make sure cable 5115417 slot 12R J31 is connected as marked. it should be routed and connect to Coax Power Splitter "1".
- ④ Connect cable 5115417-2 as marked. Route and connect to Coax Power Splitter "2". Remove 50 Ohm terminator if installed at "2".
- ⑤ Move 50 Ohm Terminator from Mux1 J6 (slots 13 & 14) to Mux2 J6 (slots 4 & 5).
- ⑥ Connect cable 5115630 as marked.
- ⑦ Connect cable 5145030-2 (HDx) or 2394870-2 (HD) to new MUX board.

2-4 INSTALL CABLES AT FRONT OF RRF MODULE

- ③ Connect cable 5121122-2 as marked to RF Slot 1F J21. Route and connect other end to MGD Slot 11F J7.
- ② Connect cable 5113182-2 as marked.
- ① Disconnect cable (2323878-3) from Slot 10F J23 and connect to Slot 1F J23.

2-5 INSTALL NEW REAR DOOR AND RATING PLATE

- ① Attach upgrade Rating Plate (46-302200P12) to top of cabinet on back side.