



System ID: _____

SR#: _____

Date: _____

PM Report Schedule A

LightSpeed VCT, Discovery CT750 HD, Revolution Discovery CT, Revolution GSI, Revolution HD, and Discovery CT

Customer Name: _____ Service Record #: _____

System ID/ Serial Number: _____ Customer Equipment ID: _____

Service Manual Document Number: _____ Software Version: _____

PM Frequency: _____ months

System Type (Select:)

PRELIMINARY TASKS

Last PM Inspection

Verify last PM Inspection

Last Inspection Date :

Last Schedule Completed:

List any Schedule A tasks that were completed earlier

PM Schedule Sys A

Review Error Logs

	Completed	Comments/Follow-up
System Application Stop and Start Errors	☐	

Initial Procedures

General Gantry Inspection

Task	Pass/Fail	Comments/Follow-up
Test Display Lights		
Test E-stops		
Test Patient Tilt Sensors		

Gantry Power Off (Covers Off)

General Tasks

Task	Pass/Fail	Comments/Follow-up
Emergency Off Switch (on wall)		

General Console Cleaning & Inspection

Checklist PM Document #: DOC0928611

Revision Number: 8.0

© 2014, General Electric Company

All Rights Reserved

Page 1 of 5



System ID:

SR#:

Date:

Task	Pass/Fail	Comments/Follow-up
Clean Component Filters, Fans & Grills		
Visual Fan Inspection		
Seismic Anchor Check (if necessary)		

Slip Ring & Brush Block Inspection and Maintenance

Task	Pass/Fail	Comments/Follow-up
Remove Slip Ring Brush Debris		
Slip Ring Tracks Inspected		
Clean Brush Blocks		
Brush Tip Inspection		

Gantry Filter Cleaning

Task	Pass/Fail	Comments/Follow-up
Clean DAS Detector Plenum Filter		
Clean Detector Face Plate		
Clean Gantry Heater Filter		
Clean Top Cover Fans		
Inspect Rotating Path Components & Cables		

NGPDU

Task	Pass/Fail	Comments/Follow-up
Vacuum Covers & Inside PDU		

Encoder

Task	Pass/Fail	Comments/Follow-up
Inspect Rotational Encoder Wheel		
Encoder Replacement		

Gantry Power On (Covers Off)

General Tasks

Task	Select	Comments/Follow-up
Grease Gantry Main Bearing		
Record Gantry Rotations		



System ID:

SR#:

Date:

Plenum Fans (5) Operational		
UPS Powered and Operational		

Tube Heat Exchanger & Pump

Task	Select	Comments/Follow-up
Inspect/Clean Tube Heat Exchanger		
Inspect/Clean JEDI Inverter Fan		
Tube Heat Exchanger Pump Inspection		

Gantry Power On (Covers On)

General Tasks

Task	Select	Comments/Follow-up
Inspect Rotational Clearances		
Scan Window Inspection/Replacement		

General Table Cleaning & Inspection

Task	Select	Comments/Follow-up
Inspect/Clean Table Pan		
Inspect Table anchors		
Inspect Drive Belts		
Clean Roller		
Grease Table Rails		
Inspect Bearing Retainer		
Grease Table Leg Bearing		
Inspect Flat Tabletop Mount		
Inspect Table Covers		
Inspect Cradle Accessory Attachment		
Check Tape Switches		

Task	QA Phantom Serial #	Location	Gantry Serial #
Scan	Scan Failures	Artifacts	Artifact Types
Series 1 (Scout/Auto Voice)			



System ID:

SR#:

Date:

Finalization

System Scanning Test

Task	QA Phantom Serial #	Location	Gantry Serial #
Series 2 (Axial/Tilt)			
Series 3 (Helical/Auto Voice)			

General Tasks

Task	Select	Comments/Follow-up
X-ray On Indicators		
Scan Control Pushbuttons Test		
	Completed	
Update Site Log	☐	
Complete PM Paper Work	☐	
Complete Customer Reports	☐	

Is follow-up work required?

Service Scheduled Date:

☐ This inspection passed and meets all GE Healthcare PM specifications

Comments:

Characters Remaining:

Tools Used:

Description	Serial Number	Bar Code	Cal Due Date	
				Add Tool
				Remove Tool

GE Representative:

Checklist PM Document #: DOC0928611
Revision Number: 8.0

© 2014, General Electric Company
All Rights Reserved



System ID:

SR#:

Date:

FSEs completing this PM

FSE Name: _____

SSO #: _____

FSE SR #: _____

FSE Signature Data...

Customer Signature: ☐ Yes ☐ No

☐ Final PM Complete

Date
Completed: _____