GE Healthcare


# High Order Shim Option

## Service Methods

OPERATING DOCUMENTATION

5314688  
Revision 1

# Important Information

## LANGUAGE

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**WARNING****(EN)**

- THIS SERVICE MANUAL IS AVAILABLE IN ENGLISH ONLY.
- IF A CUSTOMER'S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER'S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR, OR PATIENT FROM ELECTRIC SHOCK, OR FROM MECHANICAL OR OTHER HAZARDS.

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This page left intentionally blank.

## Revision History

| Rev. | Date         | Primary Reason(s) for Change |
|------|--------------|------------------------------|
| 1    | May 28, 2008 | Initial Release              |

This page left intentionally blank.

# Table of Contents

**CHAPTER 1 PREINSTALLATION . . . . . 17**

1 Unpacking the High Order Shim Hardware . . . . . 17

1.1 Personnel Requirements . . . . . 17

1.2 Overview . . . . . 17

1.3 Product Locator . . . . . 18

**CHAPTER 2 SAFETY . . . . . 19**

1 Lock Out/Tag Out . . . . . 19

1.1 Safety . . . . . 19

1.2 Required Contitions . . . . . 19

**CHAPTER 3 INSTALLATION . . . . . 21**

1 Shim Filter Installation . . . . . 21

2 Shim Power Supply Installation . . . . . 23

3 Cable Installation . . . . . 25

**CHAPTER 4 PARTS . . . . . 29**

1 Replacement Parts . . . . . 29

This page left intentionally blank.

# Chapter 1 Preinstallation

## 1 Unpacking the High Order Shim Hardware

### 1.1 Personnel Requirements

| Required Persons | Preliminary Reqs | Procedure          | Finalization   |
|------------------|------------------|--------------------|----------------|
| 1                | Not Applicable   | 2 hours to install | Not Applicable |

**NOTE:** It will take an additional 2 hours to run the functional check and to build a shim calibration file, so plan accordingly.

### 1.2 Overview

The commercial catalog M3000PP, High Order Shim Option, is compatible with and can be installed on the following product configurations:

- Signa Excite 1.5T
- Signa HD 1.5T
- Signa HDx 1.5T and 3.0T

Verify that all parts are present in the shipment (see [Table 1-1](#) ).

**Table 1-1: Catalog M3000PP Contents List**

| Item | Qty | Part Number | Description                                                                  |
|------|-----|-------------|------------------------------------------------------------------------------|
| 1    | 1   | 2352973     | Shim Filter Box Assembly                                                     |
| 2    | 1   | 2352972     | Shim Power Supply Enclosure Assembly (Y)                                     |
| 3    | 1   | 2303317     | CAN Termination (Male) (Y)                                                   |
| 4    | 1   | 2363363-4   | Cable; From Shim PS J4 to TAC I/O J119                                       |
| 5    | 1   | 2374156-2   | Cable; From Shim PS J3 To TAC I/O J126                                       |
| 6    | 1   | 2384104-2   | Cable; Run# 1152; From PP1-A5-P1 To MG2A12A3                                 |
| 7    | 1   | 5113502-5   | Cable; Run# 1233 From TAC I/O J116 To PP1-PHPS-J155                          |
| 8    | 1   | 5117890     | Cable; From Shim PS J5 to TAC I/O 116                                        |
| 9    | 1   | 5118283     | Cable; Run# 611; From TAC I/O J126 to PP1-A5-P1                              |
| 10   | 1   | 2322427     | Cable; Power Cord NEMA Male To IEC320F 7ft long                              |
| 11   | 4   | 2279742-2   | ID 0.2040 OD 0.5630 THK .0478STL F70B5A (for securing the Shim Power Supply) |
| 12   | 4   | 2279742-2   | Screw; Pan Head 5 mm x 16 mm (for securing the Shim Power Supply)            |

**Table 1-1: Catalog M3000PP Contents List (cont'd)**

|    |   |                   |                                                                                                                                                        |
|----|---|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| 13 | 4 | 46-170012P35      | 6-32, 1/4 inch Hex Keps Nut, 3/32 Thick Overall Height 7/64; Washer OD 0/32 Washer Material Thk .013; CAD PL Keps (for mounting connector to I/O J126) |
| 14 | 4 | 46-208560P28      | Screw; Pan Head; 6-32 X .500 inch; Zinc Plated (for mounting connector to I/O J126)                                                                    |
| 15 | 6 | 46-221417P2       | Female Screw lock for Sub-d Connectors 4-40 Thread; 0.503 inch Overall Length; Individual Piece Without Nut Or Washer                                  |
| 16 | 1 | Direction 5314688 | Service Manual                                                                                                                                         |

### 1.3 Product Locator

The Global Install Base Database (also known as Product Locator System) tracks shipment, trans-shipment, and field location of the serialized models. There are now two methods for submitting Product Locator information.

At this time, for U.S. ONLY, the preferred method of submitting information is the FE Site Verification Web Site at: <http://gein2.med.ge.com/gib>

The FE Site Verification consists of three components that are available on the web from the main menu. They are:

- Install/deinstall product locator model and serial numbers
- Add/modify ship to address information
- Update CARES FE data for primary/secondary FEs

To obtain a copy of the FE training tutorial for using this website, a downloadable copy is available at the following: Product Locator Support Central Page – <http://supportcentral.ge.com/15563>

One “Shipping Card” is filled out and submitted when shipped (extra cards are supplied for trans-shipments between storage and distribution points), and the “Installation Card” and extra shipment cards are attached.

1. Verify that serial and model number on each rating plate matches the installation card numbers before removing installation card.

**NOTE:** There may be one or more shipment cards and bar code labels with the installation card. These shipment cards are used to trace the transfer of serialized units between various inventory storage and distribution points until the product reaches its final installation destination. Process just the installation card and discard any extra shipment cards and labels.

2. Enter the model and serial number information into the GIB database.

## Chapter 2 Safety

### 1 Lock Out/Tag Out

#### 1.1 Safety

#### **⚠ DANGER**

**ELECTROCUTION HAZARD.  
HIGH VOLTAGE PRESENT.  
USE PROPER LOCKOUT / TAGOUT PROCEDURES BEFORE INSTALLING THE  
COMPONENTS.**

#### **NOTICE**

**Follow ALL required safety and PPE procedures customary for your organization when working on this product.**

#### 1.2 Required Contitions

Before starting any hardware installation, the MRI system must be powered down and LOTO procedures executed.

| Condition                            | Reference                                                                                                             | Eff |
|--------------------------------------|-----------------------------------------------------------------------------------------------------------------------|-----|
| System PDU locked out and tagged out | Refer to the Lockout / Tagout for Gradient Amplifier / PDU Prcedure for the specific Signa system you are working on. |     |

This page left intentionally blank.

## Chapter 3 Installation

### 1 Shim Filter Installation

#### 1.1 Personnel Requirements

| Required Persons | Preliminary Reqs | Procedure | Finalization   |
|------------------|------------------|-----------|----------------|
| 1                | Not Applicable   | 30 mins   | Not Applicable |

#### 1.2 Preliminary Requirements

##### 1.2.1 Tools and Test Equipment

| Item                                           | Qty | Effectivity | Part#          | Manufacturer |
|------------------------------------------------|-----|-------------|----------------|--------------|
| 9/32 inch socket and driver or open-end wrench | 1   | -           | Field supplied |              |

##### 1.2.2 Consumables


| Item         | Qty | Effectivity | Part#          | Manufacturer |
|--------------|-----|-------------|----------------|--------------|
| Masking tape | 1   | -           | Field supplied |              |

#### 1.3 Procedure


1. Remove the blank panel located below the Gradient Filter on the Penetration Panel.
2. Remove the nuts and washers attached to the studs on the Shim Filter Assembly.
3. Install the Shim Filter Assembly, p/n 2352973, into the opening on the Penetration Panel from the Equipment Room side of the panel. See [Illustration 3-1](#) , Equipment Room View.

**Illustration 3-1: Penetration Panel Views**

Equipment Room View

Magnet Room View

4. Apply strips of masking tape across the corners of the Shim Filter assembly to temporarily hold it in place.
5. Secure the Shim Filter assembly using the nuts and washers, which were previously removed in step 2. See [Illustration 3-1](#) , Magnet Room View.
6. Make sure all the nuts are tightened to enable a good RF seal on the panel.

## **1.4 Finalization**

No finalization steps.

## 2 Shim Power Supply Installation

### 2.1 Personnel Requirements

| Required Persons | Preliminary Reqs | Procedure | Finalization   |
|------------------|------------------|-----------|----------------|
| 1                | Not Applicable   | 30 mins   | Not Applicable |

### 2.2 Preliminary Requirements

#### 2.2.1 Tools and Test Equipment

| Item                   | Qty | Effectivity | Part#                     | Manufacturer |
|------------------------|-----|-------------|---------------------------|--------------|
| #2 Philips Screwdriver | 1   | -           | 5109894-3 or equivalent   |              |
| Allen key 5/32 inch    | 1   | -           | 5109894\2-8 or equivalent |              |

#### 2.2.2 Safety

**NOTE:** The Shim Power Supply weighs 36 lbs (16.4 Kgs). Use proper lifting techniques when handling the Shim Power Supply.

### 2.3 Procedure

1. Remove the Front Cover from the TAC Cabinet using a 5/32" Allen key to open the cover lock.
2. Install the Shim Power Supply onto the Slide Rails located above the Gradient Switch Assembly. See [Illustration 3-2](#).

**Illustration 3-2: Shim Power Supply Front View**

3. Secure the Shim Power Supply to the vertical rails using 4 flat washers and 4 pan head Philips screws.

## **2.4 Finalization**

No finalization steps.

### 3 Cable Installation

#### 3.1 Personnel Requirements

| Required Persons | Preliminary Reqs | Procedure | Finalization |
|------------------|------------------|-----------|--------------|
| 1                | Not Applicable   | 1 hr      | 2 hrs        |

#### 3.2 Preliminary Requirements

##### 3.2.1 Tools and Test Equipment

| Item                             | Qty | Effectivity | Part#                   | Manufacturer |
|----------------------------------|-----|-------------|-------------------------|--------------|
| 1/ inch Flat Blade Screwdriver   | 1   | -           | 5109895-4 or equivalent |              |
| 3/16 inch Flat Blade Screwdriver |     | -           | 5109895-7 or equivalent |              |
| 3/16 inch Nut Driver             | 1   | -           | Field Supplied          |              |


##### 3.2.2 Consumables

| Item       | Qty | Effectivity | Part#          | Manufacturer |
|------------|-----|-------------|----------------|--------------|
| Cable Ties | 20  | -           | Field supplied | -            |

#### 3.3 Procedure

1. Install the internal cabinet cables between the Shim Power Supply and the I/O Panel as shown in wiring diagram [Illustration 3-3](#) .

**Illustration 3-3: TAC Cabinet Wiring Diagram**

2. Attach the bulkhead connector of each cable to the I/O Panel as defined in [Illustration 3-3](#) .

**NOTE:** For mounting connector to J126 use the following mounting hardware:

- 4 each; 46-208560P28; Screw; Pan Head; 6-32 X .500 inch
- 4 each; 46-170012P35; 6-32, 1/4 inch Hex Keys Nut

**NOTE:** For mounting connectors to J116 and J119 use the following mounting hardware:

2 each; 46-221417P2 – Female screw lock for Sub-d Connectors 4-40 Thread


3. Route and secure the cables to the vertical rails inside the cabinet as shown in [Illustration 3-4](#) using cable ties.

**Illustration 3-4: Internal Cabinet Cable Routing**

4. Install the external cabinet cables between the TAC I/O Panel and the Penetration Panel as defined in [Illustration 3-5](#).

**Illustration 3-5: External Cabinet Cable Diagram**

### **3.4 Finalization**

1. Refer to the service methods CD for the Signa system you are working on and perform the following tasks:
  - Remove LOTO for Gradient Amplifier / PDU
  - Power Up the System
  - Perform High Order Shim Supply Functional Check
  - Perform HOShim Calibration Reference Map Creation and Procedures
  - Perform HOShim Calibration Reference Map Creation and Procedures

This page left intentionally blank.

# Chapter 4 Parts

## 1 Replacement Parts

| Item                       | Qty | Effectivity | Part Number | Manufacturer |
|----------------------------|-----|-------------|-------------|--------------|
| Shim Filter Assembly       | 1   |             | 2352973     |              |
| Shim Power Supply Assembly | 1   |             | 2352972     |              |
| CAN Terminator             | 1   |             | 2303317     |              |

This page left intentionally blank.

© 2008 General Electric Company.  
GE Medical Systems, a General Electric Company, going to market as GE Healthcare.

[www.gehealthcare.com](http://www.gehealthcare.com)