GE Healthcare Coils

GE Signa[®] 1.5T 8-Channel Cardiac Array Service Manual

GE Catalog Part Number: M3335JD

2304255-5
Rev 6
© 2009, General Electric Company,
All Rights Reserved

Revision History

Date	Rev	ECO	Change Description
10/03	1	N/A	First Issue
10/07	2	N/A	Added, "Please see the locations of the preamps shown in Figure 9A,B below. The wrist strap MUST be worn all times when handling Preamps" on page 16. Added figure 9A & B to page 17. Re-numbered figures.
11/07	3	2007460	Applied New Template, Added preamplifier to table 6.1 and revised Manual, Replaced the word "coil" to "image" on section 1-5 /line 2, Added the line "Note: The wrist strap MUST be worn all times when handling Coils" to sections 3-3, 5-1 and 5-2, Removed 2 channels from the table on section 3-5
08/08	4	2008215	Table 1-1: Revised p/n's listed in table. Table 6-1: Removed Preamplifier. Table 6-2: Added GE p/n for Pad Kit. Section 1-3: Changed Op manual p/n. Section 3-2-3: Revised step 5. Section 3-5: Revised Symptom #2 table. Removed Figure 9A nd 9B.
05/09	5	23717797	Replace U1-150478 with 2423518, U1-150497 with 2423519 and U1-150477 with 2423523 in Shipping List – Table 1-1. Add Row for Patient comfort pad U1-155164 and remove U1-155164 from pad kit row on Table 6-2. Add Croatian Service Disclaimer. Updated Figure 1.
11/09	6		Updated damage in transportation, language policy, shipping contents, and new MCQA procedure for Unified phantom set per the Phantom Elimination Program.

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "**damage in shipment**" written on **all** copies of the freight or express bill **before** delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier **immediately** upon discovery, or in any event, within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this **14** day period.

File a report with

- Call 1–800–548–3366 and use option 8.
- Fill out a report on <http://egems.med.ge.com/edq/home.jsp>
- Contact your local service coordinator for more information on this process.

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TABLE OF CONTENTS

SECTION 1 – INTRODUCTION	12
1-1 Product Identification and Shipping List	12
Note: Coils do not ship with phantoms. Phantoms come in a unified phantom set with the MR System.	13
1-2 Compatibility	13
1-3 Related Documentation	13
1-4 Environmental Requirements	13
1-5 Theory of Operation	14
SECTION 2 – SETUP AND CALIBRATION	17
2-1 Coil Installation	17
2-1-1 Special Install Notes	17
2-1-2 Configuration for 11x Software	17
2-1-3 Configuration for 10x Software	17
2-2 Installation Functional Checks	17
2-3 Periodic Quality Assurance Check.....	18
SECTION 3 – FUNCTIONAL CHECKS	19
3-1 Scanner Verification	19
3-2 Coil Imaging Performance Verification	19
3-2-1 Tools Required.....	19
3-2-2 Explanation of Procedure	19
3-2-3 SNR Setup for 11.0 M3A and Greater Software Versions with Legacy Phantom	20
3-2-4 SNR Setup for 11.0 M3A and Greater Software Versions with Unified Phantom Set	21
3-2-5 MCQA / SNR Test.....	23
3-3 External Cable Check	23
3-4 PIN Diode Check.....	24
3-6 Coil Element Information	26
SECTION 4 – MAINTENANCE	27
4-1 Coil Care	27
4-2 Special Care Requirements.....	27
SECTION 5 –REPLACEMENT	28
5-1 Anterior Cable Replacement	28
5-2 Posterior Cable Replacement	29
SECTION 6 – RENEWAL PARTS.....	31

6-1 Field Replaceable Units.....31

6-2 Other Replaceable Accessories.....31

7-1 SNR Data Sheet.....32

7-2 Schematic.....33

SECTION 1 – INTRODUCTION

1-1 Product Identification and Shipping List

To identify the 8-Channel Cardiac Array, refer to the coil label (see *Figure 1*).


Figure 1: Coil label.

Figure 2 shows a picture of the 8-Channel Cardiac Array.

Figure 2: The 8-Channel Cardiac Array.

SHIPPING LIST – TABLE 1-1

Box #	Part Name	GE Part #	USAI Part #	Qty
1	Coil	N/A	U1-100391	1
2	Patient Comfort Pad	2423518	N/A	1
2	Patient Headrest Pad	2423519	N/A	1
2	Patient Lower Body Pad	2423523	N/A	1
1	Operator Guide Packet	2422147-500	N/A	1
1	Service Manual	2304255-5	N/A	1

Note: Coils do not ship with phantoms. Phantoms come in a unified phantom set with the MR System.

1-2 Compatibility

The 8-Channel Cardiac Array is compatible with the Excite 1.5T system and requires 10.0 M4 or greater software. 11.0 M4 software is recommended.

The automated QA Tool is only available for software 11.0 M3 and greater.

1-3 Related Documentation

Operator's Guide, GE Part Number 2422147.

Signa Excite Service Methods CD, GE Part Number 2333500-3.

1-4 Environmental Requirements

Storage Requirements

Coil should be stored in the scanner room.

Dimensions

19.4" L x 16.5" W x 4.7" H (49.4cm x 42.0cm x 11.9cm)

Weight

12 lbs (5.4 kg)

1-5 Theory of Operation

The block diagram of 8-Channel Cardiac Array is shown in *Figure 3*. The coil is designed to give a high-resolution image for the cardiovascular system. It works at 1.5T field strength and is capable of doing Array Spatial Sensitivity Encoding Technique (ASSET). This design reduces the breath-hold time, which will aid the patient to successfully hold their breath for the designated time while scanning. The coil is composed of two separate pieces, anterior and posterior. The posterior piece is rigid. The anterior piece is flexible. The flexible design of the coil can accommodate different patient habits. Thus the optimum Signal to Noise Ratio (SNR) can be achieved.

The 8-Channel Cardiac Array is a receive-only coil, as indicated in the schematic shown in the Appendix. The coil is actively decoupled from the RF transmit coil during transmit by means of a RF choking circuit. The choking circuit elements are switched on and off by pin diodes or by small-signal diodes. During body coil transmission, the pin diodes are turned on, as shorted circuits, with a forward bias DC source provided by the MR system (active decoupling). The small-signal diodes are turned on by the induced RF voltage in the chokes, coupled from the transmit field (passive decoupling). When either the pin diodes or the small-signal diodes are turned on, the RF choking elements become very high in impedance (typically above 2 kilo-ohm), compared with the other circuit elements (less than 70 ohm). These high impedance elements in the coil segregate the coil circuitry into several isolated electrical segments, preventing any current flow in the coil circuit.


Figure 3: Block diagram of the coil.


Figure 4: Block Diagram of the coil to connector interface.


Figure 5: Block diagram of the circuit layout.

SECTION 2 – SETUP AND CALIBRATION

2-1 Coil Installation

2-1-1 Special Install Notes

None.

2-1-2 Configuration for 11x Software

Use the Configuration File Manager to install the coil configurations listed in Table 2-1-3.

2-1-3 Configuration for 10x Software

The 8-Channel Cardiac Array requires 10.0 M4 or greater software. 11.0 M4 software is recommended.

Updated coil configurations in the system are required for software versions prior to 11.0 M4. Updated coil configurations are recommended for 11.0 M4 and greater software revisions.

To update the coil configuration files:

1. Using the configuration manager, remove all currently installed (if any) configurations for the 8-Channel Cardiac Array.
2. Update the coil configurations to the latest version (Refer to Section 9 of “Host Software Install Release 10.X”, direction OW2SCA8G. as found on the Service Methods CD, GE Part Number 2333500-3).
3. Install the new configurations listed in table 2-1-3.

CONFIGURATION MODES – TABLE 2-1-3

8CARDIAC	8CARDIACANT	8CARDIACPOST
----------	-------------	--------------

2-2 Installation Functional Checks

1. From the Scan Desktop, start new scan by selecting [**New Pt**]; set **Patient ID** to “geservice” and **Patient Weight** to “111” pounds. Click [**Patient Position**] to open protocols window.
2. Plug the coil into the patient table.

3. At the console, verify the Coil has been properly identified by the system: correct name in **Coil** field and correct picture on the screen.
4. Perform Section 3-2 Coil Imaging Performance Verification.

2-3 Periodic Quality Assurance Check

On a periodic basis, such as during planned maintenance, perform the quality assurance checks as outlined below to ensure the coil is operating properly.

1. Check external cable for cracks or cuts.
2. Perform Section 3-2 Coil Imaging Performance Verification and record data values in Data Sheet (See section 7).

SECTION 3 – FUNCTIONAL CHECKS

3-1 Scanner Verification

Perform system level Signal to Noise Check. Refer to Service Methods CD; System Level Procedures; Functional Checks; Signal to Noise Check.

3-2 Coil Imaging Performance Verification

3-2-1 Tools Required

TOOLS REQUIRED USING LEGACY PHANTOM – TABLE 3-2-1-1

Description	GE Part #	USAI Part #	Qty
Phantom Kit	2304255-3	U1-150511	1

TOOLS REQUIRED USING UNIFIED PHANTOM – TABLE 3-2-1-2

Description	GE Part #	Qty
TL Unified Phantom	5343347	1

3-2-2 Explanation of Procedure

The image quality check uses two consecutive acquisitions for the same slice. The signal is measured at the center of the phantom with a ROI of about 60% of the phantom size. The noise is measured on the subtracted image of the two consecutive signal images using the same ROI. SNR is calculated by dividing the Signal mean by the noise standard deviation. All measurements are performed on composite images. Interim images are acquired for reference only.

3-2-3 SNR Setup for 11.0 M3A and Greater Software Versions with Legacy Phantom

1. Remove any other surface coils (if present) from the cradle.
2. Place the coil on the patient table. Position the phantom set on the coil without the pad, lining up the pegs on the phantom to the holes in the posterior section of the coil (see *Figure 6*).

Figure 6: Positioning the phantom on the coil.

3. Using the raised notches on top of the phantom, align the anterior section on the phantom and fasten the Velcro straps.
4. Landmark the coil at the indicated line and advance to scan (see *Figure 7*).


Figure 7: Landmarking the coil.

3-2-4 SNR Setup for 11.0 M3A and Greater Software Versions with Unified Phantom Set


Visual	Description
 Figure 3.2.4-1	1. Place the posterior section of the coil on the table as shown in Figure 3.2.4-1.
 Figure 3.2.4-2  Figure 3.2.4-3	2. Place the TL unified phantoms over the posterior section in such a way that the phantom set is well centered in both S/I and R/L directions as shown in Figure 3.2.4-2 and Figure 3.2.4-3.

Figure 3.2.4-4

3. Place the anterior section on the phantom set such that the coil is well centered in both S/I and R/L directions with respect to the phantom set as shown in the Figure 3.2.4-4

Figure 3.2.4-5

4. Fasten the Velcro straps so that the coil firmly snugs the phantom as shown in Figure 3.2.4-5.
5. Connect the coil to respective port of the system LPCA.


Figure 3.2.4-6

6. Landmark the coil at the cross mark as shown in Figure 3.2.4-6 and advance to scan.

3-2-5 MCQA / SNR Test

1. Move back to the operator workspace and end the previous exam, if required.
2. For 12.x or earlier systems: From the Common Service Desktop, select **Image Quality**, and then **Coil SNR**. Then select **Click here to start this tool**.

For 14.x or later systems: From the Common Service Desktop, select **Image Quality**, and then **MCQA Tool**. Then select **Click here to start this tool**.

3. Select the coil type you want to run, enter the serial number and select **Start Scan**.

Note: Additional documentation can be accessed by selecting **File, Help, [coil type]**.

3-3 External Cable Check

Diode Check from the Connector:

Note: The wrist strap **MUST** be worn all times when handling Coils.

Step 1 – Using a multi-meter set to the diode check function, check the coil diodes in the following manner.

Step 2 – Orient the connector as shown in *Figure 8*.

Step 3 – Follow Table 3-3 for the correct pins and lead locations to verify the function of the diodes.

EXTERNAL CABLE EXPECTED READINGS – TABLE 3-3

PIN Number	Positive Lead	Negative Lead	DMM Reading
1	B12	A12	0.5 +/- 0.2V
2	B14	A14	0.5 +/- 0.2V
3	B16	A16	0.5 +/- 0.2V
4	B18	A18	0.5 +/- 0.2V
5	B20	A20	0.5 +/- 0.2V
6	B22	A22	0.5 +/- 0.2V
7	B24	A24	0.5 +/- 0.2V
8	B26	A26	0.5 +/- 0.2V


Figure 8: 60 pin Bendix connector.

Note: This connector is system side. Coil side has rows A and B reversed.

3-4 PIN Diode Check

Same as External Cable Check.

The posterior section of the coil must be plugged into the system connector for the pin diode check.

3-5 Troubleshooting Tips

Symptom #1: The system reports a coil fault during prescan or does not recognize the coil connection to the system when selected in the software.

Probable Cause	Suggested Actions	Resolution
The coil connector has become disconnected from the system interface.	Check to make sure the coil connector is fully engaged.	Engage connector and try the scan again. Make sure the indicator light on the patient table turns green with the proper coil selected on the console.
One or more PIN diodes have failed.	From inside the coil, check the PIN diodes using a DVM. Use diode mode for forward-bias measurements and resistance mode for reverse-bias measurements. Note that there should be no reverse-bias leakage. Bad PIN diodes are usually shorted or exhibit reverse-bias leakage current.	Replace the coil.
There is a DC short or open somewhere within the DC path inside the coil.	Check the coil and associated DC paths with a DVM, in both diode check and resistance modes.	If a short or open is observed, check the output cable by itself (see below). If the cable and PIN diodes are okay, then replace the coil.
The output cable has a short or an open.	Disconnect the cable at the coil and check each coaxial line with a DVM (resistance). Try moving and twisting the cable as you watch the meter to find any intermittent connections. Inspect the center pins of the SMB and Bendix connections for wear and proper engagement into the mating connection.	Replace the coil.
There is a problem with the DC bias from the system interface or the coil ID circuitry.	Check to ensure the correct DC bias is being provided to the coil during transmit and receive per GEMS System Service Manual. Try a similar coil to see if the same problem occurs.	Correct MRI system problem.
There is a problem with the coil ID; <i>"The prescribed RF coil and detected RF coil do not match."</i>	Check online instructions by accessing the service browser at Diagnostics; Peripherals; SRI Functional Tests; Coil ID	Start the diagnostic by selecting SRI Coil ID and Run . The diagnostic will display the contents of the coil ID chip whenever the connector is inserted. Follow the troubleshooting instructions if the correct ID is not displayed. When the troubleshooting is complete and the ID is still not displayed, replace the cable.

Symptom #2: The coil does not pass SNR tests or exhibits poor image quality on patient scans.

Probable Cause	Suggested Actions	Resolution
One or more of the coil channels has a high noise level.	Look at the uncombined images to determine which channel has the problem. Compare the noise standard deviation measurements between channels and system performance logs; they should be approximately the same.	Replace the coil.
One or more of the coil channels has low signal.	Look at the uncombined images to determine which channel has the problem. Compare signal mean measurements between channels and system performance logs.	Replace the coil.
There is a problem with the DC bias from the system interface or the coil ID circuitry.	Check to ensure the correct DC bias is being provided to the coils during transmit and receive per GEMS System Service Manual.	Correct MRI system DC bias problem.

3-6 Coil Element Information

Mode	Rcv 1	Rcv 2	Rcv 3	Rcv 4	Rcv 5	Rcv 6	Rcv 7	Rcv 8	PIN 1	PIN 2	PIN 3	PIN 4	PIN 5	PIN 6	PIN 7	PIN 8
8CARDIAC	ON	ON	ON	ON	ON	ON	ON	ON	-5	-5	-5	-5	-5	-5	-5	-5
8CARDIAC POST					ON	ON	ON	ON	+7	+7	+7	+7	-5	-5	-5	-5
8CARDIAC ANT	ON	ON	ON	ON					-5	-5	-5	-5	N/A	N/A	N/A	N/A

SECTION 4 – MAINTENANCE

4-1 Coil Care

Detach the coil connector from the scanner before attempting to clean the coil. Do not touch the connectors with bare fingers. Never press a sharp object against the surface of the connector. Do not reattach the connector after cleaning the coil until the coil has dried completely. Electrical shock may result if the coil is attached to the system during cleaning or when it is wet.

Do not spray or pour cleaning solution directly on coil. Do not submerge coil in solution. The coil contains sensitive electronic components that could be damaged by the solution. The coil cannot be sterilized and should be cleaned only according to the procedure outlined in this section.

The following solutions are recommended for the coil and pad surfaces: (1) a ten percent bleach solution (some discoloration may occur), (2) one ounce commercial dishwashing liquid mixed with one gallon of water or (3) warm water. Apply cleaning solution to a soft cotton cloth and proceed to clean. To prevent soiling of the coil, the user should place a cotton sheet over the coil before positioning the patient. If the coil is soiled, clean the coil as described above.

4-2 Special Care Requirements

Prior to returning a coil for service, use a ten percent bleach solution (as described above) to eliminate risk of exposure to potentially infectious materials.

SECTION 5 –REPLACEMENT

Simple removals that are clearly obvious are not described here.

Unless otherwise noted, the steps for re-assembly are simply the reverse order of the steps described for disassembly.

5-1 Anterior Cable Replacement

Note: The wrist strap **MUST** be worn all times when handling Coils.

Step 1 – Remove coil cover.

Step 2 – Remove the two screws that hold down the cable assembly.

Step 3 – Disconnect SMB plug and the DC Plug on the cable end from the interface board.

Step 4 – Remove the old cable assembly.

Step 5 – Fix a new cable assembly on the anterior section of the coil and reattach the cable lock with two screws. Connect all four SMB plugs and the DC Plug to the corresponding receptacles on the interface board (see *Figure 12*).

Step 6 – Close the coil cover.


Figure 10: Anterior cable assembly.


Figure 11: Anterior DC plug assembly.


Figure 12: Anterior SMB connector assembly.

5-2 Posterior Cable Replacement

Note: The wrist strap **MUST** be worn all times when handling Coils.

Step 1 – Remove coil cover.

Step 2 – Remove the two screws that hold down the cable assembly.

Step 3 – Disconnect SMB plug and the DC Plug on the cable end from the interface board.

Step 4 – Remove the old cable assembly.

Step 5 – Fix a new cable assembly on the posterior section of the coil and reattach the cable lock with two screws. Connect all four SMB connectors and the DC Plug to the corresponding receptacles on the interface board (see *Figure 15*).

Step 6 – Close the coil cover.


Figure 13: Posterior cable assembly.


Figure 14: Posterior DC plug assembly.


Figure 15: Posterior SMB connector assembly.

SECTION 6 – RENEWAL PARTS

6-1 Field Replaceable Units

FIELD REPLACEABLE UNITS – TABLE 6-1

Part Name	GE Part #	USAI Part #
Coil	2304255-2	U1-165148
Cable Assembly	2304255-7	U1-155157
Phantom Kit (Phantom, Rings and Positioners)	2304255-3	U1-155158
TL Unified Phantom	5343347	N/A

6-2 Other Replaceable Accessories

OTHER REPLACEABLE ACCESSORIES LIST – TABLE 6-2

Part Name	GE Part #	USAI Part #
Patient Comfort Pad Kit	2304255-6	N/A
Patient Comfort Pad	N/A	U1-155164
Patient Headrest Pad	N/A	U1-155165
Patient Liner Pad	N/A	U1-155166
Velcro Straps (Set of 2)	N/A	U1-155198


Figure 16: Parts of the coil.

SECTION 7 – APPENDIX

7-1 SNR Data Sheet

Use the table provided below to record the calculated signal to noise ratio (SNR) data obtained from the Functional Checks (see section 3).

Date	Configuration	Slice Location	Corresponding Images in Series	Signal Mean	Noise Std. Dev.	SNR

Comments	

Date	Configuration	Slice Location	Corresponding Images in Series	Signal Mean	Noise Std. Dev.	SNR

Comments	

The diagram illustrates a 4-channel audio amplifier system, divided into two identical sections labeled 'Anterior' and 'Posterior'. Each section contains a pre-amplifier stage (P1), a power amplifier stage (P2), and a speaker output stage. The pre-amplifier stage uses a 2N3904 transistor and a 2N3906 transistor. The power amplifier stage uses a 2N3904 transistor and a 2N3906 transistor. The speaker output stage uses a 2N3904 transistor and a 2N3906 transistor. The system is connected to a 'Connector board' and a 'System Connector'.

GE Healthcare Coils
USA INSTRUMENTS

1515 Danner Drive
Aurora, Ohio 44202
USA

www.gehealthcare.com

Imagination at work