GE Healthcare

Technical Publications

Direction 5158952-100
Revision 6

Functool 4 Installation/Service Manual

Copyright © 2007 by the General Electric Company
All Rights Reserved

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation *Damage in shipment* written on **all** copies of the freight or express bill **before** delivery is accepted or signed for by a GE Healthcare representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier immediately upon discovery or within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14-day period.

To file a report:

- Call 1-800-548-3366, and use option 8.
- Fill out a report at <http://egems.med.ge.com/edq/home.jsp>.

Contact your local service coordinator for more information on this process.

Rev. 09/17/2006

OMMISSIONS & ERRORS

Customers, contact your GE Healthcare Sales or Service representatives. GE personnel, use the GEMS iTrak Process to report all omissions, errors, and defects in this publication.

IMPORTANT PRECAUTIONS

WARNING

(EN)

- This Service Manual is available in English only.
- If a customer's service provider requires a language other than English, it is the customer's responsibility to provide translation services.
- Do not attempt to service the equipment unless this service manual has been consulted and is understood.
- Failure to heed this warning may result in injury to the service provider, operator, or patient from electric shock, or from mechanical or other hazards.

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REVISION HISTORY

Rev	Date	Primary Reasons for Change
1	09/20/05	Initial release.
2	11/21/05	Add FRU Part Number to Section 7
3	06/09/06	Add Installation finish tasks as Section 6 and Update FRU Part Number in Section 8
4	01/15/07	Updated FRU part number in Section 8.
5	06/14/07	Added Korean warning above, changed "cd/" to "cd /" in Section 2-2-2; updated FRU part number in Section 8.
6	02/29/08	Added Omissions & Errors. Updated the following: • Damage in Transportation • Important Precautions Changed the part number in Section 8, FRU List.

TABLE OF CONTENTS

DAMAGE IN TRANSPORTATION	2
OMMISSIONS & ERRORS.....	2
IMPORTANT PRECAUTIONS	3
1- INTRODUCTION	11
1-1 Functool 4 Advanced Application “Clinical” Summary	11
2- PRELIMINARY REQUIREMENTS	11
2-1 License Key Number	11
2-2 Save Configuration	12
2-2-1 Reading Config Media → PC.....	12
2-2-2 Reading Config Media → AW	12
3- INSTALLATION PROCEDURE	12
4- AW AND FUNCTOOL 4 ERRORS.....	13
5- ADVANCED SOFTWARE APPLICATION OPERATION AND DEMO EXAM INFORMATION	13
6- INSTALLATION FINISH TASKS	13
7- SPECIAL CONSIDERATIONS.....	14
8- FRU LIST.....	14

1- INTRODUCTION

The following manual gives a brief overview of the installation process for the Functool 4 advanced application software option. The Functool 4 application is designed to run on the AW4.1 and forward AW platforms only. This guide makes references to all the appropriate standard Advantage Workstation service manuals for more detailed information about installing an advanced application software option on the AW. The current revision of the appropriate AW Service Manual should be used when installing this software. Unique installation processes for Functool 4 and NOT found in the AW service manual can be found in this document in the section, SPECIAL CONSIDERATIONS.

1-1 Functool 4 Advanced Application “Clinical” Summary

Functool 4 is an application that is used to analyze Magnetic Resonance Imaging (MRI) data sets that contain multiple images for each scan location, representing changes in intensity over time (“time series” or “dynamic” data sets), b-value (diffusion imaging), gradient orientation (diffusion tensor imaging) or frequency (spectroscopy).

Functool 4 provides the user with two main tools for analysis:

- Graphs - These show the changes in pixel value over time, b-value or frequency at a given point on the images, either at the cursor position or within a region of interest (ROI) defined by the user.
- Functional Maps - The pixel values from the data set are used to compute a characteristic parameter for every pixel location by means of a function. This can be a basic algorithm or a function derived from such an algorithm. A functional map is then constructed by displaying the value of this parameter for every pixel location.

Functool 4 is supplied in the form of protocols. Each consists of one or more protocol panels that contain the necessary information and the controls required to set up an analysis.

2- PRELIMINARY REQUIREMENTS

2-1 License Key Number

The Functool 4 option requires a license key # based upon the workstation’s Host ID. This key number may be provided by:

- Paperwork shipped with the option that will include one or more 8 to 10-digit license key number(s), created based on the workstation’s/application’s FDO and the workstation’s Host ID/License ID.
- The software package may include a portable media (i.e. floppy disk) with the license key information on it.

The license key information will also be available online at the eLicense website link: <http://elicense.gehealthcare.com/> by searching under the order number or the workstation’s Host ID/License ID.

All of these elements may exist (hard copy license #(s), media, and eLicense data), and at least one of them is necessary to load the software.

2-2 Save Configuration

Create an AW 'CONFIG' media before and after the application software installation for use as needed throughout this and other AW service activities.

2-2-1 Reading Config Media → PC

Configuration information can be read on a PC by viewing the *config.txt* file on the media using text readers like WordPad, Internet Explorer or MS Word. (MS NotePad may not work.)

2-2-2 Reading Config Media → AW

To read a Save Configuration using a floppy disk insert the floppy and type the following in a C shell:

- **su**
- Password: **operator**
- **mount /mnt/floppy**
- **cd /mnt/floppy**
- **more config.txt**
- **cd /**
- **umount /mnt/floppy**

To read a Save Configuration using a CDROM disk insert the CD, and type the following in a C shell:

- **su**
- Password: **operator**
- **mount /mnt/cdrom**
- **cd /mnt/cdrom**
- **more <backup directory name>**
- **cd /**
- **umount /mnt/cdrom**

3- INSTALLATION PROCEDURE

The Standard AW Software Options Installation – Easy Install process applies. For ongoing process revision viability, the details will not be included here, but can be found in the AW Service Manual.

The current revisions for this manual can be viewed online at: <http://aw-ib.euro.med.ge.com>.

1. Open the appropriate Advantage Workstation (AW) Service Manual.
2. Proceed to Job Card titled, *Software Options Installation*.
3. Follow the procedure using the Easy Install method.
4. When complete, perform a *save configuration*.

4- AW AND FUNCTOOL 4 ERRORS

Functool 4 and AW System errors can be examined and investigated in either the usual AW error logs – such as `/var/adm/messages`, and/or `~sdc/logfiles/EasyInstallLog`. Errors specific to Functool 4 can be found in `/export/home/sdc/logfiles/functool.log`.

5- ADVANCED SOFTWARE APPLICATION OPERATION AND DEMO EXAM INFORMATION

To test the functionality of the software, load the DEMO exam CD shipped with the product. The 'Standard' AW Software Demo Exam Installation process applies. For ongoing process revision viability, the details will not be printed here, but can be found in the AW Service Manual.

The current revisions for this manual can be viewed online at: <http://aw-ib.euro.med.ge.com>.

1. Open the appropriate Advantage Workstation (AW) Service Manual.
2. Proceed to Job Card titled, *Advantage Workstation Demo Exam Reinstallation*.
3. Follow the procedure to load the Demo Exams from the CD provided.

To test this application after it is installed, log out and log back into the AW (as suggested by the messages after the installation). Select any data set and click on **Functool** in the *Applications* menu. The *Functool* window will open (on the left side of the screen). This ensures that Functool is installed correctly. Click on the **[Quit]** button to exit Functool.

6- INSTALLATION FINISH TASKS

GIB Update - The Functool4 application must be entered into Global Installed Base (GIB) database after the initial successful installation in order to fulfill product-tracking requirements.

Failure to enter the system and all software applications into GIB will result in GE Healthcare being in violation of FDA tracking compliance, and this system will not be in the affectivity list for future upgrades.

If not familiar with how to get the GIB updated, contact the OLC AW Support. Here is the global GIB link: http://egems.gemedicalsystems.com/gib/gib_entry.jsp

7- SPECIAL CONSIDERATIONS

In addition to Functool, Functool CSI, Functool DTI and Functool FiberTrak will show up at the package installation time as a separate license key entry.

If Functool CSI, Functool DTI, and/or Functool FiberTrak was purchased and they have license keys, they must be selected and the license key numbers entered at the installation prompt for it. Otherwise, leave it unselected.

There can be up to four license keys with this option:

- o Functool **MUST** exist
- o Functool CSI may exist
- o Functool DTI may exist
- o Functool FiberTrak may exist

8- FRU LIST

The replacement part for this option are:

Part #5160880-5 – Functool4 Software and Documentation CD