

Cryogen Fill Line Vacuum Report

1290 Elm st Suite B Ontario, Ca. 91761 USA
Main: 951.734.3737, Fax: 951-735-3373



Start Date	_____	Repair person	_____
SMR barcode #	_____	Kit ID	_____
Part Description	_____		
Any Issues	_____		

Submit by Email

Print Form

Part Repair Status/ Bench Test

Fill out all fields that apply to the complete repair of the part

- ☐ Setup vacuum kit and vacuum down to 10 microns or less

Pull plunger out and watch micron meter for highest reading and note reading

- ☐ Re-vacuum line down to read under 10 microns and then close vacuum valve for 5 mins and watch to see if reading climbs

Note micron reading after the 5 mins and make sure it stops climbing, check for leaks if bad

- ☐ Open vacuum valve to get back under 10 microns, push plunger in and wait 10 mins.

After the 10 mins Pull plunger out and watch micron meter for highest reading and note reading

- ☐ Final vacuum down of line under 10 microns and note reading

Final micron reading

Part repair notes

N
o
t
e
s

Verify Part is good and Complete

- ☐ visual inspection of part making sure everything is intact and complete
- ☐ After done put back in kit, return to parts or send for repair

Status of Part

- ☐ Part is good and going back in kit or on shelf
- ☐ Part is bad and going out for repair
- ☐ Part is un-repairable

N
o
t
e
s

Verification signatures

Repair Person

Supervisor

Date complete

Date complete