

Standard List Of Components And Reference Designators

Note

Area Designators have a space between each level for ease of interpreting the sub-levels. However, on block diagrams where space is limited, it is permissible to omit the spaces between levels.

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION	
BT1	BOOT TERMINAL (Signa 3.X, 4.X)	
CC1	MAIN COMPUTER CABINET (Signa 3.X, 4.X)	
	<u>Designator</u>	<u>Description</u>
	CC1 A1	Expansion Chassis (for MV/4000) or MV/7800 Chassis
	CC1 A2	PDU Mode Control
	CC1 A3	MV/4000 Main Chassis
	CC1 A4	Array Processor
	CC1 A5	5-1/4" Disk Drive Module
	CC1 A5 A1	5-1/4" Disk Drive #1
	CC1 A5 A2	5-1/4" Disk Drive #2
	CC1 A5 A3	5-1/4" Disk Drive #3
	CC1 A5 A4	Power Supply #1
	CC1 A5 A5	Power Supply #2
	CC1 A6	Image Processor Power Supply
	CC1 A7	Image Processor 2
	CC1 A8	IP Fan Module
	CC1 A9	Mag Tape Drive
	CC1 A10	Cabinet Blower
	CC1 A11	Cabinet Blower
	CC1 A12	Temperature Sensor Control
	CC1 A13	Hum Eliminator Panel
	CC1 A14	Nitrogen Monitor
	CC1 A15	Helium Monitor
	CC1 A16	EthernetR Cluster Module
	CC1 A17	Image Processor 1 Power Supply (4.0 forward production)
	CC1 A18	Image Processor 1 (4.0 forward production)
	CC1 A19	Image Processor 1 Fan Module (4.0 forward production)
	CC1 A20	Power Panel
	CC1 A21	Interconnect Panel
	CC1 A22	Temperature Sensor
	CC1 A23	Temperature Sensor
	CC1 A24	Image Processor 1 Hum Eliminator (4.0 forward production)
	CC1 A25	50 Hz Transformer (4.0 forward production)
	CC1 A26	Ethernet Transceiver

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
CC1	COMPUTER CABINET (for Signa 5.X Genesis Operator Console)
	<u>Designator</u> <u>Description</u>
	CC1 A1 15-Slot Card Cage
	CC1 A1 A1 thru A15 (boards as provided)
	CC1 A2 I/O Panel
	CC1 A3 SCSI Tray
	CC1 A3 A1 5-1/4" Disk Drive #1
	CC1 A3 A2 5-1/4" Disk Drive #2
	CC1 A3 A3 DAT Drive (for Release 5.x)
	CC1 A3 A4 Optical Disk Drive #1
	CC1 A3 A5 Optical Disk Drive #2
	CC1 A3 A11 Optical Disk Drive Controller Board
	CC1 A4 Card Cage Fan Tray
	CC1 PS1 System Power Supply
	CC1 CB1 Power Switch Circuit Breaker
CC2	AUXILIARY COMPUTER CABINET (Signa 3.X, 4.X)
	<u>Designator</u> <u>Description</u>
	CC2 A1 Magnetic Tape Drive
	CC2 A2 Image Processor Power Supply
	CC2 A3 Image Processor 3
	CC2 A4 IP Fan Module
	CC2 A5 Disk Drive #1
	CC2 A6 Disk Drive #2
	CC2 A7 Image Processor Power Supply
	CC2 A8 Image Processor 1
	CC2 A9 IP Fan Module
	CC2 A10 Cabinet Blower
	CC2 A11 Cabinet Blower
	CC2 A12 50 Hz Transformer (Mag Tape Power-50 Hz Only)
	CC2 A13 Hum Eliminator Panel
	CC2 A14 Hum Eliminator Panel
	CC2 A20 Power Panel
	CC2 A21 Interconnect Panel

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
CII	COMPUTER CABINET (for Genesis 5.X Independent Console)
	<u>Designator</u> <u>Description</u>
	CII A1 15-Slot Card Cage
	CII A1 A1 thru A15 (boards as provided)
	CII A2 I/O Panel
	CII A3 SCSI Tray
	CII A3 A1 5-1/4" Disk Drive #1
	CII A3 A2 5-1/4" Disk Drive #2
	CII A3 A3 DAT Drive (for IC software release 3.0 and later) 1/4" Cartridge Tape Drive
	CII A3 A4 Optical Disk Drive #1
	CII A3 A5 Optical Disk Drive #2
	CII A3 A10 5-1/4" Disk Drive Controller Board
	CII A3 A11 Optical Disk Drive Controller Board
	CII A4 Card Cage Fan Tray
	CII PS1 System Power Supply
	CC1 PS2 SCSI Power Supply (on early ICs only)
CM1	CARDIAC MONITOR (Signa 3.X, 4.X, 5.X)
	<u>Designator</u> <u>Description</u>
	CM1 A1 Cardiac Lockout Board
DR1	AMPLIFIER CABINET (Distant Remote) (Signa 3.X, 4.X)
	<u>Designator</u> <u>Description</u>
	DR1 A1 Fiber Optic Multiplexer Board
	DR1 A2 Hum Eliminator
	DR1 A3 Hum Eliminator
	DR1 A4 Video Amplifier
	DR1 A5 Video Amplifier
	DR1 A6 Power Supply Connector Board
	DR1 A7 Video Amplifier Power Supply
	DR1 A8 Fiber Optic Power Supply
DR2	CONNECTOR BOX (Distant Remote) (Signa 5.X)
	<u>Designator</u> <u>Description</u>
	DR2 A1 Fiber Optic Multiplexer Board
DR3	HITCH POST

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION												
DR4	TRANSMIT BOX (Genesis Distant Remote) (Signa 5.X)												
	<table> <tr> <th><u>Designator</u></th><th><u>Description</u></th></tr> <tr> <td>DR4 A1</td><td>Mux Board</td></tr> <tr> <td>DR4 A2</td><td>Mux Power Supply</td></tr> <tr> <td>DR4 A3</td><td>Video Transceiver Power Supply</td></tr> <tr> <td>DR4 A4</td><td>Video Fiber Optic Receiver</td></tr> <tr> <td>DR4 A5</td><td>Ethernet Transceiver</td></tr> </table>	<u>Designator</u>	<u>Description</u>	DR4 A1	Mux Board	DR4 A2	Mux Power Supply	DR4 A3	Video Transceiver Power Supply	DR4 A4	Video Fiber Optic Receiver	DR4 A5	Ethernet Transceiver
<u>Designator</u>	<u>Description</u>												
DR4 A1	Mux Board												
DR4 A2	Mux Power Supply												
DR4 A3	Video Transceiver Power Supply												
DR4 A4	Video Fiber Optic Receiver												
DR4 A5	Ethernet Transceiver												
DR5	Receive Box (Genesis Distant Remote) (Signa 5.X)												
	<table> <tr> <th><u>Designator</u></th><th><u>Description</u></th></tr> <tr> <td>DR5 A1</td><td>Mux Board</td></tr> <tr> <td>DR5 A2</td><td>Mux Power Supply</td></tr> <tr> <td>DR5 A3</td><td>Video Transceiver Power Supply</td></tr> <tr> <td>DR5 A4</td><td>Video Fiber Optic Receiver</td></tr> <tr> <td>DR5 A5</td><td>Ethernet Transceiver</td></tr> </table>	<u>Designator</u>	<u>Description</u>	DR5 A1	Mux Board	DR5 A2	Mux Power Supply	DR5 A3	Video Transceiver Power Supply	DR5 A4	Video Fiber Optic Receiver	DR5 A5	Ethernet Transceiver
<u>Designator</u>	<u>Description</u>												
DR5 A1	Mux Board												
DR5 A2	Mux Power Supply												
DR5 A3	Video Transceiver Power Supply												
DR5 A4	Video Fiber Optic Receiver												
DR5 A5	Ethernet Transceiver												
DR6	Hitchpost (Genesis Distant Remote)												
EO1	EMERGENCY OFF BUTTON												
ER1	OXFORD EMERGENCY RUNDOWN UNIT												
HC1	HARD COPIER												
HE1	HEAT EXCHANGER (For Gradient Coil) (Signa 5.X, 8.X)												
HT1	PM												
	<table> <tr> <th><u>Designator</u></th><th><u>Description</u></th></tr> <tr> <td>HT1A1</td><td>PM MUX</td></tr> <tr> <td>HT2A1</td><td>PM DMUX</td></tr> </table>	<u>Designator</u>	<u>Description</u>	HT1A1	PM MUX	HT2A1	PM DMUX						
<u>Designator</u>	<u>Description</u>												
HT1A1	PM MUX												
HT2A1	PM DMUX												

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
IC1	INDEPENDENT CONSOLE (Signa 5.X Genesis)
	<u>Designator</u> <u>Description</u>
	IC1 S1 Console On/Off Switch
	IC1 A1 Plasma Display/Touch Screen/GCP/Autovoice Assembly
	IC1 A1 A1 Plasma Display/Touch Screen Assembly
	IC1 A1 A2 Genesis Console Processor (GCP) Board
	IC1 A2 Image Monitor
	IC1 A3 A/N Keyboard
	IC1 A4 System Keypad
	IC1 A5 IP Keypad
	IC1 A6 Encoder #1 (Leftmost)
	IC1 A7 Encoder #2 (Center)
	IC1 A8 Encoder #3 (Rightmost)
	IC1 A9 Trackball
	IC1 A10 Console I/F Module
	IC1 A10 A1 Video Relay
	IC1 A11 Console Swing Panel
	IC1 A11 A1 Image Monitor Remote Board
	IC1 A11 A2 Security Board
	IC1 A11 A3 Console Power Supply
	IC1 A12 Console Power Distribution Module
	IC1 CB1 Main Circuit Breaker
LC1	LASER CAMERA
LP1	LINE PRINTER

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
MF1	MULTIFORMAT UNIT (Signa 3.X, 4.X, 5.X)
	<u>Designator</u> <u>Description</u>
	MF1 A1 Monitor Assembly Chassis
	MF1 A1 A1 Horizontal Deflection Board
	MF1 A1 A2 CRT Board
	MF1 A1 A3 Vertical Deflection Board
	MF1 A1 A4 Video Amplifier Board
	MF1 A1 A5 Motherboard
	MF1 A1 A5 A1 Interlace Board
	MF1 A1 A6 High Voltage Power Supply
	MF1 A1 A7 Output Board
	MF1 A1 A8 CRT Yoke
	MF1 A1 A9 CRT (cathode ray tube)
	MF1 A1 A10 Monitor Rotation Board
	MF1 A1 A11 Light Sensing Circuit Board
	MF1 A1 A12 Photocell Board
	MF1 A2 Electronics Drawer
	MF1 A2 A1 Processor Board
	MF1 A2 A2 I/O Board
	MF1 A2 A3 Motor And Shutter Driver Board
	MF1 A2 A4 Sensor Interconnect Board
	MF1 A3 Power Supply Drawer
	MF1 A6 Interconnect Panel
	MF1 A6 A1 Video Buffer Board
	MF1 A6 A2 Remote Interconnect Board
	MF1 A7 Autoloader
	MF1 A7 A1 Autoloader Board
	MF1 A8 Control Console
MG1	OXFORD SUPERCONDUCTING MAGNET

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
MG2	MAGNET ENCLOSURE
	<u>Designator</u> <u>Description</u>
	MG2 A1 Status/Voltage Reference Module
	MG2 A1 A1 Status Indicator Board #1
	MG2 A1 A2 Status Indicator Board #2
	MG2 A1 A3 Status Indicator Board #3
	MG2 A1 A4 Voltage Reference Board
	MG2 A1 A5 Status Indicator Interface
	MG2 A2 LRIO Board (Left)
	MG2 A3 LRIO Board (Right)
	MG2 A4 Control Switch Board (Left)
	MG2 A4 Left Operator Control/Display (for Release 5.5 only)
	MG2 A4 Operator Display (MR/i and HFO/i only)
	MG2 A5 Control Switch Board (Right)
	MG2 A5L Operator Control – Left (MR/i and HFO/i only)
	MG2 A5R Operator Control – Right (MR/i and HFO/i only)
	MG2 A5 Right Operator Control/Display (for Release 5.5 only)
	MG2 A6 Limit Switch Board (Up Limit)
	MG2 A7 Limit Switch Board (Dock)
	MG2 A8 Limit Switch Board (Left-home)
	MG2 A9 Limit Switch Board (Right-home)
	MG2 A10 Limit Switch Board (Longitudinal-home)
	MG2 A11 PAC Assembly (Physiological Acquisition Controller) (for Release 5.x)
	Gated Studies Assembly (for Release 3.x & 4.x)
	MG2 A11 A1 PAC Module (Physiological Acquisition Controller) (for Release 5.x)
	Gated Studies Transmitter Module for Release 3.x & 4.x
	MG2 A11 A1 A1 PAC Board (Physiological Acquisition Controller) (for Release 5.x)
	Gated Studies Transmitter Board (for Release 3.x & 4.x)
	MG2 A11 A1 A2 Oscillator Board
	MG2 A11 A1 A4 Gated Studies Transmitter Communications Board
	MG2 A11 A2 Respiratory Bar Graph (for Release 5.x)
	Peripheral Gating Inductor Box (for Release 3.x & 4.x)
	MG2 A11 A3 Peripheral Pulse Bar Graph
	MG2 A12 RF/Gradient Coil Assembly
	MG2 A12 A1 RF/Gradient Coil
	MG2 A12 A2 I Drive Loop Disable Switch
	MG2 A12 A3 Q Drive Loop Disable Switch
	MG2 A12 A4 Dynamic Disable Switch Box 1
	MG2 A12 A5 Dynamic Disable Switch Box 2
	MG2A12 A5 A1 Diode Board (DDSB)
	MG2 A12 A6 Dynamic Disable Switch Box 3
	MG2 A12 A6 A1 Diode Board (DDSB)
	MG2 A12 A7 Dynamic Disable Switch Box 4
	MG2 A12 A8 Dynamic Disable Switch Box 5
	MG2 A12 A9 Dynamic Disable Switch Box 6
	MG2 A12 A10 Dynamic Disable Switch Box 7
	MG2 A12 A11 Dynamic Disable Switch Box 8
	MG2 A12 A12 Front Dynamic Disable Bias Splitter

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
MG2	MAGNET ENCLOSURE (continued)
	<u>Designator</u> <u>Description</u>
MG2 A12 A13	Rear Dynamic Disable Bias Splitter
MG2 A13	Magnet Splitter Box
MG2 A14	Patient Comfort Module
MG2 A15	Bore Light
MG2 A16	Carriage Cover Assembly
MG2 A16 A1	Head Coil Lower Half
MG2 A16 A1 A1	Isolation Network Board
MG2 A16 A1 A2	Head Power Splitter
MG2 A16 A1 A3	RF Coil
MG2 A16 A2	Head Coil Upper Half (for Release 3.x & 4.x)
MG2 A16 A3	Head Preamplifier Assembly
MG2 A16 A3 A1	Preamplifier Board
MG2 A16 A4	Head Preamplifier T/R Switch (3.2/Perf Plus: Quick Disconnect)
MG2 A16 A4 S1	Head-coil-in-position Switch
MG2 A16 A5	Head TR Assembly (for spectroscopy: Spectroscopy/Head Coil Switch)
MG2 A16 A6	Coil Adapter
MG2 A16 A7	Spectroscopy T/R Switch Module (Tardis spectroscopy)
MG2 A16 A7 A1	Spectroscopy T/R Switch Board
MG2 A16 A7 A2	Preamplifier Board
MG2 A16 A7 A3	Spectroscopy Coil Adapter
MG2 A16 A8	Multi-coil Preamp Protection
MG2 A16 A9	Multi-coil Preamp #1
MG2 A16 A10	Multi-coil Preamp #2
MG2 A16 A11	Multi-coil Preamp #3
MG2 A16 A12	Multi-coil Preamp #4
MG2 A16 A13	Multi-coil Adapter
MG2 A16 A14	Surface Coils
MG2 A16 A15 A1	Prostate/Uterine Multicoil Array-top
MG2 A16 A15 A2	Prostate/Uterine Multicoil Array-bottom
MG2 A16 A15 A3	Cable Interface Board (Pelvic Coil)
MG2 A16 A15 A4	0.5T Multicoil Connector Board
MG2 A16 A16 A1	0.5T TMJ Combiner
MG2 A16 A16 A1 A1	0.5T Extremity Coil Matching Board
MG2 A16 A16 A1 A2	0.5T Splitter-Extremity Coil
MG2 A16 A17 A1	T/L Multicoil
MG2 A16 A17 A2	T/L Multicoil Phase Adjustment Board
MG2 A16 A18 A1	Cervical Multicoil Input Board
MG2 A16 A19	Multicoil Preamp #5
MG2 A16 A20	Multicoil Preamp #6
MG2 A16 A21	Multicoil Preamp #7
MG2 A16 A22	Multicoil Preamp #8
MG2 A17	Head Coil Cable Take-up Interface
MG2 A18	Bias Bulkhead
MG2 A19	(Not Used)
MG2 A20	Left Lasers
MG2 A20 A1	Left Axial Laser

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
MG2	MAGNET ENCLOSURE (continued)
	<u>Designator</u> <u>Description</u>
	MG2 A20 A2 Sagittal Laser
	MG2 A20 A3 Left Coronal Laser
	MG2 A21 Right Lasers
	MG2 A21 A1 Right Axial Laser
	MG2 A21 A2 Right Coronal Laser
	MG2 A22 Left Alignment Encoder
	MG2 A23 Left Hydraulic Alignment Cylinder
	MG2 A24 Right Alignment Encoder
	MG2 A25 Right Hydraulic Alignment Cylinder
	MG2 A26 Cradle Latch Cylinder
	MG2 A27 Stop Cylinder
	MG2 A28 Vertical Hydraulic Motor
	MG2 A29 Magnet Interface #2 Or Non-hydraulic Dock Assembly
	MG2 A30 Patient Alignment Lights (Halogen)
	MG2 A31 Dock Limit Module
	MG2 A31 A1 Dock Limit Switch Board
	MG2 A32 (Not Used)
	MG2 A33 Scan Room Interface (SRI) (for Release 5.x)
	Display System Interface Module (DSIM) (for Release 4.x)
	MG2 A34 Peripheral Gating Module
	MG2 A35 Front Cover Microphone
	MG2 TS1 Front Bias Terminal Strip
MG3	REAR PEDESTAL
	<u>Designator</u> <u>Description</u>
	MG3 A1 Magnet Interface #1
	MG3 A1 A1 Limit Switch Board
	MG3 A2 Patient Communication Box
	MG3 A2 A1 Patient Communication Board (for Release 4.x & 3.x)
	Remote Intercom Board (for Release 5.x)
	MG3 A3 Body Quadrature Hybrid Assembly (Body Hybrid Assembly)
	MG3 A3 Body Coil Tr Switch Assembly Panel (linear drive only)
	MG3 A3 A1 Body T/R Switch & Hybrid (Body T/R Switch)
	MG3 A3 A2 Body Preamplifier T/R Switch (Body Splitter)
	MG3 A3 A3 Body Coil T/R Switch (IDD Switch)
	MG3 A3 A5 Body Coil Preamplifier (Quad Li Coil)
	MG3 A4 (Not Used)
	MG3 A5 (Not Used)
	MG3 A6 Longitudinal Encoder
	MG3 A7 Longitudinal Motor Drive
	MG3 A8 (Future)
	MG3 A9 (Future)
	MG3 A10 (Future)

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION																														
MG3	REAR PEDESTAL (continued)																														
	<table> <tr> <th><u>Designator</u></th><th><u>Description</u></th></tr> <tr> <td>MG3 A11</td><td>Body Quadrature Attenuator Assembly</td></tr> <tr> <td>MG3 A12</td><td>Limit Switch Board (Longitudinal-End Of Table)</td></tr> <tr> <td>MG3 A13</td><td>Gradient Cooling Assembly</td></tr> <tr> <td>MG3 A14</td><td>(Future)</td></tr> <tr> <td>MG3 A15</td><td>(Future)</td></tr> <tr> <td>MG3 A16</td><td>(Future)</td></tr> <tr> <td>MG3 A17</td><td>Coil Cable Take-up Interface</td></tr> <tr> <td>MG3 A32</td><td>Bore Venting Assembly</td></tr> <tr> <td>MG3 A33</td><td>I Inductive Drive Disable Switch (1.5T only)</td></tr> <tr> <td>MG3 A34</td><td>Q Inductive Drive Disable Switch (1.5T only)</td></tr> <tr> <td>MG3 TS2</td><td>Rear Bias Terminal Strip</td></tr> <tr> <td>MG3 A35</td><td>Drive Disable Assembly (0.5T only)</td></tr> <tr> <td>MG3 A36</td><td>Multicoil Select Switch</td></tr> <tr> <td>MG3 TS1</td><td>Gradient Terminal Strip</td></tr> </table>	<u>Designator</u>	<u>Description</u>	MG3 A11	Body Quadrature Attenuator Assembly	MG3 A12	Limit Switch Board (Longitudinal-End Of Table)	MG3 A13	Gradient Cooling Assembly	MG3 A14	(Future)	MG3 A15	(Future)	MG3 A16	(Future)	MG3 A17	Coil Cable Take-up Interface	MG3 A32	Bore Venting Assembly	MG3 A33	I Inductive Drive Disable Switch (1.5T only)	MG3 A34	Q Inductive Drive Disable Switch (1.5T only)	MG3 TS2	Rear Bias Terminal Strip	MG3 A35	Drive Disable Assembly (0.5T only)	MG3 A36	Multicoil Select Switch	MG3 TS1	Gradient Terminal Strip
<u>Designator</u>	<u>Description</u>																														
MG3 A11	Body Quadrature Attenuator Assembly																														
MG3 A12	Limit Switch Board (Longitudinal-End Of Table)																														
MG3 A13	Gradient Cooling Assembly																														
MG3 A14	(Future)																														
MG3 A15	(Future)																														
MG3 A16	(Future)																														
MG3 A17	Coil Cable Take-up Interface																														
MG3 A32	Bore Venting Assembly																														
MG3 A33	I Inductive Drive Disable Switch (1.5T only)																														
MG3 A34	Q Inductive Drive Disable Switch (1.5T only)																														
MG3 TS2	Rear Bias Terminal Strip																														
MG3 A35	Drive Disable Assembly (0.5T only)																														
MG3 A36	Multicoil Select Switch																														
MG3 TS1	Gradient Terminal Strip																														
MG4	OXFORD SHIM POWER SUPPLY-Model 2314																														
MG5	OXFORD MAGNET POWER SUPPLY-Model 2137																														
	<table> <tr> <th><u>Designator</u></th><th><u>Description</u></th></tr> <tr> <td>MG5 A3</td><td>Security Board</td></tr> <tr> <td>MG5 A4</td><td>15V, +10V Supply Board</td></tr> <tr> <td>MG5 A5</td><td>Switch Heater Supply</td></tr> <tr> <td>MG5 A6</td><td>Switch Heater Rectifier</td></tr> </table>	<u>Designator</u>	<u>Description</u>	MG5 A3	Security Board	MG5 A4	15V, +10V Supply Board	MG5 A5	Switch Heater Supply	MG5 A6	Switch Heater Rectifier																				
<u>Designator</u>	<u>Description</u>																														
MG5 A3	Security Board																														
MG5 A4	15V, +10V Supply Board																														
MG5 A5	Switch Heater Supply																														
MG5 A6	Switch Heater Rectifier																														
MG6	BLOWER BOX (Signa 5.X, 8.X)																														

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
MM1	MOBILE MAGNET
	<u>Designator</u> <u>Description</u>
	MM1 A1 Service Turret
	MM1 A1 A1 Bursting Disc Port
	MM1 A1 A2 Probe Monitoring Port
	MM1 A1 A3 Liquid Helium Fill Port
	MM1 A1 A4 Current Probe Port
	MM1 A1 A4 A1 Current Probe
	MM1 A1 A4 A2 Current Probe Exhaust Valve
	MM1 A1 A5 Diagnostic Probe Port
	MM1 A1 A6 Helium Exhaust Port
	MM1 A1 A7 Switch And Quench Heater Port
	MM1 A1 A8 Diagnostic Meter Connection
	MM1 A2 Vacuum Turret Cryostat
	MM1 A3 Cryostat
	MM1 A4 Resistive Shim Coil Assembly
MM2	(Not Used)
MM3	(Not Used)
MM4	OXFORD SHIM POWER SUPPLY-Model 2330/2331
MM5	OXFORD MAGNET POWER SUPPLY-Model 2140/2141
MM6	OXFORD DIAGNOSTIC UNIT

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
MR1	RF CABINET/RF Pen/RF Pen II (Signa 3.X, 4.X, 5.X, 8.X)
	<u>Designator</u> <u>Description</u>
	MR1 A1 ETO Amplifier Chassis
	MR1 A2 ETO Manual Driver
	MR1 A2 A1 Manual Driver Control Board
	MR1 A3 RF Power Monitor
	MR1 A3 A1 RF Power Monitor Control Board #1
	MR1 A3 A2 RF Power Monitor Control Board #2
	MR1 A3 A3 RF Sense Analog Board #1
	MR1 A3 A4 RF Sense Analog Board #2
	MR1 A3 A5 T/R Driver
	MR1 A3 A6 Coil Switch Driver (Quad only)
	MR1 A3 A7 Power Monitor I/F Board (spectro only)
	MR1 A3 PS1 +5V/±15V Power Supply For A1 & A3
	MR1 A3 PS2 +5V/±15V Power Supply For A1 & A3
	MR1 A3 PS3 ±15V Power Supply For A2 & A4
	MR1 A3 PS4 +5V Power Supply For A2 & A4
	MR1 A3 PS5 +5V Power Supply For A5
	MR1 A3 PS6 +15V Power Supply For A6
	MR1 A3 PS7 -5V Power Supply For A6 (Quad only)
	MR1 A4 ETO Controller Chassis
	MR1 A4 A1 Control Processor
	MR1 A4 A1 A1 Control Processor
	MR1 A4 A1 A2 Tuning Memory
	MR1 A4 A1 A3 Processor Display
	MR1 A4 A2 SCM Interface
	MR1 A4 A3 A/D Converter
	MR1 A4 A4 Digital Input
	MR1 A4 A5 Digital Output
	MR1 A4 A6 Extender Card
	MR1 A4 A7 RF Attenuator
	MR1 A4 A8 25W Amplifier
	MR1 A4 A9 (Not Used)
	MR1 A4 A10 (Not Used)
	MR1 A5 ETO Power Supply Chassis
	MR1 A5 A1 Power Supply Control Board
	MR1 A5 A2 Rectifier Board
	MR1 A5 A3 Filter Board
	MR1 A6 Sensor Assembly (Envelope Feedback: Envelope Feedback Sensor Assembly)
	MR1 A6 A1 T/R Filter-body (Envelope Feedback: Dual Directional Coupler)
	MR1 A6 A2 T/R Filter-head (Envelope Feedback: Dual Log Amp Board)
	MR1 A6 A3 T/R Filter-surface (Envelope Feedback: Voltage Controlled Attenuator)
	MR1 A6 A4 Power Monitor Samplers-body Elements (Envelope Feedback: Body Sensor Assembly)

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
MR1	RF CABINET/RF Pen/RF Pen II (Signa 3.X, 4.X, 5.X, 8.X) (continued)
	<u>Designator</u> <u>Description</u>
	MR1 A6 A4 A1 Forward Sensor #1
	MR1 A6 A4 A2 Forward Sensor #2
	MR1 A6 A5 Power Monitor Samplers-Head Elements (Envelope Feedback: Head Power Monitor Sensor)
	MR1 A6 A5 A1 FWD (Envelope Feedback: Forward #1)
	MR1 A6 A5 A2 (Not used) (Envelope Feedback: Forward #2)
	MR1 A6 A6 Power Monitor Samplers-Surface Elements (Envelope Feedback: EFB Enable Interface Board)
	MR1 A6 A6 A1 FWD
	MR1 A6 A6 A2 (Not Used)
	MR1 A6 A7 Power Monitor Samplers-Body Elements Manual Attenuator (for Release 4.x)
	MR1 A6 A7 A1 FWD
	MR1 A6 A7 A2 (Not Used)
	MR1 A6 A8 Power Monitor Samplers-Head Elements
	MR1 A6 A8 A1 FWD
	MR1 A6 A8 A2 (Not Used)
	MR1 A6 A9 Power Monitor Samplers-Surface Elements
	MR1 A6 A9 A1 FWD
	MR1 A6 A9 A2 (Not Used)
	MR1 A6 A10 EFB Filter
	MR1 A7 I/F Panel
	MR1 A8 Power Monitor AC Power Strip
	MR1 A9 Dynamic Disable Module
	MR1 A9 A1 Dynamic Disable Driver
	MR1 A9 A2 T/R Driver
	MR1 A9 A3 Dynamic Disable T/R Driver Board
	MR1 A9 A4 Multi-coil T/R Driver Board
	MR1 A9 A5 Spectroscopy T/R Driver Board
	MR1 A9 A6 Serial Peripheral Controller (SPC) Board
	MR1 A9 PS1 -12V Power Supply
	MR1 A9 PS2 +15V Power Supply
	MR1 A9 PS3 -15V Power Supply
	MR1 A9 PS4 +5V Power Supply 400V Power Supply (for Release 4.5 and later Dynamic
	MR1 A10 Erbttec 20kW RF Amplifier
	MR1 A11 Power Module
	MR1 A12 High Speed Blanker Board (0.5T)
	MR1 A13 RF Coupler (1.5T only)
	MR1 A14 Magnet Enclosure Power Supply (MEPS)
	MR1 A14 A1 Transformer
	MR1 A14 A2 High Current Power Supply
	MR1 A14 A3 Low Current Power Supply
	MR1 A14 A4 Longitudinal Drive Amplifier
	MR1 A14 A5 Magnet Enclosure Interface Panel Board
	MR1 A14 K1 Longitudinal Drive Amp Output Relay
	MR1 A14 K2 Bore Vent Relay
	MR1 A15 RF System Controller (RFSC)
	MR1 A15 A1 Analog Processor Board

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION																								
MR1	RF CABINET/RF Pen/RF Pen II (Signa 3.X, 4.X, 5.X, 8.X) (continued)																								
	<table> <tr> <th><u>Designator</u></th><th><u>Description</u></th></tr> <tr> <td>MR1 A15 A2</td><td>Communication Manager/Power Monitor Board</td></tr> <tr> <td>MR1 A15 A3</td><td>RFSC Power Supply Board</td></tr> <tr> <td>MR1 A15 A4</td><td>PIN Switch Driver Board</td></tr> <tr> <td>MR1 A15 A5</td><td>RFSC Interface Panel Board</td></tr> <tr> <td>MR1 A16</td><td>MNS ENI/Analogic RF Amplifier</td></tr> <tr> <td>MR1 A17</td><td>Spectro Module Assembly</td></tr> <tr> <td>MR1 A17 A3</td><td>Relay Assembly</td></tr> <tr> <td>MR1 A17 A2</td><td>3 Band Low-Pass Filter</td></tr> <tr> <td>MR1 A17 A1</td><td>Spectro Directional Coupler</td></tr> <tr> <td>MR1 A18</td><td>Adjustable 10dB Attenuator</td></tr> <tr> <td>MR1 A20</td><td>System Support Module (SSM)</td></tr> </table>	<u>Designator</u>	<u>Description</u>	MR1 A15 A2	Communication Manager/Power Monitor Board	MR1 A15 A3	RFSC Power Supply Board	MR1 A15 A4	PIN Switch Driver Board	MR1 A15 A5	RFSC Interface Panel Board	MR1 A16	MNS ENI/Analogic RF Amplifier	MR1 A17	Spectro Module Assembly	MR1 A17 A3	Relay Assembly	MR1 A17 A2	3 Band Low-Pass Filter	MR1 A17 A1	Spectro Directional Coupler	MR1 A18	Adjustable 10dB Attenuator	MR1 A20	System Support Module (SSM)
<u>Designator</u>	<u>Description</u>																								
MR1 A15 A2	Communication Manager/Power Monitor Board																								
MR1 A15 A3	RFSC Power Supply Board																								
MR1 A15 A4	PIN Switch Driver Board																								
MR1 A15 A5	RFSC Interface Panel Board																								
MR1 A16	MNS ENI/Analogic RF Amplifier																								
MR1 A17	Spectro Module Assembly																								
MR1 A17 A3	Relay Assembly																								
MR1 A17 A2	3 Band Low-Pass Filter																								
MR1 A17 A1	Spectro Directional Coupler																								
MR1 A18	Adjustable 10dB Attenuator																								
MR1 A20	System Support Module (SSM)																								
MR1	RF/PDU (SRFD)																								
	<table> <tr> <th><u>Designator</u></th><th><u>Description</u></th></tr> <tr> <td>MR1 A7</td><td>System Support Module (SSM)</td></tr> <tr> <td>MR1 A21</td><td>Analogic AN8102 RF Amplifier</td></tr> <tr> <td>MR1 A22</td><td>Analogic AN8102 RF Amplifier</td></tr> <tr> <td>MR1 A23</td><td>RF Interface (RFI)</td></tr> <tr> <td>MR1 PD1</td><td>Power Distribution Module</td></tr> </table>	<u>Designator</u>	<u>Description</u>	MR1 A7	System Support Module (SSM)	MR1 A21	Analogic AN8102 RF Amplifier	MR1 A22	Analogic AN8102 RF Amplifier	MR1 A23	RF Interface (RFI)	MR1 PD1	Power Distribution Module												
<u>Designator</u>	<u>Description</u>																								
MR1 A7	System Support Module (SSM)																								
MR1 A21	Analogic AN8102 RF Amplifier																								
MR1 A22	Analogic AN8102 RF Amplifier																								
MR1 A23	RF Interface (RFI)																								
MR1 PD1	Power Distribution Module																								
GRF	GRADIENTS/RF/PDU (GRFD)																								
	<table> <tr> <th><u>Designator</u></th><th><u>Description</u></th></tr> <tr> <td>GRF A7</td><td>System Support Module (SSM)</td></tr> <tr> <td>GRF A21</td><td>Analogic AN8102 RF Amplifier</td></tr> <tr> <td>GRF A23</td><td>RF Interface (RFI)</td></tr> <tr> <td>GRF SCA</td><td>SCA Module</td></tr> <tr> <td>GRF GRAD AMP1</td><td>Analogic AN8280 Gradient Amplifier</td></tr> <tr> <td>GRF GRAD AMP2</td><td>Analogic AN8280 Gradient Amplifier</td></tr> <tr> <td>GRF GRAD AMP3</td><td>Analogic AN8280 Gradient Amplifier</td></tr> <tr> <td>GRF PD1</td><td>Power Distribution Module</td></tr> </table>	<u>Designator</u>	<u>Description</u>	GRF A7	System Support Module (SSM)	GRF A21	Analogic AN8102 RF Amplifier	GRF A23	RF Interface (RFI)	GRF SCA	SCA Module	GRF GRAD AMP1	Analogic AN8280 Gradient Amplifier	GRF GRAD AMP2	Analogic AN8280 Gradient Amplifier	GRF GRAD AMP3	Analogic AN8280 Gradient Amplifier	GRF PD1	Power Distribution Module						
<u>Designator</u>	<u>Description</u>																								
GRF A7	System Support Module (SSM)																								
GRF A21	Analogic AN8102 RF Amplifier																								
GRF A23	RF Interface (RFI)																								
GRF SCA	SCA Module																								
GRF GRAD AMP1	Analogic AN8280 Gradient Amplifier																								
GRF GRAD AMP2	Analogic AN8280 Gradient Amplifier																								
GRF GRAD AMP3	Analogic AN8280 Gradient Amplifier																								
GRF PD1	Power Distribution Module																								

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
MR2	SYSTEM CABINET (3.x/4.x/5.x)
	<u>Designator</u> <u>Description</u>
	MR2 A1 PTS #1 (Transmit)
	MR2 A2 PTS #2 (Receive)
	MR2 A3 Nitrogen Monitor
	MR2 A4 Ethernet Cluster Module/(Thin-net for 4.x)
	MR2 A5 RF Transceiver
	MR2 A5 A1 RF Output
	MR2 A5 A2 Dual Modulator
	MR2 A5 A3 Dual DAC
	MR2 A5 A4 Transceiver I/O
	MR2 A5 A5 Quad Phase Detector (QPD)
	MR2 A5 A6 Base Band Processor
	MR2 A5 A7 A/D And Digital I/O
	MR2 A5 A8 Power Supply Bin
	MR2 A5 A9 RF Patch Panel
	MR2 A5 A13 Gated Studies Receiver
	MR2 A6 Helium Monitor
	MR2 A7 PCM or (CPS) New/changed slot locations/functions for CPS shown in ().
	MR2 A7 A3 Digital Gating Board
	MR2 A7 A4 PSG-1 Frequency Encoding
	MR2 A7 A5 PSG-1 Phase Encoding
	MR2 A7 A6 PSG-1 Slice Selection
	MR2 A7 A7 Oblique Plane Imaging
	MR2 A7 A8 PSG-1 I Envelope (PSG-1 Rho)
	MR2 A7 A9 PSG-1 Q Envelope (PSG-1 Theta)
	MR2 A7 A10 PSG-3
	MR2 A7 A11 (A10) PSG-2
	MR2 A7 (A11) SCM I/F Board
	MR2 A7 (A12) SCM CPU
	MR2 A7 (A12 A1) SBX Gating
	MR2 A7 (A12 A2) ISBX 586 Ethernet
	MR2 A7 A14 Ethernet Controller
	MR2 A7 A15 PCM CPU
	MR2 A7 A16 Splitter # 1
	MR2 A7 A17 Splitter # 2
	MR2 A7 A18 Splitter # 3
	MR2 A7 A19 Splitter # 4
	MR2 A8 SCM
	MR2 A8 A1 Vertical Logic
	MR2 A8 A2 Longitudinal Logic, Translator, DC Motor Controller
	MR2 A8 A3 (Empty)
	MR2 A8 A4 Alignment Logic
	MR2 A8 A5 (Empty)
	MR2 A8 A6 Status I/F
	MR2 A8 A7 Display I/F

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
MR2	SYSTEM CABINET (3.x/4.x/5.x) (continued)
	<u>Designator</u> <u>Description</u>
	MR2 A8 A8 Switch I/F
	MR2 A8 A9 (Empty)
	MR2 A8 A10 Multi-bus Extension
	MR2 A8 A11 Digital I/O
	MR2 A8 A12 512kb RAM
	MR2 A8 A13 (Empty)
	MR2 A8 A14 Ethernet Controller
	MR2 A8 A15 CPU And Diagnostic Display Board
	MR2 A9 Blower
	MR2 A10 Blower
	MR2 A11 I/F Panel
	MR2 A12 Power Panel
	MR2 A13 Gated Studies Receiver (GSR)
	MR2 A14 TPS Power Supply
	MR2 A15 TPS Chassis [New/changed slot locations for Release 5.x shown in ()]
	MR2 A15 A1 CPU
	MR2 A15 A2 Array Processor (AP)
	MR2 A15 A3 VME Interface
	MR2 A15 A4 (Empty)
	MR2 A15 A5 HIB (Memory #0)
	MR2 A15 A6 Memory #4 (Memory #1)
	MR2 A15 A7 Memory #3 (Memory #2)
	MR2 A15 A8 Memory #2 (Memory #3)
	MR2 A15 A9 Memory #1 (DAB)
	MR2 A15 A10 Memory #0 (Receiver #3)
	MR2 A15 A11 DAB
	MR2 A15 A12 Receiver #3 (Receiver #2)
	MR2 A15 A14 Receiver #2 (Receiver #1)
	MR2 A15 A16 Receiver #1 (Receiver #0)
	MR2 A15 A16 A1 (Receiver A/D)
	MR2 A15 A16 A2 (NB Receiver)
	MR2 A15 A16 A3 (BB Receiver)
	MR2 A15 A16 A4 $\pm$ 18V PS Decoupling Board
	MR2 A15 A18 Receiver #0 (Exciter)
	MR2 A15 A18 A1 Receiver A/D (Exciter DAC)
	MR2 A15 A18 A2 NB Receiver (NB Exciter Module)
	MR2 A15 A18 A3 BB Receiver (BB Exciter)
	MR2 A15 A18 A4 ($\pm$ 18V PS Decoupling Board)
	MR2 A15 A20 Exciter (IPG Board; Integrated Pulse Generator)
	MR2 A15 A20 A1 Exciter DAC
	MR2 A15 A20 A2 NB Exciter Module
	MR2 A15 A20 A3 BB Exciter
	MR2 A15 A22 PCM I/F Board
	MR2 A15 A23 TYME Board
	MR2 A16 TPS Fan Module

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
MR2	SYSTEM CABINET (3.x/4.x/5.x) (continued)
	<u>Designator</u> <u>Description</u>
	MR2 A17 RF SPitter Module
	MR2 A18 Temperature Sensor Control
	MR2 A19 PDU Remote Control
	MR2 A20 GSR Power Supply
	MR2 A21 Temperature Sensor (Bay 1)
	MR2 A22 Temperature Sensor (Bay 2)
	MR2 A23 A1 Interface Bracket & Harness
	MR2 A23 A2 SCB (Signal Collecting Board)
	MR2 A24 TNS Assembly
	MR2 A24 A1 TNS Module
	MR2 A24 A2 TNS Expansion Module
	MR2 A24 A3 Temperature Sensor Control
	MR2 A24 A4 Incident Tracking & Communication
	MR2 A25 Fast Receiver 0
	MR2 A26 Fast Receiver 1
	MR2 A27 Fast Receiver 2
	MR2 A28 Fast Receiver 3
MR2	SYSTEM CABINET (8.x)
	<u>Designator</u> <u>Description</u>
	MR2 A3 Nitrogen Monitor
	MR2 A6 Helium Monitor
	MR2 A11 I/F Panel
	MR2 A12 Power Panel
	MR2 A24 TNF Assembly
	MR2 A24 A1 TNF Module
	MR2 A24 A2 TNF Expansion Module
	MR2 A24 A3 Temperature Sensor Control
	MR2 A24 A4 Incident Tracking & Communication
	MR2 TNF PS TNF Power Supply
	MR2 A25 Fast Receiver 0
	MR2 A26 Fast Receiver 1
	MR2 A27 Fast Receiver 2
	MR2 A28 Fast Receiver 3
	MR2 ISE Integrated System Electronics (ISE) Chassis
	MR2 CPU ISE TPS CPU Board
	MR2 BIT3 A2 ISE Bit3 Interface Board
	MR2 BIT3 A3 ISE Bit3 Interface Board (MRT)
	MR2 AP ISE AP Board
	MR2 IPG ISE IPG Board
	MR2 DIG I/F ISE Digital Interface Board (Fast Rcvr and/or Spectroscopy Rcvr)
	MR2 SPECTRO EX ISE Spectroscopy Exciter Board
	MR2 CERD ISE CERD Board
	MR2 CERD PS ISE CERD Power Supply
	MR2 NAV/TNF Navigator/TNF module
	MR2 TEMP CNTL Temperature Controller
	MR2 SHIM PS Z2 Shim Power Supply

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
MR2	SYSTEM CABINET (8.x) (continued)
	<u>Designator</u> <u>Description</u>
	MR2 TYME ISE TYME II Board
	MR2 BKPL PS ISE Backplane Power Supply
	MR2 SPECTRO PS ISE Spectroscopy Power Supply
	MR2 FAN/PWR ISE Fan & AC Power Module
	MR2 FO I/O Bit3 Fiber Optic I/O Module
	MR2 FO I/O MRT Bit3 Fiber Optic I/O (MRT) Module
MR3	GRADIENT AMPLIFIER CABINET (Signa 3.X, 4.X, 5.X)
	<u>Designator</u> <u>Description</u>
	MR3 A1 X Amplifier (Master)
	MR3 A1 A1 Grafidy2 Board
	MR3 A1 A2 Gradient Piggyback Board
	MR3 A2 X Amplifier (Slave)
	MR3 A3 Y Amplifier (Master)
	MR3 A3 A1 Grafidy2 Board
	MR3 A3 A2 Gradient Piggyback Board
	MR3 A4 Y Amplifier (Slave)
	MR3 A5 Z Amplifier (Master)
	MR3 A5 A1 Grafidy2 Board
	MR3 A5 A2 Gradient Piggyback Board
	MR3 A6 Z Amplifier (Slave)
	MR3 A7 Gradient I/F Module
	MR3 A7 A1 Gradient I/F Board
	MR3 A7 A2 SPG (Serial to Peripheral Gradient Board)
	MR3 A9 I/F Panel
	MR3 A10 Power Panel
	MR3 A11 Cabinet Blower (Left Side)
	MR3 A12 Cabinet Blower (Right Side)
MR3	GRADIENT DRIVER (Signa 5.X, 8.X)
	<u>Designator</u> <u>Description</u>
	MR3 A1 GAP
	MR3 A2 X MIF
	MR3 A2 A1 X ASM
	MR3 A3 Y MIF
	MR3 A3 A1 Y ASM
	MR3 A4 Z MIF
	MR3 A4 A1 Z ASM
	MR3 A5 X High Power Module
	MR3 A5 A1 8645 Main Board
	MR3 A5 A2 Output Shelf
	MR3 A5 A4 A1 Toroid Power Supply
	MR3 A5 A4 A2 Toroid Power Supply

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
MR3	GRADIENT DRIVER (Signa 5.X, 8.X) (continued)
	<u>Designator</u> <u>Description</u>
	MR3 A5 A4 A3 Toroid Power Supply
	MR3 A5 A5 140 Amp External Shunt
	MR3 A6 Y High Power Module
	MR3 A7 Z High Power Module
	MR3 A8 X Low Power Module
	MR3 A9 Y Low Power Module
	MR3 A7 Z Low Power Module
	MR3 A13 A1 Fan
	MR3 A13 A2 Fan Control Board
	MR7 A1 X GRAM
	MR7 A1 A1 X GASM
	MR7 A1 A2 X Analog Tuning Board
	MR7 A1 A3 X Digital Tuning Board
	MR7 A2 Y GRAM
	MR7 A2 A1 Y GASM
	MR7 A2 A2 Y Analog Tuning Board
	MR7 A2 A3 Y Digital Tuning Board
	MR7 A3 Z GRAM
	MR7 A3 A1 Z GASM
	MR7 A3 A2 Z Analog Tuning Board
	MR7 A3 A3 Z Digital Tuning Board
MR4	GRADIENT RESISTOR BOX (Signa 4.X)
MR5	ENVIRONMENT TEMPERATURE BOX (Mobile only)
MR6	BROADBAND RF CABINET (Signa 4.X, 5.X, 8.X)
	<u>Designator</u> <u>Description</u>
	MR6 A1 ENI RF Amplifier
	MR6 A2 RF I/F Assembly
	MR6 A2 A1 Power Supply
	MR6 A2 A2 Three-band Low Pass Filter
	MR6 A2 A3 Relay Assembly
	MR6 A2 A4 RF Interface Board
	MR6 A2 A5 Sensor Driver Board
	MR6 A2 A6 Power Monitor Sensor
	MR6 A2 A7 RF Input Board
	MR6 A2 A8 Interconnect Panel
	MR6 A3 Power I/F Panel

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION																																																		
MNS	8.X MNS SPECTROSCOPY CABINET																																																		
	<table> <tr> <th><u>Designator</u></th><th><u>Description</u></th></tr> <tr> <td>MNS AMP</td><td>Analogic RF Amplifier</td></tr> <tr> <td>MNS I/F BOX</td><td>RF I/F Assembly</td></tr> <tr> <td>MNS RELAYS</td><td>Relay Assembly</td></tr> <tr> <td>MNS FILTER</td><td>Three-band Low Pass Filter</td></tr> <tr> <td>MNS DIR CPLR</td><td>Spectro Directional Coupler</td></tr> <tr> <td>MNS ATTEN</td><td>Adjustable 10dB Attenuator</td></tr> </table>	<u>Designator</u>	<u>Description</u>	MNS AMP	Analogic RF Amplifier	MNS I/F BOX	RF I/F Assembly	MNS RELAYS	Relay Assembly	MNS FILTER	Three-band Low Pass Filter	MNS DIR CPLR	Spectro Directional Coupler	MNS ATTEN	Adjustable 10dB Attenuator																																				
<u>Designator</u>	<u>Description</u>																																																		
MNS AMP	Analogic RF Amplifier																																																		
MNS I/F BOX	RF I/F Assembly																																																		
MNS RELAYS	Relay Assembly																																																		
MNS FILTER	Three-band Low Pass Filter																																																		
MNS DIR CPLR	Spectro Directional Coupler																																																		
MNS ATTEN	Adjustable 10dB Attenuator																																																		
MS1	GE SUPERCONDUCTING MAGNET																																																		
	<table> <tr> <th><u>Designator</u></th><th><u>Description</u></th></tr> <tr> <td>MS1 A1</td><td>Cryostat</td></tr> <tr> <td>MS1 A1 A1</td><td>Table End Flange</td></tr> <tr> <td>MS1 A1 A2</td><td>Vacuum Pump Out Port</td></tr> <tr> <td>MS1 A1 A3</td><td>Vacuum Monitoring Port</td></tr> <tr> <td>MS1 A1 A4</td><td>Helium Fill Valve (V1)</td></tr> <tr> <td>MS1 A1 A5</td><td>Helium Vent Valve (V2)</td></tr> <tr> <td>MS1 A1 A6</td><td>Shim Lead Vent Valve (V3)</td></tr> <tr> <td>MS1 A1 A7</td><td>Instrument Lead Vent Valve (V4)</td></tr> <tr> <td>MS1 A2</td><td>Shield Cooler Cold Head</td></tr> <tr> <td>MS1 A3</td><td>Magnet Electrical Harness</td></tr> <tr> <td>MS1 A3 A1</td><td>Terminal Box</td></tr> <tr> <td>MS1 A3 A2</td><td>Harness Lead</td></tr> <tr> <td>MS1 A3 A3</td><td>Service Turret</td></tr> <tr> <td>MS1 A4</td><td>Helium Venting</td></tr> <tr> <td>MS1 A4 A1</td><td>Burst Disc</td></tr> <tr> <td>MS1 A4 A2</td><td>Vent Adaptor</td></tr> <tr> <td>MS1 A4 A3</td><td>Vent Plumbing</td></tr> <tr> <td>MS1 A4 A4</td><td>Exhaust Plumbing</td></tr> <tr> <td>MS1 A5</td><td>Service Equipment/Instrumentation</td></tr> <tr> <td>MS1 A5 A1</td><td>Pressure Gauge</td></tr> <tr> <td>MS1 A5 A2</td><td>Flow Meters</td></tr> <tr> <td>MS1 A5 A3</td><td>Service Platform</td></tr> <tr> <td>MS1 A5 A4</td><td>Lifting Shackles</td></tr> <tr> <td>MS1 A6</td><td>Heater/RTD Distribution I/F (HFO)</td></tr> </table>	<u>Designator</u>	<u>Description</u>	MS1 A1	Cryostat	MS1 A1 A1	Table End Flange	MS1 A1 A2	Vacuum Pump Out Port	MS1 A1 A3	Vacuum Monitoring Port	MS1 A1 A4	Helium Fill Valve (V1)	MS1 A1 A5	Helium Vent Valve (V2)	MS1 A1 A6	Shim Lead Vent Valve (V3)	MS1 A1 A7	Instrument Lead Vent Valve (V4)	MS1 A2	Shield Cooler Cold Head	MS1 A3	Magnet Electrical Harness	MS1 A3 A1	Terminal Box	MS1 A3 A2	Harness Lead	MS1 A3 A3	Service Turret	MS1 A4	Helium Venting	MS1 A4 A1	Burst Disc	MS1 A4 A2	Vent Adaptor	MS1 A4 A3	Vent Plumbing	MS1 A4 A4	Exhaust Plumbing	MS1 A5	Service Equipment/Instrumentation	MS1 A5 A1	Pressure Gauge	MS1 A5 A2	Flow Meters	MS1 A5 A3	Service Platform	MS1 A5 A4	Lifting Shackles	MS1 A6	Heater/RTD Distribution I/F (HFO)
<u>Designator</u>	<u>Description</u>																																																		
MS1 A1	Cryostat																																																		
MS1 A1 A1	Table End Flange																																																		
MS1 A1 A2	Vacuum Pump Out Port																																																		
MS1 A1 A3	Vacuum Monitoring Port																																																		
MS1 A1 A4	Helium Fill Valve (V1)																																																		
MS1 A1 A5	Helium Vent Valve (V2)																																																		
MS1 A1 A6	Shim Lead Vent Valve (V3)																																																		
MS1 A1 A7	Instrument Lead Vent Valve (V4)																																																		
MS1 A2	Shield Cooler Cold Head																																																		
MS1 A3	Magnet Electrical Harness																																																		
MS1 A3 A1	Terminal Box																																																		
MS1 A3 A2	Harness Lead																																																		
MS1 A3 A3	Service Turret																																																		
MS1 A4	Helium Venting																																																		
MS1 A4 A1	Burst Disc																																																		
MS1 A4 A2	Vent Adaptor																																																		
MS1 A4 A3	Vent Plumbing																																																		
MS1 A4 A4	Exhaust Plumbing																																																		
MS1 A5	Service Equipment/Instrumentation																																																		
MS1 A5 A1	Pressure Gauge																																																		
MS1 A5 A2	Flow Meters																																																		
MS1 A5 A3	Service Platform																																																		
MS1 A5 A4	Lifting Shackles																																																		
MS1 A6	Heater/RTD Distribution I/F (HFO)																																																		
MS2	GE SUPERCONDUCTING SHIM POWER SUPPLY-PHASE II																																																		
MS3	GE RESISTIVE SHIM POWER SUPPLY-PHASE II																																																		
MS4 A1	GE EMERGENCY RUNDOWN UNIT																																																		
	<table> <tr> <th><u>Designator</u></th><th><u>Description</u></th></tr> <tr> <td>MS4 A2</td><td>Magnet Rundown Unit Cable</td></tr> </table>	<u>Designator</u>	<u>Description</u>	MS4 A2	Magnet Rundown Unit Cable																																														
<u>Designator</u>	<u>Description</u>																																																		
MS4 A2	Magnet Rundown Unit Cable																																																		
MS5	COLDHEAD COMPRESSOR (MR/i)																																																		
	<table> <tr> <th><u>Designator</u></th><th><u>Description</u></th></tr> <tr> <td>MS5 A2</td><td>He Supply Line</td></tr> <tr> <td>MS5 A3</td><td>He Return Line</td></tr> <tr> <td>MS5 A4</td><td>Coldhead Cable</td></tr> </table>	<u>Designator</u>	<u>Description</u>	MS5 A2	He Supply Line	MS5 A3	He Return Line	MS5 A4	Coldhead Cable																																										
<u>Designator</u>	<u>Description</u>																																																		
MS5 A2	He Supply Line																																																		
MS5 A3	He Return Line																																																		
MS5 A4	Coldhead Cable																																																		

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
MS5	COLDHEAD COMPRESSOR (MR/i) (continued)
	<u>Designator</u> <u>Description</u>
	MS5 A5 Water Flow
	MS5 BV Compressor Bypass Valve for HFO
MS6 A1	Magnet Power Supply
	MS6 A2 Main Coil Power Leads
	MS6 A3 Input Power Cord
MS7	S/C Shim Power Supply
MS8	COLDHEAD COMPRESSOR (HFO/i single stage coldhead)
MSM1	Magnet Monitor
MSM2	Magnet Monitor Remote Alarm
MSM3	Magnet Monitor Universal Modem
MSM4	Magnet Monitor UPS (0.7T)
MUX	Phone Line Multiplexer
OC1	OPERATOR CONSOLE (for Release 3.x & 4.x)
	<u>Designator</u> <u>Description</u>
	OC1 A1 Display Touch Screen
	OC1 A2 Image Display Area
	OC1 A2 A1 Image Monitor
	OC1 A2 A2 Volume Pots
	OC1 A3 Scan Touch Screen
	OC1 A3 A4 Scan Time Remaining Display
	OC1 A4 A1 SCM/Intercom/IF Board
	OC1 A4 A2 Speaker
	OC1 A4 S1 Mute Switch
	OC1 A5 Operator Controls
	OC1 A5 A1 A/N Keyboard
	OC1 A5 A2 Width Control/Encoder
	OC1 A5 A3 Track ball
	OC1 A5 A4 Level Control/Encoder
	OC1 A5 A5 IP Keypad
	OC1 A5 A6 A/N Keyboard
	OC1 A6 Underside Of Operator Controls
	OC1 A6 S1 Power Switch
	OC1 A7 Pedestal Components
	OC1 A7 A1 Console Power Supply
	OC1 A7 A2 Line Filter
	OC1 A7 A3 Display Touch Screen Power Supply
	OC1 A7 A4 Scan Touch Screen Power Supply
	OC1 A7 A5 Console I/O

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
OC1	OPERATOR CONSOLE (for Release 3.x & 4.x) (continued)
	<u>Designator</u> <u>Description</u>
	OC1 A8 Operator Console Interconnect Panel
	OC1 A9 Video Port (Rear Panel)
	OC1 A7 A6 Console Control
OC1	OPERATOR CONSOLE (Signa 5.X Genesis)
	<u>Designator</u> <u>Description</u>
	OC1 S1 Console On/Off Switch
	OC1 S2 Emergency Stop Switch
	OC1 SP1 Patient Intercom Speaker
	OC1 MK1 Patient Intercom Microphone
	OC1 R1 Patient Volume Control
	OC1 R2 Console Volume Control
	OC1 A1 Touch Screen/GCP/Autovoice Assembly
	OC1 A1 A1 Touch Screen Assembly
	OC1 A1 A2 Genesis Console Processor (GCP) Board
	OC1 A1 A3 Autovoice Board
	OC1 A2 Image Monitor
	OC1 A3 A/N Keyboard
	OC1 A4 System Keypad
	OC1 A5 IP Keypad
	OC1 A6 Encoder #1 (Leftmost)
	OC1 A7 Encoder #2 (Center)
	OC1 A8 Encoder #3 (Rightmost)
	OC1 A9 Trackball
	OC1 A10 Console I/F Module
	OC1 A10 A1 Video Relay
	OC1 A11 Console Swing Panel
	OC1 A11 A1 Image Monitor Remote Board
	OC1 A11 A2 Security Board
	OC1 A11 A3 Console Power Supply
	OC1 A12 Console Power Module
	OC1 CB1 Main Circuit Breaker
OM1	OXYGEN MONITOR
OM2	OXYGEN MONITOR (2nd unit in mobile)
OM3	REMOTE OXYGEN SENSOR MODULE
OW1	OPERATOR WORKSPACE (Signa 8.X Indigo)
	<u>Designator</u> <u>Description</u>
	OW1 A1 Operator Workspace Table
	OW1 A1 A1 Power Module
	OW1 A1 A2 Workspace Interface Board
	OW1 A1 A3 Intercom Board
	OW1 A1 A4 Workspace Interface Module
	OW1 A2 Operator Workspace Cabinet
	OW1 A2 A1 IBM Compatible Personal Computer

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
OW1	OPERATOR WORKSPACE (Signa 8.X Indigo) (continued)
	<u>Designator</u> <u>Description</u>
	OW1 A2 A3 SGI Host Computer (Base or Performance Model)
	OW1 A2 A2 Specialix Expansion Adapter RS232
	OW1 A2 A3 A1 Mardi Gras Power Supply
	OW1 A2 A3 A2 Internal Peripheral Power Supply
	OW1 A2 A3 A3 CD-ROM Drive
	OW1 A2 A3 A4 Hard Drive
	OW1 A2 A3 A5 Floptical Drive
	OW1 A2 A3 A6 Mid Plane Extender
	OW1 A2 A3 A7 SGI IP22 Mother Board
	OW1 A2 A3 A8 32 MB SIMM (8)
	OW1 A2 A3 A9 64 MB SIMM (4)
	OW1 A2 A3 A10 SCSI Interface Board (Host)
	OW1 A2 A3 A11 MTA/D Mdir Xacvr Assy Sub D
	OW1 A2 A3 A12 Mardi Gras Graphics (3 Board Assembly)
	OW1 A2 A3 A13 Mardi Gras Video Card
	OW1 A2 A3 A14 Max Optic 1.3 GB MOD
	OW1 A2 A3 A15 Bit3 Board
	OW1 A2 A3 A16 EISA Board (Specialix Board)
	OW1 A2 A3 A17 Full House I/O (AUI I/F)
	OW1 A2 A3 A18 SCSI Interface Board
	OW1 A2 A4 Model 800 Ethernet HUB
	OW1 A2 A5 Power Connection Module
	OW1 A2 A6 Bit 3 Fiber Optic Interface
	OW1 A2 A11 Security Key Connection
	OW1 A3 LCD Display
	OW1 A4 Host Monitor (B/W)
	OW1 A5 DASM
	OW1 A6 Host Monitor (Color)
	OW1 A7 Keyboard Assembly
	OW1 A8 Three Button Mouse
	OW1 A9 Remote Microphone
	OW1 A10 Multi Media Speakers
	OW1 A2 A12 Laptop Ethernet Connection
	OW1 A13 Modem
	OW1 A14 Legacy Media Enclosure
OW1	OPERATOR WORKSPACE (Signa 8.X LX2 Octane)
	<u>Designator</u> <u>Description</u>
	OW1 A1 Operator Workspace Table
	OW1 A1 A1 Power Module
	OW1 A1 A2 Workspace Interface Board
	OW1 A1 A3 Intercom Board
	OW1 A1 A4 Workspace Interface Module
	OW1 A2 Duplex Assembly
	OW1 A2 A1 LX PC in Duplex
	OW1 A2 A2 Specialix Expansion Adapter RS232
	OW1 A2 A3 SGI Host Computer (Indigo)

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
OW1	OPERATOR WORKSPACE (Signa 8.X LX2 Octane) (continued)
	<u>Designator</u> <u>Description</u>
	OW1 A2 A3 A1 Mardi Gras Power Supply
	OW1 A2 A3 A2 Internal Peripheral Power Supply
	OW1 A2 A3 A3 CD-ROM Drive
	OW1 A2 A3 A4 Hard Drive
	OW1 A2 A3 A5 Floptical Drive
	OW1 A2 A3 A6 Mid Plane Extender
	OW1 A2 A3 A7 SGI IP22 Mother Board
	OW1 A2 A3 A8 32 MB SIMM (8)
	OW1 A2 A3 A9 64 MB SIMM (4)
	OW1 A2 A3 A10 SCSI Interface Board (Host)
	OW1 A2 A3 A11 MTA/D Mdir Xacvr Assy Sub D
	OW1 A2 A3 A12 Mardi Gras Graphics (3 Board Assembly)
	OW1 A2 A3 A13 Mardi Gras Video Card
	OW1 A2 A3 A14 Max Optic 1.3 GB MOD
	OW1 A2 A3 A15 Bit3 Board
	OW1 A2 A3 A16 EISA Board (Specialix Board)
	OW1 A2 A3 A17 Full House I/O (AUI I/F)
	OW1 A2 A3 A18 SCSI Interface Board
	OW1 A2 A4 LAN HUB
	OW1 A2 A5 AC Power Module (Octane)
	OW1 A2 A6 Bit 3 Fiber Optic Interface
	OW1 A2 A7 SCSI Expansion Box (Octane)
	OW1 A2 A7 A1 Max Optic Drive # 1
	OW1 A2 A7 A2 Host CD-ROM
	OW1 A2 A7 A3 Host Floppy
	OW1 A2 A7 A4 Max Optic Drive # 2 (Optional)
	OW1 A2 A7 A5 Pioneer MOD
	OW1 A2 A7 A6 DAT Tape Drive
	OW1 A2 A11 Security Key
	OW1 A2 A12 Laptop Ethernet Connection
	OW1 A3 LCD Display
	OW1 A4 Host Monitor (B/W)
	OW1 A5 DASM
	OW1 A6 Host Monitor (Color)
	OW1 A7 Keyboard Assembly
	OW1 A8 Three Button Mouse
	OW1 A9 Remote Microphone
	OW1 A10 Multi Media Speakers (Stereo)
	OW1 A13 Modem
	OW1 A14 Legacy Media Enclosure
	OW1 A15 SGI - Octane
	OW1 A16 SCSI Serial Expansion
PD1	POWER DISTRIBUTION UNIT (Signa 3.X, 4.X, 5.X, 8.X)
	<u>Designator</u> <u>Description</u>
	PD1 A1 Large Door Assembly
	PD1 A2 Connector Panel

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
PD1	POWER DISTRIBUTION UNIT (continued)
	<u>Designator</u> <u>Description</u>
	PD1 A3 Contactor Compartment
	PD1 A3 A1 Contactor Panel
	PD1 A3 A2 Contactor Panel
	PD1 A3 A3 Contactor Panel
	PD1 A4 DC Power Panel
	PD1 A4 A1 24Vdc Power Supply
	PD1 A5 Load Center Compartment
	PD1 A5 A1 Circuit Breakers
	PD1 A5 A2 Temperature Exceeded
	PD1 A5 A3 Phases
	PD1 A5 A4 Service 120 Vac 15A Max
	PD1 A6 Transformer Compartment
	PD1 A6 A1 Transformer
	PD1 A6 A1 A1 Ground
	PD1 A6 A2 Primary Connections-Primary Inputs
	PD1 A6 A3 Secondary Outputs
	PD1 A6 A4 MOV Panel
	PD1 A7 Lower Panel
	PD1 A7 A1 Lower Fans
PIF	“PUCK” INTERFACE (NAVIGATOR FOR HFO)
PP1	PENETRATION PANEL
	<u>Designator</u> <u>Description</u>
	PP1 A1 Penetration Panel-Feedthrough Area
	PP1 A2 X Gradient Filters
	PP1 A3 Y Gradient Filters
	PP1 A4 Z Gradient Filters
	PP1 A5 Shim Interconnect Box
	PP1 A6 Magnet Ramp Port
	PP1 A7 Gradient Filter Box
	PP1 A8 Quad II Inductive Drive Filter
	PP1 A9 Dynamic Disable Drive Filter
	PP1 A10 Fiber Optic Interface
	PP1 A11 Dynamic Disable And Inductive Drive Filter Module
	PP1 A12 Spectroscopy Filter
	PP1 A13 20 dB Gain Block
	PP1 A14 Multi-coil Receive Select Switch
	PP1 A15 Multi-coil Switching Filter
	PP1 A16 Multicoil Bias Filter (G2)
	PP1 A17 FO Repeater Panel
	PP1 A17 A1 FO Repeater Board

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
PP2	PENETRATION CABINET
	<u>Designator</u> <u>Description</u>
	PP2 A1 TPA Assembly (Tunnel Lights, Power Supplies, Air Compressors)
	PP2 A1 A1 Bore Ventilation Compressor
	PP2 A1 A2 Coil Select Compressor
	PP2 A1 A3 Longitudinal Drive Power Amplifier
	PP2 A1 A4 +12 Vdc Tunnel Light Power Supply
	PP2 A1 A5 +28 Vdc Longitudinal Drive Power Supply
	PP2 A1 A6 Longitudinal Drive Relay
	PP2 A4 Laser Power Supplies Module
	PP2 A5 Hydraulic Manifold I/F
	PP2 A5 A1 Manifold Interface Board #1
	PP2 A5 A2 Manifold Interface Board #2
	PP2 A6 Power Module
	PP2 A7 Hydraulic Manifold
	PP2 A7 A1 Hydraulic Motor
	PP2 A7 A2 Hydraulic Pump/Reservoir
	PP2 A7 A3 Double Solenoid For Left Laser
	PP2 A7 A4 Double Solenoid For Right Laser
	PP2 A7 A5 Double Solenoid For Vertical Hydraulic Motor
	PP2 A7 A6 Single Solenoid For Stop Cylinder
	PP2 A7 A7 Single Solenoid For Cradle Latch Cylinder
	PP2 A11 Longitudinal Drive/Bore Light Power Supply
	PP2 A11 PA1 Longitudinal Drive Power Amplifier
	PP2 A11 PS1 +12 Vdc Tunnel Light Power Supply
	PP2 A11 PS2 +28 Vdc Longitudinal Drive Power Supply
	PP2 A11 K3 Longitudinal Drive Relay
	PP2 A12 Patient Alignment Power Supply
	PP2 A13 Air Compressor Assembly
	PP2 A13 A1 Bore Ventilation Compressor
	PP2 A13 A2 Coil Select Compressor
	PP2 A14 Dock Power Supply
	PP2 A15 Power Distribution and Interface Panel
	PP2 A16 Grad Cool/Bore Vent Power Supply
	PP2 A17 Power Supply Module
	PP2 A17 A1 Manifold Interface Board
	PP2 A17 A2 Long Drive Power Amp Board
	PP2 A17 K1 Longitudinal Drive Relay
	PP2 A17 PS1 Bore Light Power Supply
	PP2 A17 PS2 Patient Alignment Lights (PAL) Power Supply
	PP2 A17 PS3 Bore Vent Power Supply
	PP2 A17 PS4 24V SRI Power Supply
	PP2 A17 PS5 12V SRI Power Supply
	PP2 A17 PS6 ± 12V PAC Power Supply
	PP2 A17 PS7 ± 15V PAC Power Supply
PP3	RELIQUIFICATION PENETRATION PANEL

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
PT1	PATIENT TRANSPORT
	<u>Designator</u> <u>Description</u>
	PT1 A1 Pump Assembly And Coupling
	PT1 A2 Reservoir
	PT1 A3 Table Lift Cylinder
	PT1 A4 Latch Cylinder
	PT1 A5 Accumulator Monitor Valve
	PT1 A6 Sequence Valve
	PT1 A7 Accumulator
	PT1 A8 Undocking Valve
	PT1 A9 Table Down Valve
	PT1 A10 Foot Pump
	PT1 A11 Manifold
	PT1 A12 Cradle Latch Release Cylinder
	PT1 A13 Cradle Position Indicating Valve
RC1/RC2	REMOTE CONSOLE (Signa 3.X, 4.X)
	<u>Designator</u> <u>Description</u>
	RC1/RC2 A1 Touch Screen
	RC1/RC2 A2 CRT Assembly
	RC1/RC2 A3 Fans (B1, B2)
	RC1/RC2 A4 S1 Power Switch
	RC1/RC2 A5 Keyboard Housing Assembly
	RC1/RC2 A5 A1 A/N Keyboard
	RC1/RC2 A5 A2 Width Control/Encoder
	RC1/RC2 A5 A3 Level Control/Encoder
	RC1/RC2 A5 A4 Trackball
	RC1/RC2 A5 A5 IP Keypad
	RC1/RC2 A6 Video Panel
	RC1/RC2 A7 Panel Assembly
	RC1/RC2 A7 A1 Console Power Supply
	RC1/RC2 A7 A2 Line Filter
	RC1/RC2 A7 A3 Plasma Power Supply
	RC1/RC2 A7 A4 (Not Used)
	RC1/RC2 A7 A5 Console I/O Board
	RC1/RC2 A7 A6 Video Switch
	RC1/RC2 A7 S1 Primer Cell Switch
	RC1/RC2 A8 Interconnect Panel
RG1	OXFORD REFRIGERATION/LIQUIFICATION
	<u>Designator</u> <u>Description</u>
	RG1 A1 Compressor
	RG1 A2 Electric Control Module (for compressor)
	RG2 A1 Cold Head

SUBSYSTEM REFERENCE DESIGNATOR	SUBSYSTEM DESCRIPTION
RM1	REMOTE VIDEO MONITOR
RS1	COLD BOX
RS2	COMPRESSOR
RS3	MOTOR STARTER
RS4	CONTROL PANEL
RS5	COLD BOX MOTOR CONTROLLER
RC1/RC2	REMOTE CONSOLE
RS6	OIL COOLER
SC1	SPECTROSCOPY 1280 COMPUTER
SCC	SYSTEM COOLING CABINET (COMBINED CHILLER)
SCC2	RF COOLING ASSEMBLY
SD1	SPECTROSCOPY HARD DISK DRIVE
SK1	SPECTROSCOPY A/N KEYBOARD
SM1	SPECTROSCOPY COLOR MONITOR
SP1	SPECTROSCOPY DIGITAL PLOTTER
ST1	STRIP CHART RECORDER
TD1	REEL TO REEL TAPE DRIVE
WC1	WATER CHILLER

REVISION HISTORY

REV	DATE	AUTHOR	PRIMARY REASONS FOR CHANGE
0	2/8/99	Senski	Initial version Reference Designator list in Word format.
1	Jul 19, 1999	Laughlin-Parker	Updated for nMR
2	May 10, 2000	Laughlin-Parker	Updated for 0.7T SCC and MSM2-4
3	June 2, 2000	Laughlin-Parker	Updated for Navigator components
4	July 14, 2000	Don Thome'	Updated System Cabinet for 0.7T, updated for GRFD, and updated for 8X Broadband hardware. Added RFSC, MEPS and sub-components to MR1 list. Added MNS components to MR1 list. Added SSM to MR1 list.