

TABLE OF CONTENTS

TABLE OF CONTENTS	1
1- REMOVE, INSPECT, & REPLACE AIR FILTERS.....	2
2- APPLICATION OF FILTER COAT ADHESIVE	2
REVISION HISTORY.....	4

Description - The System Cabinet uses forced air cabinet cooling. The cabinet air passes through viscous impingement air filters that must be replaced as necessary to assure proper cabinet cooling.

Note

Filters should be replaced during initial installation of the site before turning system over to customer.

1- REMOVE, INSPECT, & REPLACE AIR FILTERS

The System Cabinet has replaceable air filters. When replacing air filters, be sure that new filters are the same as the old ones, since there are more than one size. Refer to *Renewal Parts* for correct part number.

System Cabinet model 2138300, -2 through -14:

Replace filter(s) as follows:

- a. Turn recessed screw at top center of cover 1/4-turn counterclockwise to release cover latch.
- b. Tilt top of cover forward, and lift up and out to remove cover for access to filter(s).
- c. Lift bottom of filter frame up and out to remove filter from cover.
- d. Inspect filter for dust build-up, and replace as required. If replacement is not necessary proceed to step g.
- e. Coat new filter with Filter Adhesive. Refer to Section 2- Application Of Filter Coat Adhesive.
- f. Orient replacement filter with inlet side toward grill opening on cover, and place top edge into opening under retaining lip.
- g. Slide filter upward and under retaining lip until bottom clears lower retaining lip. Lower filter into slot.
- h. Replace cover on cabinet by inserting lower cover lip behind trim strip at bottom of cabinet, and swinging top of cover into spring loaded "slam-latch".

2- APPLICATION OF FILTER COAT ADHESIVE

Air filter adhesive is low viscosity oil, held in gel-like structure that traps dust particles on surface of aluminum mesh filter. Particles recede, or are absorbed into interior of adhesive. Circulation of low viscosity oil establishes renewed impurity absorbing surface. Effectiveness of filter is a result of continuous "fly-paper" action of absorbing surface layer.

See Table 1, Required Materials.

TABLE 1
REQUIRED MATERIALS

Item	Description	Part Number	Quantity
1	RP® Super Filter Coat air filter adhesive, 1 pint can	46-251301P160	1 can/6 filters

Procedure

1. Note instructions on spray can and coat aluminum mesh with RP® Super Filter Coat. Since cabinet blower filters are not thicker than 1/2", spray only from intake side.
2. Apply several coats until aluminum mesh is uniformly coated.

Aside: Contents of can should adequately coat about 6 filters. Overspray and through-spray can also contribute to coating if air filters are grouped and stacked during spraying.

3. After propellant evaporates, this adhesive will have grease-like consistency and filter is ready for replacement.

REVISION HISTORY

REV	DATE	AUTHOR	PRIMARY REASONS FOR CHANGE
0	June 2, 1998	J. Saperstein	Initial Conversion from Toolbook to Word.