

# **Chapter 5**

## **CONSOLE**

---

Section 1.0 [GRE CONSOLE,16slice GDAS with XW8200, 1IG and Ethernet Switch \(2377708-56\), Page 370](#)

Section 2.0 [GRE CONSOLE,16slice MDAS and GDAS with XW8200 and 1IG, JARRELL \(2377708-556\), Page 386](#)

Section 3.0 [Associated Items, Page 394](#)

Section 4.0 [Options, Page 404](#)

**Section 1.0 — GRE CONSOLE,16slice GDAS with XW8200, 1IG and Ethernet Switch (2377708-56)**

| MARK | NUMBER     | QTY | FRU | DESCRIPTION                                                                                                                                                                                                    |
|------|------------|-----|-----|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1    | 2378862-56 | 1.0 | No  | Assy, GOC4 - 16 SLICE MDAS and GDAS with XW8200, 1IG and Ethernet Switch.<br>{For further detail, see <a href="#">Section 1.1, beginning on page 371</a> , including <a href="#">Table 5-2, on page 375</a> .} |
| 2    | 5172474    | 1.0 | No  | Rating plate of GOC4 for LightSpeed Xtra and RT <sup>16</sup>                                                                                                                                                  |


**Table 5-1 ASSY,GRE CONSOLE,16slice GDAS for GOC4 with XW8200, 1IG and Ethernet Switch (2377708-56 BOM Rev2)**

For information on Associated Items, See [Section 3.0](#).

For information on Options, See [Section 4.0](#).

### 1.1 GOC4 w/xw8200 (2378862-56)

$FRU = NO$



**Figure 5-1 GOC4 (2378862-56, sheet 1, rev 1)**







Figure 5-3 GOC4 (2378862-56, sheet 3, rev 1)



**Figure 5-4 GOC4 (2378862-56, sheet 4, rev 1)**

| MARK | NUMBER       | QTY | FRU | DESCRIPTION                                                                                                                              |
|------|--------------|-----|-----|------------------------------------------------------------------------------------------------------------------------------------------|
| 2    | 5117866-28   | 1.0 | Yes | Host Computer, CT USA VCT XW8200                                                                                                         |
| 3    | 5114772      | 1.0 | Yes | DARC Computer - Dual Disk H16 MDAS, GDAS                                                                                                 |
| 4    | 2362872      | 1.0 | Yes | IMAGE GENERATION COMPUTER - GRE                                                                                                          |
| 5    | 2404320      | 1.0 | Yes | CONSOLE INTERCOM PRESCRIBED-TILT HUB<br>{For further detail, see <a href="#">Section 1.3, beginning on page 380</a> }                    |
| 6    | 5167410      | 1.0 | Yes | Peripheral Tower, dual-bay, MOD-DVD Ram                                                                                                  |
| 7    | 5115334      | 1.0 | Yes | Switch, Ethernet 8-port GBe                                                                                                              |
| 8    | 5132092      | 2.0 | Yes | Three Button Optical Mouse                                                                                                               |
| 9    | 2399753-2    | 1.0 | Yes | TRACKBALL, RoHS - 3 BTTN USB for GRE                                                                                                     |
| 11   | 5128181      | 1.0 | Yes | Cover Assy - Console Front<br>{For further detail, see <a href="#">Section 1.4, beginning on page 382</a> }                              |
| 12   | 2337534      | 1.0 | Yes | SIDE COVERS GOC                                                                                                                          |
| 13   | 2304737-2    | 2.0 | Yes | CABLE, DVI VIDEO CABLE - OCTANE II                                                                                                       |
| 17   | 2324697-2    | 1.0 | Yes | SCSI CABLE, VHD68 TO CENT 50, LINUX                                                                                                      |
| 19   | 2374965-2    | 3.0 | Yes | CABLE - ETHERNET CAT5E X 610MM LG                                                                                                        |
| 25   | 46-136246P14 | 1.0 | Yes | 2~ WIDE, WHT, ADH BACK, VELCRO LOOP 25 YD                                                                                                |
| 26   | 46-136246P21 | 1.0 | Yes | 2~ WIDE X 25 YD WHITE VELCRO HOOKS                                                                                                       |
| 27   | 46-208781P2  | 4.0 | Yes | WRITE-ON LABEL FOR MARKING CABLES WHITE ADHESIVE BACKED 0.75~ X 1.90~.<br>UNMARKED, 2,500 PC FAN-FOLD SHEETS (FULL SIZE = 1.9~ X 3.15~.) |
| 29   | 2117390-2    | 1.0 | Yes | CE Compliance and WEEE Labels                                                                                                            |
| 30   | 2132620      | 1.0 | Yes | EMC LABEL CLASS A GROUP 1                                                                                                                |
| 31   | 46-279769P1  | 1.0 | Yes | IEC TYPE B SYMBOL                                                                                                                        |
| 35   | 2214164      | 1.0 | Yes | 5.25~ 2.3GB 512 BYTES/SECTOR MOD                                                                                                         |
| 36   | 2337414      | 2.0 | Yes | DVD RAM MEDIA, 9.4 GB                                                                                                                    |
| 37   | 2227726      | 1.0 | Yes | MOUSE PAD FOR LIGHTSPEED                                                                                                                 |
| 39   | 46-208758P3  | 6.0 | Yes | 7.31~ X .184~ SELF-LOCKING CABLE TIE 1/16~ TO 1-3/4~ WIRE BUNDLE RANGE. 50LB<br>TENSILE STRENGTH, NYLON MATR~L.                          |
| 42   | 2103585-2    | 4.0 | Yes | SCREW PAN HEAD 4 MM 10 MM SPRING WASHER, ZINC PLATED, STEEL                                                                              |


Table 5-2 Assy, GOC4 - 16 SLICE MDAS and GDAS with XW8200, 1IG and Ethernet Switch (2378862-56 BOM, Rev 8)

| MARK | NUMBER       | QTY  | FRU | DESCRIPTION                                                                                                  |
|------|--------------|------|-----|--------------------------------------------------------------------------------------------------------------|
| 43   | 46-208758P2  | 1.0  | Yes | 5.5~ X .140~ SELF-LOCKING CABLE TIE. 1/16~ TO 1-1/4~ WIRE BUNDLE RANGE. 30LB TENSILE STRENGTH, NYLON MATR~L. |
| 44   | 46-208747P1  | 1.0  | Yes | ADHESIVE MOUNT - CABLE TIE PANDUIT CORP, ABMM-A(STYLE A)                                                     |
| 1    | 5162580      | 1.0  | No  | New GOC4 Chassis<br>{For further detail, see <a href="#">Section 1.2, beginning on page 377</a> }            |
| 14   | 2366442-27   | 2.0  | No  | BRACKET, RACK MOUNT HOST COMPUTER                                                                            |
| 32   | 46-269337P6  | 1.0  | No  | LABEL OLD                                                                                                    |
| 33   | 2240540      | 1.0  | No  | NETSCAPE LICENSE                                                                                             |
| 41   | 46-312358P51 | 14.0 | No  | M5 X 0.8 X 12MM METRIC FREEDRIV ZINC PLATED CROSS RECESSED PAN HEAD MACHINE SCREW                            |
| 45   | 5133331      | 1.0  | No  | ETL Label for OC                                                                                             |

Table 5-2 Assy, GOC4 - 16 SLICE MDAS and GDAS with XW8200, 1IG and Ethernet Switch (2378862-56 BOM, Rev 8) (Continued)

## 1.2 New GOC4 Chassis (5162580)

FRU = NO



**Figure 5-5 New GOC4 Chassis (5162580 Sheet 1, Rev 1)**

| MARK | NUMBER       | QTY  | FRU | DESCRIPTION                                                                                                                                                                |
|------|--------------|------|-----|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 2    | 5160898      | 1.0  | Yes | SubAssy, GRE Console Rear Bulkhead                                                                                                                                         |
| 3    | 2315117      | 1.0  | Yes | CABLE, PS/2 MOUSE CABLE, LINUX                                                                                                                                             |
| 4    | 2315117-2    | 1.0  | Yes | CABLE, PS/2 KEYBOARD CABLE, LINUX                                                                                                                                          |
| 5    | 5123478      | 1.0  | Yes | Cable, USB Bulkhead, USB-Af to USB-Am x 3ft Lg                                                                                                                             |
| 6    | 46-328421P3  | 6.0  | Yes | SCREW PAN HEAD 2.5 MM 5 MM ZINC PLATED, STEEL                                                                                                                              |
| 7    | 2142225      | 4.0  | Yes | MONITOR POWER CORD                                                                                                                                                         |
| 8    | 2142226      | 5.0  | Yes | POWER CORD FOR EMC AREA                                                                                                                                                    |
| 9    | 5124071      | 4.0  | Yes | Fan, Tube Axial - 115VAC, 50-60Hz, 119mm x 119mm x 25mm x 95-115cfm, ball bearing motor - Sinwan S125AP-11-1                                                               |
| 13   | 46-208758P3  | 22.0 | Yes | 7.31~ X .184~ SELF-LOCKING CABLE TIE 1/16~ TO 1-3/4~ WIRE BUNDLE RANGE. 50LB TENSILE STRENGTH, NYLON MATR~L.                                                               |
| 14   | 2103585      | 19.0 | Yes | SCREW PAN HEAD 4 MM 8 MM SPRING WASHER, ZINC PLATED, STEEL                                                                                                                 |
| 15   | 46-328431P2  | 8.0  | Yes | WASHER PLAIN - LARGE 6.4 MM 18 MM ZINC PLATED, STEEL                                                                                                                       |
| 16   | 2109881-3    | 8.0  | Yes | WASHER LOCK - INTERNAL 6.4 MM 11 MM                                                                                                                                        |
| 20   | 46-221417P4  | 2.0  | Yes | FEMALE SCREWLOCKS FOR HD CONNECTORS. 4-40 THREAD, 0.470~ OVERALL LENGTH YELLOW CHROMATE, TWO PER KIT. 0.312~ BOTTOM, 0.094~ PANEL THK. MAX (KIT INCLUDES NUTS AND WASHERS) |
| 21   | 46-208747P1  | 1.0  | Yes | ADHESIVE MOUNT - CABLE TIE PANDUIT CORP, ABMM-A(STYLE A)                                                                                                                   |
| 23   | 46-328431P1  | 10.0 | Yes | WASHER PLAIN - LARGE 4.3 MM 12 MM ZINC PLATED, STEEL                                                                                                                       |
| 28   | 46-328417P5  | 2.0  | Yes | SCREW HEXAGON SOCKET 4 MM 16 MM ZINC PLATED, STEEL                                                                                                                         |
| 29   | 46-272301P11 | 4.0  | Yes | FAN ACCESSORIES FINGER GUARD PAPST LZ30                                                                                                                                    |
| 30   | 5127629      | 1.0  | Yes | Cable Assy - Exhaust Fan Power                                                                                                                                             |
| 31   | 2381089-3    | 16.0 | Yes | SCREW, SS NYLOCK HEX NUT, M4                                                                                                                                               |
| 32   | 46-328430P2  | 16.0 | Yes | WASHER PLAIN - NORMAL 4.3 MM 9 MM ZINC PLATED, STEEL                                                                                                                       |
| 33   | 46-312358P44 | 16.0 | Yes | M4 X 0.7 X 35MM METRIC FREEDRIV ZINC PLATED CROSS RECESSED PAN HEAD MACHINE SCREW                                                                                          |
| 34   | 5128406      | 1.0  | Yes | Assy, ICOM cable harness for GOC4                                                                                                                                          |
| 37   | 46-170686P3  | 0.0  | Yes | LOCTITE 242, 50CC 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER.                                                                                         |

Table 5-3 New GOC4 Chassis (5162580 BOM, Rev 2)

| MARK | NUMBER    | QTY | FRU | DESCRIPTION                         |
|------|-----------|-----|-----|-------------------------------------|
| 1    | 2388619-3 | 1.0 | No  | GOC4 chassis from supplier          |
| 10   | 5162696   | 1.0 | No  | GOC4 Outlet Box Assy                |
| 11   | 2385780   | 1.0 | No  | TABLE TOP ASSEMBLY, MONITOR         |
| 12   | 2385783   | 1.0 | No  | TABLE TOP ASSEMBLY, KEYBOARD        |
| 17   | 2337083-2 | 8.0 | No  | SCREW, M6 X 20 SOC HEAD, 12.9 CLASS |
| 22   | 2150839   | 1.0 | No  | CAUTION RECEPTACLES NOT FOR GEN USE |
| 25   | 5114923-2 | 1.0 | No  | Bracket, ICOM Lower                 |
| 26   | 5114924-2 | 1.0 | No  | Bracket, ICOM Upper                 |
| 27   | 5122710   | 2.0 | No  | M4 Tapped Hole Speed Nut            |
| 38   | 5166489   | 1.0 | No  | GOC4 chassis label                  |

Table 5-3 New GOC4 Chassis (5162580 BOM, Rev 2) (Continued)

## 1.3 Console Intercom Prescribed-Tilt Hub (2404320)

FRU = YES




Figure 5-6 Console Intercom Prescribed-Tilt Hub (2404320 Sheet 1, Rev 12)

| MARK | NUMBER       | QTY  | FRU | DESCRIPTION                                                                                                                                                                 |
|------|--------------|------|-----|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 3    | 2374694      | 2.0  | Yes | 8A TIME DELAY FUSES                                                                                                                                                         |
| 10   | 46-312577P14 | 22.0 | Yes | SCREW HEAD CAP SCREW STEEL ZINC PLATED (FULL THREAD) PROPERTY CLASS 8.8<br>***THIS PART IS REPLACED BY 46-328417P3***                                                       |
| 11   | 46-328430P2  | 16.0 | Yes | WASHER PLAIN - NORMAL 4.3 MM 9 MM ZINC PLATED, STEEL                                                                                                                        |
| 13   | 46-312346P1  | 2.0  | Yes | M3 X 0.5 ZINC PLATED METRIC HEX NUT WITH ATTACHED TOOTH LOCK WASHER                                                                                                         |
| 15   | 46-214941P2  | 6.0  | Yes | FLAT WASHER, NYLON, M4                                                                                                                                                      |
| 19   | 46-221417P1  | 3.0  | Yes | FEMALE SCREWLOCKS, TWO PER KIT FOR HD CONNECTORS 4-40 THREAD 0.5~<br>OVERALL LENGTH, YELLOW CHROMATE 0.312~ BOTTOM, 0.062~ PANEL THK. MAX (KIT<br>INCLUDES NUTS AND WASHERS |
| 1    | 2370567      | 1.0  | No  | POWER ENTRY MODULE 2 POLE FUSING                                                                                                                                            |
| 2    | 2374174      | 1.0  | No  | POWER ENTRY FUSE DRAWER                                                                                                                                                     |
| 4    | 2399750      | 1.0  | No  | GRE INTERCOM-TILT                                                                                                                                                           |
| 5    | 5114920-2    | 1.0  | No  | Chassis Base, ICOM                                                                                                                                                          |
| 6    | 5114921      | 1.0  | No  | ct Intercom Box Chassis Cover                                                                                                                                               |
| 7    | 5115027      | 1.0  | No  | Power supply triple output 90-264V 47-63 HZ AC IN 5V 7A +12V 3A -12V 2.5A OUT                                                                                               |
| 8    | 5115179      | 1.0  | No  | cable harness power entry to power supply                                                                                                                                   |
| 9    | 5115180      | 1.0  | No  | cable harness power supply to pwa                                                                                                                                           |
| 12   | 2276462      | 6.0  | No  | METRIC PLAIN FLAT WASHER 3.2MM ID X 7MM OD                                                                                                                                  |
| 14   | 46-312577P4  | 2.0  | No  | SCREW HEAD CAP SCREW STEEL ZINC PLATED (FULL THREAD) PROPERTY CLASS 8.8                                                                                                     |
| 16   | 2209282-3    | 5.0  | No  | HEX STANDOFF MALE/FEMALE M4 X 25MM                                                                                                                                          |
| 17   | 5128595      | 1.0  | No  | M4 x 25mm male female nylon hex stand off.                                                                                                                                  |
| 18   | 46-328417P73 | 4.0  | No  | SCREW HEXAGON SOCKET 3 MM 6 MM                                                                                                                                              |

Table 5-4 Console Intercom Prescribed-Tilt Hub (2404320 BOM, Rev 15)

1.4 Cover Assy - Console Front (5128181)

FRU = YES




Figure 5-7 Cover Assy - Console Front (5128181 Sheet 1, Rev 1)

| MARK | NUMBER       | QTY  | FRU | DESCRIPTION                                      |
|------|--------------|------|-----|--------------------------------------------------|
| 1    | 2369619      | 1.0  | Yes | COVER, CONSOLE FRONT - GRE                       |
| 2    | 5124020      | 1.0  | Yes | Foam, Sound Absorbing - 54in Wd x 1in Thick      |
| 3    | 46-136246P18 | 3.21 | Yes | 1~ WIDE, WHT, ADH BACK, VELCRO HOOK 75 FT PER RL |
| 4    | 5128180      | 1.0  | Yes | Filter, 25 PPI Quadrafoam UL94 HF-1              |

Table 5-5 Cover Assy - Console Front (5128181 BOM, Rev 1)

1.5 Side Covers GOC (2337534)


FRU = YES

| MARK | NUMBER    | QTY | FRU | DESCRIPTION                      |
|------|-----------|-----|-----|----------------------------------|
| 1    | 2315656   | 1.0 | Yes | RIGHT SIDE COVER, GLOBAL CONSOLE |
| 2    | 2315656-2 | 1.0 | Yes | LEFT SIDE COVER, GLOBAL CONSOLE  |


Table 5-6 Side Covers GOC (2337534 BOM, Rev 0)



*FRU = NO*



**Figure 5-8 GRE CONSOLE, 16slice MDAS and GDAS with XW8200 and 1IG, JARRELL (2377708-556, sheet 1, rev 1)**



**Figure 5-9 GRE CONSOLE,16slice MDAS and GDAS with XW8200 and 1IG, JARRELL (2377708-556, sheet 2, rev 1)**




Figure 5-10 GRE CONSOLE,16slice MDAS and GDAS with XW8200 and 1IG, JARRELL (2377708-556, sheet 3, rev 1)



## Chapter 5 - CONSOLE

| MARK | NUMBER       | QTY | FRU | DESCRIPTION                                                                                                                              |
|------|--------------|-----|-----|------------------------------------------------------------------------------------------------------------------------------------------|
| 2    | 5183547-38   | 1.0 | Yes | GEHC MI and CT Integrated Buy Level PC Workstation. HP8400 with 4GB Memory                                                               |
| 3    | 5147442      | 1.0 | Yes | Computer - Data Acquisition Recon Control, 2nd Generation                                                                                |
| 4    | 5147443-2    | 1.0 | Yes | Jarrell GRE IG with RAC rev.D or Higher                                                                                                  |
| 5    | 2404320-2    | 1.0 | Yes | CONSOLE INTERCOM PRESCRIBED-TILT HUB<br>{For further detail, see <a href="#">Section 2.1, beginning on page 392</a> }                    |
| 6    | 5167410      | 1.0 | Yes | Peripheral Tower, dual-bay, MOD-DVD Ram                                                                                                  |
| 7    | 5115334      | 1.0 | Yes | Switch, Ethernet 8-port GBe                                                                                                              |
| 8    | 5132092      | 2.0 | Yes | Three Button Optical Mouse                                                                                                               |
| 9    | 2399753-2    | 1.0 | Yes | TRACKBALL, RoHS - 3 BTTN USB for GRE                                                                                                     |
| 11   | 5128181      | 1.0 | Yes | Cover Assy - Console Front<br>{For further detail, see <a href="#">Section 1.4, beginning on page 382</a> }                              |
| 12   | 2337534      | 1.0 | Yes | SIDE COVERS GOC                                                                                                                          |
| 13   | 2304737-2    | 2.0 | Yes | CABLE, DVI VIDEO CABLE - OCTANE II                                                                                                       |
| 17   | 2324697-2    | 1.0 | Yes | SCSI CABLE, VHD68 TO CENT 50, LINUX                                                                                                      |
| 19   | 5127222-2    | 3.0 | Yes | CABLE, ETHERNET, 2FT Lg                                                                                                                  |
| 25   | 46-136246P14 | 1.0 | Yes | 2~ WIDE, WHT, ADH BACK, VELCRO LOOP 25 YD                                                                                                |
| 26   | 46-136246P21 | 1.0 | Yes | 2~ WIDE X 25 YD WHITE VELCRO HOOKS                                                                                                       |
| 27   | 46-208781P2  | 4.0 | Yes | WRITE-ON LABEL FOR MARKING CABLES WHITE ADHESIVE BACKED 0.75~ X 1.90~.<br>UNMARKED, 2,500 PC FAN-FOLD SHEETS (FULL SIZE = 1.9~ X 3.15~.) |
| 29   | 2117390-2    | 1.0 | Yes | CE Compliance and WEEE Labels                                                                                                            |
| 30   | 2132620      | 1.0 | Yes | EMC LABEL CLASS A GROUP 1                                                                                                                |
| 35   | 2214164      | 1.0 | Yes | 5.25~ 2.3GB 512 BYTES/SECTOR MOD                                                                                                         |
| 36   | 2337414      | 2.0 | Yes | DVD RAM MEDIA, 9.4 GB                                                                                                                    |
| 37   | 2227726      | 1.0 | Yes | MOUSE PAD FOR LIGHTSPEED                                                                                                                 |
| 39   | 46-208758P3  | 6.0 | Yes | 7.31~ X .184~ SELF-LOCKING CABLE TIE 1/16~ TO 1-3/4~ WIRE BUNDLE RANGE. 50LB<br>TENSILE STRENGTH, NYLON MATR~L.                          |
| 42   | 2103585-2    | 4.0 | Yes | SCREW PAN HEAD 4 MM 10 MM SPRING WASHER, ZINC PLATED, STEEL                                                                              |


Table 5-7 GRE CONSOLE, 16 slice MDAS and GDAS with XW8200 and 1IG, JARRELL(2377708-556 BOM, Rev 5)

| MARK | NUMBER       | QTY  | FRU | DESCRIPTION                                                                                                  |
|------|--------------|------|-----|--------------------------------------------------------------------------------------------------------------|
| 43   | 46-208758P2  | 1.0  | Yes | 5.5~ X .140~ SELF-LOCKING CABLE TIE. 1/16~ TO 1-1/4~ WIRE BUNDLE RANGE. 30LB TENSILE STRENGTH, NYLON MATR~L. |
| 44   | 46-208747P1  | 1.0  | Yes | ADHESIVE MOUNT - CABLE TIE PANDUIT CORP, ABMM-A(STYLE A)                                                     |
| 1    | 5162580      | 1.0  | No  | New GOC4 Chassis<br>{For further detail, see <a href="#">Section 1.2, beginning on page 377</a> }            |
| 14   | 2366442-27   | 2.0  | No  | BRACKET, RACK MOUNT HOST COMPUTER                                                                            |
| 23   | 5199454      | 1.0  | No  | Jarrell GOC4 for RT16 Rating Plate Assembly                                                                  |
| 28   | 5133331      | 1.0  | No  | ETL Label for OC                                                                                             |
| 33   | 2240540      | 1.0  | No  | NETSCAPE LICENSE                                                                                             |
| 41   | 46-312358P51 | 14.0 | No  | M5 X 0.8 X 12MM METRIC FREEDRIV ZINC PLATED CROSS RECESSED PAN HEAD MACHINE SCREW                            |

Table 5-7 GRE CONSOLE, 16 slice MDAS and GDAS with XW8200 and 1IG, JARRELL(2377708-556 BOM, Rev 5)

## 2.1 Console Intercom Prescribed-Tilt Hub (2404320-2)

FRU = YES



**Figure 5-12 Console Intercom Prescribed-Tilt Hub (2404320 Sheet 1, Rev 12)**

| MARK | NUMBER       | QTY  | FRU | DESCRIPTION                                                                                                                                                                 |
|------|--------------|------|-----|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 3    | 2374694      | 2.0  | Yes | 8A TIME DELAY FUSES                                                                                                                                                         |
| 10   | 46-312577P14 | 22.0 | Yes | SCREW HEAD CAP SCREW STEEL ZINC PLATED (FULL THREAD) PROPERTY CLASS 8.8<br>***THIS PART IS REPLACED BY 46-328417P3***                                                       |
| 11   | 46-328430P2  | 16.0 | Yes | WASHER PLAIN - NORMAL 4.3 MM 9 MM ZINC PLATED, STEEL                                                                                                                        |
| 13   | 46-312346P1  | 2.0  | Yes | M3 X 0.5 ZINC PLATED METRIC HEX NUT WITH ATTACHED TOOTH LOCK WASHER                                                                                                         |
| 15   | 46-214941P2  | 6.0  | Yes | FLAT WASHER, NYLON, M4                                                                                                                                                      |
| 19   | 46-221417P1  | 3.0  | Yes | FEMALE SCREWLOCKS, TWO PER KIT FOR HD CONNECTORS 4-40 THREAD 0.5~<br>OVERALL LENGTH, YELLOW CHROMATE 0.312~ BOTTOM, 0.062~ PANEL THK. MAX (KIT<br>INCLUDES NUTS AND WASHERS |
| 1    | 2370567      | 1.0  | No  | POWER ENTRY MODULE 2 POLE FUSING                                                                                                                                            |
| 2    | 2374174      | 1.0  | No  | POWER ENTRY FUSE DRAWER                                                                                                                                                     |
| 4    | 5233006      | 1.0  | No  | GRE INTERCOM PWA                                                                                                                                                            |
| 5    | 5114920-2    | 1.0  | No  | Chassis Base, ICOM                                                                                                                                                          |
| 6    | 5114921      | 1.0  | No  | ct Intercom Box Chassis Cover                                                                                                                                               |
| 7    | 5304856      | 1.0  | No  | Power supply triple output 90-264V 47-63 HZ AC IN 5V 7A +12V 3A -12V 2.5A OUT                                                                                               |
| 8    | 5115179      | 1.0  | No  | cable harness power entry to power supply                                                                                                                                   |
| 9    | 5115180      | 1.0  | No  | cable harness power supply to pwa                                                                                                                                           |
| 12   | 2276462      | 6.0  | No  | METRIC PLAIN FLAT WASHER 3.2MM ID X 7MM OD                                                                                                                                  |
| 14   | 46-312577P4  | 2.0  | No  | SCREW HEAD CAP SCREW STEEL ZINC PLATED (FULL THREAD) PROPERTY CLASS 8.8                                                                                                     |
| 16   | 2209282-3    | 5.0  | No  | HEX STANDOFF MALE/FEMALE M4 X 25MM                                                                                                                                          |
| 17   | 5128595      | 1.0  | No  | M4 x 25mm male female nylon hex stand off.                                                                                                                                  |
| 18   | 46-328417P73 | 4.0  | No  | SCREW HEXAGON SOCKET 3 MM 6 MM                                                                                                                                              |

Table 5-8 Console Intercom Prescribed-Tilt Hub (2404320-2 BOM, Rev 2)

## Section 3.0 — Associated Items

- [LCD Monitor - 19 Inch \(B71292CA/B7858LC\), Page 396](#)
- [Keyboard w/SCIM, Collectors, Page 397](#)
- [Console Slide Service Kit \(2379135\), Page 402](#)
- [Power Switch & Cables \(2379134\), Page 402](#)
- [Modem Kit \(B7700MG\), Page 402](#)

## 3.1 Monitors




Figure 5-13 Monitor Selection Flowchart

**3.1.1 LCD Monitor - 19 Inch (B71292CA/B7858LC)**

| MARK | NUMBER    | QTY | FRU | DESCRIPTION/NAME                |
|------|-----------|-----|-----|---------------------------------|
| 1    | 5128434-5 | 2.0 | No  | NEC 19 inch LCD monitor - white |

Table 5-9 LCD Monitor &amp; Cable (B71292CA BOM, rev 2)

| MARK | NUMBER  | QTY | FRU | DESCRIPTION/NAME                                 |
|------|---------|-----|-----|--------------------------------------------------|
| 1    | 5169069 | 2.0 | Yes | OEM 19 INCHES COLOR MONITOR LCD 1990SXi - SILVER |

Table 5-10 LCD Monitor &amp; Cable (B7858LC BOM, rev 3)

## 3.2 Keyboard w/SCIM, Collectors



Figure 5-14 Keyboard & SCIM Assembly (Rev 11/15/00)

**3.2.1 Keyboard w/SCIM Collector - English (2275758)**

FRU = YES

| MARK | NUMBER  | QTY | FRU | DESCRIPTION/NAME                                                                                                       |
|------|---------|-----|-----|------------------------------------------------------------------------------------------------------------------------|
| 1    | 2275756 | 1.0 | Yes | KEYBOARD, NORTH AMERICAN                                                                                               |
| 3    | 2275752 | 1.0 | Yes | SCIM COLLECTOR - ENGLISH<br>{For further detail, see <a href="#">SCIM Collector - English (2275752) on page 398.</a> } |
| 4    | 2323942 | 1.0 | Yes | CABLE, HD50 TO HD50                                                                                                    |
| 2    | 2275757 | 1.0 | No  | FUNCTION KEY TEMPLATE - E                                                                                              |

Table 5-11 Keyboard w/SCIM Collector - English (2275758 BOM, rev 5)

**3.2.1.1 SCIM Collector - English (2275752)**

FRU = YES

| MARK | NUMBER     | QTY | FRU | DESCRIPTION/NAME          |
|------|------------|-----|-----|---------------------------|
| 1    | 2262158    | 1.0 | No  | SCAN CONTROL INTERCOM MOD |
| 2    | 2262158-11 | 1.0 | No  | OVERLAY, SCIM CT, TILT -  |

Table 5-12 SCIM Collector - English (2275752 BOM, rev 2)

**3.2.2 Chinese (5167383)**

| MARK | NUMBER     | QTY | FRU | DESCRIPTION/NAME                                  |
|------|------------|-----|-----|---------------------------------------------------|
| 1    | 2275758-10 | 1.0 | No  | KEYBOARD and SCIM COLLECTOR - ENGLISH and CHINESE |

Table 5-13 Keyboard w/SCIM Collector - French (5167383 BOM, rev 1)

**3.2.2.1 SCIM Collector - Chinese (2275758-10)**

| MARK | NUMBER    | QTY | FRU | DESCRIPTION/NAME                 |
|------|-----------|-----|-----|----------------------------------|
| 1    | 2275756   | 1.0 | Yes | KEYBOARD, NORTH AMERICAN PC PS/2 |
| 2    | 2275757   | 1.0 | Yes | FUNCTION KEY TEMPLATE - ENGLISH  |
| 3    | 2275757-8 | 1.0 | Yes | FUNCTION KEY TEMPLATE - MANDARIN |

SCIM Collector - Chinese (2275758-10 BOM, rev 1)

| MARK | NUMBER     | QTY | FRU | DESCRIPTION/NAME                     |
|------|------------|-----|-----|--------------------------------------|
| 4    | 2275752-10 | 1.0 | Yes | SCIM COLLECTOR - ENGLISH and CHINESE |
| 5    | 5149890    | 1.0 | No  | Warning Label                        |
| 6    | 5149891    | 1.0 | No  | Warning Label                        |

SCIM Collector - Chinese (2275758-10 BOM, rev 1)

### 3.2.3 French (2275758-2)

| MARK | NUMBER    | QTY | FRU | DESCRIPTION/NAME                                                                                                        |
|------|-----------|-----|-----|-------------------------------------------------------------------------------------------------------------------------|
| 1    | 2275756-2 | 1.0 | Yes | KEYBOARD, EUROPEAN FRENCH                                                                                               |
| 3    | 2275752-2 | 1.0 | Yes | SCIM COLLECTOR - FRENCH<br>{For further detail, see <a href="#">SCIM Collector - Chinese (2275758-10)</a> on page 398.} |
| 4    | 2323942   | 1.0 | Yes | CABLE, HD50 TO HD50                                                                                                     |
| 2    | 2275757-2 | 1.0 | No  | FUNCTION KEY TEMPLATE - F                                                                                               |

Table 5-14 Keyboard w/SCIM Collector - French (2275758-2 BOM, rev 3)

#### 3.2.3.1 SCIM Collector - French (2275752-2)

| MARK | NUMBER     | QTY | FRU | DESCRIPTION/NAME          |
|------|------------|-----|-----|---------------------------|
| 1    | 2262158    | 1.0 | No  | SCAN CONTROL INTERCOM MOD |
| 2    | 2262158-12 | 1.0 | No  | OVERLAY, SCIM CT, TILT -  |

SCIM Collector - French (2275752-2 BOM, rev 1)

### 3.2.4 German (2275758-3)

| MARK | NUMBER    | QTY | FRU | DESCRIPTION/NAME                                                                                                       |
|------|-----------|-----|-----|------------------------------------------------------------------------------------------------------------------------|
| 1    | 2275756-3 | 1.0 | Yes | KEYBOARD, GERMAN PC PS/2                                                                                               |
| 3    | 2275752-3 | 1.0 | Yes | SCIM COLLECTOR - GERMAN<br>{For further detail, see <a href="#">SCIM Collector - German (2275752-3) on page 400.</a> } |
| 4    | 2323942   | 1.0 | Yes | CABLE, HD50 TO HD50                                                                                                    |
| 2    | 2275757-3 | 1.0 | No  | FUNCTION KEY TEMPLATE - GERMAN                                                                                         |

Table 5-15 Keyboard w/SCIM Collector - German (2275758-3 BOM, rev 3)

#### 3.2.4.1 SCIM Collector - German (2275752-3)

| MARK | NUMBER     | QTY | FRU | DESCRIPTION/NAME          |
|------|------------|-----|-----|---------------------------|
| 1    | 2262158    | 1.0 | No  | SCAN CONTROL INTERCOM MOD |
| 2    | 2262158-13 | 1.0 | No  | OVERLAY, SCIM CT, TILT -  |

Table 5-16 SCIM Collector - German (2275752-3 BOM, rev 1)

### 3.2.5 Swedish (2275758-4)

| MARK | NUMBER    | QTY | FRU | DESCRIPTION/NAME                                                                                                       |
|------|-----------|-----|-----|------------------------------------------------------------------------------------------------------------------------|
| 1    | 2275756-4 | 1.0 | Yes | KEYBOARD, SWEDISH PC PS/2                                                                                              |
| 3    | 2275752   | 1.0 | Yes | SCIM COLLECTOR - ENGLISH<br>{For further detail, see <a href="#">SCIM Collector - English (2275752) on page 398.</a> } |
| 4    | 2323942   | 1.0 | Yes | CABLE, HD50 TO HD50                                                                                                    |
| 2    | 2275757   | 1.0 | No  | FUNCTION KEY TEMPLATE - E                                                                                              |

Table 5-17 Keyboard w/SCIM Collector - Swedish (2275758-4 BOM, rev 3)

### 3.2.6 International (2275758-5)

| MARK | NUMBER    | QTY | FRU | DESCRIPTION/NAME                                                                                                                    |
|------|-----------|-----|-----|-------------------------------------------------------------------------------------------------------------------------------------|
| 1    | 2275756   | 1.0 | Yes | KEYBOARD, NORTH AMERICAN                                                                                                            |
| 3    | 2275752-4 | 1.0 | Yes | SCIM COLLECTOR - INTERNATIONAL<br>{For further detail, see <a href="#">SCIM Collector - International (2275752-4)</a> on page 401.} |
| 4    | 2323942   | 1.0 | Yes | CABLE, HD50 TO HD50                                                                                                                 |
| 2    | 2283830   | 1.0 | No  | FUNCTION KEY TEMPLATE COL                                                                                                           |

Table 5-18 Keyboard w/SCIM Collector - International (2275758-5 BOM, rev 3)

#### 3.2.6.1 SCIM Collector - International (2275752-4)

| MARK | NUMBER     | QTY | FRU | DESCRIPTION/NAME                  |
|------|------------|-----|-----|-----------------------------------|
| 1    | 2262158    | 1.0 | No  | SCAN CONTROL INTERCOM MODULE      |
| 2    | 2262158-10 | 1.0 | No  | OVERLAY, SCIM CT, TILT - BLANK    |
| 3    | 2262158-14 | 1.0 | No  | OVERLAY, SCIM CT, TILT - JAPANESE |
| 4    | 2262158-15 | 1.0 | No  | OVERLAY, SCIM CT, TILT - MANDARIN |

Table 5-19 SCIM Collector - International (2275752-4 BOM, rev 2)

### 3.3 Console Slide Service Kit (2379135)

| MARK | NUMBER    | QTY   | FRU | DESCRIPTION/NAME          |
|------|-----------|-------|-----|---------------------------|
| 3    | 2317595   | 1.000 | No  | SLIDE RAIL FOR DRAWER ASS |
| 4    | 2317501   | 1.000 | No  | BRACKET, DRAWER BALL LATC |
| 5    | 2317502   | 1.000 | No  | LATCH, DRAWER ASSEMBLY    |
| 6    | 2366442-7 | 4.000 | No  | BRACKET, RACK MOUNT       |

Table 5-20 Console Slide Service Kit (2379135 BOM, rev 0)

### 3.4 Power Switch & Cables (2379134)

| MARK | NUMBER    | QTY   | FRU | DESCRIPTION/NAME                  |
|------|-----------|-------|-----|-----------------------------------|
| 1    | 2150833   | 1.000 | No  | BREAKER SWITCH                    |
| 2    | 2108900   | 1.000 | No  | CONNECTOR 8 CONTACTS 150 VOLTS    |
| 5    | 2318308   | 1.000 | No  | ETHERNET CABLE, CHASSIS SVC PANEL |
| 7    | 2341434-2 | 1.000 | No  | CABLE, KEYBOARD - OCTANE 2        |
| 8    | 2341434   | 1.000 | No  | CABLE, MOUSE - OCTANE 2           |

Table 5-21 Power Switch &amp; Cables (2379134 BOM, rev 0)

### 3.5 Modem Kit (B7700MG)

| MARK | NUMBER  | QTY   | FRU | DESCRIPTION   |
|------|---------|-------|-----|---------------|
| 1    | 2382276 | 1.000 | No  | GRE MODEM KIT |

Table 5-22 GLOBAL MODEM FOR XTREAM W/ GRE (B7700MG BOM, rev 0)



## Section 4.0 — Options

- [Video Upgrade For Linux Console \(B7868ME\), Page 405](#)
- [LCD & SUSP \(Boom In Room\) LNX \(B7710LN\), Page 407](#)

## 4.1 Video Upgrade For Linux Console (B7868ME)

| MARK | NUMBER  | QTY | FRU | DESCRIPTION/NAME                                                                             |
|------|---------|-----|-----|----------------------------------------------------------------------------------------------|
| 1    | 5115182 | 1.0 | No  | VIDEO UPGRADE FOR LINUX CONSOLE<br>{For further detail, see <a href="#">Section 4.1.1.</a> } |

Table 5-23 Video Upgrade For Linux Console (B7868ME BOM, rev 0)

### 4.1.1 Video Upgrade (5115182)

| MARK | NUMBER       | QTY  | FRU | DESCRIPTION/NAME                                             |
|------|--------------|------|-----|--------------------------------------------------------------|
| 2    | 2403438-3    | 1.0  | Yes | 5 BNC MALE TO HD 15 MALE 75 FEET                             |
| 6    | 5161129      | 1.0  | Yes | LCD Bracket Short Arm                                        |
| 7    | 2212717      | 1.0  | Yes | CABLE COVER FOR LCD SUSPENSION 3.6 METERS LONG               |
| 8    | 5116219      | 1.0  | Yes | Grounding Cable For LCD Console To LCD                       |
| 13   | 2376772      | 1.0  | Yes | 1-IN 4-OUT RGBHV DISTRIBUTION AMP                            |
| 14   | 2403438      | 2.0  | Yes | 5 BNC MALE TO HD 15 MALE 6 FEET                              |
| 16   | 2403438-4    | 1.0  | Yes | HD 15 FEMALE TO HD 15 MALE 15 FEET                           |
| 27   | 2384323      | 10.0 | Yes | 75 OHM TERMINATOR                                            |
| 28   | 2384324      | 1.0  | YES | DVI TO HD15 CONVERTER                                        |
| 30   | 46-328421P17 | 4.0  | Yes | SCREW PAN HEAD 4 MM 20 MM ZINC PLATED, STEEL                 |
| 31   | 46-328425P1  | 2.0  | Yes | NUT HEXAGON 4 MM                                             |
| 32   | 46-328432P2  | 6.0  | Yes | WASHER LOCK - SPRING 4.1 MM 7.6 MM ZINC PLATED, SPRING STEEL |
| 33   | 46-328430P2  | 6.0  | Yes | WASHER PLAIN - NORMAL 4.3 MM 9 MM ZINC PLATED, STEEL         |
| 34   | 46-208758P5  | 10.0 | Yes | SELF LOCKING CABLE TILE                                      |
| 1    | 46-302200P10 | 1.0  | No  | 1.72 X 1.10, MODEL#, SN,                                     |
| 17   | 5125372-100  | 1.0  | No  | NEC 1980SXi LDC Service Manual                               |
| 18   | 2390072-100  | 1.0  | No  | LS LINUX CRT AND LCD REMOTE MONITOR                          |
| 25   | 2398987      | 1.0  | No  | CEILING PLATE ASSEMBLY                                       |
| 26   | 5115171      | 1.0  | No  | Handle LCD Monitor                                           |

Table 5-24 Video Upgrade (5115182 BOM, rev 6)

| MARK | NUMBER  | QTY | FRU | DESCRIPTION/NAME          |
|------|---------|-----|-----|---------------------------|
| 29   | 5115172 | 1.0 | No  | Ground Wire Ceiling Plate |

Table 5-24 Video Upgrade (5115182 BOM, rev 6)

## 4.2 LCD & SUSP (Boom In Room) LNX (B7710LN)

| MARK | NUMBER    | QTY   | FRU | DESCRIPTION/NAME                                                                                          |
|------|-----------|-------|-----|-----------------------------------------------------------------------------------------------------------|
| 1    | 5115174-2 | 1.000 | No  | Boom-in-a-Room for Linux with new suspension<br>{For further detail, see <a href="#">Section 4.2.1</a> .} |

Table 5-25 LCD &amp; SUSP (Boom In Room) LNX (B7710LN BOM, rev 1)

### 4.2.1 Boom-in-a-Room for Linux (5115174-2)

| MARK | NUMBER       | QTY  | FRU | DESCRIPTION                                                                        |
|------|--------------|------|-----|------------------------------------------------------------------------------------|
| 2    | 2403438-3    | 1.0  | Yes | 5 BNC MALE TO HD 15 MALE 75 FEET                                                   |
| 3    | 5160969-3    | 1.0  | Yes | Complete Suspension containing ceiling column 5160969 and suspension arm 5160969-2 |
| 5    | 5169069      | 1.0  | Yes | OEM 19 INCHES COLOR MONITOR LCD 1990SXi - SILVER                                   |
| 6    | 5161129      | 1.0  | Yes | LCD BRACKET - SHORT                                                                |
| 7    | 2212717      | 1.0  | Yes | CABLE COVER FOR LCD SUSPENSION 3.6 METERS LONG                                     |
| 8    | 5116219      | 1.0  | Yes | Grounding Cable for LCD CONSOLE TO LCD                                             |
| 13   | 2376772      | 1.0  | Yes | 1-IN 4-OUT RGBHV DISTRIBUTION AMP                                                  |
| 14   | 2403438      | 2.0  | Yes | 5 BNC MALE TO HD 15 MALE 6 FEET                                                    |
| 16   | 2403438-4    | 1.0  | Yes | HD 15 FEMALE TO HD 15 MALE 15 FEET                                                 |
| 27   | 2384323      | 10.0 | Yes | 75 OHM TERMINATOR                                                                  |
| 28   | 2384324      | 1.0  | Yes | DVI TO HD15 CONVERTER                                                              |
| 1    | 46-302200P10 | 1.0  | No  | 1.72 X 1.10, MODEL#, SN, DATE, DESC.                                               |
| 10   | 2213219      | 1.0  | No  | POWER CABLE FOR LCD-CONSOLE TO LCD                                                 |
| 18   | 2390072-100  | 1.0  | No  | LS LINUX CRT AND LCD REMOTE MONITOR                                                |
| 25   | 2398987      | 1.0  | No  | CEILING PLATE ASSEMBLY                                                             |
| 26   | 5115171      | 1.0  | No  | HANDLE LCD MONITOR                                                                 |
| 29   | 5115172      | 1.0  | No  | Ground Wire Ceiling Plate                                                          |
| 30   | 5161926      | 1.0  | No  | Hardware collecotr for Boom-in-a room                                              |

Table 5-26 Boom-in-a-Room for Linux (5115174-2 BOM, Rev 2)

